CHLORIDE ELECTRODE

ION-SELECTIVE ELECTRODE FOR MEASUREMENT OF CHLORIDE IONS

CHLORIDE ELECTRODE

Cat. No E-708

PACKAGE

1 CHLORIDE ELECTRODE

INTENDED USE

The CHLORIDE ELECTRODE is a membrane electrode designed for in vitro chloride ions measurements in blood, serum, plasma and urine samples. Electrode is intended for use in electrolyte analysers based on ISE systems.

STABILITY AND STORAGE

- Store at room temperature (18-25°C).
- The electrode must be installed before the expiry date indicated on the product label, but can remain in the instrument after the date.

INSTALLATION

1. Check that the electrode is $\frac{3}{4}$ filled with fill solution.
2. Check for air bubbles in the fill chamber near sample path. If there are any, tap the electrode body to dislodge air bubbles.
3. Verify that the O-rings are in place on the electrode body.
4. Carefully install electrode in your system.
5. Allow the electrode to equilibrate for at least 15 minutes.
6. Perform a two-point calibration.
7. Analyse a control material sample.

NOTES

- Product for in vitro diagnostic use only.
- Electrodes should be used by suitably qualified laboratory personnel only in accordance with intended purpose.

WASTE MANAGEMENT

Please refer to local legal requirements.

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PZ CORMAY S.A.