## ALP

	(PL)	HC-ALP	OS-ALP	B50-ALP	
<b>Nazwa zestawu</b>	<b>Nr kat.</b>	1-REAGENT	6 x 76,5 ml	3 x 44 ml	2 x 58,5 ml
Liquick Cor-ALP 30	1-218	2-REAGENT	6 x 19,5 ml	3 x 14 ml	2 x 17,5 ml
Liquick Cor-ALP 60	1-212				
Liquick Cor-ALP 120	3-327				
HC-ALP	4-512				
OS-ALP	9-420				
B50-ALP	5-528				

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### INTENDED USE


ALP reagent is intended to determine quantitatively alkaline phosphatase activity in serum and plasma. It is used to monitoring and as an aid to diagnosis of clinical conditions associated with abnormal ALP activity. ALP reagent is intended to use on automatic analyzers and manually. It is only for *in vitro* diagnosis, for healthcare professional users.

### SUMMARY<sup>1-6</sup>

Determinations of alkaline phosphatase activity in serum and plasma are mainly used as an aid in diagnosis, differentiation and monitoring of hepatobiliary disease and bone diseases. Elevated levels of the enzyme in the blood are observed, e.g. in cholestatic liver diseases, osteoporosis or Paget's disease. Decreased ALP levels are associated with Wilson's disease, Hypophosphatasia, and deficiency of magnesium, zinc and vitamin B12. Measurement of ALP activity is also used as a predictive marker in outcomes (death or liver transplantation) of patients with primary sclerosing cholangitis (PSC) and primary biliary cholangitis (PBC). The determination of the level of alkaline phosphatase is also used to monitoring the effectiveness of the PBC and IgG4-related cholangiopathy treatment.

### METHOD PRINCIPLE<sup>7,8</sup>

Kinetic method recommended by International Federation of Clinical Chemistry (IFCC).

The rate of 4-nitrophenol formation is directly proportional to the ALP activity.

### REAGENTS

Package	Liquick Cor-ALP 30	Liquick Cor-ALP 60	Liquick Cor-ALP 120
1-REAGENT	5 x 24 ml	5 x 48 ml	5 x 96 ml
2-REAGENT	1 x 30 ml	1 x 60 ml	1 x 120 ml
	HC-ALP	OS-ALP	B50-ALP
1-REAGENT	6 x 76.5 ml	3 x 44 ml	2 x 58.5 ml
2-REAGENT	6 x 19.5 ml	3 x 14 ml	2 x 17.5 ml

### CONCENTRATIONS OF THE ACTIVE INGREDIENTS IN THE REAGENT

1-REAGENT	
2-amino-2-methyl-1-propanol (AMP)	425 mmol/l
Mg <sup>2+</sup>	2.5 mmol/l
Zn <sup>2+</sup>	1.25 mmol/l
metal ion chelator	
pH adjuster	
preservative	
2-REAGENT	
p-nitrophenylphosphate	83.9 mmol/l
2-amino-2-methyl-1-propanol (AMP)	50 mmol/l
pH adjuster	
preservative	

### REAGENT STABILITY

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyzer at 2-10°C are stable for 12 weeks (Biolis 30i).

### Working reagent preparation and stability

Assay can be performed with use of separate 1-REAGENT and 2-REAGENT reagents or with use of working reagent. For working reagent preparation mix gently 4 parts of 1-REAGENT with 1 part of 2-REAGENT. Avoid foaming!

Stability of working reagent:	4 weeks at 2-8°C
	5 days at 15-25°C

### WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- The reagents are usable when the absorbance of the working reagent is less than 1.300 (read against distilled water, wavelength  $\lambda=405$  nm, cuvette l=1 cm, at temp. 25°C).
- During the reaction p-nitrophenol is produced. Do not swallow or inhale, avoid contact with skin.
- Slightly yellow colour of working reagent is normal and does not influence the result.
- Do not use reagent beyond expiry date printed on the package.
- Do not mix reagents from different kits or lots.
- Use personal protective equipment to prevent contact with samples, reagents and controls.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT and 2-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

### Ingredients:

1-REAGENT and 2-REAGENT contain 2-methyl-2H-isothiazol-3-one.

### Warning

H317 May cause an allergic skin reaction.  
P280 Wear protective gloves, protective clothing, eye protection or face protection.  
P302 + P352 IF ON SKIN: Wash with plenty of soap and water.

### EQUIPMENT REQUIRED BUT NOT PROVIDED

- automatic analyzer or photometer able to read at 410 nm (405/412 nm);
- thermostat at 37°C;
- general laboratory equipment;

### SPECIMEN<sup>1</sup>

Serum, heparinized plasma free from hemolysis.

Do not use EDTA, citrate and oxalate as anticoagulants because of ALP activity inhibition!

ALP activity remains stable in specimen up to 4 hours at 15-25°C. Freezing of sample causes a decrease of enzyme activity. Frozen specimens should be thawed and kept at room temperature for 18 to 24 hours before measurement to achieve full enzyme reactivation.

Nevertheless it is recommended to perform the assay with freshly collected samples!

Follow tube manufacturers' instructions carefully when using collection tubes.

Human-origin material should be handled as potentially infectious. Standard precautions in normal laboratory work are required.

### ASSAY PROCEDURE


Reagents are ready to use.

Carry-over effect may occur when the reagent is used on automated chemistry analyzers. It is based on interference between specific test combinations and results in either underestimating or inflating patient sample results. To minimize this effect, the preventive actions should be taken, such as ordering an extra washing program or performing tests in a separate order.

### General procedure for automatic determinations:

R1: 200  $\mu$ l

Sample: 7  $\mu$ l

**Main Wavelength:** 405 nm

**2nd wavelength:** 700 nm

**Method:** Kinetic

**Calibration Type:** Linear

**Direction:** Increasing

**R1 Volume:** 200  $\mu$ l

**Sample Volume:** 7  $\mu$ l

**R2 Volume:** 50  $\mu$ l

Application guides for analyzers are available on request.

The application guide should be verified and validated by the user prior to testing patient samples.

### Manual procedure

wavelength	410 nm (405/412 nm)
temperature	37°C
cuvette	1 cm

### Sample Start method

Pipette into the cuvette:

working reagent	1000 $\mu$ l
sample	18 $\mu$ l

Bring up to the temperature of determination. Then add:  
Mix and incubate at adequate temperature. After about 1 min. read the absorbance against air or water. Repeat the reading after exactly 1, 2 and 3 minutes. Calculate the mean absorbance change per minute ( $\Delta A/\text{min.}$ ).

### Reagent Start method

The determination can be also performed with use of separate 1-REAGENT and 2-REAGENT reagents.

Pipette into the cuvette:

1-REAGENT	1000 $\mu$ l
sample	17 $\mu$ l
2-REAGENT	250 $\mu$ l

Mix well; perform measurement as described for Sample Start method.

### Calculation

ALP activity [U/l] =  $\Delta A/\text{min.}$  x calibrator concentration or

ALP activity [U/l] =  $\Delta A/\text{min.}$  x F

Sample Start method	Reagent Start method
F = 3038	F = 3442

### QUALITY CONTROL

Pay attention to preparation of calibrator and control before ALP determination. For internal quality control, it is recommended to use with each batch of samples the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of manual assay the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or LEVEL 2 (Cat. No 5-175; 5-177) or appropriate assigned factor (F) is recommended.

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every week (Biolis 30i).

Calibration is recommended in the following cases:

- after each change of lot,
- after instrument service,
- if controls lie outside the expected range,
- each time a new reagent kit is used.

If the quality control results do not fall within the expected values or within the range determined in the laboratory, despite a successful calibration procedure, do not report results. In this case, please take the following actions:

- verify reagents are not out of expiration date.

- verify that the required maintenance has been carried out.
- verify that the procedure has been performed in accordance with the instructions for use.
- contact the Service Department or distributor for assistance.

#### REFERENCE VALUES <sup>9-11</sup>

Gender	Age	U/l (37°C)	µkat/l (37°C)
female	1 – 30 days	48 – 406	0.80 – 6.77
	31 days – 1 year	124 – 341	2.07 – 5.68
	1 year – 3 years	108 – 317	1.80 – 5.28
	4 – 15 years	54 – 369	0.91 – 6.23
	16 – 18 years	35 – 124	0.58 – 2.07
	19 – 20 years	39 – 118	0.65 – 1.97
	20 – 50 years	42 – 98	0.71 – 1.67
	50 – 60 years	39 – 118	0.65 – 1.97
male	≥ 60 years	53 – 141	0.90 – 2.40
	1 – 30 days	75 – 316	1.25 – 5.27
	31 days – 1 year	82 – 383	1.37 – 6.38
	1 year – 3 years	104 – 345	1.73 – 5.75
	4 – 15 years	54 – 369	0.91 – 6.23
	16 – 18 years	58 – 331	0.97 – 5.52
	19 – 20 years	41 – 137	0.68 – 2.28
	20 – 50 years	53 – 128	0.90 – 2.18
50 – 60 years	41 – 137	0.68 – 2.28	
≥ 60 years	56 – 119	0.95 – 2.02	

It is recommended for each laboratory to establish its own reference ranges for local population and method principle. The assay results should be used in conjunction with other data, such as symptoms, results of other tests and clinical history to make clinical decisions. It is not recommended to make clinical diagnosis based on a single result.

#### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser Biolis 30i. Results may vary if a different instrument is used.

- LoB (Limit of Blank):**  
1.7 U/l (0.028 µkat/l)
- LoD (Limit of Detection):**  
3.9 U/l (0.065 µkat/l)
- LoQ (Limit of Quantitation):**  
7 U/l (0.12 µkat/l)
- Linearity:**  
up to 900 U/l (µkat/l)

For higher activity, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

- Measurement range:**  
7 U/l (0.12 µkat/l) – 900 U/l (15 µkat/l)

#### Specificity / Interferences

Haemoglobin up to 0.625 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

#### Precision

Repeatability (run to run) n = 20	Mean [U/l]	SD [U/l]	CV [%]
level 1	103.9	0.80	0.77
level 2	445.8	1.65	0.37
Reproducibility (day to day) n = 80	Mean [U/l]	SD [U/l]	CV [%]
level 1	106.7	4.41	4.1
level 2	484.5	15.86	3.3

#### Method comparison

A comparison between ALP values determined at **Biolis 30i** (y) and at **ADVIA SIEMENS 1800** (x) using 63 serum samples gave following results:  
y = 0.9829 x – 2.1291 U/l;  
R = 0.999 (R – correlation coefficient)

#### WASTE MANAGEMENT <sup>12</sup>

After use, the reagents should be handled as potentially infectious and disposed of in accordance with local legal requirements.

- Reagents residues: 18 01 06\*
- Empty packages: 15 01 10\*
- Wastewater from the analyzer: 18 01 03\*

#### INCIDENTS <sup>13</sup>

Any serious incident that has occurred in relations to the device shall be reported to the manufacturer (website address: incidents@cormay.pl) and the competent authority of the Member State in which the user and/or the patient is established.

**Serious incident** means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- a serious public health threat.

#### LITERATURE

- Rifai N., Horvath A.R., Wittwer C., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 6th ed., St. Louis, Missouri, Elsevier, 417, 1393 (2018).
- European Association for the Study of the Liver, EASL Clinical Practice Guidelines: The diagnosis and management of patients with primary biliary cholangitis, Journal of Hepatology, 67, 145-172, (2017).
- Hartleb M. i wsp., Choroby cholestatyczne u dorosłych - wytyczne postępowania Sekcji Hepatologicznej Polskiego Towarzystwa Gastroenterologii, Gastroenterologia Praktyczna, vol. 2, 5-31 (2018).
- Alvarez L, et al., Usefulness of Biochemical Markers of Bone Turnover in Assessing Response to the Treatment of Paget's Disease, Bone, 29 (5), 441-452 (2001).
- Kravets I., Paget's Disease of Bone: Diagnosis and Treatment, The American Journal of Medicine, 131 (11), 1298-1303, (2018)
- Lammers W.J., et al., Levels of Alkaline Phosphatase and Bilirubin Are Surrogate End Points of Outcomes of Patients With Primary Biliary Cirrhosis: An International Follow-up Study, Gastroenterology, 1:47, 1338-1349, (2014)
- Bessey O.A., Lowry O.H. And Brock M.J., A method for the rapid determination of alkaline phosphates with five cubic millimeters of serum, Biol. Chem., 164, 321-9 (1946).
- Bowers G.N.Jr. and McComb R.B., A continuous spectrophotometric method for measuring the activity of serum alkaline phosphatase, Clin Chem., 12(2), 70-89 (1966)
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 4th ed., PA: WB Saunders, 2290 (2006)
- Soldin S. J., Brugnara C., Wong E. C.: Pediatric reference Ranges, 4th ed, AACC Press, 10-11 (2003).
- Hay, W.W., et al, Current pediatric diagnosis and treatment (15th ed.). New York: Lange MedicalBooks/McGraw Hill. (2000).
- European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
- Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices.

#### LIST OF CHANGES

Previous version: 05	Current version: 06
Sections updated: <i>REAGENTS; WARNING AND NOTES; WASTE MANAGEMENT</i>	

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- Rifai N., Horvath A.R., Wittwer C., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 6th ed., St. Louis, Missouri, Elsevier, 417, 1393 (2018).
- European Association for the Study of the Liver, EASL Clinical Practice Guidelines: The diagnosis and management of patients with primary biliary cholangitis, Journal of Hepatology, 67, 145-172, (2017).
- Hartleb M. i wsp., Choroby cholestatyczne u doroslych - wytyczne postepowania Sekcji Hepatologicznej Polskiego Towarzystwa Gastroenterologii, Gastroenterologia Praktyczna, vol. 2, 5-31 (2018).
- Alvarez L, et al., Usefulness of Biochemical Markers of Bone Turnover in Assessing Response to the Treatment of Paget's Disease, Bone, 29 (5), 441-452 (2001).
- Kravets I., Paget's Disease of Bone: Diagnosis and Treatment, The American Journal of Medicine, 131 (11), 1298-1303, (2018)
- Lammers W.J., et al., Levels of Alkaline Phosphatase and Bilirubin Are Surrogate End Points of Outcomes of Patients With Primary Biliary Cirrhosis: An International Follow-up Study, Gastroenterology, 1:47, 1338-1349, (2014)
- Bessey O.A., Lowry O.H. And Brock M.J., A method for the rapid determination of alkaline phosphates with five cubic millimeters of serum, Biol. Chem., 164, 321-9 (1946).
- Bowers G.N.Jr. and McComb R.B., A continuous spectrophotometric method for measuring the activity of serum alkaline phosphatase, Clin Chem., 12(2), 70-89 (1966)
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 4th ed., PA: WB Saunders, 2290 (2006)
- Soldin S. J., Brugnara C., Wong E. C.: Pediatric reference Ranges, 4th ed, AACC Press, 10-11 (2003).
- Hay, W.W., et al, Current pediatric diagnosis and treatment (15th ed.). New York: Lange MedicalBooks/McGraw Hill. (2000).
- European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
- Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices.

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