CK

| | (PL) |
|----------------------|----------------|
| Nazwa zestawu | Nr kat. |
| Liquick Cor-CK 30 | 1-219 |
| Liquick Cor-CK 60 | 1-220 |
| Liquick Cor-CK 120 | 3-330 |
| HC-CK | 4-520 |
| OS-CK | 9-421 |
| B50-CK | 5-531 |

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Data wydania: 08. 2023.

CK

(EN) Stability of working reagent: 4 days at 2-8°C

| Kit name | Cat. No |
|--------------------|---------|
| Liquick Cor-CK 30 | 1-219 |
| Liquick Cor-CK 60 | 1-220 |
| Liquick Cor-CK 120 | 3-330 |
| HC-CK | 4-520 |
| OS-CK | 9-421 |
| B50-CK | 5-531 |

INTENDED USE


Diagnostic kit for determination of creatine kinase activity intended to use both for manual assay and in several automatic analysers. The reagents must be used only for in vitro diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Creatine kinase (CK) catalyzes the transfer of phosphate group between creatine phosphate and adenosine diphosphate (ADP). The product of this reaction is adenosine triphosphate (ATP) – molecular source of energy. CK is a dimer, composed of two different subunits called M and B. Three different isoenzymes formed from these subunits are found in brain and smooth muscle (BB), skeletal muscle (MM) and cardiac muscle (MM and MB). Increased level of CK is usually the result of muscle injury, myocardial or pulmonary infarction.

METHOD PRINCIPLE

Optimized kinetic method according to International Federation of Clinical Chemistry (IFCC).

The rate of NADPH formation, measured as the change in absorbance at $\lambda=340$ nm, is directly proportional to the activity of creatine kinase.

REAGENTS

| Package | Liquick Cor-CK 30 | Liquick Cor-CK 60 | Liquick Cor-CK 120 |
|-----------|-------------------|-------------------|--------------------|
| 1-REAGENT | 5 x 25 ml | 5 x 50 ml | 5 x 100 ml |
| 2-REAGENT | 1 x 25 ml | 1 x 50 ml | 1 x 100 ml |
| | HC-CK | OS-CK | B50-CK |
| 1-REAGENT | 6 x 87.5 ml | 3 x 44.5 ml | 3 x 58.5 ml |
| 2-REAGENT | 6 x 18.5 ml | 3 x 14 ml | 3 x 14.4 ml |

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 4 weeks on board the analyser (Prestige 24i) at 2-10°C.

Working reagent preparation and stability

Assay can be performed with use of separate 1-REAGENT and 2-REAGENT reagents or with use of working reagent. For working reagent preparation mix gently 5 parts of 1-REAGENT with 1 part of 2-REAGENT. Avoid foaming!

Concentrations in the test

| | |
|--|------------|
| 1-REAGENT | |
| imidazole buffer | 100 mmol/l |
| glucose | 20 mmol/l |
| N-acetylcysteine | 20 mmol/l |
| magnesium acetate | 10 mmol/l |
| EDTA | 2 mmol/l |
| NADP | 2 mmol/l |
| ADP | 2 mmol/l |
| AMP | 5 mmol/l |
| HK | > 2.5 U/ml |
| 2-REAGENT | |
| diadenosinepentaphosphate | 10 µmol/l |
| glucose-6-phosphate-dehydrogenase (G6P-DH) | > 1.5 U/ml |
| creatine phosphate preservative | 30 mmol/l |

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Do not freeze reagents.
- Do not interchange caps among reagents.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

1-REAGENT contains imidazole.

Danger

H360D May damage the unborn child.

P201 Obtain special instructions before use.

P308+P313 IF exposed or concerned: Get medical advice.

P405 Store locked up.

ADDITIONAL EQUIPMENT

- automatic analyzer or photometer able to read at 340 nm (334/365 nm); with resolving power of absorbance 0.0001;
- thermostat at 37°C;
- general laboratory equipment;

SPECIMEN

Serum, free from hemolysis.

CK activity is unstable and is rapidly lost during storage. Probes should be stored tightly closed and protected from light.

Specimens can be stored up to 4-8 hours at 15-25 °C or 1-2 days at 2-8°C or 1 month at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-REAGENT and 2-REAGENT are ready to use.

Applications for analyzers are available on request.

Manual procedure

| | |
|-------------|--------------------------|
| wavelength | 340 nm (365 nm / 334 nm) |
| temperature | 37°C |
| cuvette | 1 cm |

Sample Start method

Pipette into the cuvettes:

| | reagent blank (RB) | standard (S) | test (T) |
|-----------------|--------------------|--------------|----------|
| working reagent | 1000 µl | 1000 µl | 1000 µl |

Bring up to the temperature of determination. Then add:

| | | | |
|------------|---|-------|-------|
| sample | - | - | 40 µl |
| calibrator | - | 40 µl | - |

Mix and incubate at adequate temperature (37 °C). After about 2 min. read the absorbance A of standard sample A(S) and test sample A(T) against reagent blank (RB). Repeat the reading after exactly 1, 2 and 3 minutes. Calculate the mean absorbance change per minute for the standard sample $\Delta A/\text{min.}(S)$ and the test sample $\Delta A/\text{min.}(T)$.

Calculation

$$\text{CK activity [U/l]} = \frac{\Delta A/\text{min.}(T)}{\Delta A/\text{min.}(S)} \times \text{calibrator concentration [U/l]}$$

Reagent Start method

The determination can be also performed with use of separate 1-REAGENT and 2-REAGENT reagents.

Pipette into the cuvettes:

| | reagent blank (RB) | standard (S) | test (T) |
|------------|--------------------|--------------|----------|
| 1-REAGENT | 1000 µl | 1000 µl | 1000 µl |
| sample | - | - | 40 µl |
| calibrator | - | 40 µl | - |

Mix gently, incubate for 5 min. Then add:

| | | | |
|-----------|--------|--------|--------|
| 2-REAGENT | 200 µl | 200 µl | 200 µl |
|-----------|--------|--------|--------|

Mix and incubate at adequate temperature (37 °C). After about 2 min. read the absorbance A of standard sample A(S) and test sample A(T) against reagent blank (RB). Repeat the reading after exactly 1, 2, 3 and 4 minutes. Calculate the mean absorbance change per minute for the standard sample $\Delta A/\text{min.}(S)$ and the test sample $\Delta A/\text{min.}(T)$.

Calculation

$$\text{CK activity [U/l]} = \frac{\Delta A/\text{min.}(T)}{\Delta A/\text{min.}(S)} \times \text{calibrator concentration [U/l]}$$

REFERENCE VALUES *

| serum | 37°C | |
|--------|-----------|---------------|
| female | < 167 U/l | < 2.78 µkat/l |
| male | < 190 U/l | < 3.17 µkat/l |

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) with each batch of samples.

For calibration the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or LEVEL 2 (Cat. No 5-175; 5-177) are recommended.

The calibration curve should be prepared every 4 weeks (Prestige 24i), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

These metrological characteristics have been obtained using automatic analysers Biolis 24i Premium or/and Prestige 24i. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity:** 4.4 U/l (0.072 µkat/l).
- Linearity:** up to 1600 U/l (26.7 µkat/l).

For higher activities dilute sample with 0.9% NaCl and repeat the assay. Multiply the result by the dilution factor.

Specificity / Interferences

Haemoglobin up to 0.156 g/dl, bilirubin up to 20 mg/dl, ascorbate up to 62 mg/l and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision

| Repeatability (run to run) n = 10 | Mean [U/l] | SD [U/l] | CV [%] |
|--|------------|----------|--------|
| level 1 | 137.22 | 0.78 | 0.57 |
| level 2 | 509.97 | 1.14 | 0.22 |
| Reproducibility (day to day) n = 20 | Mean [U/l] | SD [U/l] | CV [%] |
| level 1 | 140.43 | 2.33 | 1.66 |
| level 2 | 521.19 | 5.06 | 0.97 |

Method comparison

A comparison between CK values determined at **Biolis 24i Premium** (y) and at **OLYMPUS AU400** (x) using 24 samples gave following results:

$$y = 0.9355x + 2.3019 \text{ U/l;}$$

$$R = 1.0 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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