

ACP

For use on Diatron Pictus® series analyzers

Method: α -Naphthylphosphate – Fast Red

Product code: 1519-0234

Package: 6 x 6 ml (R1) + 1 x 5 ml (R2) + 1 x 5 ml (R3)

Store at: 2 – 8°C

For *in vitro* use only

INTENDED USE


Reagent for the quantitative determination of total and prostatic acid phosphatase in human serum specifically for use with Diatron Pictus® series analyzers. For *in vitro* diagnostic use only.

CLINICAL SIGNIFICANCE

The α -naphthylphosphate – Fast Red method is applied. Acid phosphatase activity is increased in prostate cancer and especially but not always with metastases, multiple myeloma, Paget's disease, sickle cell crisis, Gaucher disease, cirrhosis, hyperparathyroidism, thrombocytosis.

METHOD PRINCIPLE

The determination of acid phosphatase is based on the following reactions:



The α -Naphthol / Fast Red TR formed complex absorbs light at 405/750 nm.

For the determination of prostatic acid phosphatase, the prostatic fraction is specifically inhibited by L-Tartrate, leaving the non prostatic fraction active. The activity of the prostatic fraction is calculated by extracting the activity of non prostatic from total acid phosphatase.

METHOD LIMITATIONS

Refer to the book "Effects of Preanalytical Variables on Clinical Laboratory Tests" for possible interference of other pharmaceutical agents in this particular test. Interference of other agents is described in the "Clinical Guide to Laboratory Tests".

The reagent is designed especially for use with the Diatron Pictus® series of chemistry analyzers. For chemistry protocols and further information please contact the customer support unit at Diatron.

REAGENT COMPOSITION

Reagent 1 (ACP):

Citrate buffer (pH 5.3±0.1): 60 mM

α -naphthyl phosphate: 3 mM

Reagent 2 (Buffer):

Acetate buffer (pH 5.0): 5 M

Reagent 3 (L-Tartrate):

Sodium L-Tartrate: 2 mM

Citric Acid: 70 mM

Sodium Citrate: 10 mM

WARNINGS-PRECAUTIONS

- This reagent is designed for *in vitro* diagnostic use. *In vitro* diagnostic reagents can be hazardous. They should be handled according to good laboratory techniques. Avoid inhalation and contact with eyes and skin.
- Samples should be considered as potentially infectious. Handle with special caution.
- Avoid swallowing and contact with skin and mucous membranes.
- Dispose all waste according to national laws.
- MSDS is available by Diatron or MEDICON HELLAS (manufacturer) upon request.

REAGENT PREPARATION

Reconstitute the ACP reagent R1 by adding exactly 6 ml deionized H₂O. Shake gently. Avoid foaming. Reconstitute reagent R3 (L-Tartrate) with exactly 5 ml deionized H₂O. If needed, warm the reagent at 40 – 50°C to completely dissolve the material.

Reagent R2 (Acetate Buffer) is ready to use. Transfer contents of all vials to analyzer-specific reagent containers supplied with barcodes for recognition by Pictus® series analyzers.

- For the determination of ACP use reconstituted reagent R1 (ACP).
- For the determination of NP-ACP use a vial of reagent R1 (ACP), in which 60 μ l of reagent R3 (L-Tartrate) have been added.
- Reagent R2 (Acetate Buffer) is used for the stabilization of ACP in serum.

REAGENT DETERIORATION

The reagent should not be used

- When the blank measurement of the reconstituted reagent of ACP is greater than 0.4 when measured against H₂O at 405 nm.
- When the L-Tartrate has precipitated. Warm the L-Tartrate solution at 40° – 50°C to redissolve it.

SHELF LIFE

Unopened, the reagents are stable up to the expiry date stated on the label when stored at 2° – 8°C. The reconstituted ACP reagent is stable for 1 day at room temperature and for 7 days stored in the cooled reagent tray of Pictus® series analyzers. The reconstituted L-Tartrate reagent is stable up to the expiry date stated on the label when stored at 2° – 8°C. If recrystallization occurs warm the solution at 40° – 50°C to redissolve the content.

The acetate buffer is stable up to the expiry date stated on the label when stored at 2° – 8°C.

SAMPLE

Specimens of non-hemolyzed serum must be used. Acid phosphatase in serum is particularly unstable. It is necessary to stabilize the ACP in the sample by adding 20 μ l Acetate buffer for every 1 mL serum. ACP in control sera is more unstable. Do not omit, whenever ACP/NP-ACP is to be measured, to add Acetate Buffer in control sera at the same volume ratio as for serum. NOTE: The strongly acidic environment this will create may affect the determination of other analytes, (urea, enzymes, etc.) therefore the samples for ACP/NP-ACP determination should be aliquoted immediately after centrifugation and prior to Acetate buffer addition.

CALIBRATION

Diatron provides MEDI-CAL (1578-0891) for calibration of ACP. Calibrate the assay when a new lot of reagent is installed. The analyzer will automatically perform a Reagent Blank measurement every 14 days. Calibration should be repeated when a new lot of reagent is used, after major maintenance is performed on the analyzer, after a critical part is replaced, or when a significant shift in control values occurs. The calibration factor of NP-ACP is set automatically from the ACP calibration factor.

QUALITY CONTROL

Diatron provides Clinical Chemistry Control Level 1 & 2 (1578-0901-12 & 1578-0902-12 respectively) for serum quality control. Control recovery should lie within the acceptable range. Results outside the acceptable range even after recalibration could be due to reagent deterioration, unsuitable storage conditions or control deterioration, instrument malfunction, or error during test procedure

MATERIALS REQUIRED BUT NOT SUPPLIED WITH THE KIT

- ACP Calibrator
- Quality control materials
- Diatron Pictus® P400/P700/P500
- Common laboratory equipment

CALCULATION

Total ACP (U/L) = $\Delta A/\text{min} \times 853$

Non Prostatic ACP (U/L) = $\Delta A/\text{min} \times 860$

(Prostatic ACP) = (Total ACP) – (Non prostatic ACP)

$\Delta A/\text{min}$: Change of Optical Density / min

Results are automatically calculated by all automated biochemistry analyzers, based on the Calibration Factor and $\Delta A/\text{min}$.

REFERENCE INTERVALS

Total ACP up to 9 U/L

Prostatic ACP up to 3 U/L

SPECIFIC PERFORMANCE CHARACTERISTICS

The following values are representative of the reagent performance on Diatron Pictus® series analyzers. The reagent performance has been evaluated on other types of analyzers, covering all requirements of the 98/79 IVD Directive. A list of analyzers with the corresponding performance characteristics is available in the special leaflet accompanying the insert. The results taken in your laboratory may differ from these values

Linearity	Pictus® P400			Pictus® P700/P500		
	Total:	Prostatic:	CV%	Total:	Prostatic:	CV%
Lowest Detection Limit:	0.4 – 65 U/L	0.4 – 40 U/L	0.3 U/L	0.4 – 65 U/L	0.4 – 40 U/L	0.3 U/L
	0.3 U/L	0.3 U/L	0.3 U/L	0.4 U/L	0.4 U/L	0.4 U/L

The lowest detection limit (LDL) is defined as the lowest concentration of analyte that is distinguishable from zero. A sample free of analyte is assayed 20 times in one run and the LDL is calculated as the absolute mean plus three standard deviations.

Precision: Precision is estimated on two concentration levels of analyte according to CLSI protocol EP5-A (20 consecutive days, 2 runs per day, 2 repeats per run).

	Pictus 400			Pictus 700 / P500		
	Level (U/L)	CV%	Total CV%	Level (U/L)	CV%	Total CV%
Total ACP	10.2	1.55	3.75	8.6	2.79	4.29
	29.1	2.30	3.38	20.0	2.13	4.68
Prostatic ACP	6.3	2.41	3.52	4.6	5.14	5.35
	14.9	1.72	3.02	9.5	3.83	4.53

Interferences Criterion: recovery within $\pm 20\%$ from target value

	Pictus® P400	Pictus® P700/P500
Lipemia	Insignificant up to Intralipid® 1000 mg/dl	Insignificant up to Intralipid® 1000 mg/dl
Haemoglobin	Insignificant up to 75 mg/dL	Insignificant up to 200 mg/dL
Non conj. Bilirubin	Insignificant up to 20 mg/dL	Insignificant up to 20 mg/dL
Conj. Bilirubin	Insignificant up to 2.5 mg/dL	Insignificant up to 2.5 mg/dL
Ascorbate	Insignificant up to 3 mg/dL	Insignificant up to 8 mg/dL
Prostatic ACP		
Lipemia	Insignificant up to Intralipid® 1000 mg/dl	Insignificant up to Intralipid® 1000 mg/dl
Haemoglobin	Insignificant up to 75 mg/dL	Insignificant up to 150 mg/dL
Non conj. Bilirubin	Insignificant up to 20 mg/dL	Insignificant up to 20 mg/dL
Conj. Bilirubin	Insignificant up to 2.5 mg/dL	Insignificant up to 2.5 mg/dL
Ascorbate	Insignificant up to 3 mg/dL	Insignificant up to 8 mg/dL

Correlation A comparison was performed between this reagent on a Pictus® series analyzer, and a BECKMAN COULTER AU-series system. The results were as follows:

Pictus® P400

Total Y = 0.975X + 0.32 R=0.969 N=90 Sample range: 0.55 – 52.40U/L

Prostatic Y = 0.968X + 0.21 R=0.975 N=90 Sample range: 0.15 – 3.49U/L

Pictus® P700/P500

Total Y = 0.983X + 0.22 R=0.974 N=90 Sample range: 0.55 – 52.40U/L

Prostatic Y = 0.975X + 0.14 R=0.982 N=90 Sample range: 0.15 – 3.49U/L

BIBLIOGRAPHY

- Bergmeyer, H.V., Methods of Enzymatic analysis. Weinheim, Verlag chemie, 3rd p. 92 (1984)
- Tietz, N.W., Fundamentals of Clinical Chemistry, Philadelphia: W.B. Saunders Company Ltd., p. 614 (1976)
- Amador, E et al, Am. J. Clin. Path. 51:202
- Hillman, G.Z., Klin. Chem. Klin. Biochem 3:323 (1971)
- Shaw, L.M., et al, Am. J. Clin. Path. 68:57 (1977)

SYMBOLS



Temperature Limits (L/H)



Manufacturer



Read the Instructions



Catalog Number (ISO 15223 / rev. EN980)



Batch Code (ISO 15223 / rev. EN980)



For *in vitro* use (ISO 15223 / rev. EN980)



Date of Expiry (ISO 15223 / rev. EN980)

