


AccuDiagTM Total Human IgA ELISA Kit

Cat# 1802-9



Test	Total Human IgA
Method	Enzyme Linked Immunosorbent Assay
Principle	Sandwich Complex
Detection Range	0.031 mg/mL - 2 mg/mL
Sample	10 μL serum/plasma
Specificity	94 %
Sensitivity	0.031 mg/mL
Total Time	~150 min
Shelf Life	12-14 Months from the manufacturing date

INTENDED USE

To quantitate total human Immunoglobulin A (IgA)

TEST PRINCIPLE

Solid phase capture sandwich ELISA assay using a microwell format.

Patient and Standard Dilutions:

Dilute each serum or plasma specimen to be tested initially 1:100 with phosphate buffered saline (PBS), then dilute 1:100 in the IgA specimen dilute provided. The final dilution factor will be 1:10,000.

Prepare serial two fold dilutions of the human IgA standard: Neat, 1:2, 1:4, 1:8 etc. with the specimen diluent provided. Use the specimen diluent alone as the blank control well.

MATERIALS AND COMPONENTS

Materials provided with the test kits

- 1. Anti-Human IgA coated microwell strips 12x8 with plastic frame
- 2. HRP conjugated goat anti-human IgA -12mL
- 3. IgA standard (pre-diluted)- 1 mL (Store at -20 $^{\circ}$ C)
- 4. TMB/peroxide substrate color developer -12mL
- 5. IgA specimen diluent (Specimen Diluent Green II) -60mL
- 6. Sulfuric acid termination reagent (0.5N) -12mL
- 7. 15 X Wash buffer concentrate 60mL

ASSAY PROCEDURE

* Caution: All human fluids should be treated as infectious agents that could carry HIV.

*Allow each reagent to reach room temperature before use.

- 1. Add 100uL of *diluted* specimen or standard to each microwell.
- 2. Incubate at room temperature for 60 minutes.
- 3. Decant and wash each microwell four times with wash buffer (dilute buffer 1:15 with reagent grade water).
- 4. Add 100uL of HRP conjugated goat anti-human IgA to each well.
- 5. Incubate at room temperature for 60 minutes.
- Decant and wash as in step 3.
- Add 100uL of TMB/peroxide substrate and incubate at room temperature for 30 minutes.
- 8. Terminate the reaction with 100uL of 0.5N sulfuric acid.
- Zero the microwell reader at 450nm using the specimen diluent zero control well.
- 10. Determine the optical density (O.D.) of the remaining wells.
- Construct a standard curve using the O.D. values obtained for each of the standards.
- 12. Interpolate the unknowns from the standard curve.

*Interpolated concentrations greater than 2 mg/dL should be sub diluted 1:4 and re-assayed then corrected mathematically.

LIMITATIONS OF PROCEDURE

No single assay should be used as the only basis for arriving at a diagnostic conclusion. For research use only.

Dynamic Range:


0.031 mg/mL- 2.0 mg/mL

Reproducibility:

C.V. 6%-10% depending upon the region of the standard curve.

Shelf Life

The expiration date for the package and each component is stated on the label(s). Store all components at $2^{\circ}-8^{\circ}$ degrees C. Store all components at $2^{-8^{\circ}}$ C with the exception of the standard, which should be stored frozen.



Diagnostic Automation/Cortez Diagnostics, Inc. 23961 Craftsman Road, Suite D/E/F, Calabasas, CA 91302 USA Phone: 818.591.3030 Fax 818.591.8383 Email:onestep@rapidtest.com Website: www.rapidtest.com

DAI CODE #9 Page 1 of 1