PRESTIGE 24i CYSTATIN C

Nr kat. 4-200, 4-400 (PL)

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Cat. No **4-200, 4-400** (EN)

INTENDED USE

Diagnostic kit for determination of cystatin C concentration used in automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Cystatin C is a low molecular weight protein (13kD), one of the cysteine proteinases inhibitors. Cystatin C is produced in all nucleated cells and secreted into extracellular space at constant rate. Cystatin C molecule stability and dependence of its concentration solely on GFR (Glomerular Filtration Rate) decide about high diagnostic efficiency of cystatin C determination. Cystatin C level is not affected by muscle mass or diet and its increase is observed even at slight GFR reduction. Clinical applications of cystatin C are for monitoring GFR in children, elderly patients, patients with potentially nephrotoxic drug therapy, for assessment of renal transplantation status, for kidney function monitoring in acute and chronic kidney diseases including diabetic nephropathy.

METHOD PRINCIPLE

Turbidimetric method. An antigen-antibody reaction occurs between cystatin C and antibodies coated on polystyrene particles and immuno-complexes are formed. The change of turbidity is related to the quantity of cystatin C in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

REAGENTS

Package

	Cat. No 4-200 (24-TRAY)	Cat. No 4-400 (36-TRAY)
1-Reagent	1 x 29 ml	2 x 18 ml
2-Reagent	1 x 7 ml	2 x 4.4 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 9 weeks on board the analyser at 2-10°C.

Concentrations in the test

suspension of polystyrene particles coated with anti-cystatin C antibodies	8.5 g/l
MOPS buffer [3-(N-morpholino)-propanesulfonic acid]	8 g/l
gentamycin	12.5 mg/l
amphotericin B	1.25 mg/l
preservative	

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- The reagents contain antibiotics and must be handled with due cautions.
- Mix samples well before analysing.

- The samples can be shipped without special cooling and must then be analysed within 14 days after shipment.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum or EDTA/heparinised plasma. Sample stability testing showed that cystatin C in serum and plasma samples are stable for 14 days at room temperature (8-25°C) and for 21 days if stored at 2-8°C, otherwise store frozen at below -20°C. Nevertheless it is recommended to perform the assay with freshly collected samples.

PROCEDURE

1-Reagent and 2-Reagent are ready to use.
1-Reagent put on basic position in reagent tray.
2-Reagent put on start position in reagent tray.
For reagent blank 0.9% NaCl is recommended.

REFERENCE VALUES ⁹

serum	mg/l
newborns < 1 month	1.49 – 2.85
children < 5 months	1.01 – 1.92
< 12 months	0.75 – 1.53
12-24 months (F)	0.60 – 1.20
12-24 months (M)	0.77 – 1.85
2-19 years	0.62 – 1.11
adults (F)	0.61 – 1.05
adults (M)	0.71 – 1.21

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY CYSTATIN C CONTROLS (Cat. No 4-460).

For the calibration of automatic analysers the CORMAY CYSTATIN C CALIBRATORS (Cat. No 5-185) is recommended. The calibration curve should be prepared every 9 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using an automatic analyser Biolis 30i. Results may vary if a different instrument is used.

- LoB (Limit of Blank):** 0.07 mg/l
- LoD (Limit of Detection):** 0.11 mg/l
- LoQ (Limit of Quantitation):** 0.25 mg/l
- Linearity:** up to 7.8 mg/l

Specificity / Interferences

Haemoglobin up to 0.7 g/dl, bilirubin up to 80 mg/dl, triglycerides up to 1400 mg/dl and ascorbate up to 300 mg/l do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean [mg/l]	SD [mg/l]	CV [%]
level 1	0.91	0.01	1.06
level 2	3.62	0.09	2.57
Reproducibility (day to day) n = 80	Mean [mg/l]	SD [mg/l]	CV [%]
level 1	1.0	0.05	4.9
level 2	3.5	0.08	2.4

Method comparison

A comparison between cystatin C values determined at **Biolis 30i** (y) and at **BS-800** (x) using 60 serum samples gave following results:

$$y = 1.0505x - 0.0921 \text{ mg/l}$$

$$R = 0.998 \quad (R - \text{correlation coefficient})$$

A comparison between cystatin C values determined at **Biolis 30i** (y) and at **BS-800** (x) using 30 plasma samples gave following results:

$$y = 0.9146x + 0.143 \text{ mg/l}$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	68	CYS C
Auto Rerun SW		
OFF		
Auto Rerun Range (Result)		
	OFF	OFF
	Lower	Higher
Serum		
Urine		
Plasma		
CSF		
Dialysis		
Other		
Auto Rerun Condition (Absorbance)		
Absorbance Range	Lower	OFF
	Higher	OFF
Auto Rerun Condition (Prozone)		
Prozone Range	OFF	

• **Prestige 24i Premium, Biolis 24i Premium**

Item No.	68	Item Name	CYS C	Optical
Data information				
Units	mg/l			
Decimals	2			
Analysis				
Type	END method			
Main Wave Length	546nm			
Sub Wave Length				
Method	Immunoturbidymetric			
Correlation				
Slope	Intercept			
Y=	1	X+	0	
Calibration				
Type	Spline1			
Std sample conc.				
Blank	0	#1 *	#2 *	*
#3	*	#4 *	#5 *	*
#6	*			

Item No.	68	Item Name	CYS C	Optical
Aspiration				
Kind	Double			
Vol.				
Kind	Vol.	Add	Units	
Sample	2	5	µl	
Reagent 1	160	10	µl	
Reagent 2	30	10	µl	
Data Process				
Read	Start	End		
	Main	51	52	
	Sub	33	34	
Abs.Limit				
Low	-3	High	3	
Blank value				
Water Blank				
Reaction Monitor				
0 Level Point	1			
Span	3			
Correction value				
Blank correction				1
End Point Limit				2
Linear Check (%)				90
Prozone Check				
	Start	End	Limit (%)	
First				
Second			Low	

Item No.	68	Item Name	CYS C	Optical
Normal Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				
Panic Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	68	Item Name	CYS C	Optical			
Auto Rerun SW							
OFF							
Auto Rerun Range (Conc.)							
	First Dil	Low		High			
		Re	Value	Dil	Re	Value	Dil
Serum							
Urine							
Plasma							
CSF							
Dialysis							
Other							
Auto Rerun Condition (Absorbance)							
Lower		OFF					
Higher		OFF					
Auto Rerun Condition (Prozone)							
OFF							
Dilution							
99:Dil1							

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