## PRESTIGE 24i LQ TG mono

Nr kat. 4-226, 4-426 (PL)

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Cat. No **4-226, 4-426** (EN)

### INTENDED USE

Diagnostic kit for determination of triglycerides concentration used in automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.


The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Triglycerides are built of glycerol molecule esterified with three fatty acids molecules. Triglycerides are delivered with food or are synthesized endogenously in liver. Triglycerides stored in adipose tissue constitute a reserve of energy. Elevated triglycerides serum level is a risk factor of atherosclerosis. Triglycerides measurement is useful for hyperlipidemia diagnosis and treatment or for estimation of atherosclerosis progression.

### METHOD PRINCIPLE

Colorimetric, enzymatic method with glycerophosphate oxidase.

The colour intensity is proportional to the triglycerides concentration.

### REAGENTS

#### Package

	Cat. No <b>4-226</b> <b>(24-TRAY)</b>	Cat. No <b>4-426</b> <b>(36-TRAY)</b>
1-Reagent	6 x 58 ml	8 x 33 ml

The reagent when stored at 2-8°C are stable up to expiry date printed on the package. The reagent is stable for 8 weeks on board the analyser at 2-10°C (Prestige 24i).

### Concentrations in the test

buffer TRIS (PH 8.0)	200 mmol/l
4-aminoantipyrine (4-AA)	< 0.4 mmol/l
ATP	< 1.5 mmol/l
Mg <sup>2+</sup>	< 1.6 mmol/l
4-chlorophenol	< 2.5 mmol/l
chlorophenicol	1.6 mmol/l
potassium hexacyanoferrate (II)	< 1 mmol/l
FAD-2Na	< 1 mmol/l
glycerol kinase (GK)	~2500 U/l
glycerol phosphate oxidase (GPO)	~2500 U/l
peroxidase (POD)	~1900 U/l
lipoprotein lipase (LPL)	~2000 U/l
detergents, preservatives	

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

### SPECIMEN

Serum, EDTA or heparinized plasma (recommended: heparine lithium, sodium or ammonium salt) **free from hemolysis**.

Blood should be collected only if the patient has been fasting for minimum of 12 hours. Before blood collection patient should stay in rest position for about 30 minutes. Venous blood is recommended for triglycerides measurement.

Plasma triglycerides values have been reported to be 2% to 4% lower than serum triglycerides values.

Serum and plasma can be stored up to 3 days at 2-8°C or up to 3 months at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent is ready to use.

1-Reagent put on basic position in reagent tray.

For reagent blank deionized water is recommended.

#### Actions required:

- Biolis 24i Premium:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: GLUCOSE – TG MONO, HDL DIRECT – TG MONO, LDL DIRECT – TG MONO, LIPASE – TG MONO, LIPOPROTEIN(a) – TG MONO. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_008\_BIOLIS\_24i\_PREMIUM\_CARRYOVER.
- Biolis 30i:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: TG MONO – CALCIUM ARSENAZO, HDL DIRECT – TG MONO. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_009\_BIOLIS\_30i\_CARRYOVER.

### REFERENCE VALUES <sup>7</sup>

serum, plasma	< 150 mg/dl < 1.7 mmol/l
---------------	-----------------------------

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control, it is recommended to use, with each batch of samples the CORMAY SERUM HN (Cat. No 5-172), CORMAY SERUM HP (Cat. No 5-173) and CORMAY LIPID CONTROL LEVEL 1 (Cat. No 5-179), CORMAY LIPID CONTROL LEVEL 2 (Cat. No 5-180).

For the calibration of automatic analysers Prestige 24i the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

For the calibration of automatic analysers: Biolis 24i Premium, Biolis 30i the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 8 weeks (Prestige 24i), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**  
0.2 mg/dl (0.0023 mmol/l)
- LoD (Limit of Detection):**  
1.4 mg/dl (0.016 mmol/l)
- LoQ (Limit of Quantitation):**  
3.5 mg/dl (0.04 mmol/l)
- Linearity:**  
up to 1165 mg/dl (13.2 mmol/l)

For higher triglycerides concentrations dilute the sample with 0.9% NaCl in the ratio of 1 to 4 and repeat the assay. Multiply the result by 5.

#### Specificity / Interferences

Haemoglobin up to 0.31 g/dl, bilirubin up to 8.6 mg/dl and ascorbate up to 31 mg/l do not interfere with the test.

#### Precision

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	85.1	0.62	0.73
level 2	192.7	1.07	0.55
Reproducibility (day to day) n = 80	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	91.7	2.47	2.7
level 2	178.4	2.26	1.3

#### Method comparison

A comparison between triglycerides values determined at **Biolis 30i** (y) and at **ADVIA SIEMENS 1800** (x) using 69 serum samples, gave following results:

$$y = 1.0129x - 1.3173 \text{ mg/dl}$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	11	TG mono
<b>Auto Rerun SW</b>		
ON		
<b>Auto Rerun Range (Result)</b>		
	ON	ON
	Lower	Higher
Serum	5.3	1100
Urine		
Plasma		
CSF		
Dialysis		
Other		
<b>Auto Rerun Condition (Absorbance)</b>		
Absorbance Range	Lower	OFF
	Higher	OFF
Prozone Range	OFF	

▪ Prestige 24i Premium, Biolis 24i Premium

Item No.	11	Item Name	TG mono	Optical	
<b>Data information</b>					
Units	mg/dl				
Decimals	1				
<b>Calibration</b>					
Type	Linear2				
Std sample conc.					
Blank	0	#1	*	#2	*
#3		#4		#5	
#6					
<b>Analysis</b>					
Type	END method				
Main Wave Length	546nm				
Sub Wave Length	800nm				
Method	GPO-POD				
<b>Correlation</b>					
Y=	Slope	1	X+	Intercept	0

Item No.	11	Item Name	TG mono	Optical
<b>Aspiration</b>				
Kind	Single			
<b>Data Process</b>				
Read	Main	Start	End	
	Sub	9	10	
<b>Abs.Limit</b>				
Low	-3.0	High	3.0	
<b>Blank value</b>				
Water Blank				
<b>Reaction Monitor</b>				
0 Level Point	1			
Span	3			
<b>Correction value</b>				
Blank correction	1			
End Point Limit	2			
Linear Check (%)	0			
<b>Prozone Check</b>				
	Start	End	Limit (%)	
First				
Second				Low
<b>Third mixing</b>				
OFF				

Item No.	11	Item Name	TG mono	Optical
<b>Normal Range</b>				
	Male		Female	
	Low	High	Low	High
Serum	0	150	0	150
Urine				
Plasma				
CSF				
Dialysis				
Other				
<b>Panic Range</b>				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	11	Item Name	TG mono	Optical	
<b>Auto Rerun SW</b>					
ON					
<b>Auto Rerun Range (Conc.)</b>					
	First Dil	Low Value	Dil	High Value	Dil
Serum		2.6		1000	
Urine					
Plasma					
CSF					
Dialysis					
Other					
<b>Auto Rerun Condition (Absorbance)</b>					
	Lower	OFF			
	Higher	OFF			
<b>Auto Rerun Condition (Prozone)</b>					
OFF					
<b>Dilution</b>					
100:Di12					

## PRESTIGE 24i LQ TG mono

▪ **Biolis 30i**

Item no	11	Item name	TG MONO	Specimen	SERUM/ PLASMA	OPTICAL
<b>Data information</b>			<b>Aspiration volume</b>			
UNITS	mg/dL		TYPE	Single		
DECIMALS	1					
<b>Analysis</b>			<b>Data processing read</b>			
METHOD	END method					
Main Wave Length	546 nm		VOL. (µL)	SAMPLE	REAGENT 1	REAGENT 2
Sub Wave Length	800 nm		BOTTLE (ml)	6	300	
<b>CORRELATION (Y= AX + B)</b>			FIRST DIL.			
A =	1					
B =	0		<b>Blank value</b>			
<b>Calibration</b>			ABS LIMIT	-3 TO 3		
TYPE	Linear 2		END POINT	2.5		
STABILITY			LINEARITY CHECK (%)	0		
<b>Prozone check</b>						
			START	END	LIMIT (%)	
			FIRST			
			SECOND			
			MINIMUM ABS.			
			°HIGH	MEAN		
			°LOW	VARIATE		

Item No	11	Item Name	TG MONO	Specimen	SERUM/ PLASMA	OPTICAL				
<b>Reference intervals</b>			<b>Auto rerun</b>							
MALE		FEMALE		°ON °OFF						
LOW	HIGH	LOW	HIGH							
0	150	0	150	<b>Auto rerun range (conc.)</b>						
<b>Panic range</b>			Re	Value	Dil.	Re	Value	Dil.		
MALE		FEMALE		3.5			1165			
LOW	HIGH	LOW	HIGH	<b>Auto rerun condition (abs.)</b>						
									DIL.	
<b>Decision limit</b>			LOWER	°ON		°OFF				
			HIGH	°ON		°OFF				
			<b>Auto rerun condition (prozone)</b>							
			°ON		°OFF					
			SAMPLE VOL.							
			<b>Dilution</b>							
			°DIL 1		° DIL 2					
<b>Reaction check</b>										
			°ON		°OFF					
			CHECK							
			LOW							
			HIGH							
<b>VL CHECK</b>			°ON		°OFF					
			°OFF							
<b>VH CHECK</b>			°ON		°OFF					
			°OFF							

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