## PRESTIGE 24i LQ LDH

Nr kat. **4-239, 4-439** (PL)

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Cat. No **4-239, 4-439** (EN)

### INTENDED USE

Diagnostic kit for determination of lactate dehydrogenase activity intended to use in automatic analyzers: automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.


The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Lactate dehydrogenase (LDH, LD) is intracellular enzyme occurred in all tissues. LDH catalyzes the reversible conversion of lactate to pyruvate using NAD<sup>+</sup> as a cofactor. LD is a tetramer containing two possible forms of subunits: H and M. The result is five isoenzymes termed LD-1 (H<sub>4</sub>) through LD-5 (M<sub>4</sub>). The isoenzymes are present in different proportion in each tissue and have different electrophoresis mobility, what is very useful for diagnostic.

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The rate of absorbance changing at  $\lambda=340$  nm is directly proportional to lactate dehydrogenase activity.

### REAGENTS

#### Package

	Cat. No 4-239 (24-TRAY)	Cat. No 4-439 (36-TRAY)
1-Reagent	4 x 40 ml	4 x 23 ml
2-Reagent	4 x 12.5 ml	4 x 7.5 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. Stability on board of the analyser at 2-10°C: Prestige 24i - 8 weeks, Biolis 24i Premium - 8 weeks.

#### Concentrations in the test

phosphate buffer (pH 7.5)	50 mmol/l
pyruvate	0.6 mmol/l
NADH	0.25 mmol/l
preservative	

#### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

### SPECIMEN

Serum, heparinized plasma free from hemolysis.

Do not use hemolyzed blood or serum because erythrocytes contain 150 times more LDH activity than serum.

As an anticoagulant for plasma preparation use heparin lithium or ammonium salt!

LDH activity is unstable and is rapidly lost during storage.

Specimens can be stored up to 4 hours at 15-25°C or 1-2 days at 2-8°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank deionized water is recommended.

#### Actions required:

- Biolis 24i Premium:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: ALAT (II GEN) – LDH, ASAT (II GEN) - LDH, ALAT – LDH, ASAT - LDH. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_008\_BIOLIS\_24i\_PREMIUM\_CARRYOVER.
- Biolis 30i:** When performing assays in analyser, there is a probability of **cross-contamination** affecting the tests results: LDH – ALAT (II GEN), ALAT (II GEN) - LDH. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_009\_BIOLIS\_30i\_CARRYOVER.

### REFERENCE VALUES <sup>4</sup>

serum / plasma	37°C	37°C
adults	225 – 450 U/l	3.75 – 7.50 $\mu$ kat/l

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 5 weeks (Prestige 24i) or 8 weeks (Biolis 24i Premium), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):** 2.6 U/l (0.043  $\mu$ kat/l).

- LoD (Limit of Detection):** 6.1 U/l (0.10  $\mu$ kat/l).

- LoQ (Limit of Quantitation):** 11 U/l (0.18  $\mu$ kat/l).

- Linearity:** up to 1825 U/l (30.42  $\mu$ kat/l).

For higher activity dilute sample with 0.9% NaCl and repeat the assay. Multiply the result by the dilution factor.

#### Specificity / Interferences

Bilirubin up to 20 mg/dl, haemoglobin up to 5 g/dl, ascorbate up to 62 mg/l and triglycerides up to 1000 mg/dl do not interfere with the test.

#### Precision

Repeatability (run to run) n = 20	Mean [U/l]	SD [U/l]	CV [%]
level 1	337.2	5.24	1.55
level 2	802.5	4.86	0.61
Reproducibility (day to day) n = 80	Mean [U/l]	SD [U/l]	CV [%]
level 1	347.0	3.74	1.1
level 2	824.1	9.38	1.1

#### Method comparison

A comparison between LDH values determined at **Biolis 30i** (y) and at **BS-400** (x) using 60 serum samples gave following results:

$$y = 1.0134x - 0.823 \text{ U/l};$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	15	LDH
<b>Auto Rerun SW</b>		
ON		
<b>Auto Rerun Range (Result)</b>		
	ON	ON
	Lower	Higher
Serum	25	2000
Urine		
Plasma		
CSF		
Dialysis		
Other		
<b>Auto Rerun Condition (Absorbance)</b>		
Absorbance Range	Lower	OFF
	Higher	OFF
<b>Auto Rerun Condition (Prozone)</b>		
Prozone Range	OFF	

▪ **Prestige 24i Premium, Biolis 24i Premium**

Item No.	15	Item Name	LDH	Optical
<b>Data information</b>				
Units	U/L			
Decimals	1			
<b>Calibration</b>				
Type	Linear2			
Std sample conc.				
Blank	0	#1	*	#2
#3		#4		#5
#6				
<b>Analysis</b>				
Type	RATE method			
Main Wave Length	340nm			
Sub Wave Length	450nm			
Method	DGKC			
<b>Correlation</b>				
Y=	Slope	X+	Intercept	
	1	0		

Item No.	15	Item Name	LDH	Optical
<b>Aspiration</b>				
Kind	Double			
Vol.				
	Kind	Vol.	Add	Units
Sample	3	5		µl
Reagent 1	160	10		µl
Reagent 2	40	10		µl
<b>Data Process</b>				
Read	Main	Start	End	
	Sub	36	52	
<b>Abs.Limit</b>				
Low	1	~	High	2
<b>Blank value</b>				
Water Blank				
<b>Reaction Monitor</b>				
0 Level Point	1			
Span	3			
<b>Third mixing</b>				
OFF				
<b>Correction value</b>				
Blank correction				
End Point Limit				
2				
Linear Check (%)				
90				
<b>Prozone Check</b>				
First	Start	End	Limit (%)	
Second				Low
Third				Low

Item No.	15	Item Name	LDH	Optical
<b>Normal Range</b>				
	Male		Female	
	Low	High	Low	High
Serum	225	450	225	450
Urine				
Plasma	225	450	225	450
CSF				
Dialysis				
Other				
<b>Panic Range</b>				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	15	Item Name	LDH	Optical
<b>Auto Rerun SW</b>				
ON				
<b>Auto Rerun Range (Conc.)</b>				
	First Dil	Low	High	
		Re	Value	Dil
Serum			20.1	
Urine				
Plasma			2000	
CSF				
Dialysis				
Other				
<b>Auto Rerun Condition (Absorbance)</b>				
Lower	OFF			
Higher	OFF			
<b>Auto Rerun Condition (Prozone)</b>				
OFF				
<b>Dilution</b>				
100:DiI2				

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