PRESTIGE 24i LQ HBDH

Nr kat. 4-241, 4-441 (PL)

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Cat. No **4-241, 4-441** (EN)

INTENDED USE

Diagnostic kit for determination of α -hydroxybutyrate dehydrogenase activity used in automatic analysers Prestige 24i, Biolis 24i and Prestige 24i Premium, Biolis 24i Premium.


The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Lactate dehydrogenase (LDH, LD) is a tetrameric molecule containing two possible forms of subunits (H and M). The result is five isoenzymes, one of which is hydroxybutyrate dehydrogenase (HBDH, LD-1) formed by four H subunits. HBDH is present mainly in heart muscle but occur also in kidney and erythrocytes. Normal serum contains mostly LD-2 with lesser amount of LD-1. Changes in the ratio of LD-1 to LD-2 indicate myocardial infarction or hemolysis.

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The rate of absorbance changing at $\lambda=340$ nm is directly proportional to α -hydroxybutyrate dehydrogenase activity.

REAGENTS

Package

	Cat. No 4-241 (24-TRAY)	Cat. No 4-441 (36-TRAY)
1-Reagent	4 x 40 ml	4 x 23 ml
2-Reagent	4 x 12.5 ml	4 x 7.5 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. Stability on board of the analyser at 2-10°C: Prestige 24i - 8 weeks, Biolis 24i Premium - 8 weeks.

Concentrations in the test

phosphate buffer (pH 7.5)	50 mmol/l
2-oxybutyrate	3 mmol/l
NADH	0.25 mmol/l
preservative	

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum.

Do not use hemolyzed blood because erythrocytes contain very high HBDH activity. Do not chill or freeze samples. HBDH activity is unstable and is rapidly lost during storage. Specimens can be stored up to 6 hours at 15-25°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank deionized water is recommended.

Application guide for Biolis 30i analyser is available on request.

REFERENCE VALUES ⁷

serum	37°C
adults	< 182 U/l (< 3.04 μ kat/l)

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 8 weeks (Prestige 24i, Biolis 24i Premium), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers Prestige 24i and Biolis 24i Premium. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity:**
25 U/l (5.43 μ kat/l) - Prestige 24i
9.21 U/l (0.153 μ kat/l) - Biolis 24i Premium
- Linearity:**
up to 500 U/l (8.33 μ kat/l) - Prestige 24i
up to 500 U/l (8.33 μ kat/l) - Biolis 24i Premium

For higher activity dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 2.5 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision

Repeatability (run to run)		Mean [U/l]	SD [U/l]	CV [%]
Prestige 24i n=20	level 1	117.19	2.45	2.09
	level 2	343.09	3.43	1.00
Biolis 24i Premium n=20	level 1	141.63	2.12	1.50
	level 2	362.68	4.52	1.25
Reproducibility (day to day)		Mean [U/l]	SD [U/l]	CV [%]
Prestige 24i n=80	level 1	131.44	3.24	2.47
	level 2	383.35	4.78	1.25
Biolis 24i Premium n=80	level 1	148.94	2.03	1.37
	level 2	380.33	3.28	0.86

Method comparison

A comparison between HBDH values determined at **Prestige 24i** (y) and at **ACCENT-300** (x) using 63 samples gave following results:

$$y = 1.0514x - 10.717 \text{ U/l;}$$

$$R = 0.9613 \quad (R - \text{correlation coefficient})$$

A comparison between HBDH values determined at **Biolis 24i Premium** (y) and at **Prestige 24i** (x) using 101 samples gave following results:

$$y = 1.0432x - 7.6087 \text{ U/l;}$$

$$R = 0.9945 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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• **Prestige 24i Premium, Biolis 24i Premium**

Item No.	16	Item Name	HBDH	Optical	
Data information					
Units	U/L				
Decimals	1				
Calibration					
Type	Linear2				
Std sample conc.					
Blank	0	#1	*	#2	*
#3		#4		#5	
#6					
Analysis					
Type	RATE method				
Main Wave Length	340nm				
Sub Wave Length	450nm				
Method	DGKC				
Correlation					
	Slope	Intercept			
Y=	1	X+	0		
Item No.	16	Item Name	HBDH	Optical	
Aspiration					
Kind	Double				
Vol.					
	Kind	Vol.	Add	Units	
Sample	3		5	µl	
Reagent 1	200		10	µl	
Reagent 2	50		10	µl	
Data Process					
Read	Main	Start	End		
	Sub	35	52		
Abs.Limit	Low	High			
	1	2			
Correction value					
Blank correction					
End Point Limit	2				
Linear Check (%)	90				
Prozone Check					
	Start	End	Limit (%)		
First					
Second			Low		
Item No.	16	Item Name	HBDH	Optical	
Normal Range					
	Male		Female		
	Low	High	Low	High	
Serum	0	182	0	182	
Urine					
Plasma					
CSF					
Dialysis					
Other					
Panic Range					
	Male		Female		
	Low	High	Low	High	
Serum					
Urine					
Plasma					
CSF					
Dialysis					
Other					
Item No.	16	Item Name	HBDH	Optical	
Auto Rerun SW					
ON					
Auto Rerun Range (Conc.)					
	First Dil	Low Value	Dil	High Value	Dil
Serum		9.2		500	
Urine					
Plasma					
CSF					
Dialysis					
Other					
Auto Rerun Condition (Absorbance)					
	Lower	OFF			
	Higher	OFF			
Auto Rerun Condition (Prozone)					
	OFF				
Dilution					
	100:Dil2				

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