## PRESTIGE 24i LQ BIL DIRECT

Nr kat. **4-248, 4-448** (PL)

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(EN)

**2-REAGENT**

sodium metavanadate 4.1 mmol/l  
phosphate buffer  
pH adjuster  
stabilizer

### REAGENT STABILITY

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 12 weeks on board the analyser at 2-10°C (Biolis 30i).

### WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- Do not freeze reagents.
- Do not interchange caps.
- Reagent bottles should be shaken before use by gently inverting several times.
- The appearance of turbidity or control sera values outside the manufacturer's acceptable range may indicate of the reagents instability.
- Lack of significant changes in the color of the reaction mixture at the samples with low bilirubin concentration does not indicate the assay malfunction.
- Do not use the reagent beyond expiry date printed on the package.
- Do not mix reagents from different kits or lots.
- Use personal protective equipment to prevent contact with samples, reagents and controls.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT contains hydroxylammonium chloride. May produce an allergic reaction (EUH208).

### SPECIMEN<sup>6,7,8</sup>

Serum free from hemolysis.  
Serum should be separated from red blood cells as soon as possible after blood collection.  
Lipemic specimens may show falsely decreased bilirubin concentration thus fasting specimen is recommended.  
Because bilirubin is photooxidized when exposed to light, specimen should be protected from direct exposure to either artificial light or sunlight. Therefore it is essential to store specimens in the dark! Serum can be stored at 2-8°C or at room temperature (20-25°C) up to 2 days.

Nevertheless it is recommended to perform the assay with freshly collected samples!

Follow tube manufacturers' instructions carefully when using collection tubes.

Human-origin material should be handled as potentially infectious. Standard precautions in normal laboratory work are required.

### ASSAY PROCEDURE

1-REAGENT and 2-REAGENT are ready to use.  
1-REAGENT put on basic position in reagent tray.  
2-REAGENT put on start position in reagent tray.

### INTENDED USE

BIL DIRECT reagent is intended to determine quantitatively direct bilirubin level in serum. It is intended to diagnosis, monitoring and as an aid to diagnosis of clinical conditions connected with abnormal bilirubin level. BIL DIRECT reagent is intended to use on automatic analysers Biolis 24i Premium and Biolis 30i. It is only for *in vitro* diagnostics, for healthcare professional users.

### SUMMARY<sup>1,2,3</sup>

Determinations of direct bilirubin level in serum are mainly used in diagnosis of liver functions. Measurement of direct bilirubin concentration is used as an aid in diagnosis of clinical conditions associated with hyperbilirubinemia. Increased direct bilirubin levels can be observed in mechanical jaundice, Dubin-Johnson syndrome, bile ducts diseases, gallbladder diseases. Direct bilirubin is also the marker of perforation in acute sigmoid diverticulitis.

### METHOD PRINCIPLE<sup>4,5</sup>

Method is based on chemical oxidation, utilizing vanadate as an oxidizing agent.

In the presence of detergent and vanadate in an acidic solution, conjugated (direct) bilirubin is oxidized to produce biliverdin.

This oxidation reaction causes change of the yellow colour, which is specific to bilirubin to the green colour typical for biliverdin. Therefore, the direct bilirubin concentration in the sample can be obtained by measuring the absorbance before and after the vanadate oxidation.

### REAGENTS

Package	Cat. No 4-248 (24-TRAY)	Cat No 4-448 (36-TRAY)
1-REAGENT	6 x 37 ml	6 x 21 ml
2-REAGENT	6 x 11 ml	6 x 6.5 ml

### Number of tests:

Biolis 24i Premium	710	400
Biolis 30i	800	460

### CONCENTRATION OF THE ACTIVE INGREDIENTS IN THE REAGENT

<b>1-REAGENT</b>	
citrate buffer	
hydroxylammonium chloride	
stabilizers	
detergent	0.34 %
pH adjuster	
preservative	

For reagent blank deionized water is recommended.

#### Actions required:

- **Biolis 24i Premium:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: CREATININE – BIL DIRECT. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_008\_BIOLIS\_24i\_PREMIUM\_CARRYOVER.
- **Biolis 30i:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: CREATININE – BIL DIRECT. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_009\_BIOLIS\_30i\_CARRYOVER.

Application guide for Prestige 24i analyser is available on request.

The application guide should be verified and validated by the user prior to testing patient samples.

#### QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173). For the calibration of automatic analysers the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 10 weeks (Biolis 30i).

Calibration is recommended in the following cases:

- after each change of lot,
- after instrument service,
- if controls lie outside the expected range,
- each time a new reagent kit is used.

If the quality control results do not fall within the expected values or within the range determined in the laboratory, despite a successful calibration procedure, do not report results. In this case, please take the following actions:

- verify reagents are not out of expiration date.
- verify that the required maintenance has been carried out.
- verify that the procedure has been performed in accordance with the instructions for use.
- contact the Service Department or distributor for assistance.

#### REFERENCE VALUES <sup>3</sup>

serum (adults)	0.1 – 0.3 mg/dl
	1.7 – 5.1 µmol/l

It is recommended for each laboratory to establish its own reference ranges for local population and method principle. The assay results should be used in conjunction with other data, such as symptoms, results of other tests and clinical history to make clinical decisions. It is not recommended to make clinical diagnosis based on a single result.

#### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers Biolis 30i or BS-400. Results may vary if a different instrument is used.

- **LoB (Limit of Blank):**  
0.00 mg/dl (0.00 µmol/l) – BS-400
- **LoD (Limit of Detection):**  
0.01 mg/dl (0.171 µmol/l) – BS-400
- **LoQ (Limit of Quantitation):**  
0.04 mg/dl (0.684 µmol/l) – Biolis 30i
- **Linearity:**  
up to 40 mg/dl (684 µmol/l) – Biolis 30i
- **Measurement range:**  
0.04 mg/dl (0.684 µmol/l) – 40 mg/dl (684 µmol/l) – Biolis 30i
- **Specificity / Interferences**  
Ascorbic acid up to 62 mg/l and triglycerides up to 550 mg/dl (in samples containing low concentration of direct bilirubin) and up to 900 mg/dl (in samples containing high concentration of direct bilirubin) do not interfere with the test. Haemoglobin interferes even in small amount with the determination.

#### ▪ Precision (Biolis 30i)

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
poziom 1	0.35	0.01	2.9
poziom 2	1.00	0.01	1.0
Reproducibility (day to day) n = 80	Mean [mg/dl]	SD [mg/dl]	CV [%]
poziom 1	0.36	0.01	3.2
poziom 2	1.02	0.07	6.7

#### ▪ Method comparison

A comparison between direct bilirubin values determined at **Biolis 30i** (y) and at **BECKMAN COULTER AU680** (x) using 120 serum samples gave following results:

$$y = 1.0564x - 0.0196 \text{ mg/dl}$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

#### WASTE MANAGEMENT <sup>9</sup>

After use, the reagents should be handled as potentially infectious and disposed of in accordance with local legal requirements.

- Residual reagents:  
1-REAGENT – 18 01 07  
2-REAGENT – 18 01 07
- Empty packages:  
1-REAGENT – 15 01 02  
2-REAGENT – 15 01 02
- Wastewater from the analyzer:  
1-REAGENT – 18 01 03\*  
2-REAGENT – 18 01 03\*

#### INCIDENTS <sup>10</sup>

Any serious incident that has occurred in relations to the device shall be reported to the manufacturer (websites address: incidents@cormay.pl) and the competent authority of the Member State in which the user and/or the patient is established.

**Serious incident** means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or other person.
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health.
- a serious public health threat.

#### LITERATURE

1. Rifai N., Horvath A.R., Wittwer C., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 6th ed.. St. Louis, Missouri. Elsevier. 767-771 (2018).
2. Käser S.A., Fankhauser G., Glauser P.M. et al. Diagnostic Value of Inflammation Markers in Predicting Perforation in Acute Sigmoid Diverticulitis. World J Surg 34. 2717–2722 (2010).
3. Pagana K. Pagana T. Pagana T. Mosby's Diagnostic & Laboratory Test Reference 14th ed.. St. Louis, Missouri. Elsevier. 137-139 (2019).
4. Tokuda K., Tanimoto K. A new method of measuring bilirubin in serum by vanadic acid. Jpn J Clin Chem (in Japanese with English abstract) 22: 116–22 (1993).
5. Dhungana N., Morris C., Krasowski M.D. Operational impact of using a vanadate oxidase method for direct bilirubin measurements at an academic medical center clinical laboratory. Practical laboratory medicine 8: 77-85 (2017).
6. Tietz N.W., ed. Clinical Guide to Laboratory Test. 3rd ed. Philadelphia: WB Saunders 518-525 (1995).
7. Kaplan L.A., Pesce A.J., Clinical Chemistry: Theory, analysis, and correlation. 2nd ed.. The C.V. Mosby Company. St. Louis, USA. 1013 (1989).
8. World Health Organization (2002). Use of Anticoagulants in Diagnostic Laboratory Investigations (Stability of Blood, Plasma and serum samples); WHO/DIL/LAB/99.1 Rev.2: 1-64.
9. European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
10. Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices.

#### LIST OF CHANGES

Previous version: 01	Current version: 02
<i>Sections updated: INTENDED USE; SUMMARY; REAGENTS; CONCENTRATION OF THE ACTIVE INGREDIENTS; WARNING AND NOTES; SPECIMEN; ASSAY PROCEDURE; QUALITY CONTROL; REFERENCE VALUES; PERFORMANCE CHARACTERISTICS; WASTE MANAGEMENT; LITERATURE; APPLICATION</i>	
<i>Sections added: REAGENT STABILITY. INCIDENTS</i>	

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- Rifai N, Horvath A.R., Wittwer C., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 6th ed., St. Louis, Missouri. Elsevier. 767-771 (2018).
- Käser S.A., Fankhauser G., Glauser P.M. et al. Diagnostic Value of Inflammation Markers in Predicting Perforation in Acute Sigmoid Diverticulitis. World J Surg 34. 2717–2722 (2010).
- Pagana K, Pagana T, Pagana T, Mosby’s Diagnostic & Laboratory Test Reference 14th ed., St. Louis, Missouri. Elsevier. 137-139 (2019).
- Tokuda K., Tanimoto K. A new method of measuring bilirubin in serum by vanadic acid. Jpn J Clin Chem (in Japanese with English abstract) 22: 116–22 (1993).
- Dhungana N., Morris C., Krasowski M.D. Operational impact of using a vanadate oxidase method for direct bilirubin measurements at an academic medical center clinical laboratory. Practical laboratory medicine 8: 77-85 (2017).
- Tietz N.W., ed. Clinical Guide to Laboratory Test. 3rd ed. Philadelphia: WB Saunders 518-525 (1995).
- Kaplan L.A., Pesce A.J., Clinical Chemistry: Theory, analysis, and correlation. 2nd ed., The C.V. Mosby Company, St. Louis, USA. 1013 (1989).

- World Health Organization (2002). Use of Anticoagulants in Diagnostic Laboratory Investigations (Stability of Blood, Plasma and serum samples); WHO/DIL/LAB/99.1 Rev.2: 1-64.
- European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
- Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices.

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Item No.	4	Item Name	Bil D	Optical
<b>Normal Range</b>				
	Male		Female	
	Low	High	Low	High
Serum	0.1	0.3	0.1	0.3
Urine				
Plasma				
CSF				
Dialysis				
Other				
<b>Panic Range</b>				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	4	Item Name	Bil D	Optical	
<b>Auto Rerun SW</b>					
ON					
<b>Auto Rerun Range (Conc.)</b>					
	First Dil	Low		High	
		Re	Value	Dil	Re
Serum			0.04		40
Urine					
Plasma					
CSF					
Dialysis					
Other					
<b>Auto Rerun Condition (Absorbance)</b>					
	Lower	OFF			
	Higher	OFF			
<b>Auto Rerun Condition (Prozone)</b>					
OFF					
<b>Dilution</b>					
100:Dil2					

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