



PRESTIGE 24i DIGOXIN

DIAGNOSTIC KIT FOR DETERMINATION OF DIGOXIN CONCENTRATION

INTRODUCTION

Digoxin is a cardiac glycoside used to treat heart failure and atrial fibrillation with rapid ventricular rhythm. It increases the force and excitability of the heart muscle, slows down conduction in atrioventricular node and also reduces heart rate. Digoxin is a glycoside obtained from *Digitalis lanata*. It is excreted mainly unchanged via the kidneys. Digoxin concentration monitoring is recommended because of a narrow therapeutic range and serious side effects: arrhythmia, disturbance of conduction. Using digoxin under the control of concentration enables detection of dangerous side effects before clinical symptoms appear, also it allows to avoid underdosage - "sham therapy".

METHOD PRINCIPLE

Immunoturbidimetric method; inhibition of agglutination. Increase of absorbance measured at $\lambda = 700$ nm is inversely proportional to the concentration of digoxin in the sample. Digoxin which is present in the sample forms immune complexes with specific antibody contained in the first reagent. After adding a second reagent containing polystyrene latex particles coated with digoxin, agglutination is inhibited in proportion to digoxin concentration in the sample.

REAGENTS

Package

	Cat. No 4-254 (24-TRAY)	Cat. No 4-454 (36-TRAY)
1-Reagent	2 x 11.5 ml	2 x 10 ml
2-Reagent	2 x 7.5 ml	2 x 7 ml

Unopened reagents stored at 2-8°C are stable up to the expiry date printed on the package. After opening the reagents are stable for 21 days on board the analyser at 2-10°C. Protect from light!

Reagents composition

Bis-Tris buffer, monoclonal antibodies to digoxin, polystyrene latex particles coated with digoxin, sodium azide (< 0.1%).

Warnings and notes

- Products for in vitro diagnostic use only.
- The reagents should be used by suitably qualified laboratory personnel only in accordance with intended purpose.
- The reagents contain sodium azide as a preservative (< 0.1%). Avoid contact with skin and mucous membranes. Sodium azide can form high explosive metal azide combinations with lead and copper. Drains should be flushed well with a large amount of water when discarding the solution.
- Mix well Reagent 2 before first use. Avoid foam formation.
- For optimal reagents stability it is recommended to remove them from the system and store in tightly closed bottles at 2-8°C.

SPECIMEN

Serum. Samples may be stored up to 3 days at 2-8 °C. or longer at -20°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

Avoid repeated freezing and thawing. Mix well the samples before analysis.

PROCEDURE

These reagents may be used in automatic analysers Prestige 24i, Biolis 24i and Sapphire 400.

- 1-Reagent and 2-Reagent are ready to use.
 - 1-Reagent put on basic position in reagent tray.
 - 2-Reagent put on start position in reagent tray.
- For reagent blank 0.9% NaCl is recommended.

APPLICATION

Item name	DIGOX
-----------	-------

Data information

Units	ng/ml
Decimals	2

Analysis

Type	END
Main W.Length1	700
Sub W.Length2	
Method	Immunoturbidimetric

Corr

	Slope		Inter
Y=	1.000	X+	0.000

Calibration

Type	Linear
Factor	

Standard

#1	*	#4	*
#2	*	#5	*
#3	*	#6	*

Normal Range

	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Aspiration

Kind	Double	
	Volume	
Sample	10	μl
Reagent1	150	
Reagent2	100	

Third Mix.	OFF
R1 Blank	Water-B

Monitor

0 Level Point	1
Span	3.000

Data Process

Read	Absorbance Limit		
	Start	End	
Main	53	54	
Sub	34	35	
		Low	0.000
		High	3.000

Factor		Endpoint Limit	2.000
Blank correction:	1.0000	Linear Check (%)	90%

Dilution

Diluent	100 : Dil2
---------	-------------------

Prozone Check

	Start	End	Limit (%)
First			
Second			Low
Third			Low

Auto Rerun SW

ON

Auto Rerun Range (Result)

	OFF	OFF
	Lower	Higher
Serum		
Urine		
Plasma		
CSF		
Dialysis		
Other		

Auto Rerun Condition (Absorbance)

Absorbance Range		
	Lower	OFF
	Higher	OFF
Prozone Range		OFF

THERAPEUTIC RANGE

therapeutic concentration 0.8 – 2.0 ng/ml (1 – 2.6 nmol/l)

Some patients achieve the desired therapeutic response at levels outside this range, so it is recommended to consider the need to establish an individual therapeutic ranges for each patient.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY TDM CONTROLS (Cat. No 5-107) with each batch of samples. For the calibration of automatic analysers systems the CORMAY DIGOXIN CALIBRATORS (Cat. No 5-113) is recommended.

The calibration curve should be prepared every 7 days, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

Metrological characteristics may differ from values presented below with the instrument used.

- **Sensitivity / Limit of Detection:** 0.14 ng/ml (0.18 nmol/l).
- **Linearity:** up to 5.5 ng/ml (7.08 nmol/l). If the digoxin concentration exceeds 5.5 ng/ml, dilute the sample 1:5 with saline solution and repeat the assay. The dilution take into account when making the results.
- **Specificity / Interferences**
Haemoglobin up to 3 g/dl, bilirubin up to 60 mg/dl, total protein up to 10 g/dl, triglyceride up to 1.9 g/dl do not interfere with the test.

▪ Precision

Reproducibility (day to day) n = 80	Mean [ng/ml]	SD [ng/ml]	CV [%]
level 1	1.0	0.07	7.32
level 2	1.89	0.09	4.82
level 3	2.92	0.08	2.86

▪ Method comparison

A comparison between CORMAY reagent (y) and commercially available assay (x) using 50 samples gave following results:

$$y = 0.949x - 0.081 \text{ ng/ml};$$

$$R = 0.9893$$

(R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

1. Oellerich, M. Therapeutic drug monitoring. In: Thomas L, ed. Clinical Laboratory Diagnostics. Use and Assessment of Clinical Laboratory Results. 1st Edition. TH-Books, Frankfurt/Main, Germany, 1998.
2. Ruiz R, Borque L, Soria AG, Córdova MA, Asolo B. Evaluation of an immunoturbidimetric assay of serum digoxin without sample pretreatment. Eur J Clin Chem Clin Biochem 33, 171-175, 1995.
3. Scholer A, Boecker J, Engelmayer U, Feldmann K, Hannak D, Kattermann R, Oellerich M, Raith H, Schlebusch H, Wieland H, Willems D, Jarausch J, Domke I. Comparability of a new turbidimetric digoxin test with other immunochemical tests and with HPLC – a multicenter evaluation. Clin Chem 43, 92- 99, 1997.
4. Biosafety in Microbiological and Biomedical Laboratories, Richmond JY, McKinney RW, eds. US Department of Health and Human Services, 4th Edition, 1999.
5. Westgard JO, Barry PL. Cost-Effective Quality Control: Managing the Quality and Productivity of Analytical Processes, AACC Press, 1986.
6. Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th Edition, AACC Press, 2000.
7. Tietz NW. Clinical Guide to Laboratory Tests. WB Saunders Company, Philadelphia, 1990.
8. Alan H. B. Wu, Tietz Clinical Guide to Laboratory Tests, W.B. Saunders Company, 4th edition, 1456 (2006).

Date of issue: 09. 2012.

MANUFACTURER

PZ CORMAY S.A.
ul. Wiosenna 22,
05-092 Łomianki, POLAND
tel.: +48 (0) 81 749 44 00
fax: +48 (0) 81 749 44 34
<http://www.pzcormay.pl>

09/12/09/12