PRESTIGE 24i CRP ULTRA

Nr kat. 4-261, 4-480 (PL)

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Cat. No **4-261, 4-480** (EN)

INTENDED USE

Diagnostic kit for determination of C-reactive protein concentration used in automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

CRP (C-reactive protein) is an acute phase protein whose concentration is seen to increase as a result of the inflammatory process, most notably in response to pneumococcal (bacterial) infectious, histolytic disease and a variety of disease states. CRP to be used as a marker or general diagnostic indicator of infections and inflammation, in addition to serving as a monitor of patient response to therapy and surgery. Furthermore, regular measurements of CRP in infants can be a useful aid in the early diagnosis of infectious disease.

METHOD PRINCIPLE

When an antigen-antibody reaction occurs between CRP in a sample and anti-CRP antibody which has been sensitized to latex particles, agglutination results. This agglutination is detected as an absorbance change (572 nm), with the magnitude of the change being proportional to the quantity of CRP in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

REAGENTS

Package

	Cat. No 4-261 (24-TRAY)	Cat. No 4-480 (36-TRAY)
1-REAGENT	2 x 20 ml	3 x 13 ml
2-REAGENT	2 x 20 ml	3 x 13 ml

The reagents when stored at 2-10°C are stable up to expiry date printed on the package. Stability on board of the analyser at 2-10°C: Prestige 24i – 12 weeks, Biolis 24i Premium – 12 weeks.

Concentrations in the test

suspension of latex particles sensitized with anti-CRP antibodies (rabbit) (pH 7.3) 0.20 w/v%
 glycine buffer solution (pH 7.0)
 preservative

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- After measurements are taken, reagent bottles should be capped and kept at 2-10°C. Care should be taken not to interchange the caps of reagent bottles.
- Reagents with different lot numbers should not be interchanged or mixed.
- EUH210 Safety data sheet available on request.

SPECIMEN

Serum or plasma (Na-EDTA, K-EDTA, Na-Heparin, Li-Heparin, citrate).

If the test cannot be done immediately, the sample should be placed in a tightly sealable container and stored at -20°C. Repeated freezing and thawing should be avoided. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

- 1-REAGENT and 2-REAGENT are ready to use.
- 1-REAGENT put on basic position in reagent tray.
- 2-REAGENT put on start position in reagent tray.

For reagent blank 0.9% NaCl is recommended.

Actions required:

- Biolis 30i:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: LACTATE – CRP ULTRA. To avoid this effect follow the recommendations contained in the instruction 51_03_24_009_BIOLIS 30i_CARRYOVER.

REFERENCE VALUES ^{3,4}

serum, plasma	
adults	< 0.5 mg/dl (< 5 mg/l)
children (2 months – 15 years)	0.01 – 0.28 mg/dl (0.1 – 2.8 mg/l)
newborns (0 – 3 weeks)	0.01 – 0.41 mg/dl (0.1 – 4.1 mg/l)

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY IMMUNO-CONTROL I (Cat. No 4-288).

For the calibration of automatic analysers systems the CORMAY CRP ULTRA CALIBRATORS kit (Cat. No 4-276) is recommended. The calibration curve should be prepared every 12 weeks (Prestige 24i, Biolis 24i Premium), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Biolis 30i. Results may vary if a different instrument is used.

- LoB (Limit of Blank):**
0.02 mg/dl (0.2 mg/l)
- LoD (Limit of Detection):**
0.03 mg/dl (0.3 mg/l)
- LoQ (Limit of Quantitation):**
0.09 mg/dl (0.9 mg/l)
- Linearity:**
up to 35 mg/dl (350 mg/l)

For higher concentrations dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 0.5 g/dl, bilirubin up to 30 mg/dl, triglycerides up to 500 mg/dl and RF up to 500 IU/ml do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	0.96	0.02	1.64
level 2	2.79	0.05	1.80
Reproducibility (day to day) n = 80	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	0.9	0.03	2.8
level 2	2.5	0.05	2.0

Method comparison

A comparison between CRP values determined at **Biolis 30i** (y) and at **ADVIA SIEMENS 1800** (x) using 53 serum samples gave following results:

$$y = 0.994x - 0.1575 \text{ mg/dl};$$

$$R = 0.996 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	42	CRPU
Auto Rerun SW		
OFF		
Auto Rerun Condition (Absorbance)		
Absorbance Range	Lower Higher	OFF OFF
Auto Rerun Range (Result)		
	OFF Lower Higher	OFF Higher
Serum		
Urine		
Plasma		
CSF		
Dialysis		
Other		
Prozone Range		
OFF		

Prestige 24i Premium, Biolis 24i Premium

Item No.	42	Item Name	CRPU	Optical
Data information				
Units	mg/dl			
Decimals	2			
Calibration				
Type	Spline1			
Std sample conc.				
Blank	0	#1 *	#2 *	*
#3	*	#4 *	#5 *	*
#6				
Analysis				
Type	END method			
Main Wave Length	570nm			
Sub Wave Length	800nm			
Method	Immuno			
Correlation				
Y=	Slope 1	X+	Intercept 0	

Item No.	42	Item Name	CRPU	Optical
Aspiration				
Kind	Double			
Vol.				
	Kind	Vol.	Add	Units
Sample		3.2	5	µl
Reagent 1		160	10	µl
Reagent 2		160	10	µl
Data Process				
Read	Main Sub	Start End	44 45 36 37	
Abs.Limit Low High				
	-3	~	3	
Correction value				
Blank correction	1			
End Point Limit	2			
Linear Check (%)	0			
Prozone Check				
	Start End Limit (%)			
First				
Second			Low	
Reaction Monitor				
0 Level Point	1			
Span	3			
Third mixing				
ON				

Item No.	42	Item Name	CRPU	Optical
Normal Range				
	Male	Female		
	Low High	Low High	Low High	
Serum	0.00 0.50	0.00 0.50	0.00 0.50	
Urine				
Plasma	0.00 0.50	0.00 0.50	0.00 0.50	
CSF				
Dialysis				
Other				
Panic Range				
	Male	Female		
	Low High	Low High	Low High	
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	42	Item Name	CRPU	Optical
Auto Rerun SW				
OFF				
Auto Rerun Condition (Absorbance)				
	Lower Higher	OFF OFF		
Auto Rerun Range (Conc.)				
	First Dil	Low Re Value Dil	High Re Value Dil	
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				
Auto Rerun Condition (Prozone)				
OFF				
Dilution				
99:Di11				

PRESTIGE 24i CRP ULTRA

• **Biolis 30i**

Item no	42	Item name	CRP	Specimen	SERUM / PLASMA	OPTICAL	
Data information				Aspiration volume			
UNITS	mg/dL			TYPE	Double		
DECIMALS	2				SAMPLE	REAGENT 1	REAGENT 2
Analysis		END method		VOL. (µL)	3.2	160	160
METHOD				BOTTLE (ml)			
Main Wave Length	570 nm			FIRST DIL.			
Sub Wave Length	800 nm			Data processing read			
CORRELATION (Y= AX + B)					START	END	
A =	1			MAIN	44	45	
B =	0			SUB	35	36	
Blank value				ABS LIMIT			
	• WATER	° REAGENT			-3	TO	3
Calibration				Collection value			
TYPE	Spline			END POINT	2.5		
STABILITY				LINEARITY CHECK (%)	0		
				Prozone check			
		START	END	LIMIT (%)			
FIRST							
SECOND							
°HIGH				MINIMUM ABS.			
•LOW				MEAN			
				VARIATE			
Item No	42	Item Name	CRP	Specimen	SERUM / PLASMA	OPTICAL	
Reference intervals				Auto rerun			
MALE		FEMALE		•ON		°OFF	
LOW	HIGH	LOW	HIGH				
0	0.5	0	0.5				
Panic range				Auto rerun range (conc.)			
MALE		FEMALE		Re	Value	Dil.	Re
LOW	HIGH	LOW	HIGH		0.09		Value
							35
Decision limit				Auto rerun condition (abs.)			
MALE		FEMALE					DIL.
				LOWER	°ON	•OFF	
				HIGH	°ON	•OFF	
				Auto rerun condition (prozone)			
				°ON		•OFF	
				SAMPLE VOL.			
				Dilution			
				•DIL 1		° DIL 2	
Reaction check							
		°ON		•OFF			
CHECK							
LOW							
HIGH							
VL CHECK		VH CHECK					
°ON		•OFF		°ON		•OFF	

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