## PRESTIGE 24i FERRITIN

Nr kat. 4-264, 4-483 (PL)

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Cat. No **4-264, 4-483** (EN)

### INTENDED USE

Diagnostic kit for determination of ferritin concentration intended to use in automatic analyzers: Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Ferritin is an iron-containing protein with a molecular weight of approximately 450 kDa. It is found mainly in the human liver and spleen, where its function is to eliminate and store iron in the body, and is also found in small amounts in human serum. This amount varies according to the movement of iron in the body, and hepatitis and malignant tumors, may be seen to increase due to cell destruction or tumor cell production, independent of iron reserves. Consequently, the measurement of ferritin is considered to be useful in the diagnosis, treatment, assessment of disease progression, and postoperative prognosis for such disease conditions.

### METHOD PRINCIPLE

When an antigen-antibody reaction occurs between ferritin in a sample and anti-ferritin antibody which has been sensitized to latex particles, agglutination results. This agglutination is detected as an absorbance change (572 nm), with the magnitude of the change being proportional to the quantity of ferritin in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

### REAGENTS

#### Package

	Cat. No 4-264 (24-TRAY)	Cat. No 4-483 (36-TRAY)
1-Reagent	1 x 40 ml	2 x 23 ml
2-Reagent	1 x 19 ml	2 x 11.5 ml

The reagents when stored at 2-10°C are stable up to expiry date printed on the package. The reagents are stable for 11 weeks on board the analyser at 2-10°C.

### Concentrations in the test

suspension of latex particles sensitized with anti-ferritin (rabbit) antibodies (pH 7.3) 0.07 w/v%  
glycine buffer solution (pH 8.3)  
preservative

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- After measurements are taken, reagent bottles should be capped and kept at 2-10°C. Care should be taken not to interchange the caps of reagent bottles.
- Reagents with different lot numbers should not be interchanged or mixed.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

### SPECIMEN

Serum.

If the test cannot be done immediately, the sample should be placed in a tightly sealable container and stored at -20°C. Repeated freezing and thawing should be avoided. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank 0.9% NaCl is recommended.

#### Actions required:

- Biolis 24i Premium:** When performing assays in the analyzer there is a probability of **cross-contamination** affecting the tests results. To avoid this effect, the tests for determination of ferritin concentration should be performed **in the separate order** (follow the recommendations contained in the instruction 51\_03\_24\_008\_BIOLIS 24i PREMIUM\_CARRYOVER).
- Biolis 30i:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: FERRUM - FERRITIN. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_009\_BIOLIS 30i\_CARRYOVER.

### REFERENCE VALUES <sup>6</sup>

serum	ng/ml
male	20 – 250
female	10 – 120

It is recommended for each laboratory to establish its own reference ranges for local population. Diagnosis should only be made after taking clinical symptoms and the results of other tests into consideration.

### QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY IMMUNO-CONTROL II (Cat. No 4-290).

For the calibration of automatic analyzers systems the CORMAY FERRITIN CALBRATORS kit (Cat. No 4-491) is recommended.

The calibration curve should be prepared every 11 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results been obtained using an automatic analyser Biolis 30i. Results may vary if a different instrument is used.

- LoB (Limit of Blank):** 1 ng/ml
- LoD (Limit of Detection):** 3 ng/ml
- LoQ (Limit of Quantitation):** 8 ng/ml
- Linearity:** up to 750 ng/ml

For higher concentrations dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

#### Specificity / Interferences

Haemoglobin up to 0.98 g/dl, bilirubin up to 62 mg/dl, RF up to 520 IU/ml, intralipid up to 800 mg/dl do not interfere with the test.

#### Precision

Repeatability (run to run) n = 20	Mean [ng/ml]	SD [ng/ml]	CV [%]
level 1	104	1.25	1.20
level 2	390	6.99	1.79
Reproducibility (day to day) n = 80	Mean [ng/ml]	SD [ng/ml]	CV [%]
level 1	104	3.42	3.3
level 2	405.8	9.18	2.3

#### Method comparison

A comparison between ferritin values determined at **Biolis 30i** (y) and at **ADVIA SIEMENS 1800** (x) using 60 serum samples gave following results:

$$y = 0.8507x + 3.3232 \text{ ng/ml};$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	45	FERRI
<b>Auto Rerun SW</b>		
OFF		
<b>Auto Rerun Range (Result)</b>		
	OFF	OFF
	Lower	Higher
Serum		
Urine		
Plasma		
CSF		
Dialysis		
Other		
<b>Auto Rerun Condition (Absorbance)</b>		
Absorbance Range	Lower	OFF
	Higher	OFF
<b>Prozone Range</b>		
OFF		

#### • Prestige 24i Premium, Biolis 24i Premium

Item No.	45	Item Name	FERRI	Optical
<b>Data information</b>				
Units	ng/ml			
Decimals	1			
<b>Analysis</b>				
Type	RATE method			
Main Wave Length	570 nm			
Sub Wave Length				
Method	Immuno			
<b>Correlation</b>				
Y=	Slope	X+	Intercept	
	1	0		
<b>Calibration</b>				
Type	Logit 2			
Std sample conc.				
Blank	0	#1 *	#2 *	
#3	*	#4 *	#5 *	
#6				

Item No.	45	Item Name	FERRI	Optical
<b>Aspiration</b>				
Kind	Double			
<b>Data Process</b>				
Read	Main	Start	End	
	Sub	35	52	
<b>Vol.</b>				
	Kind	Vol.	Add	Units
	Sample	8	5	µl
	Reagent 1	180	10	µl
	Reagent 2	65	10	µl
<b>Blank value</b>				
Water Blank				
<b>Reaction Monitor</b>				
0 Level Point	1			
Span	3			
<b>Third mixing</b>				
ON				
<b>Correction value</b>				
Blank correction				1.00
End Point Limit				2
Linear Check (%)				90
<b>Prozone Check</b>				
	Start	End	Limit (%)	
First				
Second				Low

Item No.	45	Item Name	FERRI	Optical
<b>Normal Range</b>				
	Male		Female	
	Low	High	Low	High
Serum	20	250	10	120
Urine				
Plasma				
CSF				
Dialysis				
Other				
<b>Panic Range</b>				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	45	Item Name	FERRI	Optical
<b>Auto Rerun SW</b>				
ON				
<b>Auto Rerun Range (Conc.)</b>				
	First Dil	Low	High	
		Re	Value	Dil
Serum			3.6	980
Urine				
Plasma				
CSF				
Dialysis				
Other				
<b>Auto Rerun Condition (Absorbance)</b>				
Lower				OFF
Higher				OFF
<b>Auto Rerun Condition (Prozone)</b>				
OFF				
<b>Dilution</b>				
99:Dil 1				

PRESTIGE 24i FERRITIN

• **Biolis 30i**

Item no	45	Item name	FERRITIN	Specimen	SERUM	OPTICAL												
<b>Data information</b>																		
UNITS	ng/mL		<b>Aspiration volume</b>															
DECIMALS	1		TYPE Double															
<b>Analysis</b>			<table border="1"> <tr> <td></td> <td>SAMPLE</td> <td>REAGENT 1</td> <td>REAGENT 2</td> </tr> <tr> <td>VOL. (µL)</td> <td>8</td> <td>180</td> <td>65</td> </tr> <tr> <td>BOTTLE (ml)</td> <td></td> <td></td> <td></td> </tr> </table>					SAMPLE	REAGENT 1	REAGENT 2	VOL. (µL)	8	180	65	BOTTLE (ml)			
	SAMPLE	REAGENT 1	REAGENT 2															
VOL. (µL)	8	180	65															
BOTTLE (ml)																		
METHOD	RATE method		FIRST DIL.															
Main Wave Length	570 nm		<b>Data processing read</b>															
Sub Wave Length			<table border="1"> <tr> <td></td> <td>START</td> <td>END</td> </tr> <tr> <td>MAIN</td> <td>35</td> <td>52</td> </tr> <tr> <td>SUB</td> <td></td> <td></td> </tr> </table>					START	END	MAIN	35	52	SUB					
	START	END																
MAIN	35	52																
SUB																		
<b>CORRELATION (Y= AX + B)</b>			ABS LIMIT															
A =	1		-3 TO 3															
B =	0		<b>Collection value</b>															
<b>Blank value</b>			<table border="1"> <tr> <td>END POINT</td> <td>2.5</td> </tr> <tr> <td>LINEARITY CHECK (%)</td> <td>90</td> </tr> </table>				END POINT	2.5	LINEARITY CHECK (%)	90								
END POINT	2.5																	
LINEARITY CHECK (%)	90																	
• WATER ° REAGENT			<b>Prozone check</b>															
<b>Calibration</b>			<table border="1"> <tr> <td></td> <td>START</td> <td>END</td> <td>LIMIT (%)</td> </tr> <tr> <td>TYPE</td> <td></td> <td></td> <td></td> </tr> <tr> <td>STABILITY</td> <td>Logit 2</td> <td></td> <td></td> </tr> </table>					START	END	LIMIT (%)	TYPE				STABILITY	Logit 2		
	START	END	LIMIT (%)															
TYPE																		
STABILITY	Logit 2																	
			<table border="1"> <tr> <td></td> <td>MINIMUM ABS.</td> </tr> <tr> <td>°HIGH</td> <td>MEAN</td> </tr> <tr> <td>•LOW</td> <td>VARIATE</td> </tr> </table>					MINIMUM ABS.	°HIGH	MEAN	•LOW	VARIATE						
	MINIMUM ABS.																	
°HIGH	MEAN																	
•LOW	VARIATE																	

Item No	45	Item Name	FERRITIN	Specimen	SERUM	OPTICAL
<b>Reference intervals</b>						
MALE		FEMALE		<b>Auto rerun</b>		
LOW	HIGH	LOW	HIGH	•ON °OFF		
20	250	10	120	<b>Auto rerun range (conc.)</b>		
<b>Panic range</b>				Re	Value	Dil.
MALE		FEMALE		8	750	
LOW	HIGH	LOW	HIGH	<b>Auto rerun condition (abs.)</b>		
					DIL.	
<b>Decision limit</b>				LOWER	°ON	•OFF
MALE		FEMALE		HIGH	°ON	•OFF
				<b>Auto rerun condition (prozone)</b>		
				°ON	•OFF	
				SAMPLE VOL.		
				<b>Dilution</b>		
				•DIL 1	° DIL 2	
<b>Reaction check</b>						
		°ON		•OFF		
CHECK						
LOW						
HIGH						
<b>VL CHECK</b>			<b>VH CHECK</b>			
°ON	•OFF	°ON	•OFF			

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