



## PRESTIGE 24i ALPHA-FETOPROTEIN

### DIAGNOSTIC KIT FOR DETERMINATION OF $\alpha$ - FETOPROTEIN CONCENTRATION

#### INTRODUCTION

$\alpha$ -fetoprotein (AFP) is fetoprotein with a molecular weight of approximately 70 kD containing about 3 % sugar. While it is present in high concentrations during fetal growth, its concentration rapidly decrease after birth and is present at extremely low levels in normal human blood.

AFP shows a notable increase in primary hepatic cancer and is considered to be of great diagnostic importance. It is also thought that fluctuations in blood AFP are useful for evaluating the progress, effects of therapy, and postoperative prognosis of hepatoma.

#### METHOD PRINCIPLE

When an antigen-antibody reaction occurs between AFP in a sample and anti-AFP antibody which has been sensitized to latex particles, agglutination results. This agglutination is detected as an absorbance change, with the magnitude of the change being proportional to the quantity of AFP in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

#### REAGENTS

##### Package

	Cat. No 4-266 (24-TRAY)	Cat. No 4-485 (36-TRAY)
1-Reagent	1 x 38 ml	2 x 23 ml
2-Reagent	1 x 20 ml	2 x 12.5 ml

The reagents when stored at 2-10°C are stable up to expiry date printed on the package. Protect from light and avoid contamination!

#### Concentrations in the test

suspension of latex particles sensitized with (rabbit) anti-AFP antibodies (pH 7.3) 0.12 w/v%  
glycine buffer solution (pH 8.3)

#### Warnings and notes

- Product for in vitro diagnostic use only.
- Reagent bottles should be shaken before use by gently inverting several times.
- After measurements are taken, reagent bottles should be capped and kept at 2-10°C. Care should be taken not to interchange the caps of reagent bottles.
- Reagents with different lot numbers should not be interchanged or mixed.
- The reagents contain sodium azide (< 0.1%) as a preservative. Avoid contact with skin and mucous membranes.

#### SPECIMEN

Serum.

After blood has clotted thoroughly, the sample is centrifuged and the serum is separated from blood cells and fibrins. Samples can be stored for several weeks at 2-8°C or for 1 year at -20°C. Repeated freezing and thawing should be avoided.

Nevertheless it is recommended to perform the assay with freshly collected samples!

#### PROCEDURE

These reagents may be used in automatic analysers Prestige 24i, Biolis 24i, Sapphire 400 and Prestige 24i Premium, Biolis 24i Premium, Sapphire 400 Premium.

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank 0.9% NaCl is recommended.

#### REFERENCE VALUES <sup>5</sup>

serum	< 15 ng/ml
-------	------------

It is recommended for each laboratory to establish its own reference ranges for local population. Diagnosis should only be made after taking clinical symptoms and the results of other tests into consideration.

#### QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY IMMUNO-CONTROL I (Cat. No 4-288) with each batch of samples. For the calibration of automatic analysers systems the CORMAY AFP CALIBRATORS kit (Cat. No 4-282) is recommended.

The calibration curve should be prepared with change of reagent lot number or as required e.g. quality control findings outside the specified range.

#### PERFORMANCE CHARACTERISTICS

These metrological characteristics have been obtained using an automatic analyser HITACHI 917. Results may vary if a different instrument or manual procedure is used.

- Analytical range:** 7 – 250 ng/ml.  
For higher concentrations dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.
- Specificity / Interferences**  
Haemoglobin up to 0.3 g/dl, bilirubin up to 30 mg/dl and triglycerides up to 300 mg/dl do not interfere with the test.
- Precision**

Repeatability (run to run) n = 20	Mean [ng/ml]	SD [ng/ml]	CV [%]
level 1	9.9	0.4	4.03
level 2	26.6	0.3	1.37
level 3	96.5	0.7	0.71

- Method comparison**

A comparison between CORMAY reagent (y) and commercially available assay (x) using 78 samples gave following results:

$$y = 1.01 x + 16.73 \text{ ng/ml};$$

$$R = 0.996 \quad (R - \text{correlation coefficient})$$

#### WASTE MANAGEMENT

Please refer to local legal requirements.

#### LITERATURE

- Bergstrand C. G. et al.: Demonstration of a new protein fraction in serum from the human fetus., Scand. J. Clin. Lab. Invest., 8, 174 (1956).

- Galvin J. P. et al.: Particle enhanced photometric immunoassay systems., Clin. Lab. Assays (Pap. Annu. Clin. Lab. Assays Conf.), 4<sup>th</sup>, 73 (1983).
- Singer J. M. et al.: The latex fixation test. I. Application to the serologic diagnosis of rheumatoid arthritis, Amer. J. Med., 21, 888 (1956).
- Pesce A. J., Kaplan L.A.: Methods in Clinical Chemistry, St. Louis, Mosby, 459-465 (1987).
- Burtis CA, Ashwood ER, Bruns DE, editors. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. 4<sup>th</sup> ed, St. Louis: W. B Saunders Company; 2006, 2269.

**APPLICATION for Prestige 24i, Biolis 24i and Sapphire 400**

Item name	47	AFP		
<b>Data information</b>				
Units	ng/ml			
Decimals	1			
<b>Analysis</b>				
Type	RATE			
Main W.Length1	570			
Sub W.Length2	800			
Method	Immuno			
<b>Corr</b>				
Y=	Slope	Inter		
	1.000	0.000		
<b>Calibration</b>				
Type	Spline			
Standard				
#1	*	#4 *		
#2	*	#5		
#3	*	#6		
<b>Normal Range</b>				
	Male		Female	
	Low	High	Low	High
Serum	0	15	0	15
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item name	47	AFP
<b>Aspiration</b>		
Kind	Double	
<b>Data Process</b>		
Read	Start	End
Main	36	42
Sub		
<b>Absorbance Limit</b>		
Low	-3.000	
High	3.000	
Sample	Volume	µl
Reagent1	200	
Reagent2	100	
<b>Factor</b>		
Blank correction	1.0000	Endpoint Limit
		2.000
<b>Dilution</b>		
Diluent	99:Dil 1	
<b>Monitor</b>		
0 Level Point	1	
Span	3.000	
<b>Prozone Check</b>		
First	Start	End
		Limit (%)
Second		
		Low
Third		
		Low

Item name	47	AFP
<b>Auto Rerun SW</b>		
OFF		
<b>Auto Rerun Range (Result)</b>		
	OFF	OFF
	Lower	Higher
Serum		
Urine		
Plasma		
CSF		
Dialysis		
Other		
<b>Auto Rerun Condition (Absorbance)</b>		
Absorbance Range		
	Lower	OFF
	Higher	OFF
<b>Prozone Range</b>		
OFF		

**APPLICATION for Prestige 24i Premium, Biolis 24i Premium and Sapphire 400 Premium**

Item No.	47	Item Name	AFP	Optical
<b>Data information</b>				
Units	ng/ml			
Decimals	1			
<b>Analysis</b>				
Type	RATE			
Main Wave Length	570 nm			
Sub Wave Length	800 nm			
Method	Immuno			
<b>Calibration</b>				
Type	Spline 1			
Std sample conc.				
Blank	0	#1	*	#2 *
#3	*	#4	*	#5
#6				
<b>Correlation</b>				
Slope		Intercept		
Y=	1	X+	0	

Item No.	47	Item Name	AFP	Optical
<b>Aspiration</b>				
Kind	Double			
<b>Data Process</b>				
Read	Start	End		
Main	36	48		
Sub				
<b>Abs.Limit</b>				
Low	-3	High		
		3		
<b>Blank value</b>				
Water Blank				
<b>Reaction Monitor</b>				
0 Level Point	1			
Span	3			
<b>Third mixing</b>				
ON				
<b>Correction value</b>				
Blank correction		1		
End Point Limit		2		
Linear Check (%)		0		
<b>Prozone Check</b>				
	Start	End	Limit (%)	
First				
Second			Low	
Third			Low	

Item No.	47	Item Name	AFP	Optical
<b>Normal Range</b>				
	Male		Female	
	Low	High	Low	High
Serum	0	15	0	15
Urine				
Plasma				
CSF				
Dialysis				
Other				
<b>Panic Range</b>				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	47	Item Name	AFP	Optical	
<b>Auto Rerun SW</b>					
OFF					
<b>Auto Rerun Range (Conc.)</b>					
	First Dil	Low		High	
	Re	Value	Dil	Re	Value
Serum					
Urine					
Plasma					
CSF					
Dialysis					
Other					
<b>Auto Rerun Condition (Absorbance)</b>					
Lower		OFF			
Higher		OFF			
<b>Auto Rerun Condition (Prozone)</b>					
OFF					
<b>Dilution</b>					
99:Dil 1					

Date of issue: 01. 2013.

**MANUFACTURER**

**PZ CORMAY S.A.**  
 22 Wiosenna Street,  
 05-092 Łomianki, POLAND  
 tel.: +48 (0) 22 751 79 10  
 fax: +48 (0) 22 751 79 14  
<http://www.cormay.pl>

01/13/01/13