PRESTIGE 24i IgM

Nr kat. 4-273, 4-332 (PL)

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Cat. No 4-273, 4-332 (EN)

INTENDED USE

Diagnostic kit for determination of IgM concentration used in automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Immunoglobulins (Igs) are the instrumental proteins of immunity. Immunity is a property of the lymphoid system which is made of organs (spleen, thymus, bone marrow) and of cells (lymphocytes). Circulating immunoglobulins are secreted in the blood by B lymphocytes and they thereby export far-away the specific biological functions of humoral immunity. Immunoglobulin M (IgM) is the first Ig to appear in response to an antigenic stimulus such as an infectious agent. In many cases, the antigen-specific IgM level subsequently falls and remains low as the IgG response appears.

METHOD PRINCIPLE

The IgM present in a sample form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum measured at $\lambda=340$ nm is proportional to IgM concentration in the sample.

REAGENTS

Package

| | Cat. No 4-273 (24-TRAY) | Cat. No 4-332 (36-TRAY) |
|-----------|----------------------------|----------------------------|
| 1-Reagent | 1 x 40 ml | 2 x 23 ml |
| 2-Reagent | 1 x 10 ml | 2 x 6 ml |

Buffer (1-Reagent) stored at 2-25°C and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package.

Reagent components

1-Reagent: tricine buffer, polymer, inorganic salt, preservative.

2-Reagent: anti-human IgM antiserum, buffer, inorganic salt, preservative.

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Do not use after expiry date.
- Do not interchange caps.
- Products from human source have been tested for HBsAg and antibodies to HIV and HCV and found to be non-reactive. However this material should be handled as thought capable of transmitting infectious disease.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum. Nonhemolyzed, nonlipaemic, fresh serum is recommended.

Serum can be stored up to 3 days at 2-8°C or up to 6 months at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank 0.9% NaCl is recommended.

REFERENCE VALUES

| | |
|--------------------------------|-----------------|
| adults | 0.50 – 3.00 g/l |
| children (1 year – 12 years) | 0.45 – 2.50 g/l |
| children (1 month – 12 months) | 0.20 – 1.50 g/l |

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY IMMUNO-CONTROL III (Cat. No 4-291)

For the calibration of automatic analysers systems the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended.

The calibration curve should be prepared with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using an automatic analysers Hitachi or Biolis 30i. Results may vary if a different instrument is used.

- Measurement range:** 0.038 g/l to 7 g/l.

- Specificity / Interferences:**

Hemoglobin up to 0.32 g/dl, bilirubin up to 22 mg/dl, triglycerides up to 312 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4 g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

- Precision**

| Repeatability (run to run) n = 30 | Mean [g/l] | SD [g/l] | CV [%] |
|--|---------------|-------------|-----------|
| level 1 | 0.60 | 0.004 | 0.73 |
| level 2 | 1.25 | 0.012 | 0.87 |
| level 3 | 18.92 | 0.014 | 0.72 |
| Reproducibility (day to day) n = 60 | Mean [g/l] | SD [g/l] | CV [%] |
| level 1 | 0.56 | 0.024 | 3.99 |
| level 2 | 1.30 | 0.037 | 2.81 |
| level 3 | 19.95 | 0.050 | 2.49 |

- Method comparison**

A comparison between IgM values determined at Biolis 30i (y) and at ADVIA SIEMENS 1800 (x) using 62 serum samples gave following results:

$$y = 1.0933x - 0.0387 \text{ g/l;}$$

$$R = 0.995 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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| | | |
|--|-------|------------|
| Item name | 30 | IgM |
| Auto Rerun SW | | |
| OFF | | |
| Auto Rerun Condition (Absorbance) | | |
| Absorbance Range | | Lower OFF |
| | | Higher OFF |
| Auto Rerun Range (Result) | | |
| | OFF | OFF |
| | Lower | Higher |
| Serum | | |
| Urine | | |
| Plasma | | |
| CSF | | |
| Dialysis | | |
| Other | | |
| Prozone Range | | |
| OFF | | |

• **Prestige 24i Premium, Biolis 24i Premium**

| | | | | |
|-------------------------|------------|-----------|------|---------|
| Item No. | 30 | Item Name | IgM | Optical |
| Data information | | | | |
| Units | g/l | | | |
| Decimals | 2 | | | |
| Analysis | | | | |
| Type | END method | | | |
| Main Wave Length | 340nm | | | |
| Sub Wave Length | 700nm | | | |
| Method | Immuno | | | |
| Calibration | | | | |
| Type | Logit2 | | | |
| Std sample conc. | | | | |
| Blank | 0 | #1 * | #2 * | #3 * |
| #3 | * | #4 * | #5 * | #6 * |
| #6 | | | | |
| Correlation | | | | |
| Slope | | Intercept | | |
| Y= | 1 | X+ | 0 | |

| | | | | |
|-------------------------|--------|-----------|-----------|---------|
| Item No. | 30 | Item Name | IgM | Optical |
| Aspiration | | | | |
| Kind | Double | | | |
| Data Process | | | | |
| Read | Main | Start | End | |
| | Sub | 29 | 31 | |
| Abs.Limit | | | | |
| Low | -3 | High | 3 | |
| Blank value | | | | |
| Water Blank | | | | |
| Reaction Monitor | | | | |
| 0 Level Point | 1 | | | |
| Span | 3 | | | |
| Correction value | | | | |
| Blank correction | | | | 0.8350 |
| End Point Limit | | | | 2 |
| Linear Check (%) | | | | 0 |
| Prozone Check | | | | |
| | Start | End | Limit (%) | |
| First | | | | |
| Second | | | Low | |

| | | | | |
|---------------------|------|-----------|--------|---------|
| Item No. | 30 | Item Name | IgM | Optical |
| Normal Range | | | | |
| | Male | | Female | |
| | Low | High | Low | High |
| Serum | 0.5 | 3.0 | 0.5 | 3.0 |
| Urine | | | | |
| Plasma | | | | |
| CSF | | | | |
| Dialysis | | | | |
| Other | | | | |
| Panic Range | | | | |
| | Male | | Female | |
| | Low | High | Low | High |
| Serum | | | | |
| Urine | | | | |
| Plasma | | | | |
| CSF | | | | |
| Dialysis | | | | |
| Other | | | | |

| | | | | | | | |
|--|-----------|-----------|-------|---------|------|-------|-----|
| Item No. | 30 | Item Name | IgM | Optical | | | |
| Auto Rerun SW | | | | | | | |
| OFF | | | | | | | |
| Auto Rerun Condition (Absorbance) | | | | | | | |
| Lower | | OFF | | | | | |
| Higher | | OFF | | | | | |
| Auto Rerun Range (Conc.) | | | | | | | |
| | First Dil | Low | | | High | | |
| | | Re | Value | Dil | Re | Value | Dil |
| Serum | | | | | | | |
| Urine | | | | | | | |
| Plasma | | | | | | | |
| CSF | | | | | | | |
| Dialysis | | | | | | | |
| Other | | | | | | | |
| Auto Rerun Condition (Prozone) | | | | | | | |
| OFF | | | | | | | |
| Dilution | | | | | | | |
| 99:Dil1 | | | | | | | |

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