PRESTIGE 24i APOLIPOPROTEIN B

Nr kat. 4-275, 4-334 (PL)

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LITERATURA

- Marcovina, S.M., Albers, International Federation of Clinical Chemistry Standardization Project for Measurements of Apolipoproteins A1 and B. III.
- Tietz, N.W. Fundamentals of Clinical Chemistry. Saunders, Philadelphia 1987.
- Alan H.B. Wu, ed.: Tietz Clinical Guide to Laboratory Tests, 4th ed. W.B. Saunders Company., 146, (2006).
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders, 1802, (1999).

Data wydania: 04. 2021.

PRESTIGE 24i APOLIPOPROTEIN B

Cat. No 4-275, 4-334 (EN)

INTENDED USE

Diagnostic kit for determination of apolipoprotein B concentration used in automatic analysers Prestige 24i, Biolis 24i and Prestige 24i Premium, Biolis 24i Premium.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Lipids are transported in serum under the form of micelle known as lipoproteins. Lipoproteins are macromolecular complexes containing proteins (**apolipoproteins**), cholesterol and phospholipids in their outer layer, triglycerides and cholesterol esters in their inner structure. Lipoproteins are classified according to their increasing density respectively as chylomicrons, very low density lipoproteins (VLDL), low density lipoproteins (LDL) and high density lipoproteins (HDL). Apolipoprotein B is the major protein moiety of LDL. Apolipoprotein measurements are more discriminating than HDL and LDL cholesterol measurements in allowing to identify patients with coronary heart diseases and in assessing atherosclerotic risk.

METHOD PRINCIPLE

The apolipoprotein B presents in a sample form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum measured at $\lambda=340$ nm is proportional to apolipoprotein B concentration in the sample.

REAGENTS

Package

	Cat. No 4-275 (24-TRAY)	Cat. No 4-334 (36-TRAY)
1-REAGENT	1 x 40 ml	2 x 23 ml
2-REAGENT	1 x 10 ml	2 x 6 ml

Buffer (1-REAGENT) stored at 2-8°C and antiserum (2-REAGENT) stored at 2-8°C are stable until expiry date printed on the package.

Reagent components

1-Reagent: TRIS buffer, polymer, inorganic salt, preservative.
2-Reagent: anti-human apolipoprotein B antiserum, buffer, inorganic salt, preservative.

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Store closed.
- Do not freeze the reagents.
- Nanoparticle-based reagents can settle over time. It may be necessary to delicately mix by repeated turning.

- Human-origins products have been tested for HBsAg and antibodies to HIV 1, HIV 2, HCV and found to be non-reactive. However this material should be handled as potentially infectious.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN³

Serum or plasma.

Serum and plasma can be stored 4-7 days at 4°C or at least 6 months at -70°C.

It is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-REAGENT and 2-REAGENT are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank 0.9% NaCl is recommended.

Application guide for Biolis 30i analyser is available on request.

REFERENCE VALUES³

children (4-11 years)	0.56 – 1.13 g/l
children (12-19 years)	0.55 – 1.19 g/l
adults	0.59 – 1.73 g/l

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples the CORMAY APOLIPOPROTEIN CONTROL (Cat. No 4-293).

For the calibration of automatic analysers systems the CORMAY APOLIPOPROTEIN CALIBRATORS (Cat. No 4-289) is recommended.

The calibration curve should be prepared with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Hitachi. Results may vary if a different instrument is used.

- Measurement range:** 0.019 g/l – 4 g/l.
- Specificity / Interferences**
 Hemoglobin up to 0.32 g/dl, bilirubin up to 29,5 mg/dl, triglycerides up to 264 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4 g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

Precision

Repeatability (run to run) n = 30	Mean	SD	CV
	[g/l]	[g/l]	[%]
level 1	0.42	0.006	1.53
level 2	0.88	0.011	1.26
level 3	1.35	0.014	1.04
Reproducibility (day to day) n = 66	Mean	SD	CV
	[g/l]	[g/l]	[%]
level 1	0.33	0.001	3.44
level 2	0.95	0.003	3.04
level 3	1.51	0.004	2.71

Method comparison

A comparison between CORMAY reagent (y) and commercially available assay (x) using at least 30 samples gave following results:

$$y = 0.8998x + 0.0867 \text{ g/l};$$

$$R = 0.997 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

- Marcovina, S.M., Albers, International Federation of Clinical Chemistry Standardization Project for Measurements of Apolipoproteins A1 and B. III.
- Tietz, N.W. Fundamentals of Clinical Chemistry. Saunders, Philadelphia 1987.
- Alan H.B. Wu, ed.: Tietz Clinical Guide to Laboratory Tests, 4th ed. W.B. Saunders Company., 146, (2006).
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders, 1802, (1999).

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1. Marcovina, S.M., Albers, International Federation of Clinical Chemistry Standardization Project for Measurements of Apolipoproteins A1 and B. III.
2. Tietz, N.W. Fundamentals of Clinical Chemistry. Saunders, Philadelphia 1987.
3. Alan H.B. Wu, ed.: Tietz Clinical Guide to Laboratory Tests, 4th ed. W.B. Saunders Company., 146, (2006).
4. Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders, 1802, (1999).

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Item name	31	Apo B
Auto Rerun SW		
OFF		
Auto Rerun Range (Result)		
	OFF	OFF
	Lower	Higher
Serum		
Urine		
Plasma		
CSF		
Dialysis		
Other		
Auto Rerun Condition (Absorbance)		
Absorbance Range	Lower	OFF
	Higher	OFF
Auto Rerun Condition (Prozone)		
Prozone Range	OFF	

• **Prestige 24i Premium, Biolis 24i Premium**

Item No.	31	Item Name	Apo B	Optical	
Data information					
Units	g/L				
Decimals	2				
Analysis					
Type	END method				
Main Wave Length	340 nm				
Sub Wave Length	700 nm				
Method	Immuno				
Calibration					
Type	Logit2				
Std sample conc.					
Blank	0	#1	*	#2	*
#3	*	#4	*	#5	*
#6					
Correlation					
Slope	Intercept				
Y=	1	X+	0		

Item No.	31	Item Name	Apo B	Optical
Aspiration				
Kind	Double			
Data Process				
Read	Main	Start	End	
	Sub	30	31	
Abs.Limit				
Low	-3	High	3	
Blank value				
Water Blank				
Reaction Monitor				
0 Level Point	1			
Span	3			
Third mixing				
ON				
Correction value				
Blank correction	0.8350			
End Point Limit	2			
Linear Check (%)	0			
Prozone Check				
	Start	End	Limit (%)	
First				
Second				Low

Item No.	31	Item Name	Apo B	Optical
Normal Range				
	Male		Female	
	Low	High	Low	High
Serum	0.59	1.73	0.59	1.73
Urine				
Plasma				
CSF				
Dialysis				
Other				
Panic Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	31	Item Name	Apo B	Optical
Auto Rerun SW				
OFF				
Auto Rerun Range (Conc.)				
	First Dil	Low	High	
		Re	Value	Dil
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				
Auto Rerun Condition (Absorbance)				
Lower	OFF			
Higher	OFF			
Auto Rerun Condition (Prozone)				
OFF				
Dilution				
99:Dil 1				

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