## PRESTIGE 24i IgG

Nr kat. 4-272, 4-331 (PL)

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Cat. No 4-272, 4-331 (EN)

### INTENDED USE

Diagnostic kit for determination of IgG concentration used in automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Immunoglobulins (Igs) are the instrumental proteins of immunity. Immunity is a property of the lymphoid system which is made of organs (spleen, thymus, bone marrow) and of cells (lymphocytes). Circulating immunoglobulins are secreted in the blood by B lymphocytes and they thereby export far-away the specific biological functions of humoral immunity.

Immunoglobulin G (IgG) - the predominant serum immunoglobulin (75 % of Igs)- is of particular importance in the body's long-term defense against infection; IgG deficiency is associated with recurrent and occasionally severe pyogenic infections. IgG serum levels are increased in response to chronic or recurrent infections or autoimmune diseases.

### METHOD PRINCIPLE

The IgG present in a sample form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum is proportional to IgG concentration in the sample.

### REAGENTS

#### Package

	Cat. No 4-272 (24-TRAY)	Cat. No 4-331 (36-TRAY)
1-Reagent	1 x 40 ml	2 x 23 ml
2-Reagent	1 x 16.5 ml	2 x 9.5 ml

Buffer (1-Reagent) stored at 2-25°C and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package.

### Reagent components

**1-Reagent:** tricine buffer, polymer, inorganic salt, preservative.

**2-Reagent:** anti-human IgG antiserum, buffer, inorganic salt, preservative.

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Store closed.
- Do not freeze the reagents.
- Nanoparticle-based reagents can settle over time. It may be necessary to delicately mix by repeated turning.

- Human-origin products have been tested for HBsAg and antibodies to HIV 1, HIV 2, HCV and found to be non-reactive. However this material should be handled as potentially infectious.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

### SPECIMEN

Serum. Nonhemolyzed, nonlipaemic, fresh serum is recommended.

Serum can be stored up to 3 days at 2-8°C or up to 6 months at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank 0.9% NaCl is recommended.

### REFERENCE VALUES

adults	5.65 – 17.65 g/l
children (> 6 years)	6.50 – 16.00 g/l
children (1 month – 6 years)	2.00 – 12.40 g/l

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY IMMUNO-CONTROL III (Cat. No 4-291).

For the calibration of automatic analysers systems the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended.

The calibration curve should be prepared with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers Hitachi or Biolis 30i. Results may vary if a different instrument is used.

- Measurement range:** 0.125 g/l to 35 g/l.

#### Specificity / Interferences

Hemoglobin up to 0.32 g/dl, bilirubin up to 22 mg/dl, triglycerides up to 940 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

### Precision

Repeatability (run to run) n = 30	Mean [g/l]	SD [g/l]	CV [%]
level 1	6.83	0.01	1.43
level 2	14.03	0.16	1.12
level 3	19.86	0.36	1.80
Reproducibility (day to day) n = 66	Mean [g/l]	SD [g/l]	CV [%]
level 1	6.57	0.29	4.46
level 2	13.66	0.43	3.15
level 3	18.87	0.56	2.95

### Method comparison

A comparison between IgG values determined at **Biolis 30i** (y) and at **ADVIA SIEMENS 1800** (x) using 63 serum samples, gave following results:

$$y = 0.9221x + 0.7106 \text{ g/l};$$

$$R = 0.986 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	29	IgG
<b>Auto Rerun SW</b>		
	OFF	
<b>Auto Rerun Range (Result)</b>		
	OFF	OFF
	Lower	Higher
Serum		
Urine		
Plasma		
CSF		
Dialysis		
Other		
<b>Auto Rerun Condition (Absorbance)</b>		
Absorbance Range	Lower	OFF
	Higher	OFF
<b>Auto Rerun Condition (Prozone)</b>		
Prozone Range	OFF	

• **Prestige 24i Premium, Biolis 24i Premium**

Item No.	29	Item Name	IgG	Optical
<b>Data information</b>				
Units	g/l			
Decimals	2			
<b>Analysis</b>				
Type	END method			
Main Wave Length	600nm			
Sub Wave Length				
Method	Immuno			
<b>Calibration</b>				
Type	Logit2			
Std sample conc.				
Blank	0	#1 *	#2 *	#3 *
#3	*	#4 *	#5 *	#6 *
#6				
<b>Correlation</b>				
Slope	Intercept			
Y=	1	X+	0	

Item No.	29	Item Name	IgG	Optical
<b>Aspiration</b>				
Kind	Double			
Vol.				
Kind	Vol.	Add	Units	
Sample	3.8	5	µl	
Reagent 1	250	10	µl	
Reagent 2	95	10	µl	
<b>Data Process</b>				
Read	Main	Start	End	
	Sub	29	31	
<b>Abs.Limit</b>				
Low	-3	High	3	
<b>Blank value</b>				
<b>Water Blank</b>				
<b>Reaction Monitor</b>				
0 Level Point	1			
Span	3			
<b>Third mixing</b>				
<b>ON</b>				
<b>Correction value</b>				
Blank correction	0.7276			
End Point Limit	2			
Linear Check (%)	0			
<b>Prozone Check</b>				
	Start	End	Limit (%)	
First				
Second				Low

Item No.	29	Item Name	IgG	Optical
<b>Normal Range</b>				
	Male		Female	
	Low	High	Low	High
Serum	5.65	17.65	5.65	17.65
Urine				
Plasma				
CSF				
Dialysis				
Other				
<b>Panic Range</b>				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	29	Item Name	IgG	Optical
<b>Auto Rerun SW</b>				
	OFF			
<b>Auto Rerun Range (Conc.)</b>				
	First Dil	Low	High	
		Re	Value	Dil
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				
<b>Auto Rerun Condition (Absorbance)</b>				
	Lower	OFF		
	Higher	OFF		
<b>Auto Rerun Condition (Prozone)</b>				
	OFF			
<b>Dilution</b>				
	99:Dil1			

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