## PRESTIGE 24i LQ CK

Nr kat. 4-220, 4-420 (PL)

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Cat. No **4-220, 4-420** (EN)

### INTENDED USE


Diagnostic kit for determination of creatine kinase activity used in automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Creatine kinase (CK) catalyzes the transfer of phosphate group between creatine phosphate and adenosine diphosphate (ADP). The product of this reaction is adenosine triphosphate (ATP) – molecular source of energy. CK is a dimer, composed of two different subunits called M and B. Three different isoenzymes formed from these subunits are found in brain and smooth muscle (BB), skeletal muscle (MM) and cardiac muscle (MM and MB). Increased level of CK is usually the result of muscle injury, myocardial or pulmonary infarction.

### METHOD PRINCIPLE

Optimized kinetic method according to International Federation of Clinical Chemistry (IFCC).

The rate of absorbance changes at  $\lambda=340$  nm is directly proportional to creatine kinase activity.

### REAGENTS

#### Package

	Cat. No <b>4-220</b> <b>(24-TRAY)</b>	Cat. No <b>4-420</b> <b>(36-TRAY)</b>
1-Reagent	2 x 40 ml	2 x 22.5 ml
2-Reagent	2 x 10.5 ml	2 x 6 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. Stability on board of the analyser at 2-10°C: Prestige 24i - 12 weeks, Biolis 24i Premium – 12 weeks.

### Concentrations in the test

1-Reagent	
imidazole buffer	100 mmol/l
glucose	20 mmol/l
N-acetylcysteine	20 mmol/l
magnesium acetate	10 mmol/l
EDTA	2 mmol/l
NADP	2 mmol/l
ADP	2 mmol/l
AMP	5 mmol/l
HK	> 2.5 U/ml

### 2-Reagent

diadenosinepentaphosphate	10 μmol/l
glucose-6-phosphate-dehydrogenase (G6P-DH)	> 1.5 U/ml
creatine phosphate	30 mmol/l

### Warnings and notes

- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- Do not freeze reagents.
- Protect from direct sunlight and avoid contamination!
- Do not use reagents past the expiry date.
- Do not interchange caps among reagents.
- 1-Reagent meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

### Ingredients:

1-Reagent contains imidazole.

H360 May damage fertility or the unborn child.  
P201 Obtain special instructions before use.  
P202 Do not handle until all safety precautions have been read and understood.  
P308+P313 IF exposed or concerned: Get medical advice/attention.  
P405 Store locked up.  
P501 Dispose of the contents/containers in accordance with the current legislation on waste treatment.

### SPECIMEN

Serum, free from hemolysis.  
CK activity is unstable and is rapidly lost during storage. Probes should be stored tightly closed and protected from light.  
Specimens can be stored up to 4-8 hours at 15-25 °C or 1-2 days at 2-8°C or 1 month at -20°C.  
Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.  
1-Reagent put on basic position in reagent tray.  
2-Reagent put on start position in reagent tray.  
For reagent blank deionized water is recommended.

### Actions required:

- Biolis 24i Premium:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: CK - CREATININE. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_008\_BIOLIS\_24i\_PREMIUM\_CARRYOVER.
- Biolis 30i:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: CK – CREATININE. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_009\_BIOLIS\_30i\_CARRYOVER

### REFERENCE VALUES <sup>8</sup>

serum	37°C	
female	< 167 U/l	< 2.78 μkat/l
male	< 190 U/l	< 3.17 μkat/l

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended.

The calibration curve should be prepared every 12 weeks (Prestige 24i, Biolis 24i Premium), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**  
1.0 U/l (0.017 μkat/l)
- LoD (Limit of Detection):**  
2.0 U/l (0.033 μkat/l)
- LoQ (Limit of Quantitation):**  
8.5 U/l (0.14 μkat/l)
- Linearity:**  
up to 2800 U/l (46.67 μkat/l).

For higher activity, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by the dilution factor.

- Specificity / Interferences**  
Haemoglobin up to 0.156 g/dl, bilirubin up to 20 mg/dl, ascorbate up to 62 mg/l and triglycerides up to 1000 mg/dl do not interfere with the test.

### Precision

Repeatability (run to run) n = 20	Mean [U/l]	SD [U/l]	CV [%]
level 1	138.0	0.82	0.59
level 2	421.3	2.29	0.54
Reproducibility (day to day) n = 80	Mean [U/l]	SD [U/l]	CV [%]
level 1	135.3	6.48	4.8
level 2	434.4	5.26	1.2

### Method comparison

A comparison between CK values determined at **Biolis 30i** (y) and at **ADVIA SIEMENS 1800** (x) using 55 serum samples gave following results:

$$y = 0.922x + 2.4994 \text{ U/l};$$

$$R = 0.999$$

(R – correlation coefficient)

### WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	19	CK
<b>Auto Rerun SW</b>		
ON		
<b>Auto Rerun Range (Result)</b>		
	ON Lower	ON Higher
Serum	3.2	1200
Urine		
Plasma		
CSF		
Dialysis		
Other		
<b>Auto Rerun Condition (Absorbance)</b>		
Absorbance Range	Lower Higher	OFF OFF
<b>Prozone Range</b>		
		OFF

• Prestige 24i Premium, Biolis 24i Premium

Item No.	19	Item Name	CK	Optical	
<b>Data information</b>					
Units	U/L				
Decimals	1				
<b>Calibration</b>					
Type	Linear2				
Std sample conc.					
Blank	0	#1	*	#2	*
#3		#4		#5	
#6					
<b>Analysis</b>					
Type	RATE method				
Main Wave Length	340nm				
Sub Wave Length	660nm				
Method	IFCC				
<b>Correlation</b>					
	Slope	Intercept			
Y=	1	X+	0		

Item No.	19	Item Name	CK	Optical
<b>Aspiration</b>				
Kind	Double			
Vol.				
	Kind	Vol.	Add	Units
Sample		12	5	µl
Reagent 1		200	10	µl
Reagent 2		40	10	µl
<b>Data Process</b>				
Read	Main Sub	Start 38	End 52	
<b>Abs.Limit</b>				
Low	-3		High	3
<b>Correction value</b>				
Blank correction	1			
End Point Limit	2			
Linear Check (%)	90			
<b>Prozone Check</b>				
	Start	End	Limit (%)	
First				
Second			Low	
<b>Blank value</b>				
<b>Water Blank</b>				
<b>Reaction Monitor</b>				
0 Level Point	1			
Span	3			
<b>Third mixing</b>				
<b>OFF</b>				

Item No.	19	Item Name	CK	Optical
<b>Normal Range</b>				
	Male		Female	
	Low	High	Low	High
Serum	0	190	0	167
Urine				
Plasma	0	190	0	167
CSF				
Dialysis				
Other				
<b>Panic Range</b>				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	19	Item Name	CK	Optical	
<b>Auto Rerun SW</b>					
ON					
<b>Auto Rerun Range (Conc.)</b>					
	First Dil	Low Value	Dil	High Value	Dil
Serum		4.4		1600	
Urine					
Plasma					
CSF					
Dialysis					
Other					
<b>Auto Rerun Condition (Absorbance)</b>					
	Lower	Higher			OFF
					OFF
<b>Auto Rerun Condition (Prozone)</b>					
	OFF				
<b>Dilution</b>					
	100:DiI2				

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