PRESTIGE 24i GLUCOSE HEX

Nr kat. 4-231, 4-423 (PL)

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LITERATURA

- Bahram P., Trinder P.: *Analyst 97*, 142-145 (1972).
- Burtis C.A., Ashwood E.R., ed. *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 4th ed.*, PA: WB Saunders, 868-869, 2006.
- McPherson Richard A., Pincus Matthew R.: *Henry's Clinical Diagnosis and Management by Laboratory Methods*, 22nd ed. Philadelphia, PA: WB Saunders, 2011.
- Dujmovic and F. Deisenhammer, *Stability of cerebrospinal fluid/serum glucose ratio and cerebrospinal fluid lactate concentrations over 24 h: analysis of repeated measurements*, *Clinical Chemistry and Laboratory Medicine*, Volume 48, Issue 2, 2010, pp. 209-212.
- Zalecenia kliniczne dotyczące postępowania u chorych na cukrzycę 2014, *Diabetologia Kliniczna*, tom 3, suplement A, 2014.
- Sacks, David B., et al. *Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus*. *Clinical chemistry* 48.3 (2002): 436-472.
- Miles RR, Roberts RF, Putnam AR, Roberts WL. *Comparison of serum and heparinized plasma samples for measurement of chemistry analytes*. *Clin Chem* 2004;50:1704-1706.
- Alan H.B. Wu: *Tietz Clinical Guide to Laboratory Tests*, 4th ed. WB Saunders, 444-450 (2006).
- Young DS., *Effects of drugs on clinical laboratory tests*, 5. Vol. 2. Washington DC, USA: AACC Press (2000).

Data wydania: 04.2021

PRESTIGE 24i GLUCOSE HEX

Cat. No **4-231, 4-423** (EN)

INTENDED USE

Diagnostic kit for determination of glucose concentration used in automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.


The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Glucose is a simple six-carbon sugar. Oxidative metabolism of glucose provides the energy for most cellular processes. Glucose level in the blood is tightly controlled by several hormones. Elevated glucose level is the classic sign of diabetes mellitus. Glucose level abnormalities (hyper- or hypoglycemia) might be caused also by pancreas tumors and diseases of liver, thyroid gland or adrenal glands.

METHOD PRINCIPLE

Enzymatic method with hexokinase and glucose-6-phosphate dehydrogenase (G6P-DH).

The rate of NADH formation is directly proportional to the glucose concentration in the sample.

REAGENTS

Package

	Cat. No 4-231 (24-TRAY)	Cat No. 4-423 (36-TRAY)
1-Reagent	6 x 39.3 ml	8 x 25.2 ml
2-Reagent	6 x 10.3 ml	8 x 6.7 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. Stability on board of the analyser at 2-10°C: Prestige 24i – 12 weeks, Biolis 24i Premium – 12 weeks.

Concentrations in the reagents

1-Reagent	
PIPES buffer (pH 7.5)	80 mmol/l
Mg ²⁺	10 mmol/l
ATP	4 mmol/l
NAD	3 mmol/l
preservative	
2-Reagent	
hexokinase	≥ 4500 U/l
glucose-6-phosphate dehydrogenase (G6P-DH)	≥ 14000 U/l
preservative	

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Do not use after expiry date.
- Do not freeze reagents.
- Do not interchange caps.

- Reagents should be mixed before use by gentle inverting the bottles several times.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

EDTA or heparinized plasma / serum, free from hemolysis, cerebrospinal fluid, urine.

Plasma / Serum. Serum and plasma specimens should be separated from cells within 30 minutes after collection.

Plasma specimen which is not assayed immediately after collection should be kept in tubes containing sodium fluoride or sodium iodoacetate. These compounds adding prevent glycolysis and stabilize glucose level.

Serum and plasma can be stored up to 2 days at 4°C.³

Plasma is the specimen recommended for the glucose determination in the blood.⁵

Cerebrospinal fluid. Glucose concentration in cerebrospinal fluid should be measured directly after specimen collection. Cerebrospinal fluid must be analysed simultaneously with a blood sample.

After centrifuge CSF sample can be stored up to 24 hours at 4°C.⁴

Urine. Collect 24-hour sample in dark bottle and keep on ice. Preserve sample by adding 5 ml of glacial acetic acid to the container before starting the collection. The final pH of the sample should be between 4 and 5. Centrifuge samples with visible turbidity or precipitates before analysis.

Urine can be stored up to 24 hour at 4°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank deionized water is recommended.

Actions required:

- Biolis 30i:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: GLUCOSE HEX - CREATININE. To avoid this effect follow the recommendations contained in the instruction 51_03_24_009_BIOLIS 30i_CARRYOVER.

REFERENCE VALUES

	mg/dl	mmol/l
plasma, serum ^{5,6,7}	70 – 99	3.9 – 5.5
urine (24h) ⁸	1 – 15	0.1 – 0.8
cerebrospinal fluid ⁸	40 – 70	2.2 – 3.9

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum or CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine.

For the calibration of automatic analyser the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 12 weeks (Prestige 24i, Biolis 24i Premium), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following result have been obtained using automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**
0.6 mg/dl (0.03 mmol/l)
- LoD (Limit of Detection):**
1.2 mg/dl (0.067 mmol/l)
- LoQ (Limit of Quantitation):**
3 mg/dl (0.17 mmol/l)
- Linearity:**
up to 800 mg/dl (44,4 mmol/l)

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 1.25 g/dl, bilirubin up to 40 mg/dl, ascorbate up to 62 mg/L and triglycerides up to 1000 mg/dl do not interfere with the test. Some medicines can interfere.⁹

Precision

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	82.1	1.05	1.28
level 2	295.3	2.33	0.79
Reproducibility (day to day) n = 80	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	83	1.45	1.7
level 2	290	2.54	0.9

Method comparison

A comparison between glucose values determined at **Biolis 30i** (y) and at **BECKMAN COULTER AU680** (x) using 59 serum samples gave following results:

$$y = 1.0588x - 0.9542 \text{ mg/dl;}$$

$$R = 1.000 \quad (R - \text{correlation coefficient})$$

A comparison between glucose values determined at **Biolis 30i** (y) and at **ADVIA SIEMENS 1800** (x) using 59 plasma samples gave following results:

$$y = 1.0308x - 0.7088 \text{ mg/dl;}$$

$$R = 1.000 \quad (R - \text{correlation coefficient})$$

A comparison between glucose values determined at **Biolis 30i** (y) and at **BS-800** (x) using 60 cerebrospinal fluid samples gave following results:

$$y = 1.0069x + 0.4056 \text{ mg/dl;}$$

$$R = 0.998 \quad (R - \text{correlation coefficient})$$

A comparison between glucose values determined at **Biolis 30i** (y) and at **BS-800** (x) using 62 urine samples gave following results:

$$y = 1.034x - 1.3652 \text{ mg/dl;}$$

$$R = 1.000 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

- Barham P., Trinder P.: Analyst 97, 142-145 (1972).
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 4th ed., PA: WB Saunders, 868-869, 2006.
- McPherson Richard A., Pincus Matthew R.: Henry's Clinical Diagnosis and Management by Laboratory Methods, 22nd ed. Philadelphia, PA: WB Saunders, 2011.
- Dujmovic and F. Deisenhammer, Stability of cerebrospinal fluid/serum glucose ratio and cerebrospinal fluid lactate concentrations over 24 h: analysis of repeated measurements, Clinical Chemistry and Laboratory Medicine, Volume 48, Issue 2, 2010, pp. 209-212.
- Zalecenia kliniczne dotyczące postępowania u chorych na cukrzycę 2014, Diabetologia Kliniczna, tom 3, suplement A, 2014.
- Sacks, David B., et al. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. Clinical chemistry 48.3 (2002): 436-472.
- Miles RR, Roberts RF, Putnam AR, Roberts WL. Comparison of serum and heparinized plasma samples for measurement of chemistry analytes. Clin Chem 2004;50:1704-1706.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 444-450 (2006).
- Young DS., Effects of drugs on clinical laboratory tests, 5. Vol. 2. Washington DC, USA: AACC Press (2000).

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- Barham P., Trinder P.: Analyst 97, 142-145 (1972).
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 4th ed., PA: WB Saunders, 868-869, 2006.
- McPherson Richard A., Pincus Matthew R.: Henry's Clinical Diagnosis and Management by Laboratory Methods, 22nd ed. Philadelphia, PA: WB Saunders, 2011.
- Dujmovic and F. Deisenhammer, Stability of cerebrospinal fluid/serum glucose ratio and cerebrospinal fluid lactate concentrations over 24 h: analysis of repeated measurements, Clinical Chemistry and Laboratory Medicine, Volume 48, Issue 2, 2010, pp. 209-212.
- Zalecena kliniczne dotyczące postępowania u chorych na cukrzycę 2014, Diabetologia Kliniczna, tom 3, suplement A, 2014.
- Sacks, David B., et al. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. Clinical chemistry 48.3 (2002): 436-472.
- Miles RR, Roberts RF, Putnam AR, Roberts WL. Comparison of serum and heparinized plasma samples for measurement of chemistry analytes. Clin Chem 2004;50:1704-1706.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 444-450 (2006).
- Young DS., Effects of drugs on clinical laboratory tests, 5. Vol. 2. Washington DC, USA: AACC Press (2000).

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Item name	7	GLUC HEX
Auto Rerun SW		
ON	Auto Rerun Condition (Absorbance)	
	Absorbance Range	
	Lower	OFF
	Higher	OFF
Auto Rerun Range (Result)		
	ON	ON
	Lower	Higher
Serum	4.5	665
Urine		
Plasma		
CSF		
Dialysis		
Other		
Prozone Range		
		OFF

• Prestige 24i Premium, Biolis 24i Premium

Item No.	7	Item Name	GLUC HEX	Optical	
Data information					
Units	mg/dl				
Decimals	1				
Calibration					
Type	Linear2				
Std sample conc.					
Blank	0	#1	*	#2	*
#3		#4		#5	
#6					
Analysis					
Type	END				
Main Wave Length	340 nm				
Sub Wave Length					
Method	Hexokinase				
Correlation					
Y=	Slope	X+	Intercept		
	1		0		

Item No.	7	Item Name	GLUC HEX	Optical
Aspiration				
Kind	Double			
Vol.				
Kind	Vol.	Add	Units	
Sample	2.5	5	µl	
Reagent 1	200	10	µl	
Reagent 2	40	10	µl	
Data Process				
Read	Start	End		
Main	47	48		
Sub	30	31		
Abs.Limit				
Low	-3	High	3	
Blank value				
Water Blank				
Reaction Monitor				
0 Level Point	1			
Span	3			
Third mixing				
OFF				
Correction value				
Blank correction				
End Point Limit				
Linear Check (%)				
2				
0				
Prozone Check				
	Start	End	Limit (%)	
First				
Second				Low

Item No.	7	Item Name	GLUC HEX	Optical
Normal Range				
	Male		Female	
	Low	High	Low	High
Serum	70	99	70	99
Urine	1	15	1	15
Plasma	70	99	70	99
CSF	40	70	40	70
Dialysis				
Other				
Panic Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	7	Item Name	GLUC HEX	Optical		
Auto Rerun SW						
ON	Auto Rerun Condition (Absorbance)					
	Lower	Higher	OFF	OFF		
Auto Rerun Range (Conc.)						
	First Dil	Low	High			
	Re	Value	Dil	Re	Value	Dil
Serum		4.4			670	
Urine						
Plasma						
CSF						
Dialysis						
Other						
Auto Rerun Condition (Prozone)						
OFF						
Dilution						
100:Dil2						

PRESTIGE 24i GLUCOSE HEX

• Biolis 30i

Item no	7	Item name	GLUC HEX	Specimen	SERUM/ PLASMA/ CSF/ URINE	OPTICAL												
Data information																		
UNITS	mg/dL		Aspiration volume															
DECIMALS	1		TYPE Double															
Analysis																		
METHOD	END method		<table border="1"> <tr> <td></td> <td>SAMPLE</td> <td>REAGENT 1</td> <td>REAGENT 2</td> </tr> <tr> <td>VOL. (µL)</td> <td>2.5</td> <td>200</td> <td>40</td> </tr> <tr> <td>BOTTLE (ml)</td> <td></td> <td></td> <td></td> </tr> </table>					SAMPLE	REAGENT 1	REAGENT 2	VOL. (µL)	2.5	200	40	BOTTLE (ml)			
	SAMPLE	REAGENT 1	REAGENT 2															
VOL. (µL)	2.5	200	40															
BOTTLE (ml)																		
Main Wave Length	340 nm		FIRST DIL.															
Sub Wave Length																		
CORRELATION (Y= AX + B)																		
A =	1		Data processing read															
B =	0		<table border="1"> <tr> <td></td> <td>START</td> <td>END</td> </tr> <tr> <td>MAIN</td> <td>48</td> <td>49</td> </tr> <tr> <td>SUB</td> <td>28</td> <td>29</td> </tr> </table>					START	END	MAIN	48	49	SUB	28	29			
	START	END																
MAIN	48	49																
SUB	28	29																
Blank value																		
• WATER		° REAGENT																
Calibration																		
TYPE	Linear 2		ABS LIMIT															
STABILITY			-3 TO 3															
Collection value																		
END POINT		2.5																
LINEARITY CHECK (%)		0																
Prozone check																		
		START	END	LIMIT (%)														
FIRST																		
SECOND																		
°HIGH				MINIMUM ABS.														
•LOW				MEAN														
				VARIATE														
Item No	7	Item Name	GLUC HEX	Specimen	SERUM/ PLASMA/ CSF/ URINE	OPTICAL												
Reference intervals																		
MALE				FEMALE														
LOW	HIGH	LOW	HIGH															
70	99	70	99															
Panic range																		
MALE				FEMALE														
LOW	HIGH	LOW	HIGH															
Decision limit																		
MALE				FEMALE														
Reaction check																		
		°ON	•OFF															
CHECK																		
LOW																		
HIGH																		
VL CHECK			VH CHECK															
°ON	•OFF		°ON	•OFF														
Auto rerun																		
		•ON		°OFF														
Auto rerun range (conc.)																		
Re	Value	Dil.	Re	Value	Dil.													
	3			800														
Auto rerun condition (abs.)																		
					DIL.													
LOWER	°ON		•OFF															
HIGH	°ON		•OFF															
Auto rerun condition (prozone)																		
		°ON		•OFF														
SAMPLE VOL.																		
Dilution																		
		•DIL 1		° DIL 2														

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