PRESTIGE 24i LQ GGT

Nr kat. 4-224, 4-424 (PL)

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LITERATURA

- Szasz G., Weimann G. Suhler F., Wahlefrld A.W., Persijn J.P.: Z. Klin. Chem. Klin. Biochem., 12 (1974): 228.
- Persijn J.P., van der Silk W.: J. Clin. Chem. Clin. Biochem. 14 (1976): 421-427
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 2nd ed. Philadelphia, PA: WB Saunders, (1994), 850-1.
- Tietz N.W., ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA (1995), WB Saunders: 286
- Kaplan L.A., Pesce A.J.: Clinical Chemistry. Theory, analysis and correlation 3rd Ed. (1996), the C.V. Mosby Company, St. Louis: 1072.
- Shaw LM, Strømme JH, London JL, Theodorsen L. International Federation Clinical Chemistry, (IFCC), Scientific Committee, Analytical Section: IFCC method for gamma-glutamyltransferase J Clin ChemClin Biochem. Vol 21(1983): 633-46
- Schumann GI et al. International Federation of Clinical Chemistry and Laboratory Medicine. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 degrees C. International Federation of Clinical Chemistry and Laboratory Medicine. Part 6. Reference procedure for the measurement of gamma-glutamyltransferase. Clin Chem Lab Med. (2002), 40(7): 734-8.
- Schumann GI, Klauke R. New IFCC reference procedures for the determination of catalytic activity concentrations of five enzymes in serum: preliminary upper reference limits obtained in hospitalized subjects. Clin Chim Acta 327 (2003): 69-79

Data wydania: 10. 2023.

PRESTIGE 24i LQ GGT

Cat. No **4-224, 4-424**

(EN)

Warning

H317 May cause an allergic skin reaction.
H412 Harmful to aquatic life with long lasting effects.
P280 Wear protective gloves, protective clothing, eye protection or face protection.

P302 + P352 IF ON SKIN: Wash with plenty of soap and water.

P273 Avoid release to the environment.

INTENDED USE


Diagnostic kit for determination of γ -glutamyltransferase activity used in automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

γ -Glutamyltransferase (GGT, GGTP) is a membrane localized enzyme that catalyzes the transfer of glutamyl groups from glutathione to amino acids or peptides. Large GGT amounts are present in secretory organs: kidney, liver, bile duct, pancreas. Although the GGT activity is highest in renal tissue, serum GGT is generally elevated as a result of liver disease. Since alcohol induces GGT production, measurement of GGT activity is used for monitoring of abstinence in withdrawal treatment.

METHOD PRINCIPLE

IFCC method. Kinetic method with L- γ -glutamyl-3-carboxy-4-nitroanilide.

The rate of 5-amino-2-nitrobenzoate creation measured colorimetrically is directly proportional to γ -glutamyl-transferase activity.

REAGENTS

Package

	Cat. No 4-224 (24-TRAY)	Cat. No 4-424 (36-TRAY)
1-REAGENT	4 x 40 ml	8 x 23 ml
2-REAGENT	4 x 12.5 ml	8 x 7.5 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. Stability on board of the analyser at 2-10°C: 8 weeks (Prestige 24i) or 12 weeks (Biolis 24i Premium, Biolis 30i).

Concentrations in the test

Tris (pH 8.25)	100 mmol/l
glycylglycine	100 mmol/l
L- γ -glutamyl-3-carboxy-4-nitroanilide	4 mmol/l

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT and 2-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

1-REAGENT and 2-REAGENT contain reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1).

SPECIMEN

Serum, EDTA plasma free from hemolysis.

Do not use citrate, oxalate and fluoride as anticoagulants because of GGT activity inhibition!

Heparin causes turbidity in the reaction mixture!

GGT activity remains stable in specimen up to 2 days at 15-25°C or 1 week at 2-8°C or 1 month at -25°C. Freezing of sample causes a loss of enzyme activity. Frozen specimens should be thawed and kept at room temperature for 18 to 24 hours before measurement to achieve full enzyme reactivation.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-REAGENT and 2-REAGENT are ready to use.

1-REAGENT put on basic position in reagent tray.

2-REAGENT put on start position in reagent tray.

For reagent blank deionized water is recommended.

Actions required:

- Biolis 24i Premium:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: GGT – UA, URINE PROTEINS - GGT. To avoid this effect follow the recommendations contained in the instruction 51_03_24_008_BIOLIS_24i_PREMIUM_CARRYOVER.
- Biolis 30i:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: GGT – URINE PROTEINS. To avoid this effect follow the recommendations contained in the instruction 51_03_24_009_BIOLIS_30i_CARRYOVER.

REFERENCE VALUES ⁷

serum / plasma	37°C
female	< 38 U/l (0.633 μ kat/l)
male	< 55 U/l (0.917 μ kat/l)

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 8 weeks (Prestige 24i) or every 12 weeks (Biolis 24i Premium, Biolis 30i), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):** 0.0 U/l (0.0 μ kat/l)
- LoD (Limit of Detection):** 0.4 U/l (0.006 μ kat/l)
- LoQ (Limit of Quantitation):** 4.0 U/l (0.06 μ kat/l)
- Linearity:** up to 880 U/l (14.67 μ kat/l)

Specificity / Interferences

Haemoglobin up to 0.16 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean [U/l]	SD [U/l]	CV [%]
level 1	60.9	0.52	0.86
level 2	158.5	0.77	0.48

Method comparison

A comparison between GGT values determined at **Biolis 30i** (y) and at **BECKMAN COULTER AU680** (x) using 60 serum samples gave following results:

$$y = 0.9746x + 0.1781 \text{ U/l;}$$

$$R = 1.000 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

- Szasz G., Weimann G. Suhler F., Wahlefrld A.W., Persijn J.P.: Z. Klin. Chem. Klin. Biochem., 12 (1974): 228.
- Persijn J.P., van der Silk W.: J. Clin. Chem. Clin. Biochem. 14 (1976): 421-427
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 2nd ed. Philadelphia, PA: WB Saunders, (1994), 850-1.
- Tietz N.W., ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA (1995), WB Saunders: 286
- Kaplan L.A., Pesce A.J.: Clinical Chemistry. Theory, analysis and correlation 3rd Ed. (1996), the C.V. Mosby Company, St. Louis: 1072.
- Shaw LM, Strømme JH, London JL, Theodorsen L. International Federation Clinical Chemistry, (IFCC), Scientific Committee, Analytical Section: IFCC method for gamma-glutamyltransferase J Clin ChemClin Biochem. Vol 21(1983): 633-46
- Schumann GI et al. International Federation of Clinical Chemistry and Laboratory Medicine. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 degrees C. International Federation of Clinical Chemistry and Laboratory Medicine. Part 6. Reference procedure for the measurement of gamma-glutamyltransferase. Clin Chem Lab Med. (2002), 40(7): 734-8.
- Schumann GI, Klauke R. New IFCC reference procedures for the determination of catalytic activity concentrations of five enzymes in serum: preliminary upper reference limits obtained in hospitalized subjects. Clin Chim Acta 327 (2003): 69-79

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- Szasz G., Weimann G. Suhler F., Wahlefrld A.W., Persijn J.P.: Z. Klin. Chem. Klin. Biochem., 12 (1974): 228.
- Persijn J.P., van der Silk W.: J. Clin. Chem. Clin. Biochem. 14 (1976): 421-427
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 2nd ed. Philadelphia, PA: WB Saunders, (1994), 850-1.
- Tietz N.W., ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA (1995), WB Saunders: 286
- Kaplan L.A., Pesce A.J.: Clinical Chemistry. Theory, analysis and correlation 3rd Ed. (1996), the C.V. Mosby Company, St. Louis: 1072.
- Shaw LM, Strømme JH, London JL, Theodorsen L. International Federation Clinical Chemistry, (IFCC), Scientific Committee, Analytical Section: IFCC method for gamma-glutamyltransferase J Clin Chem Clin Biochem. Vol 21(1983): 633-46
- Schumann GI et al. International Federation of Clinical Chemistry and Laboratory Medicine. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 degrees C. International Federation of Clinical Chemistry and Laboratory Medicine. Part 6. Reference procedure for the measurement of gamma-glutamyltransferase. Clin Chem Lab Med. (2002), 40(7): 734-8.
- Schumann GI, Klauke R. New IFCC reference procedures for the determination of catalytic activity concentrations of five enzymes in serum: preliminary upper reference limits obtained in hospitalized subjects. Clin Chim Acta 327 (2003): 69-79

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Item name	18	GGT
Auto Rerun SW		
ON		
Auto Rerun Condition (Absorbance)		
Absorbance Range	Lower	OFF
	Higher	OFF
Auto Rerun Range (Result)		
	ON	ON
	Lower	Higher
Serum	11.2	580
Urine		
Plasma		
CSF		
Dialysis		
Other		
Prozone Range		
		OFF

- Prestige 24i Premium, Biolis 24i Premium

Item No.	18	Item Name	GGT	Optical	
Data information					
Units	U/L				
Decimals	1				
Calibration					
Type	Linear2				
Std sample conc.					
Blank	0	#1	*	#2	*
#3		#4		#5	
#6					
Analysis					
Type	RATE method				
Main Wave Length	405nm				
Sub Wave Length	660nm				
Method	IFCC				
Correlation					
	Slope	Intercept			
Y=	1	X+	0		

Item No.	18	Item Name	GGT	Optical
Aspiration				
Kind	Double			
Data Process				
Read	Main	Start	End	
	Sub	29	30	
Vol.				
	Kind	Vol.	Add	Units
Sample		5	5	µl
Reagent 1		200	10	µl
Reagent 2		50	10	µl
Abs.Limit				
Low	-0.1	High	2	
Blank value				
Water Blank				
Reaction Monitor				
0 Level Point	1			
Span	3			
Third mixing				
OFF				
Correction value				
Blank correction				
End Point Limit				
2				
Linear Check (%)				
10				
Prozone Check				
	Start	End	Limit (%)	
First				
Second			Low	

Item No.	18	Item Name	GGT	Optical
Normal Range				
	Male		Female	
	Low	High	Low	High
Serum	0	55	0	38
Urine				
Plasma	0	55	0	38
CSF				
Dialysis				
Other				
Panic Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	18	Item Name	GGT	Optical	
Auto Rerun SW					
ON					
Auto Rerun Condition (Absorbance)					
	Lower	OFF			
	Higher	OFF			
Auto Rerun Range (Conc.)					
	First Dil	Low Value	Dil	High Value	Dil
Serum		11.2		580	
Urine					
Plasma					
CSF					
Dialysis					
Other					
Auto Rerun Condition (Prozone)					
OFF					
Dilution					
100:DiI2					

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