PRESTIGE 24i LQ ALBUMIN

Nr kat. **4-238, 4-438** (PL)

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Cat. No **4-238, 4-438** (EN)

INTENDED USE

Diagnostic kit for determination of albumin concentration intended to use in automatic analyzers: Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium, and Biolis 30i. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Albumin is the major serum protein, but is present also in other body fluids: cerebrospinal, pleural and peritoneal. Albumin regulates blood oncotic pressure and serves as amino acids reservoir. Beyond of these functions albumin is very important transport protein – binds and keeps dispersed bilirubin, hormones, vitamins, calcium, magnesium, fatty acids and medicines. Decreased albumin blood level is caused usually by liver or kidney disease, malabsorption or malnutrition.

METHOD PRINCIPLE

Bromocresol green (BCG) forms with albumin, in succinate buffer (acid medium), a coloured complex. The absorbance of this complex is proportional to the albumin concentration in the sample. The colour intensity of the formed complex measured at 630 nm is proportional to albumin concentration in the sample.

REAGENTS

Package

	Cat. No 4-238 (24-TRAY)	Cat. No 4-438 (36-TRAY)
1-REAGENT	4 x 60 ml	6 x 25 ml

The reagent is stable up to the kit expiry date printed on the package when stored at 2-8°C. Stability on board of the analyser at 2-10°C is 8 weeks (Prestige 24i, Biolis 24i Premium) or 12 weeks (Biolis 30i).

Concentrations in the test

succinate buffer	≤ 108 mmol/l
bromocresol green (BCG)	≤ 0.35 mmol/l
sodium hydroxide	≤ 60 mmol/l
detergent	
preservative	

WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- Do not freeze the reagent.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT meets the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Warning

H319 Causes serious eye irritation.
 P280 Wear protective gloves, protective clothing, eye protection or face protection.
 P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

SPECIMEN

Serum free from hemolysis.
 Serum should be separated from red blood cells as soon as possible after blood collection.
 Serum can be stored up to 3 days at 2-8°C or 6 months at -20°C.²
 Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-REAGENT is ready to use.
 1-REAGENT put on basic position in reagent tray.
 For reagent blank deionized water is recommended.

REFERENCE VALUES 7

serum		g/dl	g/l
children	3 – 5 years	3.9 – 5.0	39 - 50
	6 – 15 years	4.1 – 5.1	41 - 51
females	16 – 54 years	3.9 – 5.0	39 – 50
males	16 – 29 years	4.6 – 5.3	46 - 53
	30 – 54 years	4.4 – 5.1	44 - 51
females and males	55 – 79 years	4.2 – 5.0	42 – 50

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173). For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended.

The calibration curve should be prepared every 7 weeks (Prestige 24i, Biolis 24i Premium) or 4 weeks (Biolis 30i), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**
0.02 g/dl (0.2 g/l)

- LoD (Limit of Detection):**
0.04 g/dl (0.4 g/l)

- LoQ (Limit of Quantitation):**
0.8 g/dl (8 g/l)

- Linearity:**
up to 6.5 g/dl (65 g/l)

For higher concentration of dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 2.5 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1200 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean	SD	CV
	[g/dl]	[g/dl]	[%]
level 1	3.08	0.02	0.61
level 2	4.39	0.02	0.46
Reproducibility (day to day) n = 80	Mean	SD	CV
	[g/dl]	[g/dl]	[%]
level 1	4.4	0.07	1.7
level 2	3.1	0.06	1.8

Method comparison

A comparison between albumin values determined at **Biolis 30i** (y) and at **ADVIA SIEMENS 1800** (x) using 68 serum samples gave following results:
 $y = 0.9798x + 0.2306$ g/dl;
 $R = 0.986$ (R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	ALB	
Auto Rerun SW		
ON		
Auto Rerun Range (Result)		
	ON	ON
	Lower	Higher
Serum	1.13	6.5
Urine		
Plasma		
CSF		
Dialysis		
Other		
Auto Rerun Condition (Absorbance)		
Absorbance Range	Lower	OFF
	Higher	OFF
Prozone Range		
Prozone Range	OFF	

▪ Prestige 24i Premium, Biolis 24i Premium

Item No.	1	Item Name	ALB	Optical	
Data information					
Units	g/dl				
Decimals	2				
Calibration					
Type	Linear2				
Std sample conc.					
Blank	0	#1	*	#2	*
#3		#4		#5	
#6					
Analysis					
Type	END method				
Main Wave Length	600nm				
Sub Wave Length	750nm				
Method	BCG				
Correlation					
Slope		Intercept			
Y=	1	X+	0		

Item No.	1	Item Name	ALB	Optical
Aspiration				
Kind	Single			
Vol.				
	Kind	Vol.	Add	Units
Sample	3	5		µl
Reagent 1	300	10		µl
Reagent 2				
Data Process				
Read	Start	End		
Main	29	30		
Sub	9	10		
Abs.Limit Low				
	-0.1	~	High	
			2.5	
Blank value				
Water Blank				
Reaction Monitor				
0 Level Point	1			
Span	3			
Third mixing				
OFF				
Correction value				
Blank correction				
End Point Limit		2		
Linear Check (%)				
Prozone Check				
	Start	End	Limit (%)	
First				
Second			Low	
Third			Low	

Item No.	1	Item Name	ALB	Optical
Normal Range				
	Male		Female	
	Low	High	Low	High
Serum	3.5	5.2	3.5	5.2
Urine				
Plasma				
CSF				
Dialysis				
Other				
Panic Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	1	Item Name	ALB	Optical			
Auto Rerun SW							
ON							
Auto Rerun Condition (Absorbance)							
Lower		OFF					
Higher		OFF					
Auto Rerun Range (Conc.)							
	First Dil	Low		High			
		Re	Value	Dil	Re	Value	Dil
Serum			1.14			6.5	
Urine							
Plasma							
CSF							
Dialysis							
Other							
Auto Rerun Condition (Prozone)							
OFF							
Dilution							
100:Dil2							

PRESTIGE 24i LQ ALBUMIN

▪ **Biolis 30i**

Item no	1	Item name	ALB	Specimen	SERUM	OPTICAL
Data information						
UNITS		g/dL		Aspiration volume		
DECIMALS		2		TYPE Single		
Analysis						
METHOD		END method				
Main Wave Length		600 nm				
Sub Wave Length		750 nm				
CORRELATION (Y= AX + B)						
A =		1				
B =		0				
Blank value						
• WATER		° REAGENT				
Calibration						
TYPE		Linear 2				
STABILITY						
Data processing read						
START		END				
MAIN		27 28				
SUB						
ABS LIMIT						
-0.1		TO			2.5	
Collection value						
END POINT		2.5				
LINEARITY CHECK (%)		0				
Prozone check						
START		END		LIMIT (%)		
FIRST						
SECOND						
°HIGH				MINIMUM ABS.		
•LOW				MEAN		
				VARIATE		

Item No	1	Item Name	ALB	Specimen	SERUM	OPTICAL
Reference intervals						
MALE		FEMALE				
LOW	HIGH	LOW	HIGH			
4.2	5.3	3.9	5			
Panic range						
MALE		FEMALE				
LOW	HIGH	LOW	HIGH			
Decision limit						
MALE		FEMALE				
Reaction check						
°ON		•OFF				
CHECK						
LOW						
HIGH						
VL CHECK						
°ON		•OFF				
•OFF		°ON				
		•OFF				
Auto rerun						
•ON		°OFF				
Auto rerun range (conc.)						
Re	Value	Dil.	Re	Value	Dil.	
	0.8			6.5		
Auto rerun condition (abs.)						
					DIL.	
LOWER	°ON		•OFF			
HIGH	°ON		•OFF			
Auto rerun condition (prozone)						
°ON		•OFF				
SAMPLE VOL.						
Dilution						
•DIL 1		° DIL 2				

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