PRESTIGE 24i LQ LDH

Nr kat. **4-239, 4-439** (PL)

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Cat. No **4-239, 4-439** (EN)

INTENDED USE

Diagnostic kit for determination of lactate dehydrogenase activity intended to use in automatic analyzers: automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.


The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Lactate dehydrogenase (LDH, LD) is intracellular enzyme occurred in all tissues. LDH catalyzes the reversible conversion of lactate to pyruvate using NAD⁺ as a cofactor. LD is a tetramer containing two possible forms of subunits: H and M. The result is five isoenzymes termed LD-1 (H₄) through LD-5 (M₄). The isoenzymes are present in different proportion in each tissue and have different electrophoresis mobility, what is very useful for diagnostic.

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The rate of absorbance changing at $\lambda=340$ nm is directly proportional to lactate dehydrogenase activity.

REAGENTS

Package

	Cat. No 4-239 (24-TRAY)	Cat. No 4-439 (36-TRAY)
1-Reagent	4 x 40 ml	4 x 23 ml
2-Reagent	4 x 12.5 ml	4 x 7.5 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. Stability on board of the analyser at 2-10°C: Prestige 24i - 8 weeks, Biolis 24i Premium - 8 weeks.

Concentrations in the test

phosphate buffer (pH 7.5)	50 mmol/l
pyruvate	0.6 mmol/l
NADH	0.25 mmol/l
preservative	

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum, heparinized plasma free from hemolysis.

Do not use hemolyzed blood or serum because erythrocytes contain 150 times more LDH activity than serum.

As an anticoagulant for plasma preparation use heparin lithium or ammonium salt!

LDH activity is unstable and is rapidly lost during storage.

Specimens can be stored up to 4 hours at 15-25°C or 1-2 days at 2-8°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank deionized water is recommended.

Actions required:

- Biolis 24i Premium:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: ALAT (II GEN) – LDH, ASAT (II GEN) - LDH, ALAT – LDH, ASAT - LDH. To avoid this effect follow the recommendations contained in the instruction 51_03_24_008_BIOLIS_24i_PREMIUM_CARRYOVER.
- Biolis 30i:** When performing assays in analyser, there is a probability of **cross-contamination** affecting the tests results: LDH – ALAT (II GEN), ALAT (II GEN) - LDH. To avoid this effect follow the recommendations contained in the instruction 51_03_24_009_BIOLIS_30i_CARRYOVER.

REFERENCE VALUES ⁴

serum / plasma	37°C	37°C
adults	225 – 450 U/l	3.75 – 7.50 μ kat/l

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 5 weeks (Prestige 24i) or 8 weeks (Biolis 24i Premium), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):** 2.6 U/l (0.043 μ kat/l).

- LoD (Limit of Detection):** 6.1 U/l (0.10 μ kat/l).

- LoQ (Limit of Quantitation):** 11 U/l (0.18 μ kat/l).

- Linearity:** up to 1825 U/l (30.42 μ kat/l).

For higher activity dilute sample with 0.9% NaCl and repeat the assay. Multiply the result by the dilution factor.

Specificity / Interferences

Bilirubin up to 20 mg/dl, haemoglobin up to 5 g/dl, ascorbate up to 62 mg/l and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean [U/l]	SD [U/l]	CV [%]
level 1	337.2	5.24	1.55
level 2	802.5	4.86	0.61
Reproducibility (day to day) n = 80	Mean [U/l]	SD [U/l]	CV [%]
level 1	347.0	3.74	1.1
level 2	824.1	9.38	1.1

Method comparison

A comparison between LDH values determined at **Biolis 30i** (y) and at **BS-400** (x) using 60 serum samples gave following results:

$$y = 1.0134x - 0.823 \text{ U/l};$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	15	LDH
Auto Rerun SW		
ON		
Auto Rerun Range (Result)		
	ON	ON
	Lower	Higher
Serum	25	2000
Urine		
Plasma		
CSF		
Dialysis		
Other		
Auto Rerun Condition (Absorbance)		
Absorbance Range	Lower	OFF
	Higher	OFF
Auto Rerun Condition (Prozone)		
Prozone Range	OFF	

■ Prestige 24i Premium, Biolis 24i Premium

Item No.	15	Item Name	LDH	Optical
Data information				
Units	U/L			
Decimals	1			
Calibration				
Type	Linear2			
Std sample conc.				
Blank	0	#1	*	#2
#3		#4		#5
#6				
Analysis				
Type	RATE method			
Main Wave Length	340nm			
Sub Wave Length	450nm			
Method	DGKC			
Correlation				
Y=	Slope	X+	Intercept	
	1		0	

Item No.	15	Item Name	LDH	Optical
Aspiration				
Kind	Double			
Vol.				
	Kind	Vol.	Add	Units
Sample	3	5		µl
Reagent 1	160	10		µl
Reagent 2	40	10		µl
Data Process				
Read	Main	Start	End	
	Sub	36	52	
Abs.Limit				
Low	1	~	High	2
Correction value				
Blank correction				
End Point Limit 2				
Linear Check (%) 90				
Prozone Check				
	Start	End	Limit (%)	
First				
Second				Low
Blank value				
Water Blank				
Reaction Monitor				
0 Level Point	1			
Span	3			
Third mixing				
OFF				

Item No.	15	Item Name	LDH	Optical
Normal Range				
	Male		Female	
	Low	High	Low	High
Serum	225	450	225	450
Urine				
Plasma	225	450	225	450
CSF				
Dialysis				
Other				
Panic Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	15	Item Name	LDH	Optical
Auto Rerun SW				
ON				
Auto Rerun Range (Conc.)				
	First Dil	Low	High	
		Re	Value	Dil
Serum			20.1	
Urine				
Plasma			2000	
CSF				
Dialysis				
Other				
Auto Rerun Condition (Absorbance)				
	Lower	OFF		
	Higher	OFF		
Auto Rerun Condition (Prozone)				
OFF				
Dilution				
100:DiI2				

PRESTIGE 24i LQ LDH

▪ Biolis 30i

Item no	15	Item name	LDH	Specimen	SERUM/ PLASMA	OPTICAL												
Data information																		
UNITS	U/L		Aspiration volume															
DECIMALS	1		TYPE															
Analysis			Double															
METHOD	RATE method		<table border="1"> <tr> <td></td> <td>SAMPLE</td> <td>REAGENT 1</td> <td>REAGENT 2</td> </tr> <tr> <td>VOL. (µL)</td> <td>3</td> <td>160</td> <td>40</td> </tr> <tr> <td>BOTTLE (ml)</td> <td></td> <td></td> <td></td> </tr> </table>					SAMPLE	REAGENT 1	REAGENT 2	VOL. (µL)	3	160	40	BOTTLE (ml)			
	SAMPLE	REAGENT 1	REAGENT 2															
VOL. (µL)	3	160	40															
BOTTLE (ml)																		
Main Wave Length	340 nm		FIRST DIL.															
Sub Wave Length	450 nm																	
CORRELATION (Y= AX + B)																		
A =	1		Data processing read															
B =	0		<table border="1"> <tr> <td></td> <td>START</td> <td>END</td> </tr> <tr> <td>MAIN</td> <td>35</td> <td>54</td> </tr> <tr> <td>SUB</td> <td></td> <td></td> </tr> </table>					START	END	MAIN	35	54	SUB					
	START	END																
MAIN	35	54																
SUB																		
Blank value																		
• WATER ° REAGENT																		
Calibration																		
TYPE	Linear 2		ABS LIMIT															
STABILITY			TO															
Collection value																		
END POINT			2.5															
LINEARITY CHECK (%)			90															
Prozone check																		
		START	END	LIMIT (%)														
FIRST																		
SECOND																		
		MINIMUM ABS.																
°HIGH		MEAN																
•LOW		VARIATE																

Item No	15	Item Name	LDH	Specimen	SERUM/ PLASMA	OPTICAL
Reference intervals						
MALE				FEMALE		
LOW	HIGH	LOW	HIGH			
225	450	225	450			
Panic range						
MALE				FEMALE		
LOW	HIGH	LOW	HIGH			
Decision limit						
		MALE	FEMALE			
Reaction check						
		°ON	•OFF			
CHECK						
LOW						
HIGH						
VL CHECK			VH CHECK			
°ON	•OFF	°ON	•OFF			
Auto rerun						
		•ON	°OFF			
Auto rerun range (conc.)						
Re	Value	Dil.	Re	Value	Dil.	
	11			1825		
Auto rerun condition (abs.)						
					DIL.	
LOWER	°ON	•OFF				
HIGH	°ON	•OFF				
Auto rerun condition (prozone)						
°ON		•OFF				
SAMPLE VOL.						
Dilution						
•DIL 1		° DIL 2				

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