PRESTIGE 24i CALCIUM ARSENAZO

Nr kat. 4-247, 4-447 (PL)

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The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Calcium plays an essential role in many cell functions: intracellularly in muscle contraction and glycogen metabolism, extracellularly, in bone mineralization, in blood coagulation and in transmission of nerve impulses. Calcium is present in plasma in three forms: free, bound to proteins or complexed with anions as phosphate, citrate and bicarbonate. Decreased total calcium levels can be associated with diseases of the bone apparatus (especially osteoporosis), kidney diseases (especially under dialysis), defective intestinal absorption and hypoparathyroidism. Increased total calcium can be measured in hyperparathyroidism, malignant diseases with metastases and sarcoidosis. Calcium measurements also help in monitoring of calcium supplementation mainly in the prevention of osteoporosis.

METHOD PRINCIPLE

Photometric test using arsenazo III. Calcium with arsenazo III at neutral pH yields a blue colored complex, whose intensity is proportional to the calcium concentration. Interference by magnesium is eliminated by addition of 8-hydroxyquinoline-5-sulfonic acid.

REAGENTS

Package	Cat. No 4-247 (24-TRAY)	Cat. No 4-447 (36-TRAY)
1-Reagent	4 x 52 ml	6 x 22 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package.

Concentrations in the test

phosphate buffer (pH 7.5)	≤ 80 mmol/l
8-Hydroxyquinoline-5-sulfonic acid	≤ 6 mmol/l
arsenazo III	≤ 120 μmol/l
detergents, preservative	

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Do not freeze the reagent.
- Contaminated glassware is the greatest source of error. The use of disposable plastic ware is recommended. Glassware should be soaked for a few hours in 2M HCl solution and then thoroughly rinsed with distilled water.

SPECIMEN

Serum. Random or 24-hours urine.
Serum can be stored up to 7 days at 20-25°C or up to 3 weeks at 4-8°C. Samples frozen at -20°C can be stored up to 8 months. Discard contaminated specimens.

24-hours urine preparation: To prevent calcium salt precipitation specimens should be collected in 10 ml of 6M HCl. In case of presence of precipitants they can be solved by lowering pH of the urine to below 2.0. Prior to determination dilute the sample with 0.9% NaCl in the ratio of 1 to 1. Multiply the result by the dilution factor. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent is ready to use.
1-Reagent put on basic position in reagent tray.
For reagent blank deionized water is recommended.

Actions required:

- Biolis 24i Premium:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: LIPASE II GEN – CALCIUM ARSENAZO, RF – CALCIUM ARSENAZO. To avoid this effect follow the recommendations contained in the instruction 51_03_24_008_BIOLIS_24i_PREMIUM_CARRYOVER.

REFERENCE VALUES 4

serum	mg/dl	mmol/l
adult	8.6 – 10.3	2.15 – 2.57
random urine	mg/dl	mmol/l
male	0.9 – 37.9	0.225 – 9.47
female	0.5 – 35.7	0.125 – 8.92
24-hours urine	mg/24h	mmol/24h
adult	100 – 300	2.5 – 7.5

It is recommended for each laboratory to establish its own reference ranges for local population.

Calcium concentration 24-hours urine - calculation

calcium concentration in 24-hours urine [mg/24h]	=	calcium concentration in sample of 24-hours urine [mg/dl]	×	urine volume of 24-hours urine [dl/24h]
--	---	---	---	---

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum or CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine.

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended.

The calibration curve should be prepared with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Biolis 24i Premium. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**
0.06 mg/dl (0.015 mmol/l)
- LoD (Limit of Detection):**
0.09 mg/dl (0.023 mmol/l)
- LoQ (Limit of Quantitation):**
0.30 mg/dl (0.075 mmol/l)
- Linearity:**
up to 23 mg/dl (5.75 mmol/l)

For higher concentration, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 5 g/dl, ascorbate up to 62 mg/l, bilirubin up to 40 mg/dl, triglycerides up to 2000 mg/dl and magnesium up to 20 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	9.67	0.09	0.91
level 2	12.8	0.10	0.79
Reproducibility (day to day) n = 80	Średnia [mg/dl]	SD [mg/dl]	CV [%]
level 1	9.5	0.14	1.4
level 2	12.7	0.13	1.0

Method comparison

A comparison between calcium values determined at **Biolis 24i Premium** (y) and at **BECKMAN COULTER AU680 0** (x) using 30 serum samples gave following results:
 $y = 0.9595x + 0.3577$ mg/dl;
 $R = 0.991$ (R – correlation coefficient)

A comparison between calcium values determined at **Biolis 24i Premium** (y) and at **BECKMAN COULTER AU680** (x) using 39 urine samples gave following results:
 $y = 1.0338x - 0.2229$ mg/dl;
 $R = 1.000$ (R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	57	CaArs			
Auto Rerun SW					
ON					
Auto Rerun Range (Result)					
	ON	ON			
	Lower	Higher			
Serum	0.62	23			
Urine					
Plasma					
CSF					
Dialysis					
Other					
Auto Rerun Condition (Absorbance)					
Absorbance Range		Lower	ON		
		Higher	OFF		
Auto Rerun Condition (Prozone)					
Prozone Range		OFF			

- Prestige 24i Premium, Biolis 24i Premium, Sapphire 400 Premium

Item No.	57	Item Name	CaArs	Optical	
Data information					
Units	mg/dl				
Decimals	2				
Calibration					
Type	Linear2				
Std sample conc.					
Blank	0	#1	*	#2	*
#3		#4		#5	
#6					
Analysis					
Type	END method				
Main Wave Length	660nm				
Sub Wave Length	700nm				
Method	Arsenazo				
Correlation					
Y=	Slope	Intercept			
	1	0			

Item No.	57	Item Name	CaArs	Optical	
Aspiration					
Kind	Single				
Vol.					
	Kind	Vol.	Add	Units	
Sample	3	5		µl	
Reagent 1	300	10		µl	
Data Process					
Read		Start	End		
	Main	30	31		
	Sub	9	10		
Abs.Limit					
Low	-3	High	3		
Correction value					
Blank correction					
End Point Limit		2			
Linear Check (%)		90			
Prozone Check					
	Start	End	Limit (%)		
First					
Second			Low		
Blank value					
Water Blank					
Reaction Monitor					
0 Level Point	1				
Span	3				
Third mixing					
OFF					

Item No.	57	Item Name	CaArs	Optical	
Normal Range					
	Male		Female		
	Low	High	Low	High	
Serum	8.6	10.3	8.6	10.3	
Urine	0.9	37.9	0.5	35.7	
Plasma					
CSF					
Dialysis					
Other					
Panic Range					
	Male		Female		
	Low	High	Low	High	
Serum					
Urine					
Plasma					
CSF					
Dialysis					
Other					

Item No.	57	Item Name	CaArs	Optical	
Auto Rerun SW					
ON					
Auto Rerun Range (Conc.)					
	First Dil	Low	High		
		Re	Value	Dil	Re
Serum			0.3		23
Urine					
Plasma					
CSF					
Dialysis					
Other					
Auto Rerun Condition (Absorbance)					
		Lower	OFF		
		Higher	OFF		
Auto Rerun Condition (Prozone)					
OFF					
Dilution					
100:Di12					