CORMAY D-DIMER CONTROLS

(PL)

CORMAY D-DIMER CONTROLS
Nr kat. 4-459 4 x 0,5 ml

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CORMAY D-DIMER CONTROLS
Cat. No 4-459 4 x 0.5 ml

INTENDED USE

CORMAY D-DIMER CONTROLS is designed for use as control serum when measuring D-dimer concentration. Kit contains lyophilised control sera containing D-dimer at two levels (low – L and high – H). The exact concentration of D-dimer ($\mu\text{g/ml}$ FEU) is printed on the label.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

RECONSTITUTION

Carefully remove the rubber cap from vial and add 0.5 ml D-DIMER DILUENT to the vial. Recap and leave at room temperature for 15 minutes, inverting gently to dissolve contents.

STABILITY AND STORAGE

1. Lyophilised serum remains stable when stored at 2-10°C until expiry date given on the product label.
2. The constituents in reconstituted control serum are stable for 2 weeks at 2-10°C.

NOTES

1. This control serum has been tested for the HIV antibody, HBsAg and HCV and found to be non-reactive. However this material should be handled as though capable of transmitting infectious disease.
2. D-DIMER DILUENT is present in kits for determination of D-dimer concentration.
3. Do not mix different lot.
4. Do not use after expiration date.
5. D-DIMER DILUENT contains sodium azide (< 1 g/l) as a preservative. Avoid contact with skin and mucous membranes.

WASTE MANAGEMENT

Please refer to local legal requirements.

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