PRESTIGE 24i ALPHA 1-MICROGLOBULIN

Nr kat. 4-268, 4-487 (PL)

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Cat. No 4-268, 4-487 (EN)

INTENDED USE

Diagnostic kit for determination of α 1-microglobulin concentration used in automatic analysers: Prestige 24i, Biolis 24i and Prestige 24i Premium, Biolis 24i Premium.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

α 1-microglobulin (α Mi) is a low molecular weight glycoprotein (24-33 kD) which was initially isolated from the urine of patients with renal tubular disorders in 1975. It is mainly synthesized in the liver and is widely distributed in various body fluids.

The measurement of α Mi in serum and urine has been considered to be useful for the diagnosis of functional renal disorders, the assessment of the progress and prognosis of diseases.

METHOD PRINCIPLE

When an antigen-antibody reaction occurs between α Mi in a sample and anti- α Mi antibody which has been sensitized to latex particles, agglutination results. This agglutination is detected as an absorbance change (572 nm), with the magnitude of the change being proportional to the quantity of α Mi in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

REAGENTS

Package

	Cat. No 4-268 (24-TRAY)	Cat. No 4-487 (36-TRAY)
1-Reagent	1 x 20 ml	2 x 13 ml
2-Reagent	1 x 20 ml	2 x 13 ml

The reagents are stable up to the kit expiry date printed on the package when stored at 2-10°C. The reagents are stable for 6 weeks on board the analyser at 2-10°C.

Concentrations in the test

suspension of latex particles sensitized with (rabbit) anti- α Mi antibodies 0.25 w/v%
 glycine buffer solution
 preservative

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Reagent bottles should be shaken before use by gently inverting several times.

- After measurements are taken, reagent bottles should be capped and kept at 2-10°C. Do not interchange the caps of reagent bottles.
- Reagents with different lot numbers should not be interchanged or mixed.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum, plasma or urine.

If the test cannot be done immediately, the sample should be placed in a tightly sealable container and stored at -20°C. Repeated freezing and thawing should be avoided. Nevertheless, it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank 0.9% NaCl is recommended.

Application guide for Biolis 30i analyser is available on request.

REFERENCE VALUES ³

serum, plasma	10.0 – 30.0 mg/l
urine	1.0 – 5.0 mg/l

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, control serum for determination of α Mi, eg.: ROCHE or BIORAD.

For the calibration of automatic analysers systems it is recommended to use: CORMAY ALPHA 1-MGLOB CALIBRATORS (S) (Cat. No 4-286) - for determination in serum and CORMAY ALPHA 1-MGLOB CALIBRATORS (U) (Cat. No 4-285) - for determination in urine.

The calibration curve should be prepared every time the test is performed, with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser HITACHI 917. Results may vary if a different instrument is used.

- Analytical range:**
 1.5 – 200 mg/l (serum)
 0.3 – 50.0 mg/l (urine)

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 0.5 g/dl, NH₄Cl do 400 mg/dl, bilirubin up to 31 mg/dl do not interfere with the test in urine.

Precision

Repeatability (run to run) n = 20	Mean [mg/l]	SD [mg/l]	CV [%]
level 1	0.5	0.0	3.97
level 2	1.6	0.0	1.81
level 3	13.9	0.1	0.42

Method comparison

A comparison between CORMAY reagent (y) and another commercially available assay (x) using 55 serum samples gave following results:

$$y = 1.00 x + 2.83 \text{ mg/l;}$$

$$R = 1.00 \quad (R - \text{correlation coefficient})$$

A comparison between CORMAY reagent (y) and another commercially available assay (x) using 55 urine samples gave following results:

$$y = 1.00 x - 0.52 \text{ mg/l;}$$

$$R = 1.00 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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Item No.	49	Item Name	AMG				Optical
Data information							
Units	mg/l						
Decimals	1						
Analysis							
Type	END method						
Main Wave Length	570 nm						
Sub Wave Length	800 nm						
Method	Immuno						
Correlation							
Y=	Slope	1		X+	Intercept	0	
Item No.	49	Item Name	AMG				Optical
Aspiration							
Kind	Double						
Vol.	Kind	Vol.	Add	Units			
Sample		3	5	µl			
Reagent 1		160	10	µl			
Reagent 2		160	10	µl			
Blank value							
Water Blank							
Reaction Monitor							
0 Level Point	1						
Span	3						
Third mixing							
ON							
Item No.	49	Item Name	AMG				Optical
Normal Range							
	Male		Female				
	Low	High	Low	High			
Serum	10	30	10	30			
Urine	1	5	1	5			
Plasma	10	30	10	30			
CSF							
Dialysis							
Other							
Panic Range							
	Male		Female				
	Low	High	Low	High			
Serum							
Urine							
Plasma							
CSF							
Dialysis							
Other							
Item No.	49	Item Name	AMG				Optical
Auto Rerun SW							
OFF							
Auto Rerun Range (Conc.)							
	First Dil	Low			High		
		Re	Value	Dil	Re	Value	Dil
Serum							
Urine							
Plasma							
CSF							
Dialysis							
Other							
Auto Rerun Condition (Absorbance)							
		Lower	OFF				
		Higher	OFF				
Auto Rerun Condition (Prozone)							
OFF							
Dilution							
99:Dil1							

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