TRANSFERRIN

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INTENDED USE

Diagnostic kit for determination of transferrin concentration intended to be used in automatic analysers according to their user manual.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Transferrin (siderophilin) is a glycoprotein synthesised in the liver. Transferrin is the major plasma transport protein for iron. Its concentration correlates with the total iron-binding capacity (TIBC).

Evaluation plasma/serum transferrin levels is useful for the differential diagnosis of anaemia and for monitoring its treatment. Causes of decreased synthesis and low plasma levels may be chronic liver disease and malnutrition. High levels of transferrin occur in pregnancy and during estrogen administration.

METHOD PRINCIPLE

Transferrin, present in a sample, form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum is proportional to transferrin concentration in the sample.

REAGENTS

	CORMAY	OS-	B50-
	TRANSFERRIN	TRANSFERRIN	TRANSFERRIN
1-REAGENT	1 x 63 ml	1 x 40.5 ml	1 x 57 ml
2-REAGENT	1 x 18 ml	1 x 12.5 ml	1 x 17 ml

Buffer (1-REAGENT) stored at 2-25°C and antiserum (2-REAGENT) stored at 2-8°C are stable until expiry date printed on the package.

Reagent components

1-REAGENT: imidazole buffer, polymer, inorganic salt, buffer, preservative.

2-REAGENT: anti-human transferrin antiserum, inorganic salt, buffer, preservative.

WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- Store closed.
- Do not freeze the reagents.
- Nanoparticle-based reagents can settle over time. It may be necessary to delicately mix by repeated turning.

- Human-origin products have been tested for HBsAg and antibodies to HIV 1, HIV 2, HCV and found to be non-reactive. However this material should be handled as potentially infectious.
- EUH210 Safety data sheet available on request.

ADDITIONAL EQUIPMENT

- automated clinical chemistry analyser capable of accommodating two-reagent assays;
- general laboratory equipment;

SPECIMEN ⁶

Serum or plasma. Nonhemolyzed, nonlipaemic, fresh serum is recommended.

Serum can be stored up to 3 days at 4°C or up to 6 months at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-REAGENT and 2-REAGENT are ready to use.

Applications for analysers are available on request.

REFERENCE VALUES ⁶

newborns	1.30 – 2.75 g/l
children (3 months – 16 years)	2.03 – 3.60 g/l
adults F	2.50 – 3.80 g/l
M	2.15 – 3.65 g/l

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples the CORMAY IMMUNO-CONTROL III (Cat. No 4-291).

For the calibration of automatic analysers the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended. **Calibrators and 0.9% NaCl** should be used for calibration.

The calibration curve should be prepared with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers: Hitachi or Biolis 24i Premium. Results may vary if a different instrument is used.

- **LOQ:** 0.076 g/l
- **Linearity:** up to 8 g/l
- **Specificity / Interferences**
 Haemoglobin up to 0.32 g/dl, bilirubin up to 29.5 mg/dl, triglycerides up to 312 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4 g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

Precision

Repeatability (run to run) n = 30	Mean [g/l]	SD [g/l]	CV [%]
level 1	1.44	0.02	1.12
level 2	2.52	0.04	1.52
level 3	3.47	0.06	1.83
Reproducibility (day to day) n = 66	Mean [g/l]	SD [g/l]	CV [%]
level 1	1.49	0.11	4.56
level 2	2.51	0.12	4.06
level 3	3.51	0.12	1.64

Method comparison

A comparison between transferrin values determined at **Biolis 24i Premium** (y) and at **ADVIA 1650** (x) using 36 samples gave following results:

$$y = 1.0074x - 0.1457 \text{ g/l};$$

$$R = 0.9780 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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