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Data wydania: 10.2019

CORMAY LIPID CONTROL 1 CORMAY LIPID CONTROL 2

(EN)

CORMAY LIPID CONTROL 1

Cat. No 5-179 3 x 1 ml

CORMAY LIPID CONTROL 2

Cat. No 5-180 3 x 1 ml

INTENDED USE

CORMAY LIPID CONTROL 1 and CORMAY LIPID CONTROL 2 are designed for routine quality control of direct HDL, direct LDL, cholesterol, and triglycerides concentrations. Components concentrations are available at 2 levels.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

RECONSTITUTION

Open the vial carefully and add exactly 1 ml of distilled water, close bottle and let stand out for 30 min. Ensure that all traces of dry material are dissolved by swirling gently. Avoid foam forming.

STABILITY AND STORAGE

1. Lyophilised serum remains stable when stored at 2-8°C until expiry date given on the product label.
2. The components in reconstituted control serum are stable for:

Room temperature	8 hours
Temp. 2-8°C	7 days
Temp. -20°C	30 days

NOTES

- This serum has been tested for HBsAg, anti- HIV-1, anti- HIV-2 antibodies, anti- HCV antibodies and found to be non-reactive. However this material should be handled as thought capable of transmitting infectious disease.

ASSIGNED VALUES

The assigned values are the means obtained by using various clinical chemistry systems.

It is recommended that each laboratory establish its own means and acceptable ranges.

WASTE MANAGEMENT

Please refer to local legal requirements.

Date of issue: 10. 2019.

CORMAY LIPID CONTROL 1 CORMAY LIPID CONTROL 2

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