CORMAY RPR CARBON

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LITERATURA

- George P. Schimid. Current Opinion in Infectious Diseases 1994; 7: 34-40.
- Sandra A Larsen et al. Clinical Microbiology Reviews 1995; 8 (1): 1-21.
- Sandra Larsen et al. A manual of Test for Syphilis American Public Health Association 1990: 1-192.
- Joseph Earle Moore et al. Gastrointestinal Haemorrhage 1952; 150(5): 467-473.
- Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACC Press, 1995.

Data wydania: 05. 2018.

CORMAY RPR CARBON

Kit name	Kit size	Cat. No
CORMAY RPR Carbon 150	150 tests	6-255
CORMAY RPR Carbon 500	500 tests	6-256

INTENDED USE

Agglutination test for detection of plasma reagins recommended for the qualitative and semi-quantitative manual assay. The reagents must be used only for in vitro diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Reagins are a group of antibodies against some components of the damage tissues from patients infected by *Treponema pallidum*, the agent which causes the syphilis. This microorganism produces some damage to the liver and heart, releasing some tissue fragments. Immunological patient system reacts producing reagins-antibodies against these fragments. Determination of plasma reagins is non-specific test for diagnosis of *Treponema pallidum* and is useful to follow up the antibiotic therapy.

METHOD PRINCIPLE

Carbon particles coated with a lipid complex (cardiolipin, lecithin, cholesterol) are agglutinated when mixed with samples containing reagins.

REAGENTS

Package

	CORMAY RPR Carbon 150	CORMAY RPR Carbon 500
RPR-Carbon	1 x 3 ml	2 x 5 ml
RPR-Control (+) (red cap)	1 x 1 ml	1 x 1 ml
RPR-Control (-) (blue cap)	1 x 1 ml	1 x 1 ml
Stirrers	3 x 25 pcs.	10 x 25 pcs.
Slides (8 circles each)	1 x 21 pcs.	3 x 21 pcs.
Dispensing vial	1 x 1 pcs.	1 x 1 pcs.
Needle	1 x 1 pcs.	1 x 1 pcs.

Reagent preparation and stability

The reagents are ready to use.

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. Do not freeze.

Concentrations in the test

carbon particles coated with a lipid complex	10 mmol/l
artificial serum solution	150 mmol/l
animal serum solution	150 mmol/l
preservatives	

Warnings and notes

- RPR-Carbon contains sodium azide (< 0.1%) as a preservative. Avoid contact with skin and mucous membranes.
- High temperature may cause test components to dry on the slide giving an agglutination aspect that can be

interpreted as false positive results. It is recommended to place the slide under a humidifying cover.

- CORMAY RPR Carbon test is non-specific for syphilis. All reactive samples should be retested with treponemic methods such as TPHA and FTA-Abs to confirm the results.
- A negative result does not exclude a diagnosis of syphilis.
- False positive results have been reported in diseases such as infectious mononucleosis, viral pneumonia, toxoplasmosis, pregnancy and autoimmune diseases.
- Diagnosis should only be made after taking clinical symptoms and the results of other tests into consideration.
- RPR-Control (+) meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients: RPR-Control (+) contains citric acid monohydrate.

Warning:

H319 Causes serious eye irritation.
 P280 Wear protective gloves/protective clothing/eye protection/face protection.
 P501 Dispose of the contents/containers in accordance with the current legislation on

waste treatment

ADDITIONAL EQUIPMENT

- mechanical rotator with adjustable speed at 80-100 r.p.m.
- general laboratory equipment.

SPECIMEN

Serum or plasma. Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged before testing. Do not use highly hemolized or lipemic samples.

It is recommended to perform the assay with freshly collected samples.

PROCEDURE

Qualitative method

- Allow the reagents and samples to reach room temperature. The sensitivity of the test may be reduced at low temperatures.
- Place 50 µl of the sample and one drop of each positive and negative controls into separate circles on the slide test.
- Swirl the RPR-Carbon reagent gently before using and add (using dispensing vial and needle) one drop (20 µl) next to the sample to be tested.
- Mix the drops with a stirrer, spreading them over the entire surface of the circle. Use different stirrers for each sample.
- Place the slide on a mechanical rotator at 80-100 r.p.m. **for 8 minutes.** False positive results could appear if the test is read later than 8 minutes.

Reading and interpretation

Examine macroscopically the presence or absence of visible agglutination immediately after removing the slide from the rotator.

The presence of agglutination indicates of plasma reagins.

Semi-quantitative method

- Make serial two fold dilutions of the sample in 0.9% NaCl solution.
- Proceed for each dilution as in the qualitative method.

Reading and interpretation

The result (titer), in the semi-quantitative method, is defined as the highest dilution showing a positive result.

QUALITY CONTROL

Positive and negative controls are recommended to monitor the performance of the procedure, as well as a comparative pattern for a better result interpretation.

PERFORMANCE CHARACTERISTICS

- Analytical sensitivity:** according to recommendations of International Reference WHO (1st Standard Human Syphilitic Serum, ref. 05/132).
- Prozone effect:** no prozone effect up to titers $\geq 1/128$.
- Diagnostic sensitivity:** 100%.
- Diagnostic specificity:** 100%.
- Interferences:** Haemoglobin up to 10 g/l, bilirubin up to 20 mg/dl, triglycerides up to 10 g/l and RF up to 300 IU/ml do not interfere with the test.

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

- George P. Schmid. Current Opinion in Infectious Diseases 1994; 7: 34-40.
- Sandra A Larsen et al. Clinical Microbiology Reviews 1995; 8 (1): 1-21.
- Sandra Larsen et al. A manual of Test for Syphilis American Public Health Association 1990: 1-192.
- Joseph Earle Moore et al. Gastrointestinal Haemorrhage 1952; 150(5): 467-473.
- Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACC Press, 1995.

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