## CORMAY VDRL

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### LITERATURA

- George P. Schmid. Current Opinion in Infectious Diseases 1994; 7: 34-40.
- Sandra A Larsen et al. Clinical Microbiology Reviews 1995; 8 (1): 1-21.
- Sandra Larsen et al. A manual of Test for Syphilis American Public Health Association 1990: 1-192.
- Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACC Press, 1995.

Data wydania: 08. 2020.

## CORMAY VDRL

Kit name	Kit size	Cat. No
CORMAY VDRL 250	250 tests	6-257
CORMAY VDRL 1500	1500 tests	6-258

### INTENDED USE

Confirming haemagglutination test for detection of plasma regains recommended for the qualitative and semi-quantitative manual assay.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Reagins are a group of antibodies against some components of the damage tissues from patients infected by *Treponema pallidum*, the agent which causes the syphilis. This microorganism produces some damage to the liver and heart, releasing some tissue fragments. Immunological patient system reacts producing reagins-antibodies against these fragments. Determination of plasma reagins is non-specific test for diagnosis of *Treponema pallidum* and is useful to follow up the antibiotic therapy. The assay is useful to follow the antibiotic therapy answer.

### METHOD PRINCIPLE

The antigen suspension, a lipid complex (cardiolipin, lecithin, cholesterol), is agglutinated when mixed with samples containing regains of patient affected by syphilis.

### REAGENTS

#### Package

	CORMAY VDRL 250	CORMAY VDRL 1500
VDRL- stabilized Antigen	1 x 5 ml	1 x 30 ml
VDRL-Control (+) (red cap)	1 x 1 ml	1 x 1 ml
VDRL-Control (-) (blue cap)	1 x 1 ml	1 x 1 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. Keep vials in vertical position. If the position is changed, gently mix to dissolve aggregates that may be present.

#### Preparation:

The reagents are ready to use.

#### Concentrations in the test

cardiolipin	0.3 g/l
lecithin	2.1 g/l
cholesterol	9.0 g/l
phosphate buffer pH 7.0	1.5 mmol/l
artificial serum solution	
animal serum solution	
preservatives	

#### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Do not freeze.
- CORMAY VDRL test is non-specific for syphilis. All reactive samples should be retested with treponemic

methods such as TPHA and FTA-Abs to confirm the results.

- A negative result does not exclude infection of *Treponema pallidum*.
- False positive results have been reported in diseases such as infectious mononucleosis, viral pneumonia, toxoplasmosis, pregnancy and autoimmune diseases.
- Diagnosis should only be made after taking clinical symptoms and the results of other tests into consideration.
- VDRL-Control (+) meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

#### Warning

H319 Causes serious eye irritation.  
 P280 Wear protective gloves/protective clothing/eye protection/face protection.

P501 Dispose of the contents/containers in accordance with the current legislation on

waste treatment

#### ADDITIONAL EQUIPMENT

- mechanical rotator with adjustable speed at 160-180 r.p.m.;
- glass slides;
- light microscope (10x);
- pipettes 50 µl
- general laboratory equipment.

#### SPECIMEN

Fresh serum or plasma. Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged before use. Do not use highly hemolized or lipemic samples.

#### PROCEDURE

##### Qualitative method

- Allow the reagents and samples to reach room temperature. The sensitivity of the test may be reduced at low temperatures.
- Place 50 µl of the sample and one drop of each positive and negative controls into separate circles on the slide test.
- Gently stir the VDRL- Reagent and add one drop (20 µl) to each test field.
- Place the slide on a mechanical rotator at 180 r.p.m. **for 4 minutes**. False positive results could appear if the test is read later than four minutes.

#### Reading and interpretation

Examine the presence or absence of agglutination immediately after rotation using the light microscope (10x). The presence of agglutination indicates of plasma reagins.

Agglutination	Reading	Report
Medium or large clumps	R	Reactive
Small clumps	W	Weakly reactive
No clumping or very slight "roughness"	N	Non reactive

#### Semi-quantitative method

- Make serial two fold dilutions of the sample in 0.9% NaCl solution.
- Proceed for each dilution as in the qualitative method.

#### Reading and interpretation

The result (titer), in the semi-quantitative method, is defined as the highest dilution showing a positive result.

#### QUALITY CONTROL

Positive and negative controls are recommended to monitor the performance of the procedure, as well as a comparative pattern for a better result interpretation. All result different from negative control result will be considered as positive.

#### PERFORMANCE CHARACTERISTICS

- Analytical sensitivity:** according to recommendations of International Reference WHO (1<sup>st</sup> Standard Human Syphilitic Serum, ref. 05/132).
- Prozone effect:** no prozone effect up to titers  $\geq 1/128$ .
- Diagnostic sensitivity:** 100 %.
- Diagnostic specificity:** 100 %.
- Interferences:**  
 Haemoglobin up to 10 g/l, bilirubin up to 20 mg/dl, triglycerides up to 10 g/l and RF up to 300 IU/ml do not interfere with the test.

#### WASTE MANAGEMENT

Please refer to local legal requirements.

#### LITERATURE

- George P. Schmid. Current Opinion in Infectious Diseases 1994; 7: 34-40.
- Sandra A Larsen et al. Clinical Microbiology Reviews 1995; 8 (1): 1-21.
- Sandra Larsen et al. A manual of Test for Syphilis American Public Health Association 1990: 1-192.
- Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACC Press, 1995.

Date of issue: 08. 2020.