CORMAY ETHANOL

| | (PL) Nr kat. |
|--------------------|-----------------|
| Nazwa zestawu | |
| CORMAY ETHANOL 30 | 2-178 |
| CORMAY ETHANOL 60 | 2-378 |
| CORMAY ETHANOL 500 | 6-378 |

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LITERATURA

- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, Toronto: WB Saunders Co., 1171-1175 (1994).
- William H., Porter Ph.D. Clinical Toxicology. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 922-923.
- Bonnichsen, R.K., Theorell, H., An Enzymatic Method for the Micro Determination of Ethanol, Scand. J. Clin. Lab. Invest. 3, 58 (1951).
- Lundquist, F., The Determination of Ethyl Alcohol in Blood and Tissues., Methods of Biochemical Analysis, Vol III, Editor, D. Glick, Interscience Publishers, New York, 1957 p. 217-251.
- CLSI Guidelines and Standards, Clinical Laboratory Tests, AACC Press, Washington, 1990.
- Preeti Dalawari: emedicine.medscape.com./article/2090019-overview#a3.
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders, 923, (1999).
- Krajowa Izba Diagnostów Laboratoryjnych, Oznaczenie etanolu w krwi i moczu pacjentów w laboratorium medycznym, Warszawa 2016. p.18-19.

Data wydania: 09.2020.

CORMAY ETHANOL

| | |
|--------------------|----------------------|
| Kit name | (EN) |
| CORMAY ETHANOL 30 | Cat. No 2-178 |
| CORMAY ETHANOL 60 | 2-378 |
| CORMAY ETHANOL 500 | 6-378 |

INTENDED USE

Diagnostic kit for determination of ethyl alcohol concentration used both for manual assay and in several automatic analysers.


The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Ethanol after consumption is absorbed in alimentary tract and transported with blood to whole organism. A little amount of ethanol is excreted with urine or exhaled with air, but majority is metabolized by liver to acetic aldehyde, then to acetate, finally to carbon dioxide and water. Ethanol consumption is often cause of various types of accidents, and ethanol poisoning may be cause of decease.

METHOD PRINCIPLE

Enzymatic method with alcohol dehydrogenase.

The rate of absorbance changes measured at $\lambda=340$ nm or $\lambda=380$ nm is proportional to the quantity of ethyl alcohol in the sample.

REAGENTS

Package

| | | | |
|-----------|-------------------|-------------------|--------------------|
| | CORMAY | CORMAY | CORMAY |
| | ETHANOL 30 | ETHANOL 60 | ETHANOL 500 |
| 1-Reagent | 6 x 30 ml | 6 x 60 ml | 4 x 500 ml |

The reagent when stored tightly capped at 2-8°C is stable up to expiry date printed on the package. On board stability of the reagents depends on type of the analyser used for analysis.

Concentrations in the test

| | |
|----------------------------|-------------|
| buffer (pH 7,3) | |
| NADanalog | 1.35 mmol/l |
| alcohol dehydrogenase | >203 kU/l |
| stabilizers, preservatives | |

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Avoid contact with skin and mucous membranes.
- Concentration of ethyl alcohol in serum or plasma is approximately 1.2 times higher than in whole blood.⁶
- Turbidity would indicate contamination or deterioration of reagent.

SPECIMEN

Serum or plasma (heparin, EDTA or citrate) free from hemolysis; urine.

Samples may be stored tightly closed in temp. 2-8°C up to 3 days. Urine samples should be stored with as little dead air space as possible.

Nevertheless it is recommended to perform the assay with freshly collected samples!

- Do not use ethanol to swab venipuncture site and to clean or sterilize glassware or other equipment used to collect the specimen or perform the assay.
- Samples should be at room temperature for testing.
- Before opening mix samples by inverting the container.

PROCEDURE

Applications for analysers are available on request.

Manual procedure

| | |
|-------------|--------|
| wavelength | 340 nm |
| temperature | 37°C |
| cuvette | 1 cm |

Pipette into the cuvettes:

| | test (T) | calibrator (C) |
|-----------|--------------|----------------|
| 1-Reagent | 1000 μ l | 1000 μ l |

Bring up to the temperature of determination. Then add:

| | | |
|------------|------------|------------|
| calibrator | - | 10 μ l |
| sample | 10 μ l | - |

Mix and incubate at adequate temperature. After 1 min. read the absorbance of the calibrator (C) and test (T) against air or water. Repeat the reading after exactly 1, 2 and 3 minutes. Calculate the mean change of absorbance per minute ($\Delta A/\text{min.}$).

Calculation

$$\text{alcohol concentration [mg/dl]} = \frac{\Delta A/\text{min}(T)}{\Delta A/\text{min}(C)} \times \text{calibrator concentration}$$

REFERENCE VALUES⁷

Ethanol consumption cause following symptoms:

| blood ethanol concentration (mg/dl) | clinical symptoms* |
|-------------------------------------|--|
| 30 – 120 | euphoria; sociability; talkativeness; diminution of attention, judgment and control; some sensory-motor impairment |
| 90 – 250 | emotional instability; loss of critical judgment; impairment of perception, memory and comprehension; increased reaction time; sensory-motor incoordination; drowsiness; |
| 180 – 300 | disorientation; mental confusion; dizziness; staggering gait; slurred speech; disturbances of vision; increased pain threshold; apathy; lethargy; |
| 250 – 400 | general inertia; markedly decreased response to stimuli; vomiting; incontinence of urine and feces; sleep or stupor; |

| | |
|-----------|---|
| 350 – 500 | complete unconsciousness; coma; depressed or abolished reflexes; impairment of circulation and respiration; possible death. |
|-----------|---|

***Effects are dependent on individual use, chronicity, tolerance and other factors.⁸**

To calculate the concentration of ethyl alcohol in blood the result received for serum or plasma should be divided by 1.2 factor. Otherwise, doctor should be informed that the result is overstated and requires a proper interpretation.⁸

Please refer to the local regulations concerning legal ethanol limits for drivers.

Conversion factor of the blood ethanol concentration: 100 mg/dl = 1 ‰.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY AMMONIA/ETHANOL CONTROLS (Cat. No 5-163).

For the calibration of manual determinations CORMAY ETHANOL CALIBRATOR 100 (Cat. No 5-106) is recommended.

For the calibration of automatic analysers the CORMAY ETHANOL CALIBRATORS (Cat. No 5-105) is recommended. Two levels of calibrator and deionised water or 0.9% NaCl (as a calibrator 0) should be used for calibration.

Calibration stability depends on type of analyser used for analysis. The calibration curve should be prepared with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers Biolis 24i Premium or Beckman CX. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**
1.2 mg/dl (0.26 mmol/l)
- LoD (Limit of Detection):**
2.9 mg/dl (0.63 mmol/l)
- LoQ (Limit of Quantitation):**
7.0 mg/dl (1.52 mmol/l)
- Linearity:**
up to 800 mg/dl (173.6 mmol/l)

Specificity / Interferences

Haemoglobin up to 0.6 g/dl, bilirubin up to 40 mg/dl, intralipid up to 1000 mg/dl, LDH up to 10 000 U/l (in the presence of 200 mg/dl lactate) do not interfere with the test.

Acetone up to 2000 mg/dl, butanol up to 200 mg/dl, ethylene glycol up to 2000 mg/dl, methanol up to 2000 mg/dl and isopropanol up to 2000 mg/dl do not interfere with the test, which confirms that assay is specific to ethanol.

Precision

| Repeatability (run to run) n = 20 | Mean [mg/dl] | SD [mg/dl] | CV [%] |
|-----------------------------------|--------------|------------|--------|
| level 1 | 32.0 | 1.17 | 3.66 |
| level 2 | 83.0 | 1.70 | 2.05 |
| level 3 | 199.9 | 4.28 | 2.14 |

| Reproducibility (day to day) n = 20 | Mean [mg/dl] | SD [mg/dl] | CV [%] |
|-------------------------------------|--------------|------------|--------|
| level 1 | 38.6 | 1.17 | 3.03 |
| level 2 | 91.05 | 1.29 | 1.41 |
| level 3 | 244.47 | 4.54 | 1.86 |

Method comparison

A comparison between CORMAY reagent (y) and commercially available assay (x) using 40 serum samples. gave following results:

$$y = 1.064 x - 2.33 \text{ mg/dl;}$$

$$R = 0.9994 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry. Toronto: WB Saunders Co., 1171-1175 (1994).
- William H. Porter Ph.D. Clinical Toxicology. In: Burtis C.A. Ashwood ER. editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 922-923.
- Bonnichsen. R.K.,Theorell. H., An Enzymatic Method for the Micro Determination of Ethanol. Scand. J. Clin. Lab. Invest. 3. 58 (1951).
- Lundquist. F., The Determination of Ethyl Alcohol in Blood and Tissues.. Methods of Biochemical Analysis. Vol III. Editor. D. Glick. Interscience Publishers. New York. 1957 p. 217-251.
- CLSI Guidelines and Standards. Clinical Laboratory Tests. AACC Press. Washington. 1990.
- Preeti Dalawari: emedicine.medscape.com/article/2090019-overview#a3.
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia. PA: WB Saunders. 923. (1999).
- Polish National Chamber of Laboratory Diagnosticians. Oznaczenie etanolu w krwi i moczu pacjentów w laboratorium medycznym. Warszawa 2016. p.18-19.

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