ACCENT-200 CYSTATIN C

Nr kat. **7-200** (PL)

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Cat. No **7-200** (EN)

INTENDED USE

Diagnostic kit for determination of cystatin C concentration intended to use in automatic analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT S120, ACCENT MC240, ACCENT M320, ACCENT 400 and ACCENT Neo200.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Cystatin C is a low molecular weight protein (13kD), one of the cysteine proteinases inhibitors. Cystatin C is produced in all nucleated cells and secreted into extracellular space at constant rate. Cystatin C molecule stability and dependence of its concentration solely on GFR (Glomerular Filtration Rate) decide about high diagnostic efficiency of cystatin C determination. Cystatin C level is not affected by muscle mass or diet and its increase is observed even at slight GFR reduction. Clinical applications of cystatin C are for monitoring GFR in children, elderly patients, patients with potentially nephrotoxic drug therapy, for assessment of renal transplantation status, for kidney function monitoring in acute and chronic kidney diseases including diabetic nephropathy.

METHOD PRINCIPLE

Turbidimetric method. An antigen-antibody reaction occurs between cystatin C and antibodies coated on polystyrene particles and immuno-complexes are formed. The change of turbidity is related to the quantity of cystatin C in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

REAGENTS

Package

1-Reagent 2 x 15 ml
2-Reagent 2 x 4 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 9 weeks on board the analyser at 2-10°C.

Concentrations in the test

suspension of polystyrene particles coated with anti-cystatin C antibodies 8.5 g/l
MOPS buffer: 8 g/l
[3-(N-morpholino)-propanesulfonic acid]
gentamycin 12.5 mg/l
amphotericin B 1.25 mg/l
preservative

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- The reagents contain antibiotics and must be handled with due cautions.
- Mix samples well before analysing.
- The samples can be shipped without special cooling and must then be analysed within 14 days after shipment.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum, plasma collected on heparin or EDTA.

Samples are stable 14 days at 8-25°C and 21 days if stored at 2-8°C. For longer storage samples should be frozen at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

For reagent blank 0.9% NaCl is recommended.

REFERENCE VALUES⁹

serum		mg/l
newborns	< 1 month	1.49 – 2.85
children	< 5 months	1.01 – 1.92
	< 12 months	0.75 – 1.53
	12-24 months (F)	0.60 – 1.20
	12-24 months (M)	0.77 – 1.85
adults (F)	2-19 years	0.62 – 1.11
		0.61 – 1.05
adults (M)		0.71 – 1.21

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY CYSTATIN C CONTROLS (Cat. No 4-460) with each batch of samples.

For the calibration of automatic analysers ACCENT-200, ACCENT-200 II GEN, ACCENT S120, ACCENT MC240, ACCENT M320 the CORMAY CYSTATIN C CALIBRATORS (Cat. No 5-185) is recommended. 0.9% NaCl should be used as a calibrator 0.

The calibration curve should be prepared every 9 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers ACCENT-200 and/or Modular P and ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity:**
0.37 mg/l – Modular P
0.27 mg/l – ACCENT MC240

Linearity:

up to 8.39 mg/l – Modular P
up to 8 mg/l – ACCENT MC240

Specificity / Interferences

Haemoglobin up to 0.7 g/dl, bilirubin up to 800 mg/l, triglycerides up to 14 g/l and ascorbate up to 300 mg/l do not interfere with the test.

Precision

Repeatability (run to run)		Mean [mg/l]	SD [mg/l]	CV [%]
Modular P n=20	level 1	1.92	0.02	1.1
ACCENT MC240 n=20	level 1	0.95	0.02	2.05
	level 2	3.76	0.03	0.79
Reproducibility (day to day)		Mean [mg/l]	SD [mg/l]	CV [%]
Modular P n=20	level 1	0.88	0.02	2.73
	level 2	5.25	0.23	4.38
ACCENT MC240 n=80	level 1	1.0	0.07	7.0
	level 2	3.5	0.11	3.1

Method comparison

A comparison between cystatin C values determined at ACCENT-200 (y) and at BN ProSpec (x) using 23 serum samples gave following results:

$$y = 0.7954x + 0.1113 \text{ mg/l}$$

$$R = 0.9986 \quad (R - \text{correlation coefficient})$$

A comparison between cystatin C values determined at ACCENT MC240 (y) and at BS-800 (x) using 65 serum samples gave following results:

$$y = 1.0858x - 0.1451 \text{ mg/l}$$

$$R = 0.992 \quad (R - \text{correlation coefficient})$$

A comparison between cystatin C values determined at ACCENT MC240 (y) and at BS-800 (x) using 60 plasma samples gave following results:

$$y = 1.0152x - 0.0385 \text{ mg/l}$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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Calibration Rule

Rule	Spline
Sensitivity	1
Replicates	1
Interval (day)	0
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT-200 II GEN

Parameters

Test Name	CysC	R1	250
Test No	69	R2	50
Full Name	Cystatin C	Sample Volume	3
Reference No	69	R1 Blank	
Analy. Type	Endpoint	Mixed Reag. Blank	
Pri. Wave.	546 nm	Concentration	
Secon. Wave.		Linearity Limit	
Trend	Increase	Substrate Limit	
Reac. Time	0 40	Factor	
Incuba. Time	4	<input type="checkbox"/> Prozone check	
Unit	mg/l	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.01	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

Calibration Rule

Rule	Spline
Sensitivity	1
Replicates	1
Interval (day)	0
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT S120

Chem	CYS C	No.	069	Sample Type	SERUM
Chemistry	CYSTATIN C	Print name	CYS C	Reaction Direction	positive
Reaction Type	Endpoint	Sec Wave		Decimal	0.01
Pri Wave	546nm	Incubation Time	18	Reaction Time	18 20
Unit	mg/L	Blank Time	-3 -1	Standard	3 μL
Sample Vol	3 μL	Aspirated	20 μL	Diluent	180 μL
Decreased	3 μL	Increased	20 μL	Reagent Vol	R1 180 μL
Sample Blank	<input type="checkbox"/>	Auto Rerun	<input type="checkbox"/>	R2	30 μL
Linearity range (Standard)	0.3 8	Linearity Limit		Substrate Depletion	
Linearity Range (Decreased)		Linearity Range (Increased)		Mixed Blank Abs	-40000 40000
R1 Blank Abs	-40000 40000	Blank Response	-40000 40000	On-board Stability	Day(s)
Twin Chemistry		Prozone Check	<input type="checkbox"/>	Reagent Alarm Limit	
Q1		Q2		Q3	
Q5		Q6		V3	
Sample Pretreatment		Control Pretreatment		Calibrator Pretreatment	
Pretreat Sample Vol	μL	Pretreat Sample Vol	μL		
CALIBRATION SETTINGS	Math model Spline	Factor		Replicates	2
AUTO CALIBRATION		Bottle Changed	<input type="checkbox"/>	Lot Changed	<input type="checkbox"/>
		Cal Time			
ACCEPTANCE LIMITS	Slope Diff	SD		Repeatability	40000
	Deter Coeff				