ACCENT-200 CHOL

Nr kat. **7-204** (PL)

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LITERATURA

- Siedel J., Hägele E.O., Zigenhorn J., Wahlefeld A.W.: Clin. Chem. 29, 6 (1983).
- Tel R.M., Berends G.T.: J Clin.Chem. Clin. Biochem. 18, 10 (1980).
- Rautela G.S., Liedtke R.J.: Clin. Chem. 24. 1 (1978)
- Schettler G., Nussel E.: Arbeitsmed. Sozialmed. Preventivmed. 10, 25 (1975).
- Richmond W.: Clin. Chem. 19, 1350 (1973).
- Roeschla P., Bernt E., Gruber W.: J. Clin. Chem. Biochem. 12, 403 (1974).
- Tietz N.W., ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA: WB Saunders, 130 (1995).
- Trinder P.: Ann. Clin. Biochem. 6, 24 (1969).
- Dembińska-Kieć A., Naskalski J.W., Solnica B.: Diagnostyka laboratoryjna z elementami biochemii klinicznej. Wyd. IV, str. 670, 875 (2017).
- Cybulska B, Szostak WB, Filipiak KJ, et al. Polish Forum for Prevention Guidelines on Dyslipidaemia: update 2016. Kardiol Pol. 2017; 75(2): 187–190, doi: 10.5603/KP.2017.0031.
- Piepoli MF, Hoes AW, Agewall S, et al. 2016 European Guidelines on cardiovascular disease prevention in clinical practice: The Sixth Joint Task Force of the European Society of Cardiology and Other Societies on Cardiovascular Disease Prevention in Clinical Practice (constituted by representatives of 10 societies and by invited experts) Developed with the special contribution of the European Association for Cardiovascular Prevention & Rehabilitation (EACPR). Eur Heart J. 2016; 37(29): 2315–2381, doi: 10.1093/eurheartj/ehw106, indexed in Pubmed: 27222591.
- NCEP Expert Panel. Arch Inter Med (148), 36-69, (1988).
- Jacobso D.S., DeMott W.R., Grady H.J., et. al., ed., Laboratory Tests Handbook, 4th ed., Hudson, Lexi-Comp, 143, (1996).

Data wydania: 06. 2023.

ACCENT-200 CHOL

Cat. No **7-204** (EN)

INTENDED USE

Diagnostic kit for determination of total cholesterol concentration intended to use in automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, BS-120, ACCENT 400 and ACCENT Neo200.


The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Cholesterol is essential structural component of cell membranes and precursor of bile acids and all steroids hormones. This is why cholesterol has enormous significance for organism normal functioning. But there is also well established association between blood cholesterol concentration and coronary heart disease. Measurement of cholesterol serum level is valuable in prevention and monitoring cardiovascular disease. This determination is useful also for evaluation of intestine absorption, liver and gallbladder function.

METHOD PRINCIPLE

Enzymatic, colorimetric method with cholesterol esterase and cholesterol oxidase.

The colour intensity is proportional to the cholesterol concentration.

REAGENTS

Package

1-REAGENT 4 x 35 ml

The reagent, stored at 2-8°C is stable up to expiry date printed on the package. The reagent stored on board of the analyzer at 2-10°C is stable for 10 weeks (ACCENT-200) or for 12 weeks (ACCENT MC240).

Concentrations in the test

Good's buffer (pH 6.4)	< 120 mmol/l
phenol	< 6 mmol/l
4-aminoantipyrine	< 0.4 mmol/l
cholesterol esterase (CHE)	< 4 µkat/l
cholesterol oxidase (CHO)	< 5 µkat/l
peroxidase (POD)	< 24 µkat/l
preservatives, detergents, stabilizers	

Warnings and notes

- Protect from direct sunlight and avoid contamination!

- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Warning

H319 Causes serious eye irritation.
P280 Wear protective gloves, protective clothing, eye protection or face protection.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing

SPECIMEN

Serum, EDTA or heparinized plasma (recommended: heparine lithium, sodium or ammonium salt) free from hemolysis. Blood should be collected only if the patient has been fasting for minimum of 12 hours. Before blood collection patient should stay in rest position for about 30 minutes. Venous blood is recommended for cholesterol measurement. Serum should be separated from red blood cells as soon as possible after blood collection.

Plasma cholesterol values have been reported to be 3% to 5% lower than serum cholesterol values.

Serum and plasma can be stored up to 3 days at 2-8°C or 6 months at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-REAGENT is ready to use.

Deionized water is recommended as a reagent blank. For analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S and BS-120 it is recommended to determine the reagent blank during each calibration. Deionized water should be used as reagent blank. When performing calibration, the task type **Calib+Rgt.Blk** should be selected.

Actions required:

When performing assays in analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S and BS-120, there is a probability of **cross-contamination** affecting the tests results: CHOL - LIPASE, CREATININE - CHOL, ALAT - CHOL, ASAT - CHOL, CK - CHOL, CHOL - MG. To avoid this effect follow the recommendations contained in the instruction 51_03_24_001_ACCENT-200_CARRYOVER.

REFERENCE VALUES

serum / plasma	mg/dl	mmol/l
children ⁹ < 1year	66.1 – 228.5	1.71 – 5.91
1 – 19 years	111.4 – 202.2	2.88 – 5.23
adults ^{9,10,11}	< 190	< 5.00

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analysers the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 10 weeks (ACCENT-200) or every 12 weeks (ACCENT MC240), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers ACCENT-200 and ACCENT MC 240. Results may vary if a different instrument or a manual procedure is used.

Sensitivity

2.85 mg/dl (0.07 mmol/l) - ACCENT-200
7.6 mg/dl (0.2 mmol/l) - ACCENT MC240

Linearity

up to 900 mg/dl (23.3 mmol/l) - ACCENT-200
up to 860 mg/dl (22.3 mmol/l) - ACCENT MC240

Specificity / Interferences

Haemoglobin up to 0.31 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 n=20	level 1	101.74	1.64	1.61
	level 2	240.78	1.50	0.62
ACCENT MC240 n=20	level 1	97.54	0.80	0.82
	level 2	243.29	2.18	0.90
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 n=80	level 1	260.25	7.88	3.03
	level 2	109.73	3.16	2.88
ACCENT MC240 n=80	level 1	242.80	4.39	1.80
	level 2	95.90	1.27	1.30

Method comparison

A comparison between cholesterol values determined at ACCENT-200 (y) and at OLYMPUS AU 400 (x) using 30 samples gave following results:

$$y = 1.047 x - 6.5774 \text{ mg/dl};$$

$$R = 0.997 \quad (R - \text{correlation coefficient})$$

A comparison between cholesterol values determined at ACCENT MC240 (y) and at ADVIA 1800 (x) using 58 serum samples gave following results:

$$y = 0.9963x - 1.8305 \text{ mg/dl};$$

$$R = 0.995 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

- Siedel J., Hägele E.O., Zigenhorn J., Wahlefeld A.W.: Clin. Chem. 29, 6 (1983).
- Tel R.M., Berends G.T.: J Clin.Chem. Clin. Biochem. 18, 10 (1980).
- Rautela G.S., Liedtke R.J.: Clin. Chem. 24. 1 (1978)
- Schettler G., Nussel E.: Arbeitsmed. Sozialmed. Praventivmed. 10, 25 (1975).
- Richmond W.: Clin. Chem. 19, 1350 (1973).
- Roeschlau P., Bernt E., Gruber W.: J. Clin. Chem. Biochem. 12, 403 (1974).
- Tietz N.W., ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA: WB Saunders, 130 (1995).
- Trinder P.: Ann. Clin. Biochem. 6, 24 (1969).
- Dembńska-Kieć A., Naskalski J.W., Solnica B.: Diagnostyka laboratoryjna z elementami biochemii klinicznej, Wyd. IV, str. 670, 875 (2017).
- Cybulska B, Szostak WB, Filipiak KJ, et al. Polish Forum for Prevention Guidelines on Dyslipidaemia: update 2016. Kardiol Pol. 2017; 75(2): 187–190, doi: 10.5603/KP.2017.0031.
- Piepoli MF, Hoes AW, Agewall S, et al. 2016 European Guidelines on cardiovascular disease prevention in clinical practice: The Sixth Joint Task Force of the European Society of Cardiology and Other Societies on Cardiovascular Disease Prevention in Clinical Practice (constituted by representatives of 10 societies and by invited experts)Developed with the special contribution of the European Association for Cardiovascular Prevention & Rehabilitation (EACPR). Eur Heart J. 2016; 37(29): 2315–2381, doi: 10.1093/eurheartj/ehw106, indexed in Pubmed: 27222591.
- NCEP Expert Panel. Arch Inter Med (148), 36-69, (1988).
- Jacobso D.S., DeMott W.R., Grady H.J., et al., ed., Laboratory Tests Handbook, 4th ed., Hudson, Lexi-Comp, 143, (1996).

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- Siedel J., Hägele E.O., Zigenhorn J., Wahlefeld A.W.: Clin. Chem. 29, 6 (1983).
- Tel R.M., Berends G.T.: J Clin.Chem. Clin. Biochem. 18, 10 (1980).
- Rautela G.S., Liedtke R.J.: Clin. Chem. 24. 1 (1978)
- Schettler G., Nussel E.: Arbeitsmed. Sozialmed. Preventivmed. 10, 25 (1975).
- Richmond W.: Clin. Chem. 19, 1350 (1973).
- Roeschlau P., Bernt E., Gruber W.: J. Clin. Chem. Biochem. 12, 403 (1974).
- Tietz N.W., ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA: WB Saunders, 130 (1995).
- Trinder P.: Ann. Clin. Biochem. 6, 24 (1969).
- Dembńska-Kieć A., Naskalski J.W., Solnica B.: Diagnostyka laboratoryjna z elementami biochemii klinicznej. Wyd. IV, str. 670, 875 (2017).
- Cybulska B, Szostak WB, Filipiak KJ, et al. Polish Forum for Prevention Guidelines on Dyslipidaemia: update 2016. Kardiol Pol. 2017; 75(2): 187–190, doi: 10.5603/KP.2017.0031.
- Piepoli MF, Hoes AW, Agewall S, et al. 2016 European Guidelines on cardiovascular disease prevention in clinical practice: The Sixth Joint Task Force of the European Society of Cardiology and Other Societies on Cardiovascular Disease Prevention in Clinical Practice (constituted by representatives of 10 societies and by invited experts) Developed with the special contribution of the European Association for Cardiovascular Prevention & Rehabilitation (EACPR). Eur Heart J. 2016; 37(29): 2315–2381, doi: 10.1093/eurheartj/ehw106, indexed in Pubmed: 27222591.
- NCEP Expert Panel. Arch Inter Med (148), 36-69, (1988).
- Jacobso D.S., DeMott W.R., Grady H.J., et. al., ed., Laboratory Tests Handbook, 4th ed., Hudson, Lexi-Comp, 143, (1996).

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Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	3
Interval (day)	70
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT-220S

Parameters

Test	CHOL	R1	250
No	2	R2	
Full Name	Cholesterol	Sample Volume	3
Standard No	2	R1 Blank	
Reac. Type	Endpoint	Mixed Rtg. Blank	
Pri. Wave.	510 nm	Linearity Range	1.95 950
Sec. Wave.	670 nm	Linearity Limit	
Direction	Increase	Substrate Limit	
Reac. Time	0 20	Factor	
Incuba. Time		<input type="checkbox"/> Prozone check	
Unit	mg/dl	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	Integer	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	3
Interval (day)	70
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• BS-120

Parameters

Test	CHOL	R1	250
No	2	R2	
Full Name	Cholesterol	Sample Volume	3
Standard No	2	R1 Blank	
Reac. Type	Endpoint	Mixed Rtg. Blank	
Pri. Wave.	510 nm	Linearity Range	1.6 800
Sec. Wave.	670 nm	Linearity Limit	
Direction	Increase	Substrate Limit	
Reac. Time	0 16	Factor	
Incuba. Time	16	<input type="checkbox"/> Prozone check	
Unit	mg/dl	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	Integer	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	3
Interval (day)	
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT S120

Chem	CHOL	No.	002	Sample Type	SERUM
Chemistry	CHOLESTEROL			Print name	CHOL
Reaction Type	Endpoint			Reaction Direction	positive
Pri Wave	510 nm			Sec Wave	670 nm
Unit	mg/dl			Decimal	0.1
Blank Time	-3	-1		Incubation Time	0
				Reaction Time	16 18
				Reagent Vol	
Standard	Sample Vol	Aspirated	Diluent	R1	200 µL
Decreased	2.5 µL	20 µL	180 µL	R2	µL
Increased	µL	µL	µL		
	<input type="checkbox"/> Sample Blank	<input checked="" type="checkbox"/> Auto Rerun			
Linearity range (Standard)	3,8	860		Linearity Limit	
Linearity Range (Decreased)				Substrate Depletion	
Linearity Range (Increased)				Mixed Blank Abs	-40000 40000
R1 Blank Abs	-40000	40000		On-board Stability	Day(s)
Blank Response	-40000	40000		Reagent Alarm Limit	
Twin Chemistry				<input type="checkbox"/> Enzyme Linear Extension	
<input type="checkbox"/> Prozone Check					
Q1 <input type="checkbox"/>	Q2 <input type="checkbox"/>	V1 <input type="checkbox"/>	Q3 <input type="checkbox"/>	Q4 <input type="checkbox"/>	V2 <input type="checkbox"/>
Q5 <input type="checkbox"/>	Q6 <input type="checkbox"/>	V3 <input type="checkbox"/>	PC1 <input type="checkbox"/>	PC2 <input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment	<input type="checkbox"/> Calibrator Pretreatment			
	Pretreat Sample Vol <input type="checkbox"/> µL	Pretreat Sample Vol <input type="checkbox"/> µL			
CALIBRATION SETTINGS	Math model <input type="checkbox"/> Multi-point linear	<input type="checkbox"/> Bottle Changed			
Factor <input type="checkbox"/>	Replicates <input type="checkbox"/> 2	<input type="checkbox"/> Lot Changed			
		<input type="checkbox"/> Cal Time			
ACCEPTANCE LIMITS	Cal Time <input type="checkbox"/> Hour				
Slope Diff <input type="checkbox"/>	SD <input type="checkbox"/>				
Sensitivity <input type="checkbox"/>	Repeatability <input type="checkbox"/> 40000				
Deter Coeff <input type="checkbox"/>					