ACCENT-200 IgM

Nr kat. 7-205

(PL)

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Cat. No **7-205** (EN)

INTENDED USE

Diagnostic kit for determination of IgM concentration intended to use in automatic analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, ACCENT 400 and ACCENT Neo200.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Immunoglobulins (Igs) are the instrumental proteins of immunity. Immunity is a property of the lymphoid system which is made of organs (spleen, thymus, bone marrow) and of cells (lymphocytes). Circulating immunoglobulins are secreted in the blood by B lymphocytes and they thereby export far-away the specific biological functions of humoral immunity. Immunoglobulin M (IgM) is the first Ig to appear in response to an antigenic stimulus such as an infectious agent. In many cases, the antigen-specific IgM level subsequently falls and remains low as the IgG response appears.

METHOD PRINCIPLE

The IgM present in a sample form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum measured at $\lambda=340$ nm is proportional to IgM concentration in the sample.

REAGENTS

Package

1-Reagent 1 x 35 ml
 2-Reagent 1 x 7 ml

Buffer (1-Reagent) stored at 2-25°C and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package.

Reagent components

1-Reagent: tricine buffer, polymer, inorganic salt, preservative.

2-Reagent: anti-human IgM antiserum, buffer, inorganic salt, preservative.

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Store closed.
- Do not freeze the reagents.
- Nanoparticle-based reagents can settle over time. It may be necessary to delicately mix by repeated turning.
- Human-origin products have been tested for HBsAg and antibodies to HIV 1, HIV 2, HCV and found to be non-reactive. However this material should be handled as potentially infectious.

- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN ⁴

Serum. Nonhemolyzed, nonlipaemic, fresh serum is recommended.

Serum can be stored up to 3 days at 2-8°C or up to 6 months at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

0.9% NaCl is recommended as a reagent blank.

REFERENCE VALUES ⁴

adults	0.50 – 3.00 g/l
children (1 – 12 years)	0.45 – 2.50 g/l
children (1 – 12 months)	0.20 – 1.50 g/l

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples the CORMAY IMMUNO-CONTROL III (Cat. No 4-291).

For the calibration of automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320 the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended. 0.9% NaCl should be used as a calibrator 0.

The calibration curve should be prepared with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers: Hitachi or ACCENT MC240. Results may vary if a different instrument is used.

- Measurement range:** 0.038 g/l to 7 g/l.

Specificity / Interferences:

Hemoglobin up to 0.32 g/dl, bilirubin up to 22 mg/dl, triglycerides up to 312 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4 g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

Precision

Repeatability (run to run) n = 30	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.60	0.004	0.73
level 2	1.25	0.012	0.87
level 3	18.92	0.014	0.72

Reproducibility (day to day) n = 60	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.56	0.024	3.99
level 2	1.30	0.037	2.81
level 3	19.95	0.050	2.49

Method comparison

A comparison between IgM values determined at **ACCENT MC240** (y) and at **ADVIA 1800** (x) using 62 samples gave following results:

$$y = 0.978x + 0.0291 \text{ g/l};$$

$$R = 0.997 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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Calibration Rule

Rule	Spline
Sensitivity	1
Replicates	1
Interval (day)	0
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT-220S

Parameters

Test	IgM	R1	250
No	42	R2	50
Full Name	Immuno-globulin M	Sample Volume	3
Standard No	42	R1 Blank	
Reac. Type	Endpoint	Mixed Rtg. Blank	
Pri. Wave	340 nm	Linearity Range	
Sec. Wave	670 nm	Linearity Limit	
Direction	Increase	Substrate Limit	
Reac. Time	-1 15	Factor	
Incuba. Time	19	<input type="checkbox"/> Prozone check	
Unit	g/l	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.01	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

Calibration Rule

Rule	Spline
Sensitivity	1
Replicates	2
Interval (day)	0
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT S120

Chem	IgM	No.	042	Sample Type	SERUM
Chemistry	IMMUNOGLOBULIN M	Print name	IgM	Reaction Direction	positive
Reaction Type	Endpoint	Sec Wave	670 nm	Decimal	0.01
Pri Wave	340 nm	Incubation Time	19	Reaction Time	14 16
Unit	g/l	Blank Time	-3 -1	Standard	2.5 μL
Sample Vol	2.5 μL	Aspirated	20 μL	Diluent	180 μL
Decreased	2.5 μL	Reagent Vol R1	200 μL	Reagent Vol R2	40 μL
Increased	μL	Sample Blank	<input type="checkbox"/>	Auto Rerun	<input checked="" type="checkbox"/>

Linearity range (Standard)	0,04	2,30	Linearity Limit	
Linearity Range (Decreased)			Substrate Depletion	
Linearity Range (Increased)			Mixed Blank Abs	-40000 40000
R1 Blank Abs	-40000	40000	On-board Stability	
Blank Response	-40000	40000	Reagent Alarm Limit	
Twin Chemistry			Enzyme Linear Extension	<input type="checkbox"/>
<input type="checkbox"/> Prozone Check				

Q1		Q2		V1		Q3		Q4		V2	
Q5		Q6		V3		PC1		PC2			
Sample Pretreatment		Control Pretreatment		Calibrator Pretreatment							
Pretreat Sample Vol		μL	Pretreat Sample Vol		μL						

CALIBRATION SETTINGS

Math model	Spline	Bottle Changed	<input type="checkbox"/>
Factor		Lot Changed	<input type="checkbox"/>
Replicates	2	Cal Time	

AUTO CALIBRATION

ACCEPTANCE LIMITS

Cal Time		Hour	
Slope Diff		SD	
Sensitivity		Repeatability	40000
Deter Coeff			

ACCENT-200 IgM

• ACCENT MC240

Chem <input type="text" value="IgM"/>	No. <input type="text" value="042"/>	Sample Type <input type="text" value="SERUM"/>
Chemistry <input type="text" value="IMMUNOGLOBULIN M"/>	Print name <input type="text" value="IgM"/>	
Reaction Type <input type="text" value="Endpoint"/>	Reaction Direction <input type="text" value="positive"/>	
Pri Wave <input type="text" value="340 nm"/>	Sec Wave <input type="text" value="700 nm"/>	
Unit <input type="text" value="g/l"/>	Decimal <input type="text" value="0.01"/>	
Blank Time <input type="text" value="-3"/> <input type="text" value="-1"/>	Incubation Time <input type="text" value="21"/>	Reaction Time <input type="text" value="10"/> <input type="text" value="11"/>
Standard <input type="text" value="2.5"/> <input type="text" value="2.5"/> <input type="text" value="2.5"/> <input type="text" value="2.5"/> <input type="text" value="2.5"/>	Aspirated <input type="text" value="20"/> <input type="text" value="20"/> <input type="text" value="20"/> <input type="text" value="20"/> <input type="text" value="20"/>	Diluent <input type="text" value="180"/> <input type="text" value="180"/> <input type="text" value="180"/> <input type="text" value="180"/> <input type="text" value="180"/>
Reagent Vol <input type="text" value="200"/> <input type="text" value="200"/> <input type="text" value="200"/> <input type="text" value="200"/> <input type="text" value="200"/>	R1 <input type="text" value="40"/> <input type="text" value="40"/> <input type="text" value="40"/> <input type="text" value="40"/> <input type="text" value="40"/>	R2 <input type="text" value="40"/> <input type="text" value="40"/> <input type="text" value="40"/> <input type="text" value="40"/> <input type="text" value="40"/>
<input type="checkbox"/> Sample Blank	<input checked="" type="checkbox"/> Auto Rerun	

Linearity range (Standard) <input type="text" value="0,04"/> <input type="text" value="2,40"/>	Linearity Limit <input type="text"/>
Linearity Range (Decreased) <input type="text"/>	Substrate Depletion <input type="text"/>
Linearity Range (Increased) <input type="text"/>	Mixed Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>
R1 Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>	On-board Stability <input type="text"/> Day(s)
Blank Response <input type="text" value="-35000"/> <input type="text" value="35000"/>	Reagent Alarm Limit <input type="text"/>
Twin Chemistry <input type="text"/>	<input type="checkbox"/> Enzyme Linear Extension
<input type="checkbox"/> Prozone Check	

Q1 <input type="text"/>	Q2 <input type="text"/>	V1 <input type="text"/>	Q3 <input type="text"/>	Q4 <input type="text"/>	V2 <input type="text"/>
Q5 <input type="text"/>	Q6 <input type="text"/>	V3 <input type="text"/>	PC1 <input type="text"/>	PC2 <input type="text"/>	
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment	<input type="checkbox"/> Calibrator Pretreatment			
Pretreat Sample Vol <input type="text"/> μ L	Pretreat Sample Vol <input type="text"/> μ L				

CALIBRATION SETTINGS	AUTO CALIBRATION
Math model <input type="text" value="Spline"/>	<input type="checkbox"/> Bottle Changed
Factor <input type="text"/>	<input type="checkbox"/> Lot Changed
Replicates <input type="text" value="2"/>	<input type="checkbox"/> Cal Time

ACCEPTANCE LIMITS
Cal Time <input type="text"/> Hour
Slope Diff <input type="text"/> SD <input type="text"/>
Sensitivity <input type="text"/> Repeatability <input type="text" value="35000"/>
Deter Coeff <input type="text"/>

• ACCENT M320

Chem <input type="text" value="IgM"/>	No. <input type="text" value="042"/>	Sample Type <input type="text" value="SERUM"/>
Chemistry <input type="text" value="IMMUNOGLOBULIN M"/>	Print name <input type="text" value="IgM"/>	
Reaction Type <input type="text" value="Endpoint"/>	Reaction Direction <input type="text" value="positive"/>	
Pri Wave <input type="text" value="340 nm"/>	Sec Wave <input type="text" value="660 nm"/>	
Unit <input type="text" value="g/l"/>	Decimal <input type="text" value="0.01"/>	
Blank Time <input type="text" value="-3"/> <input type="text" value="-1"/>	Incubation Time <input type="text" value="30"/>	Reaction Time <input type="text" value="26"/> <input type="text" value="28"/>
Standard <input type="text" value="4"/> <input type="text" value="4"/> <input type="text" value="4"/> <input type="text" value="4"/> <input type="text" value="4"/>	Aspirated <input type="text" value="20"/> <input type="text" value="20"/> <input type="text" value="20"/> <input type="text" value="20"/> <input type="text" value="20"/>	Diluent <input type="text" value="180"/> <input type="text" value="180"/> <input type="text" value="180"/> <input type="text" value="180"/> <input type="text" value="180"/>
Reagent Vol <input type="text" value="200"/> <input type="text" value="200"/> <input type="text" value="200"/> <input type="text" value="200"/> <input type="text" value="200"/>	R1 <input type="text" value="40"/> <input type="text" value="40"/> <input type="text" value="40"/> <input type="text" value="40"/> <input type="text" value="40"/>	R2 <input type="text" value="40"/> <input type="text" value="40"/> <input type="text" value="40"/> <input type="text" value="40"/> <input type="text" value="40"/>
<input type="checkbox"/> Sample Blank	<input checked="" type="checkbox"/> Auto Rerun	

Linearity range (Standard) <input type="text" value="0,07"/> <input type="text" value="2,35"/>	Linearity Limit <input type="text"/>
Linearity Range (Decreased) <input type="text"/>	Substrate Depletion <input type="text"/>
Linearity Range (Increased) <input type="text"/>	Mixed Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>
R1 Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>	On-board Stability <input type="text"/> Day(s)
Blank Response <input type="text" value="-35000"/> <input type="text" value="35000"/>	Reagent Alarm Limit <input type="text"/>
Twin Chemistry <input type="text"/>	<input type="checkbox"/> Enzyme Linear Extension
<input type="checkbox"/> Prozone Check	

Q1 <input type="text"/>	Q2 <input type="text"/>	V1 <input type="text"/>	Q3 <input type="text"/>	Q4 <input type="text"/>	V2 <input type="text"/>
Q5 <input type="text"/>	Q6 <input type="text"/>	V3 <input type="text"/>	PC1 <input type="text"/>	PC2 <input type="text"/>	
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment	<input type="checkbox"/> Calibrator Pretreatment			
Pretreat Sample Vol <input type="text"/> μ L	Pretreat Sample Vol <input type="text"/> μ L				

CALIBRATION SETTINGS	AUTO CALIBRATION
Math model <input type="text" value="Logit-Log 5P"/>	<input type="checkbox"/> Bottle Changed
Factor <input type="text"/>	<input type="checkbox"/> Lot Changed
Replicates <input type="text" value="2"/>	<input type="checkbox"/> Cal Time

ACCEPTANCE LIMITS
Cal Time <input type="text"/> Hour
Slope Diff <input type="text"/> SD <input type="text"/>
Sensitivity <input type="text"/> Repeatability <input type="text" value="35000"/>
Deter Coeff <input type="text"/>

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