## ACCENT-200 LIPASE

Nr kat. **7-209** (PL) **2-REAGENT**

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### LITERATURA

- Tietz NW et al. Lipase in serum-the elusive enzyme: An overview. Clin Chem 1993;39:746-756.
- Steinberg WM, Goldstein SS, Davies ND et al. Diagnostic assays in acute pancreatitis. (Review). Ann Intern Med 1985; 102:576-580.
- Leybold A, Junge W. Importance of colipase for the measurement of serum lipase activity. Adv Clin Enzymol 1986;4:60-67.
- Alan H. B. Wu, Tietz Clinical Guide to Laboratory Tests, W.B. Saunders Company, 4th edition, 676 (2006).

Data wydania: 07. 2023.

## ACCENT-200 LIPASE

Cat. No **7-209** (EN)

### INTENDED USE

Diagnostic kit for determination of lipase activity intended to use in automatic analyzers: ACCENT-200 (II GEN), ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, BS-120, ACCENT 400 and ACCENT Neo200.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Lipase is a digestive enzyme released into the intestine from the pancreas where it breaks down triglycerides into fatty acids and glycerol prior to absorption. Lipase measurements are used in the diagnosis and treatment of diseases of the pancreas such as acute pancreatitis, obstruction of the pancreatic duct and pancreatic tumours.

### METHOD PRINCIPLE

The colorimetric method is based on a lipase specific degradation of a chromogenic substrate. The specific lipase substrate-DGGMR [1,2-o-dilauryl-racglycero-3-glutaric acid-(6'-methylresorufin) ester] is cleaved by the catalytic action of lipase to form 1,2-o-dilauryl-racglycerol and an unstable intermediate, glutaric acid-(6-methyl resorufin) ester. This decomposes spontaneously in alkaline solution to form glutaric acid and methylresorufin. The lipase activity in the specimen is proportional to the production of methylresorufin in the reaction and can be determined photometrically.

### REAGENTS

#### Package

1-REAGENT 2 x 14 ml  
2-REAGENT 2 x 5.5 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyzer at 2-10°C are stable for 12 weeks (ACCENT-200 (II GEN), ACCENT MC240, ACCENT S120).

### Concentrations in the test

#### 1-REAGENT

TAPS [N-Tris(hydroxymethyl)methyl-3-aminopropanesulfonic acid] 100 mM  
sodium hydroxide 40 mM  
sodium deoxycholate 34 mM  
preservative

#### 2-REAGENT

tartaric acid 9,5 mM  
sodium hydroxide 19 mM  
colipase 460 IU/ml  
2-propanol 0.65 M  
DGGMR [1,2-o-dilauryl-rac-glycero-3-glutaric acid-(6'-methylresorufin)-ester] 0.4 mM

ACCENT-200 LIPASE

51\_03\_03\_050\_08

### WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- EUH210 Safety data sheet available on request.

### SPECIMEN

Serum, heparinized plasma free from hemolysis.  
Sample may be stored for up to 5 days at 2-8°C or 24 hours at 20-25°C.  
Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-REAGENT and 2-REAGENT are ready to use.  
Deionised water is recommended as a reagent blank.  
For analyzers: ACCENT-200 (II GEN), ACCENT-220S and BS-120 it is recommended to determine the reagent blank during each calibration. Deionized water should be used as reagent blank. When performing calibration, the task type **Calib+Rgt.Blk** should be selected.

### Actions required:

- ACCENT-200 (II GEN), ACCENT-220S, BS-120:** When performing assays in the analyzers, there is a probability of **cross-contamination** affecting the tests results: HDL DIRECT – LIPASE, LDL DIRECT – LIPASE, TG – LIPASE, TG mono – LIPASE, CHOL – LIPASE, CALCIUM ARSENAZO - LIPASE. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_001\_ACCENT-200\_CARRYOVER.
- ACCENT S120:** When performing assays in the analyzers, there is a probability of **cross-contamination** affecting the tests results: LDL DIRECT – LIPASE, TG – LIPASE, CHOL – LIPASE. To avoid this effect follow the recommendations contained in the instruction: 51\_03\_24\_005\_ACCENT\_S120\_CARRYOVER.
- ACCENT MC240, ACCENT M320:** When performing assays in the analyzers, there is a probability of **cross-contamination** affecting the tests results: LDL DIRECT – LIPASE, TG – LIPASE, CHOL – LIPASE. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_006\_ACCENT\_MC240\_M320\_CARRYOVER.

### REFERENCE VALUES <sup>4</sup>

Normal range	13 – 60 U/l	0.22 – 1.00 µkat/l
--------------	-------------	--------------------

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173). For the calibration of automatic analysers: ACCENT-200 (II GEN) and BS-120 the CORMAY MULTICALIBRATOR LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

For the calibration of automatic analysers: ACCENT S120, ACCENT-220S, ACCENT MC240, ACCENT M320 the CORMAY MULTICALIBRATOR LEVEL 2 (Cat. No5-175; 5-177) is recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 12 weeks (ACCENT-200 (II GEN), ACCENT MC240, ACCENT S120), with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyzers ACCENT-200 (II GEN) and ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**  
0.2 U/l (0.003 µkat/l)
- LoD (Limit of Detection):**  
0.4 U/l (0.007 µkat/l)
- LoQ (Limit of Quantitation):**  
4 U/l (0.067 µkat/l) – ACCENT-200 (II GEN)  
2 U/l (0.033 µkat/l) – ACCENT MC240
- Linearity:**  
up to 250 U/l (4.17 µkat/l) – ACCENT-200 (II GEN)  
up to 330 U/l (5.50 µkat/l) – ACCENT MC240

For higher activity dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

### Specificity / Interferences

Haemoglobin up to 0.16 g/dl, ascorbate up to 62 mg/l, bilirubin up to 15 mg/dl and triglycerides up to 750 mg/dl do not interfere with the test.

### Precision

Repeatability (run to run) n = 20		Mean [U/l]	SD [U/l]	CV [%]
ACCENT-200 (II GEN)	level 1	60.3	1.01	1.68
	level 2	119.5	1.94	1.62
ACCENT MC240	level 1	60.4	1.09	1.80
	level 2	109.4	1.20	1.10
Reproducibility (day to day) n = 80		Mean [U/l]	SD [U/l]	CV [%]
ACCENT-200 (II GEN)	level 1	61.7	2.00	3.2
	level 2	118.1	3.12	2.6
ACCENT MC240	level 1	61.4	2.47	4.0
	level 2	109.1	2.46	2.3

### Method comparison

A comparison between lipase values determined at **ACCENT-200 (II GEN)** (y) and at **ATELICA SIEMENS 930** (x), using 32 serum samples gave following results:

$$y = 1.2928x - 9.2288 \text{ U/l};$$

$$R = 0.977 \quad (R - \text{correlation coefficient})$$

A comparison between lipase values determined at **ACCENT MC240** (y) and at **ADVIA SIEMENS 1800** (x), using 31 serum samples gave following results:

$$y = 1.1127x - 5.7469 \text{ U/l};$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

### LITERATURE

- Tietz NW et al. Lipase in serum-the elusive enzyme: An overview. Clin Chem 1993;39:746-756.
- Steinberg WM, Goldstein SS, Davies ND et al. Diagnostic assays in acute pancreatitis. (Review). Ann Intern Med 1985; 102:576-580.
- Leybold A, Junge W. Importance of colipase for the measurement of serum lipase activity. Adv clin Enzymol 1986;4:60-67.
- Alan H. B. Wu, Tietz Clinical Guide to Laboratory Tests, W.B. Saunders Company, 4th edition, 676 (2006).

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Calibration Rule

Rule	One-point Linear
Sensitivity	1
Replicates	2
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT-220S

Parameters

Test	LIPASE	R1	180
No	21	R2	60
Full Name	Lipase	Sample Volume	4
Standard No	21	R1 Blank	
Reac. Type	Kinetic	Mixed Rtg. Blank	
Pri. Wave.	578 nm	Linearity Range	8 250
Sec. Wave.	670 nm	Linearity Limit	0.2
Direction	Increase	Substrate Limit	
Reac. Time	4 14	Factor	
Incuba. Time	11	<input type="checkbox"/> Prozone check	
Unit	U/l	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.1	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

Calibration Rule

Rule	One-point Linear
Sensitivity	1
Replicates	2
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• BS-120

Parameters

Test	LIPA	R1	180
No	21	R2	60
Full Name	Lipase	Sample Volume	4
Standard No	21	R1 Blank	
Reac. Type	Kinetic	Mixed Rtg. Blank	
Pri. Wave.	578 nm	Linearity Range	10 250
Sec. Wave.	670 nm	Linearity Limit	0.2
Direction	Increase	Substrate Limit	
Reac. Time	4 14	Factor	
Incuba. Time	16	<input type="checkbox"/> Prozone check	
Unit	U/l	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.1	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

Calibration Rule

Rule	One-point Linear
Sensitivity	1
Replicates	2
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT S120

Chem	LIPASE	No.	021	Sample Type	SERUM						
Chemistry	LIPASE	Print name	LIPASE	Reaction Direction	positive						
Reaction Type	Kinetic	Reaction Time	5 13	Sec Wave	670 nm						
Pri Wave	578 nm	Decimal	0.1	Incubation Time	17						
Unit	U/L	Reaction Time	5 13	Reagent Vol	R1 150 μL						
Blank Time		Standard	2 μL	Diluent							
		Aspirated	20 μL	R2	75 μL						
		Increased		Auto Rerun	<input type="checkbox"/> V						
Linearity range (Standard)	2	320	Linearity Limit	0.4							
Linearity Range (Decreased)			Substrate Depletion	40000							
Linearity Range (Increased)			Mixed Blank Abs	-40000 40000							
R1 Blank Abs	-40000	40000	On-board Stability	84 Day(s)							
Blank Response	-40000	40000	Reagent Alarm Limit								
Twin Chemistry			<input type="checkbox"/> Enzyme Linear Extension								
<input type="checkbox"/> Prozone Check											
Q1		Q2		V1		Q3		Q4		V2	
Q5		Q6		V3		PC1		PC2			
<input type="checkbox"/> Sample Pretreatment		<input type="checkbox"/> Control Pretreatment		<input type="checkbox"/> Calibrator Pretreatment							
		Pretreat Sample Vol		μL		Pretreat Sample Vol		μL			
<b>CALIBRATION SETTINGS</b>		<b>AUTO CALIBRATION</b>									
Math model	Two-point linear	<input type="checkbox"/> Bottle Changed									
Factor		Replicates	2	<input type="checkbox"/> Lot Changed							
				<input type="checkbox"/> Cal Time							
<b>ACCEPTANCE LIMITS</b>											
Cal Time	2016 Hour										
Slope Diff		SD									
Sensitivity		Repeatability	40000								
Deter Coeff											

## ACCENT-200 LIPASE

### • ACCENT MC240

Chem <b>LIPASE</b>	No. <b>021</b>	Sample Type <b>SERUM</b>
Chemistry <b>LIPASE</b>	Print name <b>LIPASE</b>	
Reaction Type <b>Kinetic</b>	Reaction Direction <b>positive</b>	
Pri Wave <b>570 nm</b>	Sec Wave <b>700 nm</b>	
Unit <b>U/L</b>	Decimal <b>0.1</b>	
Blank Time	Incubation Time <b>21</b>	Reaction Time <b>5</b> <b>17</b>
Standard <b>2.5</b> $\mu$ L	Aspirated <b>20</b> $\mu$ L	Diluent <b>180</b> $\mu$ L
Decreased <b>2.5</b> $\mu$ L	Reagent Vol R1 <b>180</b> $\mu$ L	R2 <b>90</b> $\mu$ L
Increased	Sample Blank <input type="checkbox"/>	Auto Rerun <input checked="" type="checkbox"/>

Linearity range (Standard) <b>2</b> <b>330</b>	Linearity Limit <b>0.4</b>
Linearity Range (Decreased)	Substrate Depletion <b>35000</b>
Linearity Range (Increased)	Mixed Blank Abs <b>-35000</b> <b>35000</b>
R1 Blank Abs <b>-35000</b> <b>35000</b>	On-board Stability <b>84</b> Day(s)
Blank Response <b>-35000</b> <b>35000</b>	Reagent Alarm Limit
Twin Chemistry	<input type="checkbox"/> Enzyme Linear Extension
<input type="checkbox"/> Prozone Check	
Q1 <input type="checkbox"/> Q2 <input type="checkbox"/> V1 <input type="checkbox"/> Q3 <input type="checkbox"/> Q4 <input type="checkbox"/> V2 <input type="checkbox"/>	
Q5 <input type="checkbox"/> Q6 <input type="checkbox"/> V3 <input type="checkbox"/> PC1 <input type="checkbox"/> PC2 <input type="checkbox"/>	
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment
<input type="checkbox"/> Pretreat Sample Vol <input type="checkbox"/> $\mu$ L	<input type="checkbox"/> Calibrator Pretreatment
<input type="checkbox"/> Pretreat Sample Vol <input type="checkbox"/> $\mu$ L	

<b>CALIBRATION SETTINGS</b>	<b>AUTO CALIBRATION</b>
Math model <b>Two-point linear</b>	<input type="checkbox"/> Bottle Changed
Factor <input type="checkbox"/>	<input type="checkbox"/> Lot Changed
Replicates <b>2</b>	<input type="checkbox"/> Cal Time

<b>ACCEPTANCE LIMITS</b>
Cal Time <b>2016</b> Hour
Slope Diff <input type="checkbox"/> SD <input type="checkbox"/>
Sensitivity <input type="checkbox"/> Repeatability <b>35000</b>
Deter Coeff <input type="checkbox"/>

### • ACCENT M320

Chem <b>LIPASE</b>	No. <b>021</b>	Sample Type <b>SERUM</b>
Chemistry <b>LIPASE</b>	Print name <b>LIPASE</b>	
Reaction Type <b>Kinetic</b>	Reaction Direction <b>positive</b>	
Pri Wave <b>570 nm</b>	Sec Wave <b>700 nm</b>	
Unit <b>U/L</b>	Decimal <b>0.1</b>	
Blank Time	Incubation Time <b>21</b>	Reaction Time <b>5</b> <b>17</b>
Standard <b>2.5</b> $\mu$ L	Aspirated <b>20</b> $\mu$ L	Diluent <b>180</b> $\mu$ L
Decreased <b>2.5</b> $\mu$ L	Reagent Vol R1 <b>180</b> $\mu$ L	R2 <b>90</b> $\mu$ L
Increased	Sample Blank <input type="checkbox"/>	Auto Rerun <input checked="" type="checkbox"/>

Linearity range (Standard) <b>1</b> <b>430</b>	Linearity Limit <b>0.4</b>
Linearity Range (Decreased)	Substrate Depletion <b>35000</b>
Linearity Range (Increased)	Mixed Blank Abs <b>-35000</b> <b>35000</b>
R1 Blank Abs <b>-35000</b> <b>35000</b>	On-board Stability <b>84</b> Day(s)
Blank Response <b>-35000</b> <b>35000</b>	Reagent Alarm Limit
Twin Chemistry	<input type="checkbox"/> Enzyme Linear Extension
<input type="checkbox"/> Prozone Check	
Q1 <input type="checkbox"/> Q2 <input type="checkbox"/> V1 <input type="checkbox"/> Q3 <input type="checkbox"/> Q4 <input type="checkbox"/> V2 <input type="checkbox"/>	
Q5 <input type="checkbox"/> Q6 <input type="checkbox"/> V3 <input type="checkbox"/> PC1 <input type="checkbox"/> PC2 <input type="checkbox"/>	
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment
<input type="checkbox"/> Pretreat Sample Vol <input type="checkbox"/> $\mu$ L	<input type="checkbox"/> Calibrator Pretreatment
<input type="checkbox"/> Pretreat Sample Vol <input type="checkbox"/> $\mu$ L	

<b>CALIBRATION SETTINGS</b>	<b>AUTO CALIBRATION</b>
Math model <b>Two-point linear</b>	<input type="checkbox"/> Bottle Changed
Factor <input type="checkbox"/>	<input type="checkbox"/> Lot Changed
Replicates <b>2</b>	<input type="checkbox"/> Cal Time

<b>ACCEPTANCE LIMITS</b>
Cal Time <b>2016</b> Hour
Slope Diff <input type="checkbox"/> SD <input type="checkbox"/>
Sensitivity <input type="checkbox"/> Repeatability <b>35000</b>
Deter Coeff <input type="checkbox"/>

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