ACCENT-200 ALAT

Nr kat. 7-216 (PL)

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(EN)

INTENDED USE



ALAT reagent is intended to determine quantitatively alanine aminotransferase (ALAT) activity in serum. It is used to monitoring and as an aid to diagnosis of clinical conditions associated with abnormal alanine aminotransferase activity. ALAT reagent is intended to use on automatic analysers ACCENT-200 (II GEN) / BS-200, ACCENT-220S / BS-180, ACCENT S120 / BS-230, ACCENT MC240 / BS-240Pro, ACCENT M320 / BS-360E, BS-120, ACCENT 400 and ACCENT Neo200. It is only for *in vitro* diagnostics, for healthcare professional users.

SUMMARY ^{1,2}

Determinations of elevated alanine aminotransferase activity (ALAT, ALT, GPT) in serum are mainly used to aid in the diagnosis and monitor of liver diseases e.g. acute hepatitis, hepatocyte necrosis, hepatic ischemia, cirrhosis, cholestasis, liver tumor, hepatotoxic drugs, mechanical jaundice. Increased levels of ALAT can be observed in severe burns, striated muscle injury, muscle inflammation, pancreatitis, myocardial infarction, infectious mononucleosis, trauma, patients with higher risk of COVID-19 adverse outcome.

METHOD PRINCIPLE ³

Optimized, modified method according to International Federation of Clinical Chemistry (IFCC), without pyridoxal phosphate.

The rate of absorbance changing at $\lambda=340$ nm is directly proportional to alanine aminotransferase activity.

REAGENTS

Package

1-REAGENT	2 x 35 ml
2-REAGENT	1 x 17.5 ml

Number of tests:

ACCENT-200 (II GEN)	310
ACCENT-220S	310
ACCENT S120	490
ACCENT MC240	490
ACCENT M320	370
BS-120	310

CONCENTRATION OF THE ACTIVE INGREDIENTS

1-REAGENT

L-alanine	625 mmol/l
LDH	3.3 U/ml
buffer Tris stabilizers	

ACCENT-200 ALAT (II GENERACJA/ II GENERATION/ II ПОКОЛЕНИЕ)

pH adjuster
detergent
preservative
2-REAGENT
2-oxoglutarate
NADH
buffer
pH adjuster
preservatives

62 mmol/l
1.4 mmol/l

REAGENT STABILITY

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyzer at 2-10°C are stable for 12 weeks (ACCENT-200 (II GEN), ACCENT M320, ACCENT S120).

WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- Do not use reagent beyond expiry date printed on the package.
- Do not mix reagents from different kits or lots.
- Use personal protective equipment to prevent contact with samples, reagents and controls.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN ^{4,5}

Serum free from hemolysis.

Hemolysis should be avoided, since ALAT activity in erythrocytes is 3 to 5 times higher than in normal serum. Do not freeze the samples. ALAT activity remains stable in specimen up to 3 days at 15-25°C or up to 7 days at 2-8°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

Follow tube manufacturers' instructions carefully when using collection tubes.

Human-origin material should be handled as potentially infectious. Standard precautions in normal laboratory work are required.

PROCEDURE

1-REAGENT and 2-REAGENT are ready to use.

Deionised water is recommended as a reagent blank.

For analyzers: ACCENT-200 (II GEN), ACCENT-220S and BS-120, it is recommended to determine the reagent blank during each calibration. Deionized water should be used as reagent blank. When performing calibration, the task type **Calib+Rgt.Blk** should be selected.

Actions required:

When performing assays in analyzers: ACCENT-200 (II GEN), ACCENT-220S and BS-120 there is a probability of **cross-contamination** affecting the tests results: ALAT - CHOL, ALAT - CK-MB. To avoid this effect follow the recommendations contained in the instruction 51_03_24_001_ACCENT-200_CARRYOVER.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173). For the calibration of automatic analysers the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 12 weeks (ACCENT-200 (II GEN), ACCENT M320, ACCENT S-120). Calibration is recommended in the following cases:

- after each change of lot,
- after instrument service,
- if controls lie outside the expected range,
- each time a new reagent kit is used.

If the quality control results do not fall within the expected values or within the range determined in the laboratory, despite a successful calibration procedure, do not report results. In this case, please take the following actions:

- verify reagents are not out of expiration date.
- verify that the required maintenance has been carried out.
- verify that the procedure has been performed in accordance with the instructions for use.
- contact the Service Department or distributor for assistance.

REFERENCE VALUES ⁶

serum	37°C	
women	up to 34 U/l	up to 0.58 µkat/l
men	up to 45 U/l	up to 0.77 µkat/l

It is recommended for each laboratory to establish its own reference ranges for local population and method principle.

The assay results should be used in conjunction with other data, such as symptoms, results of other tests and clinical history to make clinical decisions. It is not recommended to make clinical diagnosis based on a single result.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers ACCENT-200 (II GEN) and/or ACCENT MC240 and/or BS-400. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**
2.0 U/l (0.03 µkat/l) – BS-400
- LoD (Limit of Detection):**
3.0 U/l (0.05 µkat/l) – BS-400
- LoQ (Limit of Quantitation):**
8.0 U/l (0.13 µkat/l) – ACCENT-200 (II GEN), ACCENT MC240
- Linearity:**
up to 1000 U/l (16.7 µkat/l) – ACCENT-200 (II GEN)
up to 1100 U/l (18.3 µkat/l) – ACCENT MC240
- Measurement range:**
8.0 U/l (0.13 µkat/l) - 1000 U/l (16.7 µkat/l) – ACCENT-200 (II GEN)
8.0 U/l (0.13 µkat/l) - 1100 U/l (18.3 µkat/l) – ACCENT MC240

Specificity / Interferences

Ascorbic acid up to 62 mg/l, triglycerides up to 1000 mg/dl, bilirubin up to 20 mg/dl, haemoglobin up to 0.31 g/dl (in samples containing low activity of ALAT) and up to 2.5 g/dl (in samples containing high activity of ALAT) do not interfere with the test.

Precision

Repeatability (run to run)		Mean [U/l]	SD [U/l]	CV [%]
ACCENT-200 (II GEN) n=20	level 1	32.8	1.3	3.9
	level 2	104.5	1.5	1.4
ACCENT MC240 n=20	level 1	31.4	1.24	4.0
	level 2	99.6	1.05	1.1
Reproducibility (day to day)		Mean [U/l]	SD [U/l]	CV [%]
ACCENT-200 (II GEN) n=80	level 1	32.7	1.68	5.1
	level 2	103.3	3.53	3.4
ACCENT MC240 n=80	level 1	33.6	1.98	5.9
	level 2	102.7	3.98	3.9

Method comparison

A comparison between ALAT values determined at **ACCENT-200 (II GEN)** (y) and at **BECKMAN COULTER AU680** (x) using 60 serum samples gave following results:

$$y = 1.0368 x - 0.4275 \text{ U/l;}$$

$$R = 1.000 \quad (R - \text{correlation coefficient})$$

A comparison between ALAT values determined at **ACCENT MC240** (y) and at **BECKMAN COULTER AU680** (x) using 60 serum samples gave following results:

$$y = 0.9665 x + 1.7925 \text{ U/l;}$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT ⁷

After use, the reagents should be handled as potentially infectious and disposed of in accordance with local legal requirements.

- Reagents residues:
1-REAGENT – 18 01 07
2-REAGENT – 18 01 07
- Empty packages:
1-REAGENT – 15 01 02
2-REAGENT – 15 01 02
- Wastewater from the analyzer:
18 01 03*

INCIDENTS ⁸

Any serious incident that has occurred in relations to the device shall be reported to the manufacturer (website address: incidents@cormay.pl) and the competent authority of the Member State in which the user and/or the patient is established.

Serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- a serious public health threat.

LITERATURE

- Pagana K. Pagana T.J., Mosby's Diagnostic and Laboratory Test Reference, 10th ed., Elsevier, St. Louis, USA, 27 (2013).
- Lippi G., Plebani M., Laboratory abnormalities in patients with COVID-2019 infection, Clin Chem Lab Med, 58(7):1131-1134 (2020).
- Bergmeyer H.U., Horder M., Rej R.: J. Clin. Chem. Clin. Biochem. 24, 481-495 (1986).
- Kaplan L.A., Pesce A.J., Clinical Chemistry: Theory, analysis, and correlation, 2nd ed., The C.V. Mosby Company, St. Louis, USA, 896 (1989).
- World Health Organization (2002). Use of Anticoagulants in Diagnostic Laboratory Investigations (Stability of Blood, Plasma and serum samples); WHO/DIL/LAB/99.1 Rev.2: 1-64
- Burtis C.A., Ashwood E.R., Bruns D.E. ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 4th ed., Elsevier Saunders, St. Louis, USA, 2254 (2006).
- European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
- Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices

LIST OF CHANGES

Previous version: 06	Current version: 07
Sections updated: <i>INTENDED USE</i>	

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1. Pagana K. Pagana T.J., Mosby's Diagnostic and Laboratory Test Reference, 10th ed., Elsevier, St. Louis, USA, 27 (2013).
2. Lippi G., Plebani M., Laboratory abnormalities in patients with COVID-2019 infection, Clin Chem Lab Med, 58(7):1131-1134 (2020).
3. Bergmeyer H.U., Horder M., Rej R.: J. Clin. Chem. Clin. Biochem. 24, 481-495 (1986).
4. Kaplan L.A., Pesce A.J., Clinical Chemistry: Theory, analysis, and correlation, 2nd ed., The C.V. Mosby Company, St. Louis, USA, 896 (1989).
5. World Health Organization (2002). Use of Anticoagulants in Diagnostic Laboratory Investigations (Stability of Blood, Plasma and serum samples); WHO/DIL/LAB/99.1 Rev.2: 1-64
6. Burtis C.A., Ashwood E.R., Bruns D.E. ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 4th ed., Elsevier Saunders, St. Louis, USA, 2254 (2006).
7. European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
8. Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices

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Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	2
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• BS-120

Parameters	ALAT II GEN	R1	200
Test	ALAT II GEN	R2	50
No	15	Sample Volume	10
Full Name	ALAT	R1 Blank	
Standard No	15	Mixed Rtg. Blank	
Reac. Type	Kinetic	Linearity Range	8 1000
Pri. Wave.	340 nm	Linearity Limit	0.2
Sec. Wave.	450 nm	Substrate Limit	5600
Direction	Decrease	Factor	
Reac. Time	3 12	<input type="checkbox"/> Prozone check	
Incuba. Time	16	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Unit	U/l	PC <input type="checkbox"/> Abs <input type="checkbox"/>	
Precision	0.1		

Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	2
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT-220S

Parameters	ALAT II GEN	R1	200
Test	ALAT II GEN	R2	50
No	15	Sample Volume	12
Full Name	ALAT	R1 Blank	
Standard No	15	Mixed Rtg. Blank	
Reac. Type	Kinetic	Linearity Range	8 1000
Pri. Wave.	340 nm	Linearity Limit	0.2
Sec. Wave.	450 nm	Substrate Limit	5400
Direction	Decrease	Factor	
Reac. Time	3 12	<input type="checkbox"/> Prozone check	
Incuba. Time	11	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Unit	U/l	PC <input type="checkbox"/> Abs <input type="checkbox"/>	
Precision	0.1		

Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	2
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• **ACCENT M320**

Chem	<input type="text" value="ALAT II GEN"/>	No.	<input type="text" value="015"/>	Sample Type	<input type="text" value="SERUM"/>
Chemistry	<input type="text" value="ALAT"/>			Print name	<input type="text" value="ALAT"/>
Reaction Type	<input type="text" value="Kinetic"/>			Reaction Direction	<input type="text" value="negative"/>
Pri Wave	<input type="text" value="340 nm"/>			Sec Wave	<input type="text" value="450 nm"/>
Unit	<input type="text" value="U/l"/>			Decimal	<input type="text" value="0.1"/>
Blank Time	<input type="text"/>			Incubation Time	<input type="text" value="20"/>
				Reaction Time	<input type="text" value="10"/> <input type="text" value="25"/>
Standard	<input type="text" value="8"/> μL	Aspirated	<input type="text"/>	Diluent	<input type="text"/>
Decreased	<input type="text" value="8"/> μL	<input type="text" value="20"/> μL	<input type="text" value="180"/> μL	Reagent Vol	R1 <input type="text" value="160"/> μL
Increased	<input type="text"/>	<input type="text"/>	<input type="text"/>	R2	<input type="text" value="40"/> μL
	<input type="text"/>	Sample Blank	<input type="checkbox"/>	Auto Rerun	<input type="checkbox"/>

Linearity range (Standard)	<input type="text" value="8"/>	<input type="text" value="1200"/>	Linearity Limit	<input type="text" value="0.8"/>	
Linearity Range (Decreased)	<input type="text"/>	<input type="text"/>	Substrate Depletion	<input type="text" value="3300"/>	
Linearity Range (Increased)	<input type="text"/>	<input type="text"/>	Mixed Blank Abs	<input type="text" value="-35000"/> <input type="text" value="35000"/>	
R1 Blank Abs	<input type="text" value="-35000"/>	<input type="text" value="35000"/>	On-board Stability	<input type="text" value="84"/> Day(s)	
Blank Response	<input type="text" value="-35000"/>	<input type="text" value="35000"/>	Reagent Alarm Limit	<input type="text"/>	
Twin Chemistry	<input type="text"/>		<input type="checkbox"/> Enzyme Linear Extension		
	<input type="checkbox"/> Prozone Check				
Q1	<input type="text"/>	Q2	<input type="text"/>	V1	<input type="text"/>
Q3	<input type="text"/>	Q4	<input type="text"/>	V2	<input type="text"/>
Q5	<input type="text"/>	Q6	<input type="text"/>	V3	<input type="text"/>
PC1	<input type="text"/>	PC2	<input type="text"/>		
Sample Pretreatment	<input type="text"/>	Control Pretreatment	<input type="text"/>	Calibrator Pretreatment	<input type="text"/>
	<input type="text"/>	Pretreat Sample Vol	<input type="text"/>	μL	Pretreat Sample Vol
				μL	

CALIBRATION SETTINGS		AUTO CALIBRATION	
Math model	<input type="text" value="Multi-point linear"/>	<input type="checkbox"/> Bottle Changed	
Factor	<input type="text"/>	<input type="checkbox"/> Lot Changed	
Replicates	<input type="text" value="2"/>	<input type="checkbox"/> Cal Time	

ACCEPTANCE LIMITS	
Cal Time	<input type="text" value="2016"/> Hour
Slope Diff	<input type="text"/> SD <input type="text"/>
Sensitivity	<input type="text"/> Repeatability <input type="text" value="35000"/>
Deter Coeff	<input type="text"/>

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PZ CORMAY S.A.