ACCENT-200 TOTAL IgE

Nr kat. 7-231

(PL)

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Cat. No **7-231** (EN)

INTENDED USE

Diagnostic kit for determination of IgE levels intended to use in automatic analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT 400, ACCENT Neo200 and ACCENT M320.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

IgE is an immunoglobulin with a molecular weight of approximately 190 kD normally present in the blood in trace amounts. Continual production of IgE antibodies in response to common naturally occurring allergens, however, often results in elevated serum levels and in the development of such clinically important Type I allergic reactions as asthma, hay fever, dermatitis and food allergies. Elevated IgE levels are also seen in parasitic (helminth) diseases, IgE myeloma, and in hepatitis. The measurement of IgE in human serum is thus considered to be useful in the diagnosis, treatment, assessment of disease progression, or postoperative prognosis for such conditions.

METHOD PRINCIPLE

When an antigen-antibody reaction occurs between IgE in a sample and anti-IgE antibody which has been sensitized to latex particles, agglutination results. This agglutination is detected as an absorbance change (572 nm), with the magnitude of the change being proportional to the quantity of IgE in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

REAGENTS

Package

1-Reagent 1 x 28 ml
 2-Reagent 1 x 14 ml

The reagents, stored at 2-10°C are stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 5 weeks.

Concentrations in the test

suspension of latex particles sensitized with (mouse) anti-IgE antibodies (pH 7.3) 0.125 w/v%
 glycine buffer solution (pH 8.3)
 preservative

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- After measurements are taken, reagent bottles should be capped and kept at 2-10°C.
- Do not interchange the caps of reagent bottles.

- Reagents with different lot numbers should not be interchanged or mixed.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum, plasma collected on heparin (sodium / lithium salt), EDTA (sodium / potassium salt) or citric acid.

If the test cannot be done immediately, the sample should be stored at -20°C. Avoid repeated freezing and thawing. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.
 0.9% NaCl is recommended as a reagent blank.

REFERENCE VALUES ³

serum, plasma	< 358 IU/ml
---------------	-------------

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY IMMUNO-CONTROL II (Cat. No 4-290) with each batch of samples.

For the calibration of automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S ACCENT S120, ACCENT MC240, ACCENT M320, the CORMAY IgE CALIBRATORS kit (Cat. No 4-280) is recommended. 0.9% NaCl should be used as a calibrator 0.

The calibration curve should be prepared every 5 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers: ACCENT MC240 and TBA-30R. Results may vary if a different instrument is used.

- Measurement range:** 25 – 1000 IU/ml

For higher concentration, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 0.5 g/dl, bilirubin up to 30 mg/dl, triglycerides up to 1500 mg/dl, RF up to 500 IU/ml do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean [IU/ml]	SD [IU/ml]	CV [%]
level 1	67.14	1.50	2.23
level 2	411.70	2.80	0.68

Reproducibility (day to day) n=80	Mean [IU/ml]	SD [IU/ml]	CV [%]
level 1	67.9	4.62	6.8
level 2	404.8	12.99	3.2

Method comparison

A comparison between IgE levels determined at ACCENT MC240 (y) and at BS-800 (x) using 35 serum samples gave following results:

$$y = 0.9753 x + 3.1135 \text{ IU/ml};$$

$$R = 1.000 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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Calibration Rule

Rule	Spline
Sensitivity	1
Replicates	1
Interval (day)	0
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT-220S

Parameters			
Test	IgE	R1	160
No	34	R2	80
Full Name	Total IgE	Sample Volume	4
Standard No	34	R1 Blank	
Reac. Type	Kinetic	Mixed Rtg. Blank	
Pri. Wave.	578 nm	Linearity Range	
Sec. Wave.		Linearity Limit	
Direction	Increase	Substrate Limit	
Reac. Time	6 16	Factor	
Incuba. Time	10	<input type="checkbox"/> Prozone check	
Unit	IU/ml	q1	<input type="checkbox"/>
Precision	Integer	q2	<input type="checkbox"/>
		q3	<input type="checkbox"/>
		q4	<input type="checkbox"/>
		PC	<input type="checkbox"/>
		Abs	<input type="checkbox"/>

Calibration Rule

Rule	Spline
Sensitivity	1
Replicates	1
Interval (day)	0
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT S120

Chem	IgE	No.	034	Sample Type	SERUM
Chemistry	TOTAL IgE	Print name	IgE	Reaction Direction	positive
Reaction Type	Kinetic	Sec Wave		Decimal	0.1
Pri Wave	578 nm	Incubation Time	15	Reaction Time	4 12
Unit	IU/ml	Blank Time		Standard	Sample Vol 4 μL
				Aspirated	μL
				Diluent	μL
				Reagent Vol	R1 150 μL
					R2 75 μL
				Decreased	4 μL
				Increased	μL
				Sample Blank	<input type="checkbox"/>
				Auto Rerun	<input checked="" type="checkbox"/>
Linearity range (Standard)	25	1000	Linearity Limit	0.2	
Linearity Range (Decreased)			Substrate Depletion	40000	
Linearity Range (Increased)			Mixed Blank Abs	-40000 40000	
R1 Blank Abs	-40000	40000	On-board Stability	Day(s)	
Blank Response	-40000	40000	Reagent Alarm Limit		
Twin Chemistry			Enzyme Linear Extension	<input type="checkbox"/>	
			Prozone Check	<input type="checkbox"/>	
Q1		Q2		V1	
Q3		Q4		V2	
Q5		Q6		PC1	
Q7		Q8		PC2	
Sample Pretreatment		Control Pretreatment		Calibrator Pretreatment	
		Pretreat Sample Vol	μL	Pretreat Sample Vol	μL
CALIBRATION SETTINGS	Math model	Spline	AUTO CALIBRATION	Bottle Changed	<input type="checkbox"/>
	Factor	Replicates	2	Lot Changed	<input type="checkbox"/>
				Cal Time	<input type="checkbox"/>
ACCEPTANCE LIMITS	Cal Time	Hour	Slope Diff	SD	
	Sensitivity	Repeatability	40000	Deter Coeff	

ACCENT-200 TOTAL IgE

• ACCENT MC240

Chem IgE	No. 034	Sample Type SERUM
Chemistry TOTAL IgE	Print name IgE	
Reaction Type Kinetic	Reaction Direction positive	
Pri Wave 570 nm	Sec Wave 800 nm	
Unit IU/ml	Decimal 0.1	
Blank Time	Incubation Time 21	Reaction Time 4 14
Standard 4 μ L	Aspirated 20 μ L	Diluent 180 μ L
Decreased 4 μ L	Reagent Vol R1 150 μ L	R2 75 μ L
Increased	Sample Blank <input type="checkbox"/>	Auto Rerun <input checked="" type="checkbox"/>

Linearity range (Standard)	13	1000	Linearity Limit	0.2
Linearity Range (Decreased)			Substrate Depletion	35000
Linearity Range (Increased)			Mixed Blank Abs	-35000 35000
R1 Blank Abs	-35000	35000	On-board Stability	Day(s)
Blank Response	-35000	35000	Reagent Alarm Limit	
Twin Chemistry			Enzyme Linear Extension	<input type="checkbox"/>
Prozone Check	<input type="checkbox"/>			
Q1		Q2	V1	Q3
Q5		Q6	V3	PC1
Sample Pretreatment		Control Pretreatment		Calibrator Pretreatment
Pretreat Sample Vol	μ L	Pretreat Sample Vol	μ L	

CALIBRATION SETTINGS		AUTO CALIBRATION
Math model Spline		Bottle Changed
Factor	Replicates 2	Lot Changed
		Cal Time

ACCEPTANCE LIMITS	
Cal Time	Hour
Slope Diff	SD
Sensitivity	Repeatability 35000
Deter Coeff	

• ACCENT M320

Chem IgE	No. 034	Sample Type SERUM
Chemistry TOTAL IgE	Print name IgE	
Reaction Type Kinetic	Reaction Direction positive	
Pri Wave 570 nm	Sec Wave 800 nm	
Unit IU/ml	Decimal 0.1	
Blank Time	Incubation Time 15	Reaction Time 9 19
Standard 4.5 μ L	Aspirated 20 μ L	Diluent 180 μ L
Decreased 4.5 μ L	Reagent Vol R1 150 μ L	R2 75 μ L
Increased	Sample Blank <input type="checkbox"/>	Auto Rerun <input checked="" type="checkbox"/>

Linearity range (Standard)	13	1000	Linearity Limit	0.2
Linearity Range (Decreased)			Substrate Depletion	35000
Linearity Range (Increased)			Mixed Blank Abs	-35000 35000
R1 Blank Abs	-35000	35000	On-board Stability	Day(s)
Blank Response	-35000	35000	Reagent Alarm Limit	
Twin Chemistry			Enzyme Linear Extension	<input type="checkbox"/>
Prozone Check	<input type="checkbox"/>			
Q1		Q2	V1	Q3
Q5		Q6	V3	PC1
Sample Pretreatment		Control Pretreatment		Calibrator Pretreatment
Pretreat Sample Vol	μ L	Pretreat Sample Vol	μ L	

CALIBRATION SETTINGS		AUTO CALIBRATION
Math model Spline		Bottle Changed
Factor	Replicates 2	Lot Changed
		Cal Time

ACCEPTANCE LIMITS	
Cal Time	Hour
Slope Diff	SD
Sensitivity	Repeatability 35000
Deter Coeff	

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