## ACCENT-200 ALPHA - FETOPROTEIN

Nr kat. 7-232 (PL)

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Cat. No 7-232 (EN)

### INTENDED USE

Diagnostic kit for determination of  $\alpha$ -fetoprotein concentration used in automatic analysers ACCENT-200 and ACCENT-200 II GEN.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

$\alpha$ -fetoprotein (AFP) is fetoprotein with a molecular weight of approximately 70 kD containing about 3 % sugar. While it is present in high concentrations during fetal growth, its concentration rapidly decrease after birth and is present at extremely low levels in normal human blood.

AFP shows a notable increase in primary hepatic cancer and is considered to be of great diagnostic importance. It is also thought that fluctuations in blood AFP are useful for evaluating the progress, effects of therapy, and postoperative prognosis of hepatoma.

### METHOD PRINCIPLE

When an antigen-antibody reaction occurs between AFP in a sample and anti-AFP antibody which has been sensitized to latex particles, agglutination results. This agglutination is detected as an absorbance change, with the magnitude of the change being proportional to the quantity of AFP in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

### REAGENTS

#### Package

1-Reagent 1 x 28 ml  
2-Reagent 1 x 14 ml

The reagents when stored at 2-10°C are stable up to expiry date printed on the package.

#### Concentrations in the test

suspension of latex particles sensitized with (rabbit) anti-AFP antibodies (pH 7.3) 0.12 w/v% glycine buffer solution (pH 8.3) preservative

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Reagent bottles should be shaken before use by gently inverting several times.
- After measurements are taken, reagent bottles should be capped and kept at 2-10°C. Care should be taken not to interchange the caps of reagent bottles.
- Reagents with different lot numbers should not be interchanged or mixed.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

### SPECIMEN

Serum.

After blood has clotted thoroughly, the sample is centrifuged and the serum is separated from blood cells and fibrins. Samples can be stored for several weeks at 2-8°C or for 1 year at -20°C. Repeated freezing and thawing should be avoided.

Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.

For reagent blank 0.9% NaCl is recommended.

### REFERENCE VALUES <sup>5</sup>

serum	< 15 ng/ml
-------	------------

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples the CORMAY IMMUNO-CONTROL I (Cat. No 4-288).

For the calibration of automatic analysers systems the CORMAY AFP CALIBRATORS kit (Cat. No 4-282) is recommended. As a 0 calibrator 0.9% NaCl should be used.

The calibration curve should be prepared with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using an automatic analyser Hitachi 917. Results may vary if a different instrument is used.

- Analytical range:** 7 – 250 ng/ml.

For higher concentrations dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

- Specificity / Interferences**

Haemoglobin up to 0.3 g/dl, bilirubin up to 30 mg/dl and triglycerides up to 300 mg/dl do not interfere with the test.

- Precision**

Repeatability (run to run) n = 20	Mean [ng/ml]	SD [ng/ml]	CV [%]
level 1	9.9	0.4	4.03
level 2	22.6	0.3	1.37
level 3	96.5	0.7	0.71

- Method comparison**

A comparison between CORMAY reagent (y) and commercially available assay (x) using 78 samples gave following results:

$$y = 1.01 x + 16.73 \text{ ng/ml};$$

$$R = 0.996 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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