ACCENT-200 HBDH

Nr kat. 7-241 (PL)

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Cat. No **7-241** (EN)

INTENDED USE

Diagnostic kit for determination of α -hydroxybutyrate dehydrogenase activity used in automatic analysers ACCENT-200 and ACCENT-200 II GEN.


The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Lactate dehydrogenase (LDH, LD) is a tetrameric molecule containing two possible forms of subunits (H and M). The result is five isoenzymes, one of which is hydroxybutyrate dehydrogenase (HBDH, LD-1) formed by four H subunits. HBDH is present mainly in heart muscle but occur also in kidney and erythrocytes. Normal serum contains mostly LD-2 with lesser amount of LD-1. Changes in the ratio of LD-1 to LD-2 indicate myocardial infarction or hemolysis.

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The rate of absorbance changing at $\lambda=340$ nm is directly proportional to α -hydroxybutyrate dehydrogenase activity.

REAGENTS

Package	
1-Reagent	1 x 30 ml
2-Reagent	1 x 8 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. Stability on board of the analyser at 2-10°C: ACCENT-200 – 11 weeks, ACCENT-200 II GEN – 12 weeks.

Concentrations in the test

phosphate buffer (pH 7.5)	50 mmol/l
2-oxybutyrate	3 mmol/l
NADH	0.25 mmol/l
preservative	

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum.
 Do not use hemolyzed blood because erythrocytes contain very high HBDH activity. Do not chill or freeze samples. HBDH activity is unstable and is rapidly lost during storage. Specimens can be stored up to 6 hours at 15-25°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.
 For reagent blank deionized water is recommended.

REFERENCE VALUES ⁷

serum	37°C
adults	< 182 U/l (< 3.04 μ kat/l)

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analyser ACCENT-200 the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

For the calibration of automatic analyser ACCENT-200 II GEN the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

As a 0 calibrator deionised water should be used.

The calibration curve should be prepared every 11 weeks (ACCENT-200) or every 10 weeks (ACCENT-200 II GEN), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers ACCENT-200 and ACCENT-200 II GEN. Results may vary if a different instrument or a manual procedure is used.

• Sensitivity

25.8 U/l (0.43 μ kat/l) – ACCENT-200
 19 U/l (0.32 μ kat/l) – ACCENT-200 II GEN

• Linearity

up to 1000 U/l (16.67 μ kat/l) – ACCENT-200
 up to 1200 U/l (20.0 μ kat/l) – ACCENT-200 II GEN

For higher activity dilute the sample 10-fold with 0.9% NaCl and repeat the assay. Multiply the result by 10.

• Specificity / Interferences

Haemoglobin up to 2.5 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 10000 mg/dl do not interfere with the test.

• Precision (ACCENT-200)

Repeatability (run to run) n = 20	Mean [U/l]	SD [U/l]	CV [%]
level 1	120.05	3.35	2.79
level 2	426.73	2.27	0.53
Reproducibility (day to day) n = 80	Mean [U/l]	SD [U/l]	CV [%]
level 1	225.15	7.71	3.42
level 2	408.36	6.88	1.69

• Precision (ACCENT-200 II GEN)

Repeatability (run to run) n = 20	Mean [U/l]	SD [U/l]	CV [%]
level 1	120.05	3.35	2.79
level 2	426.73	2.27	0.53
Reproducibility (day to day) n = 10	Mean [U/l]	SD [U/l]	CV [%]
level 1	154.49	3.32	2.15
level 2	449.51	10.73	2.39

• Method comparison

A comparison between HBDH activity determined at ACCENT-200 (y) and at PRESTIGE 24i (x) using 42 samples gave following results:

$$y = 0.9387x + 8.8211 \text{ U/l};$$

$$R = 0.9920 \quad (R - \text{correlation coefficient})$$

A comparison between HBDH activity determined at ACCENT-200 II GEN (y) and at OLYMPUS AU400 (x) using 30 samples gave following results:

$$y = 0.9833x + 0.2661 \text{ U/l};$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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