ACCENT-200 URINE PROTEINS

Nr kat. 7-242 (PL)

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Cat. No **7-242** (EN)

INTENDED USE

Diagnostic kit for determination of total proteins in urine and cerebrospinal fluid intended to use in automatic analysers: ACCENT-200 (II GEN), ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, ACCENT 400, ACCENT 400, ACCENT Neo200 and BS-120.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

In healthy people with properly functioning kidneys, proteins are actively reabsorbed in the proximal tubules and only small amounts of proteins (several mg per day) are excreted in urine. The measurement of total proteins concentration in urine is used in the diagnosis and treatment of heart and thyroid diseases, which are characterized by proteinuria or albuminuria.

The measurement of total proteins concentration in cerebrospinal fluid (CSF) is especially useful in detecting increased permeability of the blood-brain barrier and in detecting increased intrathecal synthesis of immunoglobulins. Increased concentration of protein in CSF may indicate brain tumors, intracerebral hemorrhage, brain injury, bacterial and viral encephalitis and multiple sclerosis.

METHOD PRINCIPLE

Direct, colorimetric method with pyrogallol red.

At an acidic pH the protein aminoacid groups, with the pyrogallol red-molybdate complex, form a coloured compound which shows a maximum absorbance at a wavelength of 600 nm. Colour intensity is proportional to the concentration of proteins in the sample.

REAGENTS

Package 2 x 31 ml
1-REAGENT

The reagent when stored at 15-25°C is stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 12 weeks (ACCENT-200 (II GEN), ACCENT MC240, ACCENT S120).

Concentrations in the reagent

Succinate buffer ≤ 60 mmol/l
Pyrogallol red ≤ 0.07 mmol/l
Sodium molybdate ≤ 0.05 mmol/l
Stabilizers, preservative

WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- Do not freeze the reagent.
- Reagent should be mixed before use by gentle inverting the bottle several times.
- The appearance of turbidity or control urine values outside the manufacturer's acceptable range may indicate of reagent instability.
- EUH210 Safety data sheet available on request.

SPECIMEN

Urine: Urine used for analysis may come from the first morning sample, random sample or timed collection sample according to the classification by ECLM².

In order to collect and prepare the samples only dedicated tubes and containers should be used.

Do not use preservatives.

Freshly collected urine should be kept at room temperature for about an hour and then cooled to 4°C. Direct reduction of the temperature on freshly collected urine can cause precipitation of minerals.

Samples with visible turbidity should be centrifuged before analysis. Determination of uncentrifuged samples may give increased results.

Urine samples are stable for 2 days at 2-8°C⁶.

Nevertheless it is recommended to perform the assay with freshly collected samples!

Cerebrospinal fluid: Cerebrospinal fluid should be centrifuged before analysis. The presence of blood in samples of cerebrospinal fluid may result in false results of protein determination. For proper interpretation of the results, cerebrospinal fluid must be determined simultaneously with a sample of blood taken from a patient at the same time.

The stability of the cerebrospinal fluid is 3 days at 2-8°C, 6 months at the temperature -20°C⁶.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-REAGENT is ready to use.

0.9% NaCl is recommended as a reagent blank.

Actions required:

When performing assays in analysers: ACCENT-200 (II GEN), ACCENT-220S and BS-120 there is a probability of **cross-contamination** affecting the tests results: URINE PROTEINS II GEN - CREATININE, URINE PROTEINS II GEN - CREATININE ENZYMATYC, URINE PROTEINS II GEN - TG mono, URINE PROTEINS II GEN - UA / UA PLUS, GGT - URINE PROTEINS II GEN, CREATININE ENZYMATYC - URINE PROTEINS II GEN. To avoid this effect follow the recommendations contained in the instruction 51_03_24_001_ACCENT-200_CARRYOVER.

Calculation

For the calculation of proteins excreted over 24 hours, multiply the concentration (mg/dl) by the volume (dl) of the 24 hours urine.

REFERENCE VALUES^{7,10}

urine (adults)	< 15 mg/dl (0.15 g/l)	
urine 24-h (adults)	< 100 mg (0.10 g)	
cerebrospinal fluid	mg/dl	g/l
children 0 - 4 weeks	20 - 80	0.20 - 0.80
children >4 weeks, adults	15 - 45	0.15 - 0.45

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples the CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162). For the calibration of automatic analysers the CORMAY URINE PROTEINS CALIBRATORS (Cat. No 5-181) is recommended. 0.9% NaCl should be used as a calibrator 0.

The calibration curve should be prepared every 12 weeks (ACCENT-200 (II GEN), ACCENT MC240, ACCENT S120) with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers ACCENT-200 (II GEN) and ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

■ LOQ

4 mg/dl (0.04 g/l) – ACCENT-200 (II GEN)

1.85 mg/dl (0.02 g/l) – ACCENT MC240

■ Linearity

up to 130 mg/dl (1.3 g/l) - ACCENT-200 (II GEN)

up to 183 mg/dl (1.8 g/l) – ACCENT MC240

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

■ Specificity / Interferences

Hemoglobin up to 0.004 g/dl, ascorbic acid up to 20 mg/dl, creatinine up to 6 g/l, bilirubin up to 5 mg/dl, conjugated bilirubin up to 60 mg/dl, uric acid up to 85 mg/dl, glucose up to 35 g/l, citrates up to 250 mg/dl, oxalates up to 90 mg/dl, calcium ions up to 130 mg/dl, magnesium ions up to 1.8 g/l, phosphate ions up to 1.2 g/l, urea up to 50 g/l do not affect the results of the determination.

High concentration of iron (II) ions in the test sample may cause interferences¹².

Acetaminophen and some antibiotics from penicillins and aminoglycosides can interfere^{4,11-13}.

■ Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 (II GEN) n=20	level 1	25.1	0.16	0.63
	level 2	68.6	0.91	1.32
ACCENT MC240 n=20	level 1	20.9	0.11	0.5
	level 2	59.4	0.24	0.4
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 (II GEN) n=80	level 1	20.8	0.88	4.3
	level 2	59.0	1.72	2.9
ACCENT MC240 n=80	level 1	17.4	0.68	3.9
	level 2	59.2	1.85	3.1

■ Method comparison

A comparison between total proteins values determined at **ACCENT-200 (II GEN)** (y) and at **BECKMAN COULTER AU680** (x) using 60 samples of urine gave following results:

$y = 0.9357 x + 8.3994$ mg/dl;

R = 0.989 (R – correlation coefficient)

A comparison between total proteins values determined at **ACCENT-200 (II GEN)** (y) and at **BECKMAN COULTER AU680** (x) using 60 samples of CSF gave following results:

$y = 0.9521 x + 1.7825$ mg/dl;

R = 0.999 (R – correlation coefficient)

A comparison between total proteins values determined at **ACCENT MC240** (y) and at **ADVIA 1800** (x) using 65 samples of urine gave following results:

$y = 1.1266 x - 2.2936$ mg/dl;

R = 0.999 (R – correlation coefficient)

A comparison between total proteins values determined at **ACCENT MC240** (y) and at **ADVIA 1800** (x) using 66 samples of CSF gave following results:

$y = 1.1637 x - 3.8049$ mg/dl;

R = 0.999 (R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

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Calibration Rule

Rule	Two-point Linear
Sensitivity	1
Replicates	2
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT-220S

Parameters

Test	PROT U
No	55
Full Name	Urine Proteins
Standard No	55
Reac. Type	Endpoint
Pri. Wave.	630 nm
Sec. Wave.	
Direction	Increase
Reac. Time	0 11
Incuba. Time	
Unit	mg/dl
Precision	0.01

R1	200
R2	
Sample Volume	4
R1 Blank	
Mixed Rtg. Blank Concentration	
Linearity Range	3.3 145
Linearity Limit	
Substrate Limit	
Factor	
<input type="checkbox"/> Prozone check	
q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
PC <input type="checkbox"/> Abs <input type="checkbox"/>	

Calibration Rule

Rule	Two-point Linear
Sensitivity	1
Replicates	3
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• BS-120

Parameters

Test	PROT U
No	55
Full Name	Urine Proteins
Standard No	55
Reac. Type	Endpoint
Pri. Wave.	630 nm
Sec. Wave.	
Direction	Increase
Reac. Time	0 8
Incuba. Time	
Unit	mg/dl
Precision	0.1

R1	180
R2	
Sample Volume	4
R1 Blank	
Mixed Rtg. Blank Concentration	1.7 137
Linearity Limit	
Substrate Limit	
Factor	
<input type="checkbox"/> Prozone check	
q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
PC <input type="checkbox"/> Abs <input type="checkbox"/>	

Calibration Rule

Rule	Two-point Linear
Sensitivity	1
Replicates	3
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT S120

Chem	UP	No.	055	Sample Type	URINE/CSF		
Chemistry	URINE PROTEINS	Print name	UP	Reaction Direction	positive		
Reaction Type	Endpoint	Sec Wave		Decimal	0.1		
Pri Wave	578 nm	Incubation Time	0	Reaction Time	5 7		
Unit	mg/dl	Blank Time	-3 -1	Standard	Sample Vol 2 μL Aspirated μL Diluent μL Reagent Vol R1 160 μL R2 μL		
Decreased	2 μL 20 μL 180 μL	Increased	μL μL μL	Sample Blank	<input type="checkbox"/>	Auto Rerun	<input checked="" type="checkbox"/>
Linearity range (Standard)	2	274	Linearity Limit				
Linearity Range (Decreased)			Substrate Depletion				
Linearity Range (Increased)			Mixed Blank Abs	-40000 40000			
R1 Blank Abs	-40000	40000	On-board Stability	84 Day(s)			
Blank Response	-40000	40000	Reagent Alarm Limit				
Twin Chemistry		<input type="checkbox"/> Prozone Check	Enzyme Linear Extension	<input type="checkbox"/>			
Q1 <input type="checkbox"/> Q2 <input type="checkbox"/> V1 <input type="checkbox"/> Q3 <input type="checkbox"/> Q4 <input type="checkbox"/> V2 <input type="checkbox"/>		Q5 <input type="checkbox"/> Q6 <input type="checkbox"/> V3 <input type="checkbox"/> PC1 <input type="checkbox"/> PC2 <input type="checkbox"/>					
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment	<input type="checkbox"/> Calibrator Pretreatment	Pretreat Sample Vol <input type="checkbox"/> μL	Pretreat Sample Vol <input type="checkbox"/> μL			
CALIBRATION SETTINGS	Math model	Two-point linear	AUTO CALIBRATION	<input type="checkbox"/> Bottle Changed			
Factor <input type="checkbox"/>	Replicates	2	<input type="checkbox"/> Lot Changed	<input type="checkbox"/> Cal Time			
ACCEPTANCE LIMITS	Cal Time	2016 Hour	Slope Diff <input type="checkbox"/>	SD <input type="checkbox"/>			
Sensitivity <input type="checkbox"/>	Repeatability	40000	Deter Coeff <input type="checkbox"/>				

ACCENT-200 URINE PROTEINS

• ACCENT MC240

Chem	UP	No.	055	Sample Type	URINE / CSF
Chemistry	URINE PROTEINS	Print name	UP		
Reaction Type	Endpoint	Reaction Direction	positive		
Pri Wave	605 nm	Sec Wave			
Unit	mg/dl	Decimal	0.1		
Blank Time	-3	-1			
Incubation Time	0				
Reaction Time	28	30			
Standard	Sample Vol 3.2 μL	Aspirated	μL	Diluent	μL
Decreased	3.2 μL	20 μL		180 μL	
Increased	μL	μL		μL	
	<input type="checkbox"/> Sample Blank	<input checked="" type="checkbox"/> Auto Rerun			

Linearity range (Standard)	1.85	183	Linearity Limit								
Linearity Range (Decreased)			Substrate Depletion								
Linearity Range (Increased)			Mixed Blank Abs	-35000 35000							
R1 Blank Abs	-35000	35000	On-board Stability	84 Day(s)							
Blank Response	-35000	35000	Reagent Alarm Limit								
Twin Chemistry			<input type="checkbox"/> Enzyme Linear Extension								
<input type="checkbox"/> Prozone Check											
Q1		Q2		V1		Q3		Q4		V2	
Q5		Q6		V3		PC1		PC2			
<input type="checkbox"/> Sample Pretreatment		<input type="checkbox"/> Control Pretreatment		<input type="checkbox"/> Calibrator Pretreatment							
		Pretreat Sample Vol	μL	Pretreat Sample Vol	μL						

CALIBRATION SETTINGS		AUTO CALIBRATION	
Math model	Two-point linear	<input type="checkbox"/> Bottle Changed	
Factor		<input type="checkbox"/> Lot Changed	
	Replicates 2	<input type="checkbox"/> Cal Time	

ACCEPTANCE LIMITS			
Cal Time	2016 Hour		
Slope Diff		SD	
Sensitivity		Repeatability	35000
Deter Coeff			

• ACCENT M320

Chem	UP	No.	055	Sample Type	URINE / CSF
Chemistry	URINE PROTEINS	Print name	UP		
Reaction Type	Endpoint	Reaction Direction	positive		
Pri Wave	605 nm	Sec Wave			
Unit	mg/dl	Decimal	0.1		
Blank Time	-3	-1			
Incubation Time	0				
Reaction Time	38	40			
Standard	Sample Vol 2.8 μL	Aspirated	μL	Diluent	μL
Decreased	2.8 μL	20 μL		180 μL	
Increased	μL	μL		μL	
	<input type="checkbox"/> Sample Blank	<input checked="" type="checkbox"/> Auto Rerun			

Linearity range (Standard)	2	164	Linearity Limit								
Linearity Range (Decreased)			Substrate Depletion								
Linearity Range (Increased)			Mixed Blank Abs	-35000 35000							
R1 Blank Abs	-35000	35000	On-board Stability	84 Day(s)							
Blank Response	-35000	35000	Reagent Alarm Limit								
Twin Chemistry			<input type="checkbox"/> Enzyme Linear Extension								
<input type="checkbox"/> Prozone Check											
O1		O2		V1		O3		O4		V2	
O5		O6		V3		PC1		PC2			
<input type="checkbox"/> Sample Pretreatment		<input type="checkbox"/> Control Pretreatment		<input type="checkbox"/> Calibrator Pretreatment							
		Pretreat Sample Vol	μL	Pretreat Sample Vol	μL						

CALIBRATION SETTINGS		AUTO CALIBRATION	
Math model	Two-point linear	<input type="checkbox"/> Bottle Changed	
Factor		<input type="checkbox"/> Lot Changed	
	Replicates 2	<input type="checkbox"/> Cal Time	

ACCEPTANCE LIMITS			
Cal Time	2016 Hour		
Slope Diff		SD	
Sensitivity		Repeatability	35000
Deter Coeff			

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