ACCENT-200 MICROALBUMIN

Nr kat. **7-244** (PL)

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Cat. No **7-244** (EN)

INTENDED USE

Diagnostic kit for determination of albumin concentration in urine and cerebrospinal fluid intended to use in automatic analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT MC240, ACCENT M320, ACCENT 400 and ACCENT Neo200.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Albumin is a protein that is formed within the liver and it makes up approximately 60% of the serum protein. Normally only small amounts of albumin are filtered through the renal glomeruli, and that small quantity can be reabsorbed by the renal tubules. In that case there is a low albumin concentration in the urine. When renal disorders appear, level of urine albumin increase but remains still undetectable by routine screening tests (microalbuminuria). The appearance of low but abnormal levels (30-300 mg/24h) of albumin in the urine is an early clinical evidence of nephropathy (mostly diabetic) and cardiovascular disorders.

To avoid the necessity of 24-hour urine collection it is common in clinical practice to measure albumin and creatinine simultaneously and give the result as a albumin/creatinine ratio.

METHOD PRINCIPLE

Immunoturbidimetric method. Albumin in the sample forms with anti-albumin antibodies in the reagent an insoluble complex. The turbidity caused by the complexes is measured spectrophotometrically at 340 nm and is proportional to the amount of albumin in the sample.

REAGENTS

Package

1-Reagent	2 x 12 ml
2-Reagent	2 x 3.2 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyzer at 2-10°C are stable for 11 weeks (ACCENT-200) or 12 weeks (ACCENT M320, ACCENT S120).

Concentrations in the test

1-Reagent

Tris buffer (pH 7.6)	18.2 mmol/l
sodium chloride	123.2 mmol/l
PEG	< 4%

2-Reagent

sodium chloride	154 mmol/l
anti-human albumin antibodies preservatives	

Warnings and notes

- Protect from direct sunlight and avoid contamination.
- Do not freeze the reagents.
- Do not use after expiry date.

- Do not interchange caps.
- Reagents with different lot numbers should not be interchanged or mixed.
- Before use mix reagent by gently inverting each bottle.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product

SPECIMEN

Urine. Urine used for analysis may come from the first morning sample, random sample or timed collection sample [3]. Samples with visible turbidity should be centrifuged before analysis. Determination of uncentrifuged samples may give increased results.

Urine samples are stable for 2 days at room temperature, 14 days at the 8°C [7]. Nevertheless it is recommended to perform the assay with freshly collected samples!

Cerebrospinal fluid. CSF should be centrifuged before analysis. If the total protein in CSF is greater than 2000 mg/l, the CSF sample needs to be diluted 1:9 and the result multiplied by 10.

It is recommended to follow NCCLS procedures regarding specimen collecting and handling.

Samples should be stored at 2-4°C and analyzed within 2 hours after collection.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use. 0.9% NaCl is recommended as a reagent blank.

For analyzers: ACCENT-200, ACCENT-200 II GEN and ACCENT-220S it is recommended to determine the reagent blank during each calibration. Deionized water should be used as reagent blank. When performing calibration, the task type **Calib+Rgt.Blk** should be selected.

Actions required:

When performing assays in analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, there is a probability of **cross-contamination** affecting the tests results: AMYLASE - MICROALBUMIN, ALBUMIN - MICROALBUMIN, ALP - MICROALBUMIN, BIL TOTAL MALLOY-EVELYN - MICROALBUMIN, CHOLINESTERASE - MICROALBUMIN, LDL DIRECT - MICROALBUMIN, CREATININE - MICROALBUMIN, GGT - MICROALBUMIN, LDH - MICROALBUMIN, dTIBC - MICROALBUMIN, TOTAL PROTEIN - MICROALBUMIN, TG - MICROALBUMIN. To avoid this effect follow the recommendations contained in the instruction 51_03_24_001_ACCENT-200_CARRYOVER.

Calculation

For the calculation of albumin 24 hours quantity, multiply the concentration (mg/l) with the volume (l) of the 24 hours urines.

REFERENCE VALUES³

urine	mg/24h	µg/min	mg/g creatinine
normal	< 30	< 20	< 30
microalbuminuria	30 – 300	20 – 200	30 – 300
clinical albuminuria (overt nephropathy)	> 300	> 200	> 300
cerebrospinal fluid, lumbar	177 – 251 mg/l		

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY MICROALBUMIN CONTROL (Cat. No 4-461) with each batch of samples.

For the calibration of automatic analyzers the CORMAY MICROALBUMIN CALIBRATOR (Cat. No 5-193) is recommended. 0.9% NaCl should be used as a calibrator 0.

The calibration curve should be prepared every 4 weeks (ACCENT-200), 6 weeks (ACCENT S120) or 7 weeks (ACCENT M320), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyzers: ACCENT-200 and ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity** : 3.8 mg/l – ACCENT-200
- LOQ**: 2 mg/l – ACCENT MC240
- Linearity**: up to concentration of highest calibrator

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Hemoglobin up to 2.5 g/dl, ascorbate up to 200 mg/dl, creatinine up to 6 g/l, uric acid up to 100 mg/dl, glucose up to 35 g/l, urea up to 50 g/l, bilirubin conjugated up to 60 mg/dl, calcium ion up to 130 mg/dl, magnesium ion up to 1.8 g/l, do not interfere with the test.

Precision

Repeatability (run to run)		Mean [mg/l]	SD [mg/l]	CV [%]
ACCENT-200 n=10	level 1	19.06	0.54	2.84
	level 2	77.72	0.87	1.12
ACCENT MC240 n=20	level 1	21,46	0,36	1,67
	level 2	58,86	0,87	1,48
Reproducibility (day to day)		Mean [mg/l]	SD [mg/l]	CV [%]
ACCENT-200 n=20	level 1	68.02	1.87	2.74
	level 2	74.82	0.98	1.31
ACCENT MC240 n=80	level 1	28.2	1.38	4.9
	level 2	67.3	1.46	2.2

Method comparison

A comparison between microalbumin values determined at **ACCENT-200** (y) and at **COBAS INTEGRA 400** (x) using 50 samples of urine gave following results:

$$y = 0.9868x - 3.5633 \text{ mg/l;} \\ R = 0.998 \quad (R - \text{correlation coefficient})$$

A comparison between CORMAY reagent (y) and another commercially available assay (x) using 29 samples of cerebrospinal fluid, at analyzers **ACCENT-200** and **BS-400** gave following results:

$$y = 0.9785x + 8.4498 \text{ mg/dl;} \\ R = 0.998 \quad (R - \text{correlation coefficient})$$

A comparison between microalbumin values determined at **ACCENT MC240** (y) and at **BS-400** (x) using 67 samples of urine gave following results:

$$y = 0.955x + 5.8166 \text{ mg/l;} \\ R = 1.000 \quad (R - \text{correlation coefficient})$$

A comparison between microalbumin values determined at **ACCENT MC240** (y) and at **BS-400** (x) using 67 samples of cerebrospinal fluid, gave following results:

$$y = 0.955x + 5.8166 \text{ mg/l;} \\ R = 1.000 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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ACCENT-200 MICROALBUMIN

• ACCENT MC240

Chem <input type="text" value="MALB"/>	No. <input type="text" value="072"/>	Sample Type <input type="text" value="URINE/CSF"/>
Chemistry <input type="text" value="MICROALBUMIN"/>	Print name <input type="text" value="MALB"/>	
Reaction Type <input type="text" value="Endpoint"/>	Reaction Direction <input type="text" value="positive"/>	
Pri Wave <input type="text" value="340nm"/>	Sec Wave <input type="text"/>	
Unit <input type="text" value="mg/l"/>	Decimal <input type="text" value="0,1"/>	
Blank Time <input type="text" value="-3"/> <input type="text" value="-1"/>	Incubation Time <input type="text" value="21"/>	Reaction Time <input type="text" value="31"/> <input type="text" value="32"/>
Standard <input type="text" value="7.5"/> <input type="text" value="7.5"/> <input type="text" value="20"/> <input type="text" value="180"/>	Aspirated <input type="text" value="20"/> <input type="text" value="180"/>	Diluent <input type="text" value="180"/>
Decreased <input type="text" value="7.5"/> <input type="text" value="20"/> <input type="text" value="180"/>	Reagent Vol R1 <input type="text" value="150"/> <input type="text" value="30"/>	R2 <input type="text" value="30"/>
Increased <input type="text" value="7.5"/> <input type="text" value="20"/> <input type="text" value="180"/>	Sample Blank <input type="text" value="V"/>	Auto Rerun <input type="text"/>

Linearity range (Standard) <input type="text" value="2"/> <input type="text" value="320"/>	Linearity Limit <input type="text"/>
Linearity Range (Decreased) <input type="text"/>	Substrate Depletion <input type="text"/>
Linearity Range (Increased) <input type="text"/>	Mixed Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>
R1 Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>	On-board Stability <input type="text"/> Day(s)
Blank Response <input type="text" value="-35000"/> <input type="text" value="35000"/>	Reagent Alarm Limit <input type="text"/>
Twin Chemistry <input type="text"/>	Enzyme Linear Extension <input type="text"/>
<input type="text"/> Prozone Check	
Q1 <input type="text"/> Q2 <input type="text"/> V1 <input type="text"/> Q3 <input type="text"/> Q4 <input type="text"/> V2 <input type="text"/>	
Q5 <input type="text"/> Q6 <input type="text"/> V3 <input type="text"/> PC1 <input type="text"/> PC2 <input type="text"/>	
<input type="text"/> Sample Pretreatment	<input type="text"/> Control Pretreatment
<input type="text"/> Pretreat Sample Vol <input type="text"/> μ L	<input type="text"/> Calibrator Pretreatment
<input type="text"/> Pretreat Sample Vol <input type="text"/> μ L	<input type="text"/> Pretreat Sample Vol <input type="text"/> μ L

CALIBRATION SETTINGS	AUTO CALIBRATION
Math model <input type="text" value="Spline"/>	<input type="text"/> Bottle Changed
Factor <input type="text"/>	<input type="text"/> Lot Changed
Replicates <input type="text" value="2"/>	<input type="text"/> Cal Time

ACCEPTANCE LIMITS
Cal Time <input type="text"/> Hour
Slope Diff <input type="text"/> SD <input type="text"/>
Sensitivity <input type="text"/> Repeatability <input type="text" value="35000"/>
Deter Coeff <input type="text"/>

• ACCENT M320

Chem <input type="text" value="MALB"/>	No. <input type="text" value="072"/>	Sample Type <input type="text" value="URINE/CSF"/>
Chemistry <input type="text" value="MICROALBUMIN"/>	Print name <input type="text" value="MALB"/>	
Reaction Type <input type="text" value="Endpoint"/>	Reaction Direction <input type="text" value="positive"/>	
Pri Wave <input type="text" value="340nm"/>	Sec Wave <input type="text"/>	
Unit <input type="text" value="mg/l"/>	Decimal <input type="text" value="0,1"/>	
Blank Time <input type="text" value="-3"/> <input type="text" value="-1"/>	Incubation Time <input type="text" value="21"/>	Reaction Time <input type="text" value="33"/> <input type="text" value="34"/>
Standard <input type="text" value="7.5"/> <input type="text" value="7.5"/> <input type="text" value="20"/> <input type="text" value="180"/>	Aspirated <input type="text" value="20"/> <input type="text" value="180"/>	Diluent <input type="text" value="180"/>
Decreased <input type="text" value="7.5"/> <input type="text" value="20"/> <input type="text" value="180"/>	Reagent Vol R1 <input type="text" value="150"/> <input type="text" value="30"/>	R2 <input type="text" value="30"/>
Increased <input type="text" value="7.5"/> <input type="text" value="20"/> <input type="text" value="180"/>	Sample Blank <input type="text" value="V"/>	Auto Rerun <input type="text"/>

Linearity range (Standard) <input type="text" value="9.5"/> <input type="text" value="358"/>	Linearity Limit <input type="text"/>
Linearity Range (Decreased) <input type="text"/>	Substrate Depletion <input type="text"/>
Linearity Range (Increased) <input type="text"/>	Mixed Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>
R1 Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>	On-board Stability <input type="text"/> Day(s)
Blank Response <input type="text" value="-35000"/> <input type="text" value="35000"/>	Reagent Alarm Limit <input type="text"/>
Twin Chemistry <input type="text"/>	Enzyme Linear Extension <input type="text"/>
<input type="text"/> Prozone Check	
Q1 <input type="text"/> Q2 <input type="text"/> V1 <input type="text"/> Q3 <input type="text"/> Q4 <input type="text"/> V2 <input type="text"/>	
Q5 <input type="text"/> Q6 <input type="text"/> V3 <input type="text"/> PC1 <input type="text"/> PC2 <input type="text"/>	
<input type="text"/> Sample Pretreatment	<input type="text"/> Control Pretreatment
<input type="text"/> Pretreat Sample Vol <input type="text"/> μ L	<input type="text"/> Calibrator Pretreatment
<input type="text"/> Pretreat Sample Vol <input type="text"/> μ L	<input type="text"/> Pretreat Sample Vol <input type="text"/> μ L

CALIBRATION SETTINGS	AUTO CALIBRATION
Math model <input type="text" value="Spline"/>	<input type="text"/> Bottle Changed
Factor <input type="text"/>	<input type="text"/> Lot Changed
Replicates <input type="text" value="2"/>	<input type="text"/> Cal Time

ACCEPTANCE LIMITS
Cal Time <input type="text"/> Hour
Slope Diff <input type="text"/> SD <input type="text"/>
Sensitivity <input type="text"/> Repeatability <input type="text" value="35000"/>
Deter Coeff <input type="text"/>

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