## ACCENT-200 CALCIUM ARSENAZO

Nr kat. 7-247

(PL)

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Cat. No 7-247 (EN)

### INTENDED USE

Diagnostic kit for determination of calcium concentration, used in automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240 ACCENT M320 and BS-120 / BS-130.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Calcium plays an essential role in many cell functions: intracellularly in muscle contraction and glycogen metabolism, extracellularly, in bone mineralization, in blood coagulation and in transmission of nerve impulses. Calcium is present in plasma in three forms: free, bound to proteins or complexed with anions as phosphate, citrate and bicarbonate. Decreased total calcium levels can be associated with diseases of the bone apparatus (especially osteoporosis), kidney diseases (especially under dialysis), defective intestinal absorption and hypoparathyroidism. Increased total calcium can be measured in hyperparathyroidism, malignant diseases with metastases and sarcoidosis. Calcium measurements also help in monitoring of calcium supplementation mainly in the prevention of osteoporosis.

### METHOD PRINCIPLE

Photometric test using arsenazo III.

Calcium with arsenazo III at neutral pH yields a blue colored complex, whose intensity is proportional to the calcium concentration. Interference by magnesium is eliminated by addition of 8-hydroxyquinoline-5-sulfonic acid.

### REAGENTS

#### PACKAGE

1-Reagent 2 x 36 ml

The reagent, stored at 2-8°C is stable up to expiry date printed on the package.

The reagent stored on board of the analyzer at 2-10°C is stable for 12 weeks (ACCENT-200, ACCENT MC240).

#### Concentrations in the test

phosphate buffer (pH 7.5) ≤ 80 mmol/l  
 8-hydroxyquinoline-5-sulfonic acid ≤ 6 mmol/l  
 arsenazo III ≤ 120 μmol/l  
 detergents, preservative

#### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Do not freeze the reagent.
- Contaminated glassware is the greatest source of error. The use of disposable plastic ware is recommended. Glassware should be soaked for a few hours in 2M HCl solution and then thoroughly rinsed with distilled water.

### SPECIMEN

Serum. Random or 24-hours urine.

Serum can be stored up to 7 days at 20-25°C or up to 3 weeks at 4-8°C. Samples frozen at -20°C can be stored up to 8 months.

Discard contaminated specimens.

**24-hours urine preparation:** To prevent calcium salt precipitation specimens should be collected in 10 ml of 6M HCl. In case of presence of precipitants they can be solved by lowering pH of the urine to below 2.0. Prior to determination dilute the sample with 0.9% NaCl in the ratio of 1 to 1. Multiply the result by the dilution factor. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent is ready to use.

Deionized water is recommended as a reagent blank.

For analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S and BS-120 / BS-130 it is recommended to determine the reagent blank during each calibration. Deionized water should be used as reagent blank. When performing calibration, the task type **Calib+Rgt.Blk** should be selected

#### Actions required:

When performing assays in analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S and BS-120 / BS-130, there is a probability of **cross-contamination** affecting the tests results: CALCIUM ARSENAZO – PHOSPHORUS, RF - CALCIUM ARSENAZO, CALCIUM ARSENAZO - LIPASE. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_001\_ACCENT-200\_CARRYOVER.

### REFERENCE VALUES <sup>4</sup>

serum	mg/dl	mmol/l
adult	8.6 – 10.3	2.15 – 2.57
random urine	mg/dl	mmol/l
male	0.9 – 37.9	0.225 – 9.47
female	0.5 – 35.7	0.125 – 8.92
24-hours urine	mg/24h	mmol/24h
adult	100 – 300	2.5 – 7.5

It is recommended for each laboratory to establish its own reference ranges for local population.

#### Calcium concentration 24-hours urine – calculation

calcium concentration in 24-hours urine [mg/24h]	=	calcium concentration in sample of 24-hours urine [mg/dl]	×	urine volume of 24-hours urine [dl/24h]
--	---	---	---	---

### QUALITY CONTROL

For internal quality control it is recommended to use the following controls with each batch of samples: CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum; CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine. For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 12 weeks (ACCENT-200, ACCENT MC240), with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers ACCENT-200 and ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

#### LoB (Limit of Blank):

0.00 mg/dl (0.00 mmol/l) - ACCENT-200  
 0.00 mg/dl (0.00 mmol/l) – ACCENT MC240

#### LoD (Limit of Detection):

0.08 mg/dl (0.02 mmol/l) - ACCENT-200  
 0.03 mg/dl (0.0075 mmol/l) – ACCENT MC240

#### LoQ (Limit of Quantitation):

0.34 mg/dl (0.085 mmol/l) - ACCENT-200  
 0.33 mg/dl (0.0825 mmol/l) – ACCENT MC240

#### Linearity:

up to 22 mg/dl (5.5 mmol/l) – ACCENT-200  
 up to 20.1 mg/dl (5.025 mmol/l) – ACCENT MC240

For higher concentrations dilute sample with 0.9% NaCl and repeat the assay. Multiply the result by the dilution factor.

#### Specificity / Interferences

Haemoglobin up to 5 g/dl, ascorbate up to 62 mg/l, bilirubin up to 40 mg/dl, triglycerides up to 2000 mg/dl and magnesium up to 20 mg/dl do not interfere with the test.

#### Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 n=20	level 1	9.73	0.08	0.77
	level 2	12.4	0.13	1.06
ACCENT MC240 n=20	level 1	10.1	0.17	1.65
	level 2	12.5	0.06	0.45
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 n=80	level 1	9.8	0.25	2.6
	level 2	12.6	0.24	1.9
ACCENT MC240 n=80	level 1	10.3	0.20	1.9
	level 2	12.7	0.18	1.4

#### Method comparison

A comparison between calcium values determined at **ACCENT-200** (y) and at **BECKMAN COULTER AU680** (x) using 60 serum samples gave following results:

$y = 0.9846x + 0.2834$  mg/dl;  
 $R = 0.975$  (R – correlation coefficient)

A comparison between calcium values determined at **ACCENT-200** (y) and at **BECKMAN COULTER AU680** (x) using 67 urine samples gave following results:

$y = 0.956x + 0.2535$  mg/dl;  
 $R = 0.998$  (R – correlation coefficient)

A comparison between calcium values determined at **ACCENT MC240** (y) and at **BECKMAN COULTER AU680** (x) using 60 serum samples gave following results:

$y = 0.939x + 0.9582$  mg/dl;  
 $R = 0.981$  (R – correlation coefficient)

A comparison between calcium values determined at **ACCENT MC240** (y) and at **BECKMAN COULTER AU680** (x) using 70 urine samples gave following results:

$y = 0.9233x + 0.7828$  mg/dl;  
 $R = 0.997$  (R – correlation coefficient)

### WASTE MANAGEMENT

Please refer to local legal requirements.

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**Calibration Rule**

Rule	Multi-point Linear
Sensitivity	1
Replicates	2
Interval (day)	
Difference Limit	0
SD	0
Blank Response	0   50000
Error Limit	0
Coefficient	0

• ACCENT-200 II GEN

**Parameters**

Test Name	Ca ARS	R1	300
Test No	59	R2	
Full Name	Calcium Arsenazo	Sample Volume	3
Reference No	59	R1 Blank	
Analy. Type	Endpoint	Mixed Reag. Blank	
Pri. Wave.	670 nm	Concentration	0.61   22
Secon. Wave.		Linearity Limit	
Trend	Increase	Substrate Limit	
Reac. Time	0   19	Factor	
Incuba. Time		<input type="checkbox"/> Prozone check	
Unit	mg/dl	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.01	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

**Calibration Rule**

Rule	Multi-point Linear
Sensitivity	1
Replicates	2
Interval (day)	
Difference Limit	0
SD	0
Blank Response	0   50000
Error Limit	0
Coefficient	0

• ACCENT-220S

**Parameters**

Test	Ca ARS	R1	300
No	59	R2	
Full Name	Calcium Arsenazo	Sample Volume	3
Standard No	59	R1 Blank	
Reac. Type	Endpoint	Mixed Rtg. Blank	
Pri. Wave.	670 nm	Linearity Range	0.33   20.4
Sec. Wave.		Linearity Limit	
Direction	Increase	Substrate Limit	
Reac. Time	0   30	Factor	
Incuba. Time		<input type="checkbox"/> Prozone check	
Unit	mg/dl	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.01	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

**Calibration Rule**

Rule	Multi-point Linear
Sensitivity	1
Replicates	2
Interval (day)	
Difference Limit	0
SD	0
Blank Response	0   50000
Error Limit	0
Coefficient	0

• BS-120 / BS-130

**Parameters**

Test	Ca ARS	R1	300
No	67	R2	
Full Name	Calcium Arsenazo	Sample Volume	3
Standard No	67	R1 Blank	
Reac. Type	Endpoint	Mixed Rtg. Blank	
Pri. Wave.	670 nm	Linearity Range	0.31   20.9
Sec. Wave.		Linearity Limit	
Direction	Increase	Substrate Limit	
Reac. Time	0   19	Factor	
Incuba. Time	16	<input type="checkbox"/> Prozone check	
Unit	mg/dl	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.01	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

**Calibration Rule**

Rule	Multi-point Linear
Sensitivity	1
Replicates	2
Interval (day)	
Difference Limit	0
SD	0
Blank Response	0   50000
Error Limit	0
Coefficient	0

• ACCENT S120

Chem	Ca Ars	No.	059	Sample Type	SERUM / URINE
Chemistry	Calcium Arsenazo	Print name	Ca Ars	Reaction Direction	positive
Reaction Type	Endpoint	Sec Wave		Decimal	0.01
Pri Wave	670 nm	Incubation Time	0	Reaction Time	24   26
Unit	mg/dL	Blank Time	-3   -1	Reagent Vol	
Standard	Sample Vol 2 μL Aspirated μL Diluent μL	Decreased	2 μL 20 μL 180 μL	R1	200 μL
Increased		Increased		R2	
				<input type="checkbox"/> Sample Blank <input checked="" type="checkbox"/> Auto Rerun	

Linearity range (Standard)	0.31   19.6	Linearity Limit			
Linearity Range (Decreased)		Substrate Depletion			
Linearity Range (Increased)		Mixed Blank Abs	-40000   40000		
R1 Blank Abs	-40000   40000	On-board Stability		Day(s)	
Blank Response	-40000   40000	Reagent Alarm Limit		<input type="checkbox"/> Enzyme Linear Extension	
Twin Chemistry		<input type="checkbox"/> Prozone Check		Q1 <input type="checkbox"/> Q2 <input type="checkbox"/> V1 <input type="checkbox"/> Q3 <input type="checkbox"/> Q4 <input type="checkbox"/> V2 <input type="checkbox"/>	
		Q5 <input type="checkbox"/> O6 <input type="checkbox"/> V3 <input type="checkbox"/> PC1 <input type="checkbox"/> PC2 <input type="checkbox"/>		<input type="checkbox"/> Sample Pretreatment <input type="checkbox"/> Control Pretreatment <input type="checkbox"/> Calibrator Pretreatment	
		Pretreat Sample Vol <input type="checkbox"/> μL		Pretreat Sample Vol <input type="checkbox"/> μL	

<b>CALIBRATION SETTINGS</b>	<b>AUTO CALIBRATION</b>	
Math model	Multi-point linear	<input type="checkbox"/> Bottle Changed
Factor		<input type="checkbox"/> Lot Changed
Replicates	2	<input type="checkbox"/> Cal Time

<b>ACCEPTANCE LIMITS</b>			
Cal Time	<input type="checkbox"/> Hour		
Slope Diff	<input type="checkbox"/>	SD	<input type="checkbox"/>
Sensitivity	<input type="checkbox"/>	Repeatability	40000
Deter Coeff	<input type="checkbox"/>		

## ACCENT-200 CALCIUM ARSENAZO

### • ACCENT MC240

Chem	<input type="text" value="Ca Ars"/>	No.	<input type="text" value="059"/>	Sample Type	<input type="text" value="SERUM / URINE"/>
Chemistry	<input type="text" value="Calcium Arsenazo"/>	Print name	<input type="text" value="Ca Ars"/>	Reaction Direction	<input type="text" value="positive"/>
Reaction Type	<input type="text" value="Endpoint"/>	Pri Wave	<input type="text" value="660 nm"/>	Sec Wave	<input type="text" value=""/>
Pri Wave	<input type="text" value="660 nm"/>	Unit	<input type="text" value="mg/dL"/>	Decimal	<input type="text" value="0.01"/>
Unit	<input type="text" value="mg/dL"/>	Incubation Time	<input type="text" value="0"/>	Reaction Time	<input type="text" value="28"/> <input type="text" value="30"/>
Blank Time	<input type="text" value="-3"/> <input type="text" value="-1"/>	Standard	<input type="text" value="2"/> <input type="text" value="20"/>	Reagent Vol	<input type="text" value="200"/> <input type="text" value=""/>
Standard	<input type="text" value="2"/> <input type="text" value="20"/>	Aspirated	<input type="text" value="20"/> <input type="text" value="180"/>	R1	<input type="text" value="200"/> <input type="text" value=""/>
Decreased	<input type="text" value="2"/> <input type="text" value="20"/>	Diluent	<input type="text" value="180"/> <input type="text" value=""/>	R2	<input type="text" value=""/> <input type="text" value=""/>
Increased	<input type="text" value=""/> <input type="text" value=""/>	Sample Blank	<input type="checkbox"/>	Auto Rerun	<input checked="" type="checkbox"/>
Linearity range (Standard)	<input type="text" value="0.33"/> <input type="text" value="20.1"/>	Linearity Limit	<input type="text" value=""/>	Substrate Depletion	<input type="text" value=""/>
Linearity Range (Decreased)	<input type="text" value=""/>	Mixed Blank Abs	<input type="text" value="-35000"/> <input type="text" value="35000"/>	On-board Stability	<input type="text" value=""/> Day(s)
Linearity Range (Increased)	<input type="text" value=""/>	Blank Response	<input type="text" value="-35000"/> <input type="text" value="35000"/>	Reagent Alarm Limit	<input type="text" value=""/>
R1 Blank Abs	<input type="text" value="-35000"/> <input type="text" value="35000"/>	Twin Chemistry	<input type="text" value=""/>	Enzyme Linear Extension	<input type="checkbox"/>
Blank Response	<input type="text" value="-35000"/> <input type="text" value="35000"/>	Prozone Check	<input type="checkbox"/>	Q1	<input type="text" value=""/>
Twin Chemistry	<input type="text" value=""/>	Q2	<input type="text" value=""/>	V1	<input type="text" value=""/>
Prozone Check	<input type="checkbox"/>	Q3	<input type="text" value=""/>	Q4	<input type="text" value=""/>
Q1	<input type="text" value=""/>	Q5	<input type="text" value=""/>	V2	<input type="text" value=""/>
Q2	<input type="text" value=""/>	Q6	<input type="text" value=""/>	PC1	<input type="text" value=""/>
Q3	<input type="text" value=""/>	PC2	<input type="text" value=""/>	Sample Pretreatment	<input type="text" value=""/>
Q4	<input type="text" value=""/>	Sample Pretreatment	<input type="text" value=""/>	Control Pretreatment	<input type="text" value=""/>
Q5	<input type="text" value=""/>	Control Pretreatment	<input type="text" value=""/>	Calibrator Pretreatment	<input type="text" value=""/>
Q6	<input type="text" value=""/>	Calibrator Pretreatment	<input type="text" value=""/>	Pretreat Sample Vol	<input type="text" value=""/> <input type="text" value=""/> $\mu$ L
PC1	<input type="text" value=""/>	Pretreat Sample Vol	<input type="text" value=""/> <input type="text" value=""/> $\mu$ L	Pretreat Sample Vol	<input type="text" value=""/> <input type="text" value=""/> $\mu$ L
PC2	<input type="text" value=""/>	Pretreat Sample Vol	<input type="text" value=""/> <input type="text" value=""/> $\mu$ L		
Sample Pretreatment	<input type="text" value=""/>				
Control Pretreatment	<input type="text" value=""/>				
Calibrator Pretreatment	<input type="text" value=""/>				
Pretreat Sample Vol	<input type="text" value=""/> <input type="text" value=""/> $\mu$ L				
Pretreat Sample Vol	<input type="text" value=""/> <input type="text" value=""/> $\mu$ L				

### • ACCENT M320

Chem	<input type="text" value="Ca Ars"/>	No.	<input type="text" value="059"/>	Sample Type	<input type="text" value="SERUM / URINE"/>
Chemistry	<input type="text" value="Calcium Arsenazo"/>	Print name	<input type="text" value="Ca Ars"/>	Reaction Direction	<input type="text" value="positive"/>
Reaction Type	<input type="text" value="Endpoint"/>	Pri Wave	<input type="text" value="660 nm"/>	Sec Wave	<input type="text" value="700 nm"/>
Pri Wave	<input type="text" value="660 nm"/>	Unit	<input type="text" value="mg/dL"/>	Decimal	<input type="text" value="0.01"/>
Unit	<input type="text" value="mg/dL"/>	Blank Time	<input type="text" value="-3"/> <input type="text" value="-1"/>	Incubation Time	<input type="text" value="0"/>
Blank Time	<input type="text" value="-3"/> <input type="text" value="-1"/>	Standard	<input type="text" value="2"/> <input type="text" value="20"/>	Reaction Time	<input type="text" value="20"/> <input type="text" value="22"/>
Standard	<input type="text" value="2"/> <input type="text" value="20"/>	Aspirated	<input type="text" value="20"/> <input type="text" value="180"/>	Reagent Vol	<input type="text" value="200"/> <input type="text" value=""/>
Decreased	<input type="text" value="2"/> <input type="text" value="20"/>	Diluent	<input type="text" value="180"/> <input type="text" value=""/>	R1	<input type="text" value="200"/> <input type="text" value=""/>
Increased	<input type="text" value=""/> <input type="text" value=""/>	Sample Blank	<input type="checkbox"/>	R2	<input type="text" value=""/> <input type="text" value=""/>
Linearity range (Standard)	<input type="text" value="0.31"/> <input type="text" value="20.2"/>	Linearity Limit	<input type="text" value=""/>	Substrate Depletion	<input type="text" value=""/>
Linearity Range (Decreased)	<input type="text" value=""/>	Mixed Blank Abs	<input type="text" value="-35000"/> <input type="text" value="35000"/>	On-board Stability	<input type="text" value=""/> Day(s)
Linearity Range (Increased)	<input type="text" value=""/>	Blank Response	<input type="text" value="-35000"/> <input type="text" value="35000"/>	Reagent Alarm Limit	<input type="text" value=""/>
R1 Blank Abs	<input type="text" value="-35000"/> <input type="text" value="35000"/>	Twin Chemistry	<input type="text" value=""/>	Enzyme Linear Extension	<input type="checkbox"/>
Blank Response	<input type="text" value="-35000"/> <input type="text" value="35000"/>	Prozone Check	<input type="checkbox"/>	Q1	<input type="text" value=""/>
Twin Chemistry	<input type="text" value=""/>	Q2	<input type="text" value=""/>	V1	<input type="text" value=""/>
Prozone Check	<input type="checkbox"/>	Q3	<input type="text" value=""/>	Q4	<input type="text" value=""/>
Q1	<input type="text" value=""/>	Q5	<input type="text" value=""/>	V2	<input type="text" value=""/>
Q2	<input type="text" value=""/>	Q6	<input type="text" value=""/>	PC1	<input type="text" value=""/>
Q3	<input type="text" value=""/>	PC2	<input type="text" value=""/>	Sample Pretreatment	<input type="text" value=""/>
Q4	<input type="text" value=""/>	Sample Pretreatment	<input type="text" value=""/>	Control Pretreatment	<input type="text" value=""/>
Q5	<input type="text" value=""/>	Control Pretreatment	<input type="text" value=""/>	Calibrator Pretreatment	<input type="text" value=""/>
Q6	<input type="text" value=""/>	Calibrator Pretreatment	<input type="text" value=""/>	Pretreat Sample Vol	<input type="text" value=""/> <input type="text" value=""/> $\mu$ L
PC1	<input type="text" value=""/>	Pretreat Sample Vol	<input type="text" value=""/> <input type="text" value=""/> $\mu$ L	Pretreat Sample Vol	<input type="text" value=""/> <input type="text" value=""/> $\mu$ L
PC2	<input type="text" value=""/>	Pretreat Sample Vol	<input type="text" value=""/> <input type="text" value=""/> $\mu$ L		
Sample Pretreatment	<input type="text" value=""/>				
Control Pretreatment	<input type="text" value=""/>				
Calibrator Pretreatment	<input type="text" value=""/>				
Pretreat Sample Vol	<input type="text" value=""/> <input type="text" value=""/> $\mu$ L				
Pretreat Sample Vol	<input type="text" value=""/> <input type="text" value=""/> $\mu$ L				

#### CALIBRATION SETTINGS

Math model   
 Factor  Replicates

#### AUTO CALIBRATION

Bottle Changed  
 Lot Changed  
 Cal Time

#### ACCEPTANCE LIMITS

Cal Time  Hour  
 Slope Diff  SD   
 Sensitivity  Repeatability   
 Deter Coeff

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