ACCENT-200 ACP

Nr kat. 7-249

(PL)

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Cat. No **7-249** (EN)

INTENDED USE

Diagnostic kit for determination of total, non prostatic and prostatic acid phosphatase activity. These reagents may be used for determination: **ACP**: in the automatic analysers: ACCENT-200, ACCENT-200 II GEN and ACCENT-220S; **ACP-NP**: in automatic analysers: ACCENT-200, ACCENT-200 II GEN.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Acid phosphatases (ACP) are a group of enzymes that can be found in liver, spleen, erythrocytes, platelets, bone marrow and prostate gland. The ratio between total acid phosphatase and prostatic phosphatase is healthy males is about 1:1. Activity of total acid phosphatase increases in case of Paget's disease, hyperparathyroidism, bone cancer, Gaucher'a disease, Niemann-Pick disease and myelocytic leukaemia. Prostatic acid phosphatase levels increase in case of prostatic cancer.

METHOD PRINCIPLE

Acid phosphatases catalyses the hydrolysis, in acid medium, of 1-naphtylphosphate (1-NP) into 1-naphtol and phosphate. 1-naphtol reacts with diazo-2-chloro-5-toluene (Fast Red TR salt), forming an azo dye compound with absorbance increase is proportional to the total acid phosphatase activity (ACP). Prostatic acid phosphatase (ACP-P) is inhibited by tartrate and is detected with an indirect method by subtraction between ACP and non-prostatic acid phosphatase (ACP-NP).

REAGENTS

Package	
REAGENT 1a	6 x 20 ml
REAGENT 1b	6 vials (powder)
TARTRATE REAGENT	1 vial (powder)
ACETIC ACID SOLUTION	1 x 3 ml
(acetic acid 0.1 mol/l)	
small spoon	1 pcs.

The reagents, stored at 2-8°C are stable up to expiry date printed on the package.

Working reagents stability

Stability of working reagents ACP and ACP-NP after reconstitution: 7 days at 2-8°C, if contamination avoided and bottle recapped immediately after use.

Concentrations in the working reagent ACP

citrate buffer (pH 5.4)	0.1 mol/l
1-NP	10 mmol/l
Fast Red TR salt	> 1.0 mmol/l

Concentrations in the working reagent ACP-NP

citrate buffer (pH 5.4)	0.1 mol/l
1-NP	10 mmol/l
Fast Red TR salt	> 1.0 mmol/l
sodium tartrate	0.18 mol/l

Warnings and notes

- Protect from direct sunlight and avoid contamination!

ACCENT-200 ACP

- It is recommended to carry out the prostatic acid phosphatase test only on samples with total acid phosphatase activity > 3.5 U/l as it is possible to obtain not reliable with a total acid phosphatase activity lower than the above-mentioned limit.
- Controls and calibrators used in determination of ACP should be used according to the guidelines contained in their instruction for use.
- The reagents must be limpid; do not use if turbid.
- A slight pink coloration of working reagents ACP and ACP-NP does not affect the reagents performance.
- ACETIC ACID SOLUTION meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

ACETIC ACID SOLUTION contains acetic acid.

Danger

H314 Causes severe skin burns and eye damage.
 H318 Causes eye damage.
 P280: Wear protective gloves/protective clothing/eye protection/face protection
 P302+P352: IF ON SKIN: Wash with plenty

of water

P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P321: Specific treatment is urgently needed (go to see a doctor with the Safety data sheet for this product).

P501: Dispose of the contents/containers in accordance with the current legislation on waste treatment.

SPECIMEN

Fresh serum.

Do not use haemolized, lipemic and icteric samples. Since the sample activity decreases of about 50% within an hour at room temperature, acidify the specimen immediately after collection by adding 20µl ACETIC ACID SOLUTION for each ml of fresh serum.

Stability of acidified sample: 5 days at 2-8°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

Deionised water is recommended as a reagent blank.

ACP ASSAY

Working reagent ACP: Dissolve the contents of one bottle of Reagent 1b with the contents of one bottle of Reagent 1a. Let stand for about 10 minutes and gently mix. The contents pour into bottle of Reagent 1a, put in analyser and use application for ACP.

ACP-NP ASSAY

(only for ACCENT-200, ACCENT-200 II GEN)

Working reagent ACP-NP: Dissolve the contents of one bottle of Reagent 1b with the contents of one bottle of Reagent 1a. Let stand for about 10 minutes and gently mix. In order to determine ACP-NP activity add one small spoon contents (corresponding to 35-40 mg) of Tartrate Reagent for each 5 ml working reagent. Gently mix until complete solution. Such solution put in analyser and use application for ACP-NP.

CALCULATION OF ACP-P (prostatic acid phosphatase)

ACP-P [U/l] = ACP – ACP-NP

REFERENCE VALUES ⁴

ACP

serum	U/l	µkat/l
male	2.5 – 11.7	0.042 – 0.195
female	0.3 – 9.2	0.005 – 0.154

ACP-NP

serum (tartrate inhibition)	U/l	µkat/l
male	0.2 – 3.5	0.003 – 0.058
female	0.0 – 0.8	0.000 – 0.013

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For ACP internal quality control it is recommended to use with each batch of samples the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173). For the calibration of automatic analysers: ACCENT-200, ACCENT-200 II GEN, the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

For the calibration of automatic analyser ACCENT-220S, the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers: ACCENT-200 and/or ACCENT-200 II GEN and/or ACCENT-220S. Results may vary if a different instrument is used.

Sensitivity

ACP: 0.35 U/l (0.006 µkat/l)

ACP-NP: 0.39 U/l (0.007 µkat/l)

Linearity

ACP: up to 75 U/l (1.25 µkat/l)

ACP-NP: up to 75 U/l (1.25 µkat/l)

For higher activity than 75 U/l dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 0.15 g/dl, triglycerides up to 400 mg/dl, bilirubin up to 20 mg/dl and do not interfere with the test.

Precision

ACP

Repeatability (run to run) n=20	Mean [U/l]	SD [U/l]	CV [%]
level 1	6,7	0,26	3,87
level 2	13,8	0,37	2,70
level 3	28,0	0,64	2,30
Reproducibility (day to day) n=20	Mean [U/l]	SD [U/l]	CV [%]
level 1	6,64	0,27	4,06
level 2	13,84	0,37	2,64
level 3	28,10	0,62	2,22

ACP-NP

Repeatability (run to run) n=20	Mean [U/l]	SD [U/l]	CV [%]
level 1	2,8	0,19	6,73
level 2	3,9	0,18	4,50
level 3	6,0	0,23	3,86

ACCENT-200 ACP

Reproducibility (day to day) n=20	Mean [U/l]	SD [U/l]	CV [%]
level 1	2,93	0,17	5,88
level 2	3,90	0,19	4,95
level 3	5,96	0,23	3,82

Method comparison

A comparison between CORMAY reagent (y) and another commercially available assay (x) using 60 samples gave following results:

ACP

y = 1.0162 x + 0.10151 U/l;

R = 0.997 (R – correlation coefficient)

ACP-NP

y = 1.1325 x - 0.09887 U/l;

R = 0.998 (R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

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Calibration Rule

Rule	One point Linear
Sensitivity	1
Replicates	3
Interval (day)	0
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

ACCENT-220S

Parameters

Test	ACP	R1	250
No	56	R2	
Full Name	Acid phosphatase	Sample Volume	25
Standard No	56	R1 Blank	
Reac. Type	Kinetic	Mixed Rtg. Blank	
Pri. Wave.	405 nm	Linearity Range	
Sec. Wave.	670 nm	Linearity Limit	0.2
Direction	Increase	Substrate Limit	
Reac. Time	8 18	Factor	743
Incuba. Time		<input type="checkbox"/> Prozone check	
Unit	U/l	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.1	PC <input type="text"/>	Abs <input type="text"/>

Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	2
Interval (day)	7
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

ACP-NP:

ACCENT-200, ACCENT-200 II GEN

Parameters

Test Name	ACPnp	R1	250
Test No	57	R2	
Full Name	Acid phosphatase NPP	Sample Volume	25
Reference No	57	R1 Blank	
Analy. Type	Kinetic	Mixed Reag. Blank	
Pri. Wave.	405 nm	Concentration	
Secn. Wave.	670 nm	Linearity Limit	0.2
Trend	Ascending	Substrate Limit	
Reac. Time	15 22	Factor	743
Incuba. Time		<input type="checkbox"/> Prozone check	
Unit	U/l	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.1	PC <input type="text"/>	Abs <input type="text"/>

Calibration Rule

Rule	One point Linear
Sensitivity	1
Replicates	3
Interval (day)	0
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

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