ACCENT-200 TG

Nr kat. 7-253 (PL)

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Cat. No **7-253** (EN)

INTENDED USE

Diagnostic kit for determination of triglycerides concentration intended to use in automatic analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, BS-120, ACCENT 400 and ACCENT Neo200.




The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Triglycerides are built of glycerol molecule esterified with three fatty acids molecules. Triglycerides are delivered with food or are synthesized endogenously in liver. Triglycerides stored in adipose tissue constitute a reserve of energy. Elevated triglycerides serum level is a risk factor of atherosclerosis. Triglycerides measurement is useful for hyperlipidemia diagnosis and treatment or for estimation of atherosclerosis progression.

METHOD PRINCIPLE

Colorimetric, enzymatic method with glycerophosphate oxidase.

The colour intensity is proportional to the triglycerides concentration.

REAGENTS

Package

1-Reagent	4 x 30 ml
2-Reagent	2 x 15 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package.

The reagents stored on board of the analyser at 2-10°C are stable for 10 weeks (ACCENT-200) or 12 weeks (ACCENT M320).

Concentrations in the test

buffer PIPES (pH 7.0)	40 mmol/l
4-aminoantipyrine (4-AA)	0.4 mmol/l
ATP	1.5 mmol/l
Mg ²⁺	1.6 mmol/l
ADPS	0.6 mmol/l
glycerol kinase (GK)	> 66.67 µkat/l
glycerol-3-phosphate oxidase (GPO)	> 60.00 µkat/l
peroxidase (POD)	> 20.00 µkat/l
lipoprotein lipase (LPL)	> 16.67 µkat/l
preservative	

ACCENT-200 TG

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum, EDTA or heparinized plasma (recommended: heparine lithium, sodium or ammonium salt) free from hemolysis.

Blood should be collected only if the patient has been fasting for minimum of 12 hours. Before blood collection patient should stay in rest position for about 30 minutes. Venous blood is recommended for triglycerides measurement.

Plasma triglycerides values have been reported to be 2% to 4% lower than serum triglycerides values.

Serum and plasma can be stored up to 3 days at 2-8°C or 3 months at -20°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

Deionised water is recommended as a reagent blank.

For analyzers: ACCENT-200 (II GEN), ACCENT-220S and BS-120 it is recommended to determine the reagent blank during each calibration. Deionized water should be used as reagent blank. When performing calibration, the task type **Calib+Rgt.Blk** should be selected.

Actions required:

When performing assays in analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S and BS-120 there is a probability of **cross-contamination** affecting the tests results: TG - LIPASE, TG – MICROALBUMIN, HDL DIRECT II GEN - TG., To avoid this effect follow the recommendations contained in the instruction 51_03_24_001_ACCENT-200_CARRYOVER.

REFERENCE VALUES⁷:

serum, plasma	< 150 mg/dl < 1.7 mmol/l
---------------	-----------------------------

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) with each batch of samples.

For the calibration of automatic analysers: ACCENT-200, ACCENT-200 II GEN, the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

For the calibration of automatic analysers: ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, BS-120, the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 10 weeks (ACCENT-200) or 12 weeks (ACCENT M320), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers ACCENT-200 and ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

Sensitivity

5.7 mg/dl (0.06 mmol/l) - ACCENT-200

LOQ

5.8 mg/dl (0.07 mmol/l) – ACCENT MC240

Linearity

up to 2000 mg/dl (22.6 mmol/l) - ACCENT-200

up to 2350 mg/dl (26.56 mmol/l) – ACCENT MC240

For higher triglycerides concentrations dilute the sample with 0.9% NaCl in the ratio of 1 to 4 and repeat the assay. Multiply the result by 5.

Specificity / Interferences

Haemoglobin up to 2.5 g/dl, bilirubin up to 20 mg/dl and ascorbate up to 62 mg/l do not interfere with the test.

Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 n=20	level 1	62.45	0.80	1.29
	level 2	171.40	3.90	2.28
ACCENT MC240 n=20	level 1	87.45	0.64	0.73
	level 2	186.08	1.44	0.78
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 n=20	level 1	117.54	2.52	2.15
	level 2	186.06	3.76	2.02
ACCENT MC240 n=80	level 1	185.2	2.84	1.5
	level 2	91.1	1.83	2.0

Method comparison

A comparison between triglycerides values determined at **ACCENT-200** (y) and at **COBAS INTEGRA 400** (x) using 100 samples gave following results:

$$y = 0.9479x + 5.0253 \text{ mg/dl};$$

$$R = 0.994 \quad (R - \text{correlation coefficient})$$

A comparison between triglycerides values determined at **ACCENT MC240** (y) and at **ADVIA 1800** (x) using 58 serum samples gave following results:

$$y = 0.9914x + 1.3579 \text{ mg/dl};$$

$$R = 0.998 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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ACCENT-200 TG

• ACCENT MC240

Chem	TG	No.	003	Sample Type	SERUM
Chemistry	TRIGLYCERIDES	Print name	TG	Reaction Direction	positive
Reaction Type	Endpoint	Sec Wave	660 nm	Unit	mg/dL
Pri Wave	546 nm	Incubation Time	21	Reaction Time	28 30
Unit	mg/dL	Decimal	0.1	Standard	2.4 μL
Blank Time	-3 -1	Reagent Vol	R1 160 μL R2 40 μL	Aspirated	μL
Sample Vol	2.4 μL	Diluent	μL	Decreased	2.4 μL 20 μL 180 μL
Increased	μL	Auto Rerun	<input type="checkbox"/>	Blank Response	-35000 35000
Linearity range (Standard)	5.8 2350	Linearity Limit		On-board Stability	Day(s)
Linearity Range (Decreased)		Substrate Depletion		Reagent Alarm Limit	
Linearity Range (Increased)		Mixed Blank Abs	-35000 35000	Twin Chemistry	<input type="checkbox"/>
R1 Blank Abs	-35000 35000	On-board Stability	Day(s)	Prozone Check	<input type="checkbox"/>
Blank Response	-35000 35000	Reagent Alarm Limit		Q1 Q2 V1 Q3 Q4 V2	
Twin Chemistry	<input type="checkbox"/>	Enzyme Linear Extension	<input type="checkbox"/>	Q5 Q6 V3 PC1 PC2	
Sample Pretreatment	<input type="checkbox"/>	Calibrator Pretreatment	<input type="checkbox"/>	Control Pretreatment	<input type="checkbox"/>
Factor		Replicates	2	Pretreat Sample Vol	μL
Cal Time	Hour	Math model	Multi-point linear	Pretreat Sample Vol	μL
Slope Diff		SD		Factor	
Sensitivity		Repeatability	35000	Replicates	2
Deter Coeff		Cal Time		Math model	Multi-point linear

• ACCENT M320

Chem	TG	No.	003	Sample Type	SERUM
Chemistry	TRIGLYCERIDES	Print name	TG	Reaction Direction	positive
Reaction Type	Endpoint	Sec Wave	800 nm	Unit	mg/dL
Pri Wave	505 nm	Incubation Time	18	Reaction Time	38 40
Unit	mg/dL	Blank Time	-3 -1	Standard	3 μL
Sample Vol	3 μL	Aspirated	μL	Diluent	μL
Decreased	3 μL 20 μL 180 μL	Reagent Vol	R1 160 μL R2 40 μL	Auto Rerun	<input type="checkbox"/>
Increased	μL	Sample Blank	<input type="checkbox"/>	Linearity range (Standard)	5.2 2390
Linearity Range (Decreased)		Linearity Limit		Linearity Range (Increased)	
Linearity Range (Increased)		Substrate Depletion		Mixed Blank Abs	-35000 35000
R1 Blank Abs	-35000 35000	On-board Stability	Day(s)	On-board Stability	Day(s)
Blank Response	-35000 35000	Reagent Alarm Limit		Reagent Alarm Limit	
Twin Chemistry	<input type="checkbox"/>	Enzyme Linear Extension	<input type="checkbox"/>	Q1 Q2 V1 Q3 Q4 V2	
Prozone Check	<input type="checkbox"/>	Q5 Q6 V3 PC1 PC2		Control Pretreatment	<input type="checkbox"/>
Sample Pretreatment	<input type="checkbox"/>	Calibrator Pretreatment	<input type="checkbox"/>	Pretreat Sample Vol	μL
Factor		Replicates	2	Pretreat Sample Vol	μL
Cal Time	Hour	Math model	Multi-point linear	Factor	
Slope Diff		SD		Replicates	2
Sensitivity		Repeatability	35000	Math model	Multi-point linear
Deter Coeff		Cal Time		Factor	

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