## ACCENT-200 FERRUM

Nr kat. 7-258 (PL)

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(EN)

### INTENDED USE

Diagnostic kit for determination of iron concentration used in automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, BS-120, ACCENT 400 and ACCENT Neo200.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Iron is the most abundant trace element in the organism. Most of the iron in humans is located within heme molecule which is incorporated into hemoglobin, myoglobin, catalase, peroxidase and cytochromes. Iron is stored bound to ferritin or hemosiderin and is transported by transferrin. Measurement of iron level is valuable especially in diagnosis of different types of anemia.

### METHOD PRINCIPLE

Colorimetric method with ferrozine, without deproteinization. Iron ions ( $Fe^{3+}$ ), bounded in blood to transferrin are released in acid solution and in the presence of detergents and reduced to  $Fe^{2+}$  by ascorbate.  $Fe^{2+}$  forms with 3-(2-pyridyl)-5,6-bis(2-[4-phenyl sulfonic acid])-1,2,4-triazine sodium salt (ferrozine) coloured complex,  $Cu^{2+}$  ions are bound by thiourea. The colour intensity is directly related to the iron concentration.

### REAGENTS

#### Package

1-REAGENT 3 x 31 ml  
2-REAGENT 2 x 10 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 12 weeks.

### Concentrations in the test

#### 1-REAGENT

citric acid (pH 1.9) ≤ 240 mmol/l  
thiourea ≤ 108 mmol/l  
detergent ≤ 7 %

#### 2-REAGENT

sodium ascorbate ≤ 150 mmol/l  
3-(2-pyridyl)-5,6-bis(2-[5-furyl sulfonic acid])-1,2,4-triazine sodium salt (ferrozine) ≤ 6 mmol/l  
preservative  
stabilizer

### WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- Contaminated glassware is the greatest source of error. The use of disposable plastic ware is recommended. Glassware should be soaked for a few hours in 2M HCl solution and then thoroughly rinsed with distilled water.
- 1-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

### Ingredients:

1-REAGENT contains citric acid.

### Danger

H314 Causes severe skin burns and eye damage.

P280 Wear protective gloves, protective clothing, eye protection or face protection.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER or doctor.

### SPECIMEN <sup>13</sup>

Serum free from hemolysis, collected in plastic tubes.

Serum should be separated from red blood cells as soon as possible after blood collection.

Serum can be stored up to 7 days at 15-25°C or up to 7 weeks at 2-8°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-REAGENT and 2-REAGENT are ready to use. Deionised water is recommended as a reagent blank.

#### Actions required:

When performing assays in analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S and BS-120 there is a probability of **cross-contamination** affecting the tests results: FERRUM - FERRITIN, TOTAL PROTEIN II GEN - FERRUM. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_001\_ACCENT-200\_CARRYOVER.

### REFERENCE VALUES <sup>8, 10</sup>

serum	µg/dl	µmol/l
newborns	100 – 250	17.9 – 44.8
infants	40 – 100	7.2 – 17.9
children	50 – 120	9.0 – 21.5
adult females	50 – 170	9.0 – 30.4
adult males	65 – 175	11.6 – 31.3

Samples should be taken in the morning from patients in a fasting state, since iron values decrease by 30% during the course of the day.

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) with each batch of samples.

For the calibration of automatic analysers: ACCENT-200, ACCENT-200 II GEN, BS-120, the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or CORMAY MULTICALIBRATOR LEVEL 2 (Cat. No 5-175; 5-177) is recommended depending on the calibrator lot number. Deionised water should be used as a calibrator 0.

For the calibration of automatic analyser ACCENT-220S, the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) or CORMAY MULTICALIBRATOR LEVEL 2 (Cat. No 5-175; 5-177) is recommended depending on the calibrator lot number. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 12 weeks with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers: ACCENT-200 and ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

#### Sensitivity

4.1 µg/dl (0.73 µmol/l) - ACCENT-200  
2.0 µg/dl (0.36 µmol/l) – ACCENT MC240

#### Linearity

do 1000 µg/dl (179 µmol/l) - ACCENT-200  
do 1300 µg/dl (232.7 µmol/l) – ACCENT MC240

For higher concentrations dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

#### Specificity / Interferences

Ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl, triglycerides up to 1000 mg/dl and copper up to 500 µg/dl do not interfere with the test. Haemoglobin interferes even in small amount with the determination.

#### Precision

Repeatability (run to run)		Mean [µg/dl]	SD [µg/dl]	CV [%]
ACCENT-200 n=10	level 1	66.1	0.72	1.10
	level 2	246.5	1.04	0.42
ACCENT MC240 n=20	level 1	60.4	0.37	0.6
	level 2	222.7	1.16	0.5
Reproducibility (day to day)		Mean [µg/dl]	SD [µg/dl]	CV [%]
ACCENT-200 n=10	level 1	65.8	1.23	1.87
	level 2	247.1	2.51	1.01
ACCENT MC240 n=80	level 1	62.7	1.77	2.8
	level 2	231.2	3.79	1.6

#### Method comparison

A comparison between iron values determined at **ACCENT 200** (y) and at **BS-400** (x) using 85 serum samples gave following results  
 $y = 0.956x + 4.3082$  µg/dl;  
 $R = 0.997$  (R – correlation coefficient)

A comparison between iron values determined at **ACCENT MC240** (y) and at **Advia 1800** (x) using 58 serum samples gave following results:  
 $y = 1.0372x - 0.9453$  µg/dl;  
 $R = 1.000$  (R – correlation coefficient)

### WASTE MANAGEMENT

Please refer to local legal requirements.

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**Calibration Rule**

Rule	Two-point Linear
Sensitivity	1
Replicates	3
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0   50000
Error Limit	0
Coefficient	0

• ACCENT-200 II GEN

**Parameters**

Test Name	FERRUM	R1	250
Test No	14	R2	50
Full Name	FERRUM	Sample Volume	20
Reference No	14	R1 Blank	
Analy. Type	Endpoint	Mixed Reag. Blank	
Pri. Wave.	546 nm	Concentration	3.7   1000
Secon. Wave.	670 nm	Linearity Limit	
Trend	Increase	Substrate Limit	
Reac. Time	-1   20	Factor	
Incuba. Time	19	<input type="checkbox"/> Prozone check	
Unit	µg/dl	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.1	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

**Calibration Rule**

Rule	Two-point Linear
Sensitivity	1
Replicates	3
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0   50000
Error Limit	0
Coefficient	0

• ACCENT-220S

**Parameters**

Test	FERRUM	R1	250
No	14	R2	50
Full Name	FERRUM	Sample Volume	20
Standard No	14	R1 Blank	
Reac. Type	Endpoint	Mixed Rtg. Blank	
Pri. Wave.	546 nm	Linearity Range	5.5   1000
Sec. Wave.	670 nm	Linearity Limit	
Direction	Increase	Substrate Limit	
Reac. Time	-1   20	Factor	
Incuba. Time	22	<input type="checkbox"/> Prozone check	
Unit	µg/dl	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.1	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

**Calibration Rule**

Rule	Multi-point Linear
Sensitivity	1
Replicates	3
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0   50000
Error Limit	0
Coefficient	0

• BS-120

**Parameters**

Test	FERRUM	R1	250
No	14	R2	50
Full Name	FERRUM	Sample Volume	20
Standard No	14	R1 Blank	
Reac. Type	Endpoint	Mixed Rtg. Blank	
Pri. Wave.	546 nm	Linearity Range	7.5   1180
Sec. Wave.	670 nm	Linearity Limit	
Direction	Increase	Substrate Limit	
Reac. Time	-1   16	Factor	
Incuba. Time	16	<input type="checkbox"/> Prozone check	
Unit	µg/dl	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.1	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

**Calibration Rule**

Rule	Two-point Linear
Sensitivity	1
Replicates	3
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0   50000
Error Limit	0
Coefficient	0

• ACCENT S120

Chem	FERRUM	No.	014	Sample Type	SERUM						
Chemistry	FERRUM			Print name	FERRUM						
Reaction Type	Endpoint			Reaction Direction	positive						
Pri Wave	578 nm			Sec Wave	670 nm						
Unit	ug/dL			Decimal	0.1						
Blank Time	-3	-1		Incubation Time	20						
Standard	16	µL	Aspirated	µL	Diluent	µL	Reagent Vol	R1 200	µL		
Decreased	16	µL	20	µL	180	µL	R2 40	µL			
Increased		µL		µL		µL					
			Sample Blank	<input checked="" type="checkbox"/>	Auto Rerun						
Linearity range (Standard)	4	1300		Linearity Limit							
Linearity Range (Decreased)				Substrate Depletion							
Linearity Range (Increased)				Mixed Blank Abs	-40000	40000					
R1 Blank Abs	-40000	40000		On-board Stability	84	Day(s)					
Blank Response	-40000	40000		Reagent Alarm Limit							
Twin Chemistry				<input type="checkbox"/> Enzyme Linear Extension							
<input type="checkbox"/> Prozone Check											
Q1		Q2		V1		Q3		Q4		V2	
Q5		Q6		V3		PC1		PC2			
<input type="checkbox"/> Sample Pretreatment		<input type="checkbox"/> Control Pretreatment		<input type="checkbox"/> Calibrator Pretreatment							
		Pretreat Sample Vol		µL		Pretreat Sample Vol		µL			
<b>CALIBRATION SETTINGS</b>											
Math model	Two-point linear			<b>AUTO CALIBRATION</b>							
Factor		Replicates	2	<input type="checkbox"/> Bottle Changed							
				<input type="checkbox"/> Lot Changed							
				<input type="checkbox"/> Cal Time							
<b>ACCEPTANCE LIMITS</b>											
Cal Time	2016	Hour									
Slope Diff				SD							
Sensitivity		Repeatability	40000								
Deter Coeff											

## ACCENT-200 FERRUM

### • ACCENT MC240

Chem <b>FERRUM</b>	No. <b>014</b>	Sample Type <b>SERUM</b>
Chemistry <b>FERRUM</b>	Print name <b>FERRUM</b>	
Reaction Type <b>Endpoint</b>	Reaction Direction <b>positive</b>	
Pri Wave <b>570 nm</b>	Sec Wave <b>700 nm</b>	
Unit <b>ug/dL</b>	Decimal <b>0.1</b>	
Blank Time <b>-3</b>	<b>-1</b>	Incubation Time <b>21</b>
		Reaction Time <b>26</b> <b>28</b>
Standard <b>16</b> $\mu$ L	Aspirated <b>16</b> $\mu$ L	Diluent <b>180</b> $\mu$ L
Decreased <b>16</b> $\mu$ L	<b>20</b> $\mu$ L	<b>180</b> $\mu$ L
Increased <b>16</b> $\mu$ L	<b>20</b> $\mu$ L	<b>180</b> $\mu$ L
	Sample Blank <input type="checkbox"/>	Auto Rerun <input checked="" type="checkbox"/>
Linearity range (Standard) <b>2</b>	<b>1300</b>	Linearity Limit <b>1300</b>
Linearity Range (Decreased)		Substrate Depletion
Linearity Range (Increased)		Mixed Blank Abs <b>-35000</b> <b>35000</b>
R1 Blank Abs <b>-35000</b>	<b>35000</b>	On-board Stability <b>84</b> Day(s)
Blank Response <b>-35000</b>	<b>35000</b>	Reagent Alarm Limit
Twin Chemistry		<input type="checkbox"/> Enzyme Linear Extension
<input type="checkbox"/> Prozone Check		
Q1 <b>16</b>	Q2 <b>16</b>	V1 <b>180</b>
Q3 <b>16</b>	Q4 <b>16</b>	V2 <b>180</b>
Q5 <b>16</b>	Q6 <b>16</b>	V3 <b>180</b>
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment	<input type="checkbox"/> Calibrator Pretreatment
Pretreat Sample Vol <b>16</b> $\mu$ L	Pretreat Sample Vol <b>16</b> $\mu$ L	
<b>CALIBRATION SETTINGS</b>		<b>AUTO CALIBRATION</b>
Math model <b>Two-point linear</b>		<input type="checkbox"/> Bottle Changed
Factor <b>1</b>	Replicates <b>2</b>	<input type="checkbox"/> Lot Changed
		<input type="checkbox"/> Cal Time
<b>ACCEPTANCE LIMITS</b>		
Cal Time <b>2016</b> Hour		
Slope Diff <b>1</b>	SD <b>1</b>	
Sensitivity <b>1</b>	Repeatability <b>35000</b>	
Deter Coeff <b>1</b>		

### • ACCENT M320

Chem <b>FERRUM</b>	No. <b>014</b>	Sample Type <b>SERUM</b>
Chemistry <b>FERRUM</b>	Print name <b>FERRUM</b>	
Reaction Type <b>Endpoint</b>	Reaction Direction <b>positive</b>	
Pri Wave <b>570 nm</b>	Sec Wave <b>700 nm</b>	
Unit <b>ug/dL</b>	Decimal <b>0.1</b>	
Blank Time <b>-3</b>	<b>-1</b>	Incubation Time <b>29</b>
		Reaction Time <b>28</b> <b>29</b>
Standard <b>16</b> $\mu$ L	Aspirated <b>16</b> $\mu$ L	Diluent <b>180</b> $\mu$ L
Decreased <b>16</b> $\mu$ L	<b>20</b> $\mu$ L	<b>180</b> $\mu$ L
Increased <b>16</b> $\mu$ L	<b>20</b> $\mu$ L	<b>180</b> $\mu$ L
	Sample Blank <input type="checkbox"/>	Auto Rerun <input checked="" type="checkbox"/>
Linearity range (Standard) <b>3.9</b>	<b>1300</b>	Linearity Limit <b>1300</b>
Linearity Range (Decreased)		Substrate Depletion
Linearity Range (Increased)		Mixed Blank Abs <b>-35000</b> <b>35000</b>
R1 Blank Abs <b>-35000</b>	<b>35000</b>	On-board Stability <b>84</b> Day(s)
Blank Response <b>-35000</b>	<b>35000</b>	Reagent Alarm Limit
Twin Chemistry		<input type="checkbox"/> Enzyme Linear Extension
<input type="checkbox"/> Prozone Check		
Q1 <b>16</b>	Q2 <b>16</b>	V1 <b>180</b>
Q3 <b>16</b>	Q4 <b>16</b>	V2 <b>180</b>
Q5 <b>16</b>	Q6 <b>16</b>	V3 <b>180</b>
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment	<input type="checkbox"/> Calibrator Pretreatment
Pretreat Sample Vol <b>16</b> $\mu$ L	Pretreat Sample Vol <b>16</b> $\mu$ L	Pretreat Sample Vol <b>16</b> $\mu$ L
<b>CALIBRATION SETTINGS</b>		<b>AUTO CALIBRATION</b>
Math model <b>Two-point linear</b>		<input type="checkbox"/> Bottle Changed
Factor <b>1</b>	Replicates <b>2</b>	<input type="checkbox"/> Lot Changed
		<input type="checkbox"/> Cal Time
<b>ACCEPTANCE LIMITS</b>		
Cal Time <b>2016</b> Hour		
Slope Diff <b>1</b>	SD <b>1</b>	
Sensitivity <b>1</b>	Repeatability <b>35000</b>	
Deter Coeff <b>1</b>		

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