## ACCENT-200 UIBC

Nr kat. 7-259 (PL)

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### LITERATURA

- Fairbanks VF, Klee GG. Biochemical aspects of hematology. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 1642-1710.
- Wick M, Pinger W, Lehmann P. Clinical aspects and laboratory. Iron metabolism, anemias. 5th ed. Wien, New York: Springer; 2003.
- Guder Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001. p. 46.
- Perrotta, G., Iron and Iron-Binding Capacity, In: Pesce, A.J., Kaplan, L.A. eds., Methods in Clinical Chemistry, C.V. Mosby, St. Louis, 1258 - 1261, 1987
- Tietz NW (ed). Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders; 1701-1703; 1821t (1999).
- Burtis CA, Brun DE. Tietz Fundamentals of Clinical Chemistry and Molecular Diagnostics. 7th ed. St. Louis, p. 306. 2014
- Schreiber WE. „Iron and Porphyrin Metabolism” in Pesce, A.J., Kaplan, L.A. editors. Clinical Chemistry: Theory, Analysis, Correlation. Mosby Inc., an affiliate of Elsevier Inc., St. Louis, 755-770, 2010.

Data wydania: 05.2022.

## ACCENT-200 UIBC

Cat. No **7-259** (EN)

### INTENDED USE

Diagnostic kit for determination of unsaturated iron binding capacity intended to use in automatic analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, ACCENT 400, ACCENT Neo200 and BS-120


The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

The total iron content of the body is about 3 to 3.5 g. Of this amount about 2.5 g contained in erythrocytes or their precursors in the bone marrow. Plasma contains only about 2.5 mg of iron. Iron is transported as Fe (III) bound to the plasma protein apotransferrin. The apotransferrin-Fe (III) complex is called transferrin. Normally only about one third of the iron binding sites of transferrin are occupied by Fe (III). The additional amount of iron that can be bound is the unsaturated (or latent) iron-binding capacity (UIBC). The sum of the serum iron and UIBC represents total iron binding capacity (TIBC). TIBC is a measurement for the maximum iron concentration that transferrin can bind. Serum UIBC levels vary in disorders of iron metabolism where iron capacities are often increased in iron deficiency and decreased in chronic inflammatory disorders, malignancies or in course of hemochromatosis.

### METHOD PRINCIPLE

Direct, colorimetric method with ferrozine:

In an alkaline environment known ferrous ion concentration incubated with serum, binds specifically with transferrin at unsaturated iron binding sites. Remaining unbound ferrous ions are measured with the chromogen reaction. The difference between the amount of excess iron and the total amount added to the serum is equivalent to the quantity bound to transferrin. This is the UIBC (unsaturated iron binding capacity) of the sample.

### REAGENTS

**Package**  
1-Reagent 1 x 33 ml  
2-Reagent 1 x 9 ml

The reagents, stored at 2-8°C are stable until expiry date printed on the package. Stability on board of the analyser at 2-10°C: 5 weeks (ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320), or 7 weeks (BS-120).

### Concentrations in the test

**1-Reagent**  
buffer (pH 8.4) 0.25 mol/l  
ammonium iron (II) sulfate 20 µmol/l  
thiourea 90 mmol/l  
detergent 0.1 %  
sodium azide <0.1 %

**2-Reagent**  
sodium ascorbate 150 mmol/l  
sodium chloride 75 mmol/l  
3-(2-pyridyl)-5,6-bis(2-[5-furyl sulfonic acid])-1,2,4-triazine sodium salt (ferrozine) ≥ 10 mmol/l  
preservatives 0.3 %

### Warnings and notes

- Do not freeze the reagents.
- Protect from direct sunlight and avoid contamination!
- Contaminated glassware is the greatest source of error. The use of disposable plastic ware is recommended. Glassware should be soaked for a few hours in 2M HCl solution and then thoroughly rinsed with distilled water.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- The negative UIBC value will occur when the patient's serum iron level exceeds the binding capacity of the transferrin.
- For the diagnostic purpose UIBC determination should be performed at the same time with iron determination. Obtained result have to be interpreted in relation to the result of iron concentration and the percentage saturation of transferrin with iron ions<sup>7</sup>.
- 1-Reagent contains thiourea. May produce an allergic reaction (EUH208).
- 2-Reagent contains 1-[1,3-Bis (hydroxymethyl)-2,5-dioximidazolidin-4-yl]-1,3-bis (hydroxymethyl) urea. May produce an allergic reaction (EUH208).

### SPECIMEN

Serum, heparin plasma.  
Separate serum/plasma at the latest 2 h after blood collection to avoid hemolysis. Samples should be taken in the morning from patients, since iron values decrease during the course of the day.  
Discard contaminated specimens.  
Anticoagulants such as EDTA, oxalate and citrate must not be used, as they bind iron ions and prevent reaction with chromogen<sup>4</sup>.  
Serum can be stored up to 3 days at 20-25°C, 7 days at 4-8°C or up to one month at -20°C. Plasma can be stored up to 7 days at 4-8°C or up to month at -20°C<sup>4</sup>.  
Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.  
Deionised water is recommended as a reagent blank.

### Actions required:

When performing assays in analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S and BS-120 there is a probability of **cross-contamination** affecting the tests results. To avoid this effect, the tests for determination of unsaturated iron binding capacity, using the ACCENT-200 UIBC set, should be performed **in the separate order**, (follow the recommendations contained in the instruction 51\_03\_24\_001\_ACCENT-200\_CARRYOVER).

### REFERENCE VALUES<sup>5,6</sup>

The reference values were calculated from the serum iron (SI) and TIBC ranges reported in literature, in accordance with mathematic formula:  
UIBC=TIBC-SI

Reference values for UIBC are listed in table below:

serum / plasma	µg/dl	µmol/l
Females	80 – 375	14 – 67
Males	75 – 360	13 – 64

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) with each batch of samples.  
For the calibration of automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, BS-120, the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) is recommended. Deionised water should be used as a calibrator 0.  
The calibration curve should be prepared every week (ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320) or every 5 weeks (BS-120), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers ACCENT-200 and ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

#### Sensitivity

20 µg/dl (3.6 µmol/l) - ACCENT-200  
26.5 µg/dl (4.7 µmol/l) – ACCENT MC240

#### Linearity

up to 450 µg/dl (80.6 µmol/l) - ACCENT-200  
up to 525 µg/dl (94 µmol/l) – ACCENT MC240

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

#### Specificity / Interferences

Haemoglobin interfere even in small amounts, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl, triglycerides up to 1000 mg/dl, copper up to 3.5 mg/dl and zinc up to 15 mg/dl do not interfere with the test.

#### Precision

Repeatability (run to run)		Mean [µg/dl]	SD [µg/dl]	CV [%]
ACCENT-200 n=10	level 1	86.92	2.14	2.46
	level 2	153.11	3.38	2.21
ACCENT MC240 n=20	level 1	80.33	4.80	5.98
	level 2	167.53	4.48	2.67
Reproducibility (day to day)		Mean [µg/dl]	SD [µg/dl]	CV [%]
ACCENT-200 n=10	level 1	85.42	3.45	4.04
	level 2	148.73	2.25	1.51
ACCENT MC240 n=80	level 1	79.3	3.93	5.0
	level 2	157.2	4.27	2.7

#### Method comparison

A comparison between UIBC values determined at ACCENT-200 (y) and COBAS INTEGRA 400 (x) using 63 samples gave following results:  
y = 0.957 x + 10.655 µg/dl;  
R = 0.992 (R – correlation coefficient)

A comparison between UIBC values determined at ACCENT MC240 (y) and BS-400 (x) using 57 serum samples gave following results:  
y = 0.9992 x + 5.0163 µg/dl;  
R = 0.994 (R – correlation coefficient)

### WASTE MANAGEMENT

Please refer to local legal requirements.

### LITERATURE

- Fairbanks VF, Klee GG. Biochemical aspects of hematology. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 1642-1710.
- Wick M, Pingerra W, Lehmann P. Clinical aspects and laboratory. Iron metabolism, anemias. 5th ed. Wien, New York: Springer; 2003.
- Guder Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001. p. 46.
- Perrotta, G., Iron and Iron-Binding Capacity, In: Pesce, A.J., Kaplan, L.A. eds., Methods in Clinical Chemistry, C.V. Mosby, St. Louis, 1258 - 1261, 1987
- Tietz NW (ed). Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders; 1701-1703; 1821t (1999).
- Burtis CA, Bruns DE. Tietz Fundamentals of Clinical Chemistry and Molecular Diagnostics. 7th ed. St. Louis, p. 306. 2014
- Schreiber WE. „Iron and Porphyrin Metabolism” in Pesce, A.J., Kaplan, L.A. editors. Clinical Chemistry: Theory, Analysis, Correlation. Mosby Inc., an affiliate of Elsevier Inc., St. Louis, 755-770, 2010.

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- Fairbanks VF, Klee GG. Biochemical aspects of hematology. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 1642-1710.
- Wick M, Pingerra W, Lehmann P. Clinical aspects and laboratory. Iron metabolism, anemias. 5th ed. Wien, New York: Springer; 2003.
- Guder Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001. p. 46.
- Perrotta, G., Iron and Iron-Binding Capacity, In: Pesce, A.J., Kaplan, L.A. eds., Methods in Clinical Chemistry, C.V. Mosby, St. Louis, 1258 - 1261,1987
- Tietz NW (ed). Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders; 1701-1703; 1821t (1999).
- Burtis CA, Bruns DE. Tietz Fundamentals of Clinical Chemistry and Molecular Diagnostics. 7th ed. St. Louis, p. 306. 2014
- Schreiber WE. „Iron and Porphyrin Metabolism” in Pesce, A.J., Kaplan, L.A. editors. Clinical Chemistry: Theory, Analysis, Correlation. Mosby Inc., an affiliate of Elsevier Inc., St. Louis, 755-770, 2010.

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#### Calibration Rule

Rule	Two-point Linear
Sensitivity	0
Replicates	3
Interval (day)	7
Difference Limit	0
SD	0
Blank Response	0   50000
Error Limit	0
Coefficient	0

#### • ACCENT-200 II GEN

<b>Parameters</b>			
Test Name	UIBC	R1	220
Test No	60	R2	55
Full Name	UIBC	Sample Volume	24
Reference No	60	R1 Blank	
Analy. Type	Endpoint	Mixed Reag. Blank	
Pri. Wave.	578 nm	Concentration	37   480
Secon. Wave.		Linearity Limit	
Trend	Increase	Substrate Limit	
Reac. Time	-2   25	Factor	
Incuba. Time	25	<input type="checkbox"/> Prozone check	
Unit	µg/dl	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	Integer	PC <input type="checkbox"/>	Abs <input type="checkbox"/>

#### Calibration Rule

Rule	Two-point Linear
Sensitivity	0
Replicates	3
Interval (day)	7
Difference Limit	0
SD	0
Blank Response	0   50000
Error Limit	0
Coefficient	0

#### • ACCENT-220S

<b>Parameters</b>			
Test	UIBC	R1	220
No	60	R2	55
Full Name	UIBC	Sample Volume	27
Standard No	60	R1 Blank	
Reac. Type	Endpoint	Mixed Rtg. Blank	
Pri. Wave.	578 nm	Linearity Range	31.5   460
Sec. Wave.		Linearity Limit	
Direction	Increase	Substrate Limit	
Reac. Time	-2   25	<input type="checkbox"/> Prozone check	
Incuba. Time	25	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Unit	µg/dl	PC <input type="checkbox"/>	Abs <input type="checkbox"/>
Precision	Integer		

#### Calibration Rule

Rule	Two-point Linear
Sensitivity	1
Replicates	3
Interval (day)	7
Difference Limit	0
SD	0
Blank Response	0   50000
Error Limit	0
Coefficient	0

#### • BS-120

<b>Parameters</b>			
Test	UIBC	R1	200
No	60	R2	50
Full Name	UIBC	Sample Volume	22
Standard No	60	R1 Blank	
Reac. Type	Endpoint	Mixed Rtg. Blank	
Pri. Wave.	578 nm	Linearity Range	36   475
Sec. Wave.		Linearity Limit	
Direction	Increase	Substrate Limit	
Reac. Time	-1   18	Factor	
Incuba. Time	16	<input type="checkbox"/> Prozone check	
Unit	µg/dl	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	Integer	PC <input type="checkbox"/>	Abs <input type="checkbox"/>

#### Calibration Rule

Rule	Two-point Linear
Sensitivity	1
Replicates	3
Interval (day)	21
Difference Limit	0
SD	0
Blank Response	0   50000
Error Limit	0
Coefficient	0

#### • ACCENT S120

Chem	UIBC	No.	060	Sample Type	SERUM
Chemistry	UIBC			Print name	UIBC
Reaction Type	Endpoint			Reaction Direction	Positive
Pri Wave	578nm			Sec Wave	670nm
Unit	µg/dL			Decimal	0,1
Blank Time	-3	-1		Incubation Time	17
Standard	Sample Vol 18 µL	Aspirated µL	Diluent µL	Reaction Time	30   32
Decreased	18 µL	20 µL	180 µL	Reagent Vol R1	160 µL
Increased	µL	µL	µL	R2	40 µL
	<input type="checkbox"/> Sample Blank	<input checked="" type="checkbox"/> Auto Rerun			
Linearity range (Standard)	32	500		Linearity Limit	
Linearity Range (Decreased)				Substrate Depletion	
Linearity Range (Increased)				Mixed Blank Abs	-40000   40000
R1 Blank Abs	-40000	40000		On-board Stability	Day(s)
Blank Response	-40000	40000		Reagent Alarm Limit	
Twin Chemistry				<input type="checkbox"/> Enzyme Linear Extension	
<input type="checkbox"/> Prozone Check					
Q1		Q2		V1	
Q5		Q6		V3	
<input type="checkbox"/> Sample Pretreatment		<input type="checkbox"/> Control Pretreatment		PC1	
				PC2	
		Pretreat Sample Vol	µL	Calibrator Pretreatment	
				Pretreat Sample Vol	µL
<b>CALIBRATION SETTINGS</b>	Math model	Two-point Linear		<b>AUTO CALIBRATION</b>	
Factor		Replicates	2	<input type="checkbox"/> Bottle Changed	
				<input type="checkbox"/> Lot Changed	
				<input type="checkbox"/> Cal Time	
<b>ACCEPTANCE LIMITS</b>	Cal Time		Hour		
	Slope Diff		SD		
	Sensitivity		Repeatability	40000	
	Deter Coeff				