ACCENT-200 AMYLASE EPS

Nr kat. 7-276 (PL)

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LITERATURA

- Pesce, A.J., Kaplan, L.A.: "Methods in Clinical Chemistry", Mosby Ed. (1987).
- Burtis C.A., Ashwood E.R.: "Tietz Textbook of Clinical Chemistry", W.B. Saunders Company Ed. (3rd edition, 1999).
- Lorentz K. Approved recommendation on IFCC methods for the measurement of catalytic concentration of enzymes. Part 9. IFCC method for alpha-amylase (1,4-alpha-D-glucan 4-glucohydrolase, EC 3.2.1.1). International Federation of Clinical Chemistry and Laboratory Medicine (IFCC). Committee on Enzymes. Clin Chem Lab Med. 1998 Mar;36(3): 185-203.
- Junge W, Wortmann W, Wilke B, Waldenström J, Kurlle-Weittenhiller A, Finke J, Klein G. Development and evaluation of assays for the determination of total and pancreatic amylase at 37 degrees C according to the principle recommended by the IFCC. Clin Biochem. 2001 May;34:607-15.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 100-104, (2006).
- Hohenwallner W, Hagele EO, Scholer A et al. Ber Oster Ges Klin Chem 1983;6:101-112.
- Tietz N.W., ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA: WB Saunders, 46-8 (1995).

Data wydania: 05.2022.

ACCENT-200 AMYLASE EPS

Cat. No **7-276** (EN)

INTENDED USE

Diagnostic kit for determination of α -amylase activity, intended to use in automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, BS-120, ACCENT 400 and ACCENT Neo200. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions



INTRODUCTION

α -Amylases are hydrolytic enzymes which hydrolyze 1,4- α -glucosidic bond in starch and other similar polysaccharides to maltose and other oligosaccharides. Several types of amylases can be distinguished, depending on the organ they are originating from. α -amylase is the most commonly measured in the diagnosis of acute pancreatitis, when its activity in serum is very high. Elevation of α -amylase activity in serum is also accompanied by increased excretion of enzyme in urine which can last longer than in the blood. Because of that activity in α -amylase in urine is used as a indicator of acute pancreatitis. Hyperamylasemia occurs also in chronic pancreatitis, failures of kidneys, lungs, diseases of the salivary glands, cerebral traumas, surgical interventions and macroamylasemia. To confirm pancreatic specificity it is recommended to determine also other pancreas specific enzyme like lipase.

METHOD PRINCIPLE

Enzymatic colorimetric method, with EPS substrate, in accordance to recommendations of IFCC – International Federation of Clinical Chemistry and Laboratory Medicine (modified IFCC method).

α -Amylase catalyzes hydrolysis of substrate 4,6-ethylidene-(G7)-p-nitrophenyl-(G1)- α -D-maltoheptaoside (EPS, Ethylidene Protected Substrate). Ethylidene group prevents the substrate from breaking down because of exo-enzymes activity, therefore in absence of α -amylase no increase of absorbance is observed. α -Amylase hydrolyses the substrate into smaller fragments which are acted upon by α -glucosidase, causing the ultimate release of chromophore p-nitrophenol (pNP) and glucose.

Increase of absorbance related to formation of p-nitrophenol is proportional to the α -amylase activity in sample and is measured spectrophotometrically at 405 nm wavelength.

REAGENTS

Package

1-Reagent 1 x 25 ml
2-Reagent 1 x 7 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 12 weeks (ACCENT-200, ACCENT M320).

Concentrations in the test

HEPES buffer, pH 7.2	52.5 mmol/l
sodium chloride	87 mmol/l
magnesium chloride	12.6 mmol/l
calcium chloride	0.075 mmol/l
α -glucosidase	$\geq 4\text{ kU/l}$
4,6-ethylidene G7pNP(EPS)	$> 4\text{ mmol/l}$
stabilizers and preservatives	

Warnings and notes

- Prevent the reagent from microbiological contamination and from saliva and sweat α -amylase.
- Protect from direct sunlight.
- The reagents must be clear, do not use if turbid.
- A slight yellow colour of 2-Reagent does not affect the result of the determination.
- 1-Reagent and 2-Reagent meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

1-Reagent and 2-Reagent contain reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one, mixture (3:1).

Warning

H317 May cause an allergic skin reaction.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P302+P352 IF ON SKIN: Wash with plenty of soap and water.
P333+P313 If skin irritation or rash occurs: Get medical advice.

P363 Wash contaminated clothing before reuse.

SPECIMEN

Serum or plasma collected on heparin, free from hemolysis, urine.

Do not use anticoagulants: EDTA, citrates and oxalates as they inhibit amylase activity.

Serum / plasma can be stored for 7 days at 15-25°C or for one month at 2-8°C.⁷

Urine can be stored for 2 days at 15-25°C or for 10 days at 2-8°C.⁶ Amylase is very unstable in acid urine. Adjust pH to approximately 7.0 before storage.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use. Deionised water is recommended as a reagent blank.

Actions required:

When performing assays in analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S and BS-120, there is a probability of **cross-contamination** affecting the tests results: AMYLASE EPS – GLUCOSE, AMYLASE EPS - MG. To avoid this effect follow the recommendations contained in the instruction 51_03_24_001_ACCENT-200_CARRYOVER.

REFERENCE VALUES⁵

	U/l	$\mu\text{kat/l}$
serum / plasma	28 – 100	0.47 – 1.7
urine	≤ 460	≤ 7.7

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the following controls: CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum; CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine.

For the calibration of automatic analysers: ACCENT-200, ACCENT-200 II GEN, BS-120, the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

For the calibration of automatic analyser ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320 the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 12 weeks (ACCENT-200, ACCENT M320), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers: ACCENT-200 and/or ACCENT-200 II GEN and ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

Sensitivity

2.1 U/l (0.04 $\mu\text{kat/l}$) - ACCENT-200
10 U/l (0.167 $\mu\text{kat/l}$) – ACCENT MC240

Linearity

up to 2000 U/l (33.3 $\mu\text{kat/l}$) - ACCENT-200
up to 1650 U/l (27.5 $\mu\text{kat/l}$) – ACCENT MC240

For higher activity, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by the dilution factor.

Specificity / Interferences

Haemoglobin up to 0.156 g/dl, bilirubin up to 20 mg/dl, ascorbate up to 62 mg/l, triglycerides up to 1250 mg/dl and glucose up to 2000 mg/dl do not interfere with the test.

Precision

Repeatability (run to run)		Mean [U/l]	SD [U/l]	CV [%]
ACCENT 200 II GEN n=20	level 1	67.4	0.76	1.13
	level 2	385.3	3.95	1.03
ACCENT MC240 n=20	level 1	74.03	1.1	1.48
	level 2	371.19	3.79	1.02
Reproducibility (day to day)		Mean [U/l]	SD [U/l]	CV [%]
ACCENT 200 II GEN n=80	level 1	66.4	1.07	1.61
	level 2	376.4	4.11	1.09
ACCENT MC240 n=80	level 1	72.7	2.62	3.6
	level 2	367.2	11.29	3.1

Method comparison

A comparison between amylase values determined at ACCENT-200 (y) and at ADVIA 1650 (x) using 30 samples gave following results:

$y = 1.031x - 6.1311\text{ U/l}$;
 $R = 0.992$ (R – correlation coefficient)

A comparison between amylase values determined at ACCENT MC240 (y) and at ADVIA 1800 (x) using 58 serum samples gave following results:

$y = 0.9529x + 1.0293\text{ U/l}$;
 $R = 0.999$ (R – correlation coefficient)

A comparison between amylase values determined at ACCENT MC240 (y) and at BS-400 (x) using 68 urine samples gave following results:

$y = 0.9931x + 1.4537\text{ U/l}$;
 $R = 1.000$ (R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

- Pesce A.J. Kaplan. L.A.: “Methods in Clinical Chemistry”. Mosby Ed. (1987).
- Burtis C.A. Ashwood E.R.: “Tietz Textbook of Clinical Chemistry”. W.B. Saunders Company Ed. (3rd edition. 1999).
- Lorentz K. Approved recommendation on IFCC methods for the measurement of catalytic concentration of enzymes. Part 9. IFCC method for alpha-amylase (1.4-alpha-D-glucan 4-glucanohydrolase. EC 3.2.1.1). International Federation of Clinical Chemistry and Laboratory Medicine (IFCC). Committee on Enzymes. Clin Chem Lab Med. 1998 Mar;36(3): 185-203.
- Junge W. Wortmann W. Wilke B. Waldenström J. Kurrle-Weittenhiller A. Finke J. Klein G. Development and evaluation of assays for the determination of total and pancreatic amylase at 37 degrees C according to the principle recommended by the IFCC. Clin Biochem. 2001 May; 34:607-15.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests. 4th ed. WB Saunders. 100-104. (2006).
- Hohenwallner W. Hagele EO. Scholer A et al. Ber Oster Ges Klin Chem 1983; 6:101-112.
- Tietz N.W., ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA: WB Saunders, 46-8 (1995).

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- Pesce, A.J., Kaplan, L.A.: “Methods in Clinical Chemistry”, Mosby Ed. (1987).
- Burtis C.A., Ashwood E.R.: “Tietz Textbook of Clinical Chemistry”, W.B. Saunders Company Ed. (3rd edition, 1999).
- Lorentz K. Approved recommendation on IFCC methods for the measurement of catalytic concentration of enzymes. Part 9. IFCC method for alpha-amylase (1,4-alpha-D-glucanohydrolase, EC 3.2.1.1). International Federation of Clinical Chemistry and Laboratory Medicine (IFCC). Committee on Enzymes. Clin Chem Lab Med. 1998 Mar;36(3):185-203.
- Junge W, Wortmann W, Wilke B, Waldenström J, Kurlle-Weittenhiller A, Finke J, Klein G. Development and evaluation of assays for the determination of total and pancreatic amylase at 37 degrees C according to the principle recommended by the IFCC. Clin Biochem. 2001 May;34:607-15.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 100-104, (2006).
- Hohenwallner W, Hagele EO, Scholer A et al. Ber Oster Ges Klin Chem 1983;6:101-112.
- Tietz N.W., ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA: WB Saunders, 46-8 (1995).

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Calibration Rule

Rule	One-point Linear
Sensitivity	1
Replicates	3
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT-200 II GEN

Parameters			
Test Name	AMY EPS	R1	200
Test No	63	R2	50
Full Name	Amylase EPS	Sample Volume	6
Reference No	63	R1 Blank	
Analy. Type	Kinetic	Mixed Reag. Blank	
Pri. Wave.	405 nm	Concentration	2.2 2000
Secon. Wave.	670 nm	Linearity Limit	0.2
Trend	Increase	Substrate Limit	
Reac. Time	8 18	Factor	
Incuba. Time	15	<input type="checkbox"/> Prozone check	
Unit	U/l	q1	<input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>
Precision	0.1	PC	<input type="checkbox"/> Abs <input type="checkbox"/>

Calibration Rule

Rule	One-point Linear
Sensitivity	1
Replicates	3
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT-220S

Parameters			
Test	AMY EPS	R1	200
No	63	R2	50
Full Name	Amylase EPS	Sample Volume	6
Standard No	63	R1 Blank	
Reac. Type	Kinetic	Mixed Rtg. Blank	
Pri. Wave.	405 nm	Linearity Range	6 2200
Sec. Wave.	670 nm	Linearity Limit	0.2
Direction	Increase	Substrate Limit	
Reac. Time	8 18	Factor	
Incuba. Time	16	<input type="checkbox"/> Prozone check	
Unit	U/l	q1	<input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>
Precision	0.1	PC	<input type="checkbox"/> Abs <input type="checkbox"/>

Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	3
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• BS-120

Parameters			
Test	AMY EPS	R1	200
No	63	R2	50
Full Name	Amylase EPS	Sample Volume	6
Standard No	63	R1 Blank	
Reac. Type	Kinetic	Mixed Rtg. Blank	
Pri. Wave.	405 nm	Linearity Range	2.1 2000
Sec. Wave.	670 nm	Linearity Limit	0.2
Direction	Increase	Substrate Limit	
Reac. Time	8 18	Factor	
Incuba. Time	16	<input type="checkbox"/> Prozone check	
Unit	U/l	q1	<input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>
Precision	Integer	PC	<input type="checkbox"/> Abs <input type="checkbox"/>

Calibration Rule

Rule	One-point Linear
Sensitivity	1
Replicates	3
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT S120

Chem	AMY EPS	No.	063	Sample Type	SERUM
Chemistry	AMYLASE EPS	Reaction Type	Kinetic	Reaction Direction	positive
Pri Wave	405nm	Unit	U/L	Sec Wave	670nm
Blank Time		Decimals	0.1	Incubation Time	12
Standard	5	Reagent Vol	160	Reaction Time	8 18
Decreased	5	R1	160	R2	40
Increased		R2	40		
Linearity range (Standard)	5	1600	Linearity Limit	0.2	
Linearity Range (Decreased)			Substrate Depletion	40000	
Linearity Range (Increased)			Mixed Blank Abs	-40000 40000	
R1 Blank Abs	-40000	40000	On-board Stability	Day(s)	
Blank Response	-40000	40000	Reagent Alarm Limit		
Twin Chemistry			<input type="checkbox"/> Enzyme Linear Extension		
<input type="checkbox"/> Prozone Check					
Q1		Q2		V1	
Q3		Q4		V2	
Q5		Q6		PC1	
PC2					
<input type="checkbox"/> Sample Pretreatment		<input type="checkbox"/> Control Pretreatment		<input type="checkbox"/> Calibrator Pretreatment	
<input type="checkbox"/> Pretreat Sample Vol		<input type="checkbox"/> μ L		<input type="checkbox"/> Pretreat Sample Vol	<input type="checkbox"/> μ L
CALIBRATION SETTINGS	Math model	Multi-point Linear	AUTO CALIBRATION	<input type="checkbox"/> Bottle Changed	
Factor		Replicates	2	<input type="checkbox"/> Lot Changed	
				<input type="checkbox"/> Cal Time	
ACCEPTANCE LIMITS	Cal Time	Hour			
Slope Diff		SD			
Sensitivity		Repeatability	40000		
Deter Coeff					

ACCENT-200 AMYLASE EPS

• ACCENT MC240

Chem No. Sample Type

Chemistry Print name

Reaction Type Reaction Direction

Pri Wave Sec Wave

Unit Decimal

Incubation Time

Blank Time Reaction Time

Standard Sample Vol μ L Aspirated μ L Diluent μ L Reagent Vol

Decreased μ L μ L μ L R1 μ L

Increased μ L μ L μ L R2 μ L

Sample Blank Auto Rerun

Linearity range (Standard) Linearity Limit

Linearity Range (Decreased) Substrate Depletion

Linearity Range (Increased) Mixed Blank Abs

R1 Blank Abs On-board Stability Day(s)

Blank Response Reagent Alarm Limit

Twin Chemistry Enzyme Linear Extension

Prozone Check

Q1 Q2 V1 Q3 Q4 V2

Q5 Q6 V3 PC1 PC2

Sample Pretreatment Control Pretreatment Calibrator Pretreatment

Pretreat Sample Vol μ L Pretreat Sample Vol μ L

CALIBRATION SETTINGS

Math model

Factor Replicates

AUTO CALIBRATION

Bottle Changed

Lot Changed

Cal Time

ACCEPTANCE LIMITS

Cal Time Hour

Slope Diff SD

Sensitivity Repeatability

Deter Coeff

• ACCENT M320

Chem No. Sample Type

Chemistry Print name

Reaction Type Reaction Direction

Pri Wave Sec Wave

Unit Decimal

Incubation Time

Blank Time Reaction Time

Standard Sample Vol μ L Aspirated μ L Diluent μ L Reagent Vol

Decreased μ L μ L μ L R1 μ L

Increased μ L μ L μ L R2 μ L

Sample Blank Auto Rerun

Linearity range (Standard) Linearity Limit

Linearity Range (Decreased) Substrate Depletion

Linearity Range (Increased) Mixed Blank Abs

R1 Blank Abs On-board Stability Day(s)

Blank Response Reagent Alarm Limit

Twin Chemistry Enzyme Linear Extension

Prozone Check

Q1 Q2 V1 Q3 Q4 V2

Q5 Q6 V3 PC1 PC2

Sample Pretreatment Control Pretreatment Calibrator Pretreatment

Pretreat Sample Vol μ L Pretreat Sample Vol μ L

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