

ACCENT-200 CREA ENZYMATIC

PROGRAM NA ANALIZATORY / APPLICATION for / АДАПТАЦИЯ для:

• ACCENT-200, ACCENT-200 II GEN

Parameters	CREA ENZ	R1	180
Test Name	71	R2	60
Test No	Crea Enzymatic	Sample Volume	8
Full Name	71	R1 Blank	
Reference No	Endpoint	Mixed Reag. Blank	
Analy. Type	546 nm	Concentration	0.08 25.3
Pri. Wave.		Linearity Limit	
Secon. Wave.	Ascending	Substrate Limit	
Trend	-2 18	Factor	
Reac. Time	15	<input type="checkbox"/> Prozone check	
Incuba. Time	0.01	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Unit		PC <input type="checkbox"/> Abs <input type="checkbox"/>	
Precision			

Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	3
Interval (day)	56
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT-220S

Parameters	CREA ENZ	R1	180
Test	71	R2	60
No	Crea Enzymatic	Sample Volume	8
Full Name	71	R1 Blank	
Standard No	Endpoint	Mixed Rtg. Blank	
Reac. Type	546 nm	Linearity Range	0.08 30
Pri. Wave.		Linearity Limit	
Sec. Wave.	Increase	Substrate Limit	
Direction	-2 20	Factor	
Reac. Time	19	<input type="checkbox"/> Prozone check	
Incuba. Time	mg/dl	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Unit	0.01	PC <input type="checkbox"/> Abs <input type="checkbox"/>	
Precision			

Calibration Rule

Rule	Multi-point Linear
Sensitivity	1
Replicates	3
Interval (day)	56
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• BS-120

Parameters	CREA ENZ	R1	240
Test	68	R2	80
No	Crea Enzymatic	Sample Volume	14
Full Name	68	R1 Blank	
Standard No	Endpoint	Mixed Rtg. Blank	
Reac. Type	546 nm	Linearity Range	0.04 25.5
Pri. Wave.		Linearity Limit	
Sec. Wave.	670 nm	Substrate Limit	
Direction	Increase	Factor	
Reac. Time	-1 18	<input type="checkbox"/> Prozone check	
Incuba. Time	16	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Unit	mg/dl	PC <input type="checkbox"/> Abs <input type="checkbox"/>	
Precision	0.01		

Calibration Rule

Rule	Two-point / Multi-point Linear
Sensitivity	1
Replicates	3
Interval (day)	56
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT S120

Chem	CREA ENZ	No.	071	Sample Type	SERUM	
Chemistry	CREATININE ENZYMATIC	Print name	CREA ENZ			
Reaction Type	Endpoint	Reaction Direction	positive			
Pri Wave	546nm	Sec Wave				
Unit	mg/dl	Decimal	0.01			
Blank Time	-3	-1	Incubation Time	19	Reaction Time	20 22
Standard	Sample Vol	Aspirated	Diluent	Reagent Vol		
Decreased	5.3	20	180	R1 120	R2 40	
Increased						
		Sample Blank	V	Auto Rerun		
Linearity range (Standard)	0.1	32.6	Linearity Limit			
Linearity Range (Decreased)			Substrate Depletion			
Linearity Range (Increased)			Mixed Blank Abs	-40000 40000		
R1 Blank Abs	-40000	40000	On-board Stability		Day(s)	
Blank Response	-40000	40000	Reagent Alarm Limit			
Twin Chemistry			<input type="checkbox"/> Enzyme Linear Extension			
<input type="checkbox"/> Prozone Check						
Q1	Q2	V1	Q3	Q4	V2	
Q5	Q6	V3	PC1	PC2		
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment	<input type="checkbox"/> Calibrator Pretreatment				
	Pretreat Sample Vol		Pretreat Sample Vol			
CALIBRATION SETTINGS	Math model	Multi-point Linear	AUTO CALIBRATION			
Factor	Replicates	2	<input type="checkbox"/> Bottle Changed			
			<input type="checkbox"/> Lot Changed			
			<input type="checkbox"/> Cal Time			
ACCEPTANCE LIMITS	Cal Time	Hour				
Slope Diff	SD					
Sensitivity	Repeatability	40000				
Deter Coeff						

ACCENT-200 CREA ENZYMATIC

• ACCENT MC240

Chem <input type="text" value="CREA ENZ"/>	No. <input type="text" value="071"/>	Sample Type <input type="text" value="SERUM"/>
Chemistry <input type="text" value="CREATININE ENZYMATIC"/>	Print name <input type="text" value="CREA ENZ"/>	
Reaction Type <input type="text" value="Endpoint"/>	Reaction Direction <input type="text" value="positive"/>	
Pri Wave <input type="text" value="546nm"/>	Sec Wave <input type="text"/>	
Unit <input type="text" value="mg/dl"/>	Decimal <input type="text" value="0,01"/>	
Blank Time <input type="text" value="-3"/> <input type="text" value="-1"/>	Incubation Time <input type="text" value="21"/>	
	Reaction Time <input type="text" value="28"/> <input type="text" value="30"/>	
Standard <input type="text" value="5.3"/> <input type="text" value="5.3"/> <input type="text" value="5.3"/>	Aspirated <input type="text" value="20"/> <input type="text" value="20"/> <input type="text" value="20"/>	Diluent <input type="text" value="180"/> <input type="text" value="180"/> <input type="text" value="180"/>
Decreased <input type="text" value="5.3"/>	Reagent Vol <input type="text" value="120"/>	R1 <input type="text" value="40"/>
Increased <input type="text"/>	R2 <input type="text" value="40"/>	
<input type="checkbox"/> Sample Blank <input checked="" type="checkbox"/> Auto Rerun		

Linearity range (Standard) <input type="text" value="0.05"/> <input type="text" value="32"/>	Linearity Limit <input type="text"/>
Linearity Range (Decreased) <input type="text"/>	Substrate Depletion <input type="text"/>
Linearity Range (Increased) <input type="text"/>	Mixed Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>
R1 Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>	On-board Stability <input type="text"/> Day(s)
Blank Response <input type="text" value="-35000"/> <input type="text" value="35000"/>	Reagent Alarm Limit <input type="text"/>
Twin Chemistry <input type="text"/>	<input type="checkbox"/> Enzyme Linear Extension
<input type="checkbox"/> Prozone Check	

Q1 <input type="text"/>	Q2 <input type="text"/>	V1 <input type="text"/>	Q3 <input type="text"/>	Q4 <input type="text"/>	V2 <input type="text"/>
Q5 <input type="text"/>	Q6 <input type="text"/>	V3 <input type="text"/>	PC1 <input type="text"/>	PC2 <input type="text"/>	
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment	<input type="checkbox"/> Calibrator Pretreatment			
Pretreat Sample Vol <input type="text"/> μ L	Pretreat Sample Vol <input type="text"/> μ L				

CALIBRATION SETTINGS	AUTO CALIBRATION
Math model <input type="text" value="Multi-point Linear"/>	<input type="checkbox"/> Bottle Changed
Factor <input type="text"/>	<input type="checkbox"/> Lot Changed
Replicates <input type="text" value="2"/>	<input type="checkbox"/> Cal Time

ACCEPTANCE LIMITS
Cal Time <input type="text"/> Hour
Slope Diff <input type="text"/> SD <input type="text"/>
Sensitivity <input type="text"/> Repeatability <input type="text" value="35000"/>
Deter Coeff <input type="text"/>

• ACCENT M320

Chem <input type="text" value="CREA ENZ"/>	No. <input type="text" value="071"/>	Sample Type <input type="text" value="SERUM"/>
Chemistry <input type="text" value="CREATININE ENZYMATIC"/>	Print name <input type="text" value="CREA ENZ"/>	
Reaction Type <input type="text" value="Endpoint"/>	Reaction Direction <input type="text" value="positive"/>	
Pri Wave <input type="text" value="546nm"/>	Sec Wave <input type="text"/>	
Unit <input type="text" value="mg/dl"/>	Decimal <input type="text" value="0,01"/>	
Blank Time <input type="text" value="-3"/> <input type="text" value="-1"/>	Incubation Time <input type="text" value="18"/>	
	Reaction Time <input type="text" value="28"/> <input type="text" value="30"/>	
Standard <input type="text" value="5.3"/> <input type="text" value="5.3"/> <input type="text" value="5.3"/>	Aspirated <input type="text" value="20"/> <input type="text" value="20"/> <input type="text" value="20"/>	Diluent <input type="text" value="180"/> <input type="text" value="180"/> <input type="text" value="180"/>
Decreased <input type="text" value="5.3"/>	Reagent Vol <input type="text" value="120"/>	R1 <input type="text" value="40"/>
Increased <input type="text"/>	R2 <input type="text" value="40"/>	
<input type="checkbox"/> Sample Blank <input checked="" type="checkbox"/> Auto Rerun		

Linearity range (Standard) <input type="text" value="0.10"/> <input type="text" value="28"/>	Linearity Limit <input type="text"/>
Linearity Range (Decreased) <input type="text"/>	Substrate Depletion <input type="text"/>
Linearity Range (Increased) <input type="text"/>	Mixed Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>
R1 Blank Abs <input type="text" value="-35000"/> <input type="text" value="35000"/>	On-board Stability <input type="text"/> Day(s)
Blank Response <input type="text" value="-35000"/> <input type="text" value="35000"/>	Reagent Alarm Limit <input type="text"/>
Twin Chemistry <input type="text"/>	<input type="checkbox"/> Enzyme Linear Extension
<input type="checkbox"/> Prozone Check	

Q1 <input type="text"/>	Q2 <input type="text"/>	V1 <input type="text"/>	Q3 <input type="text"/>	Q4 <input type="text"/>	V2 <input type="text"/>
Q5 <input type="text"/>	Q6 <input type="text"/>	V3 <input type="text"/>	PC1 <input type="text"/>	PC2 <input type="text"/>	
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment	<input type="checkbox"/> Calibrator Pretreatment			
Pretreat Sample Vol <input type="text"/> μ L	Pretreat Sample Vol <input type="text"/> μ L				

CALIBRATION SETTINGS	AUTO CALIBRATION
Math model <input type="text" value="Multi-point Linear"/>	<input type="checkbox"/> Bottle Changed
Factor <input type="text"/>	<input type="checkbox"/> Lot Changed
Replicates <input type="text" value="2"/>	<input type="checkbox"/> Cal Time

ACCEPTANCE LIMITS
Cal Time <input type="text"/> Hour
Slope Diff <input type="text"/> SD <input type="text"/>
Sensitivity <input type="text"/> Repeatability <input type="text" value="35000"/>
Deter Coeff <input type="text"/>

Data wydania / Date of issue / Дата создания: 02. 2023.