ACCENT-200 LDL DIRECT

Nr kat. **7-280** (PL)

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Cat. No **7-280** (EN)

INTENDED USE

Diagnostic kit for determination of LDL-cholesterol concentration (direct method) used in automatic analysers ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, BS-120, ACCENT 400 and ACCENT Neo200.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Plasma lipoproteins are spherical particles containing varying amounts of cholesterol, triglycerides, phospholipids and proteins. The relative protein and lipid determine the density of these lipoproteins and provide the basis on which to begin their classification. The classes are: chylomicron, very-low-density lipoprotein (VLDL), low-density-lipoprotein (LDL) and high-density lipoprotein (HDL).

LDL are synthesized in the liver by the action of various lipolytic enzymes on triglyceride rich VLDL. LDL-cholesterol concentration is considered to be the most important clinical predictor, of all single parameters, with respect to coronary atherosclerosis.

Accurate measurement of LDL-cholesterol is of vital importance in therapies which focus on lipid reduction to prevent atherosclerosis or reduce its progress and to avoid plaque rupture.

METHOD PRINCIPLE

The assay is a homogenous method for directly measuring LDL-cholesterol concentrations in serum or plasma, without the need for any off-line pretreatment or centrifugation steps.

Liquid selective detergent method.

The method is in a two reagents format and depends on the properties of a unique detergent. This detergent (1-REAGENT) solubilizes only the non LDL particles (HDL, VLDL, CM).

The cholesterol released is consumed by cholesterol esterase and cholesterol oxidase in a non color forming reaction. A second detergent (2-REAGENT) solubilizes the remaining LDL particles and chromogenic coupler allows for color formulation. The enzyme reaction with LDL-cholesterol in the presence of the coupler produces color that is proportional to the amount of LDL cholesterol present in the sample.

REAGENTS

Package

1-REAGENT 1 x 31 ml
2-REAGENT 1 x 11 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 12 weeks (ACCENT-200, ACCENT M320).

Concentrations in the test

1-REAGENT

Buffer < 1.0 %
Detergent 1

Cholesterol esterase (<i>Pseudomonas</i> spp.)	< 1500 U/l
Cholesterol oxidase (<i>Cellulomonas</i> spp.)	< 1500 U/l
Peroxidase (horseradish)	< 1300 ppg U/l
4-aminoantipyrine	< 0.1 %
Ascorbic Acid Oxidase (<i>Curcubita</i> spp.)	< 3000 U/l
Preservative	
2-REAGENT	
Buffer	
Detergent 2	< 1.0 %
N,N-bis(sulfobutyl)-toluidine, disodium (DSBmT)	< 1.0 mM
Preservative	

Warnings and notes

- Do not freeze reagents.
- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT and 2-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

1-REAGENT and 2-REAGENT contain reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1).

Warning

 H317 May cause an allergic skin reaction.
P280 Wear protective gloves, protective clothing, eye protection or face protection.
P302 + P352 IF ON SKIN: Wash with plenty of soap and water.

SPECIMEN

Serum, sodium heparinized or EDTA plasma. Anticoagulants containing citrate should not be used. Serum: Collect whole blood by venipuncture and allow to clot. Centrifuge and remove the serum as soon as possible after collection (within 3 hours). Plasma: Specimens may be collected in EDTA or sodium heparin. Centrifuge and remove the plasma as soon as possible after collection (within 3 hours). If not analysed promptly, specimens may be stored at 2-8°C for up to 5 days. If specimens need to be stored for more than 5 days, they may be preserved at -80°C. Samples may be frozen once. Nevertheless it is recommended to perform the assay with freshly collected samples.

PROCEDURE

1-REAGENT and 2-REAGENT are ready to use. Deionized water is recommended as a reagent blank. For analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S and BS-120 it is recommended to determine the reagent blank during each calibration. Deionized water should be used as reagent blank. When performing calibration, the task type **Calib+Rgt.Blk** should be selected.

Actions required:

When performing assays in analyzers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S and BS-120 there is a probability of **cross-contamination** affecting the tests results: LDL DIRECT

II GEN - TG mono, LDL DIRECT II GEN - UA, LDL DIRECT II GEN - UA PLUS, LDL DIRECT II GEN - LIPASE, LDL DIRECT II GEN - MICROALBUMIN. To avoid this effect follow the recommendations contained in the instruction 51_03_24_001_ACCENT-200_CARRYOVER.

REFERENCE VALUES ⁷

NCEP* Classification	Adults	
	mg/dl	mmol/l
Optional	< 100	< 2.59
Near optional	< 130	< 3.37
Bordeline high	130-159	3.37-4.12
High	160-189	4.14-4.90
Very high	≥ 190	≥ 4.92

* National Cholesterol Education Program

As LDL cholesterol is affected by a number of factors such as smoking, exercise, hormones, age and sex, each laboratory should establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) with each batch of samples.

For the calibration of automatic analysers: ACCENT-200, ACCENT-200 II GEN, ACCENT-220S, ACCENT S120, ACCENT MC240, ACCENT M320, BS-120 the CORMAY HDL/LDL CALIBRATOR (Cat. No 5-178) is recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 12 weeks (ACCENT-200, ACCENT M320), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers ACCENT-200 and ACCENT MC240. Results may vary if a different instrument or a manual procedure is used.

Sensitivity

2.0 mg/dl (0.05 mmol/l) - ACCENT-200
3.6 mg/dl (0.09 mmol/l) - ACCENT MC240

Linearity

up to 700 mg/dl (18.13 mmol/l) - ACCENT-200
up to 725 mg/dl (18.78 mmol/l) - ACCENT MC240

For higher concentration of LDL cholesterol dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Triglycerides up to 1293 mg/dl, bilirubin conjugated up to 20 mg/dl, bilirubin total up to 20 mg/dl, haemoglobin up to 0.5 g/dl, ascorbic acid up to 500 mg/l, and gamma-globulins up to 5000 mg/dl do not interfere with the test.

Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 n=10	level 1	141.28	2.99	2.12
	level 2	191.66	5.79	3.02
ACCENT MC240 n=20	level 1	56.53	0.58	1.02
	level 2	132.27	0.90	0.68

Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
ACCENT-200 n=10	level 1	141.50	3.16	2.23
	level 2	189.62	3.50	1.85
ACCENT MC240 n=80	level 1	56.1	1.25	2.2
	level 2	132.4	1.78	1.3

Method comparison

A comparison between CORMAY reagent on ACCENT 200 (y) and commercially available assay on HITACHI 912 (x) using 43 samples gave following results:

y = 0.9439 x + 17.566 mg/dl;
R = 0.988 (R – correlation coefficient)

A comparison between LDL cholesterol values determined at ACCENT MC240 (y) and at ADVIA 1800 (x) using 44 serum samples gave following results:

y = 1.0986 x - 8.6484 mg/dl;
R = 0.992 (R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

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Calibration Rule

Rule	Two-point Linear
Sensitivity	1
Replicates	3
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT-200 II GEN

Parameters

Test Name	LDL D	R1	210
Test No	26	R2	70
Full Name	LDL Direct	Sample Volume	3
Reference No	26	R1 Blank	
Analy. Type	Endpoint	Mixed Reag. Blank	
Pri. Wave.	578 nm	Concentration	5 700
Secon. Wave.	670 nm	Linearity Limit	
Trend	Increase	Substrate Limit	
Reac. Time	0 40	Factor	
Incuba. Time	10	<input type="checkbox"/> Prozone check	
Unit	mg/dl	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.1	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

Calibration Rule

Rule	Two-point Linear
Sensitivity	1
Replicates	3
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT-220S

Parameters

Test	LDL D	R1	210
No	26	R2	70
Full Name	LDL Direct	Sample Volume	3
Standard No	26	R1 Blank	
Reac. Type	Endpoint	Mixed Rtg. Blank	
Pri. Wave.	578 nm	Concentration	3.9 700
Sec. Wave.	670 nm	Linearity Limit	
Direction	Increase	Substrate Limit	
Reac. Time	0 38	Factor	
Incuba. Time	11	<input type="checkbox"/> Prozone check	
Unit	mg/dl	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.1	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

Calibration Rule

Rule	Two-point Linear
Sensitivity	1
Replicates	3
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• BS-120

Parameters

Test	LDL D	R1	225
No	26	R2	75
Full Name	LDL Direct	Sample Volume	3
Standard No	26	R1 Blank	
Reac. Type	Endpoint	Mixed Rtg. Blank	
Pri. Wave.	630 nm	Concentration	1.8 680
Sec. Wave.		Linearity Limit	
Direction	Increase	Substrate Limit	
Reac. Time	4 18	Factor	
Incuba. Time	16	<input type="checkbox"/> Prozone check	
Unit	mg/dl	q1 <input type="checkbox"/> q2 <input type="checkbox"/> q3 <input type="checkbox"/> q4 <input type="checkbox"/>	
Precision	0.1	PC <input type="checkbox"/> Abs <input type="checkbox"/>	

Calibration Rule

Rule	Two-point Linear
Sensitivity	1
Replicates	3
Interval (day)	84
Difference Limit	0
SD	0
Blank Response	0 50000
Error Limit	0
Coefficient	0

• ACCENT S120

Chem	LDL D	No.	026	Sample Type	SERUM						
Chemistry	LDL DIRECT			Print name	LDL D						
Reaction Type	Endpoint			Reaction Direction	positive						
Pri Wave	578 nm			Sec Wave	670 nm						
Unit	mg/dL			Decimal	0.1						
Blank Time	-3	-1		Incubation Time	15						
				Reaction Time	22 24						
Standard	2.5	Aspirated		Diluent		Reagent Vol					
Decreased	2.5	20	180	180	180	R1	180				
Increased						R2	60				
				Sample Blank	V	Auto Rerun					
Linearity range (Standard)	2	730		Linearity Limit							
Linearity Range (Decreased)				Substrate Depletion							
Linearity Range (Increased)				Mixed Blank Abs	-40000	40000					
R1 Blank Abs	-40000	40000		On-board Stability		Day(s)					
Blank Response	-40000	40000		Reagent Alarm Limit							
Twin Chemistry				<input type="checkbox"/> Enzyme Linear Extension							
<input type="checkbox"/> Prozone Check											
Q1		Q2		V1		Q3		Q4		V2	
Q5		Q6		V3		PC1		PC2			
<input type="checkbox"/> Sample Pretreatment		<input type="checkbox"/> Control Pretreatment		<input type="checkbox"/> Calibrator Pretreatment							
		Pretreat Sample Vol		uL		Pretreat Sample Vol		uL			
CALIBRATION SETTINGS	Math model	Two-point linear		AUTO CALIBRATION		Bottle Changed					
Factor		Replicates	2			Lot Changed					
						Cal Time					
ACCEPTANCE LIMITS	Cal Time		Hour								
	Slope Diff		SD								
	Sensitivity		Repeatability	40000							
	Deter Coeff										

ACCENT-200 LDL DIRECT

• ACCENT MC240

Chem	<input type="text" value="LDL D"/>	No.	<input type="text" value="026"/>	Sample Type	<input type="text" value="SERUM"/>
Chemistry	<input type="text" value="LDL DIRECT"/>	Print name	<input type="text" value="LDL D"/>	Reaction Direction	<input type="text" value="positive"/>
Reaction Type	<input type="text" value="Endpoint"/>	Sec Wave	<input type="text" value="800 nm"/>	Incubation Time	<input type="text" value="21"/>
Pri Wave	<input type="text" value="605 nm"/>	Decimal	<input type="text" value="0.1"/>	Reaction Time	<input type="text" value="28"/> <input type="text" value="30"/>
Unit	<input type="text" value="mg/dL"/>	Blank Time	<input type="text" value="-3"/> <input type="text" value="-1"/>	Reagent Vol	<input type="text" value="180"/> <input type="text" value="60"/>
Standard	<input type="text" value="2.5"/> <input type="text" value="μL"/>	Aspirated	<input type="text" value="20"/> <input type="text" value="μL"/>	Diluent	<input type="text" value="180"/> <input type="text" value="μL"/>
Decreased	<input type="text" value="2.5"/> <input type="text" value="μL"/>	Increased	<input type="text" value="20"/> <input type="text" value="μL"/>	Sample Blank	<input type="checkbox"/>
Increased	<input type="text" value=""/> <input type="text" value="μL"/>	Auto Rerun	<input type="checkbox"/>		

Linearity range (Standard)	<input type="text" value="3.6"/> <input type="text" value="725"/>	Linearity Limit	<input type="text"/>
Linearity Range (Decreased)	<input type="text"/>	Substrate Depletion	<input type="text"/>
Linearity Range (Increased)	<input type="text"/>	Mixed Blank Abs	<input type="text" value="-35000"/> <input type="text" value="35000"/>
R1 Blank Abs	<input type="text" value="-35000"/> <input type="text" value="35000"/>	On-board Stability	<input type="text"/> Day(s)
Blank Response	<input type="text" value="-35000"/> <input type="text" value="35000"/>	Reagent Alarm Limit	<input type="text"/>
Twin Chemistry	<input type="text"/>	Enzyme Linear Extension	<input type="checkbox"/>

O1	<input type="text"/>	O2	<input type="text"/>	V1	<input type="text"/>	O3	<input type="text"/>	O4	<input type="text"/>	V2	<input type="text"/>
Q5	<input type="text"/>	Q6	<input type="text"/>	V3	<input type="text"/>	PC1	<input type="text"/>	PC2	<input type="text"/>		<input type="text"/>
Sample Pretreatment	<input type="checkbox"/>	Control Pretreatment	<input type="checkbox"/>	Calibrator Pretreatment	<input type="checkbox"/>	Pretreat Sample Vol	<input type="text"/> <input type="text" value="μL"/>	Pretreat Sample Vol	<input type="text"/> <input type="text" value="μL"/>		

CALIBRATION SETTINGS

Math model	<input type="text" value="Two-point linear"/>
Factor	<input type="text"/>
Replicates	<input type="text" value="2"/>

AUTO CALIBRATION

Bottle Changed	<input type="checkbox"/>
Lot Changed	<input type="checkbox"/>
Cal Time	<input type="text"/>

ACCEPTANCE LIMITS

Cal Time	<input type="text"/> Hour
Slope Diff	<input type="text"/>
Sensitivity	<input type="text"/>
Deter Coeff	<input type="text"/>
SD	<input type="text"/>
Repeatability	<input type="text" value="35000"/>

• ACCENT M320

Chem	<input type="text" value="LDL D"/>	No.	<input type="text" value="026"/>	Sample Type	<input type="text" value="SERUM"/>
Chemistry	<input type="text" value="LDL DIRECT"/>	Print name	<input type="text" value="LDL D"/>	Reaction Direction	<input type="text" value="positive"/>
Reaction Type	<input type="text" value="Endpoint"/>	Sec Wave	<input type="text" value="800 nm"/>	Incubation Time	<input type="text" value="30"/>
Pri Wave	<input type="text" value="570 nm"/>	Decimal	<input type="text" value="0.1"/>	Reaction Time	<input type="text" value="26"/> <input type="text" value="28"/>
Unit	<input type="text" value="mg/dL"/>	Blank Time	<input type="text" value="-3"/> <input type="text" value="-1"/>	Reagent Vol	<input type="text" value="180"/> <input type="text" value="60"/>
Standard	<input type="text" value="2.5"/> <input type="text" value="μL"/>	Aspirated	<input type="text" value="20"/> <input type="text" value="μL"/>	Diluent	<input type="text" value="180"/> <input type="text" value="μL"/>
Decreased	<input type="text" value="2.5"/> <input type="text" value="μL"/>	Increased	<input type="text" value="20"/> <input type="text" value="μL"/>	Sample Blank	<input type="checkbox"/>
Increased	<input type="text" value=""/> <input type="text" value="μL"/>	Auto Rerun	<input type="checkbox"/>		

Linearity range (Standard)	<input type="text" value="4.2"/> <input type="text" value="750"/>	Linearity Limit	<input type="text"/>
Linearity Range (Decreased)	<input type="text"/>	Substrate Depletion	<input type="text"/>
Linearity Range (Increased)	<input type="text"/>	Mixed Blank Abs	<input type="text" value="-35000"/> <input type="text" value="35000"/>
R1 Blank Abs	<input type="text" value="-35000"/> <input type="text" value="35000"/>	On-board Stability	<input type="text"/> Day(s)
Blank Response	<input type="text" value="-35000"/> <input type="text" value="35000"/>	Reagent Alarm Limit	<input type="text"/>
Twin Chemistry	<input type="text"/>	Enzyme Linear Extension	<input type="checkbox"/>

Q1	<input type="text"/>	Q2	<input type="text"/>	V1	<input type="text"/>	Q3	<input type="text"/>	Q4	<input type="text"/>	V2	<input type="text"/>
O5	<input type="text"/>	O6	<input type="text"/>	V3	<input type="text"/>	PC1	<input type="text"/>	PC2	<input type="text"/>		<input type="text"/>
Sample Pretreatment	<input type="checkbox"/>	Control Pretreatment	<input type="checkbox"/>	Calibrator Pretreatment	<input type="checkbox"/>	Pretreat Sample Vol	<input type="text"/> <input type="text" value="μL"/>	Pretreat Sample Vol	<input type="text"/> <input type="text" value="μL"/>		

CALIBRATION SETTINGS

Math model	<input type="text" value="Two-point linear"/>
Factor	<input type="text"/>
Replicates	<input type="text" value="2"/>

AUTO CALIBRATION

Bottle Changed	<input type="checkbox"/>
Lot Changed	<input type="checkbox"/>
Cal Time	<input type="text"/>

ACCEPTANCE LIMITS

Cal Time	<input type="text"/> Hour
Slope Diff	<input type="text"/>
Sensitivity	<input type="text"/>
Deter Coeff	<input type="text"/>
SD	<input type="text"/>
Repeatability	<input type="text" value="35000"/>

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