ACCENT-300 CYSTATIN C

Nr kat. 7-300 (PL)

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Cat. No **7-300** (EN)

INTENDED USE

Diagnostic kit for determination of cystatin C concentration, used in automatic analyser ACCENT-300.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Cystatin C is a low molecular weight protein (13kD), one of the cysteine proteinases inhibitors. Cystatin C is produced in all nucleated cells and secreted into extracellular space at constant rate. Cystatin C molecule stability and dependence of its concentration solely on GFR (Glomerular Filtration Rate) decide about high diagnostic efficiency of cystatin C determination. Cystatin C level is not affected by muscle mass or diet and its increase is observed even at slight GFR reduction. Clinical applications of cystatin C are for monitoring GFR in children, elderly patients, patients with potentially nephrotoxic drug therapy, for assessment of renal transplantation status, for kidney function monitoring in acute and chronic kidney diseases including diabetic nephropathy.

METHOD PRINCIPLE

Turbidimetric method. An antigen-antibody reaction occurs between cystatin C and antibodies coated on polystyrene particles and immuno-complexes are formed. The change of turbidity is related to the quantity of cystatin C in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

REAGENTS

Package

1-Reagent 2 x 13.5 ml
 2-Reagent 1 x 5.2 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 9 weeks on board the analyser at 2-10°C.

Concentrations in the test

suspension of polystyrene particles coated with anti-cystatin C antibodies 8.5 g/l
 MOPS buffer [3-(N-morpholino)-propanesulfonic acid] 8 g/l
 gentamycin 12.5 mg/l
 amphotericin B 1.25 mg/l
 preservative

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- The reagents contain antibiotics and must be handled with due cautions.
- Mix samples well before analysing.

- The samples can be shipped without special cooling and must then be analysed within 14 days after shipment.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum or EDTA/heparinised plasma.

Sample stability testing showed that cystatin C in serum and plasma samples are stable for 14 days at room temperature (8-25°C) and for 21 days if stored at 2-8°C, otherwise store frozen at below -20°C. Mix samples well before analysing. The samples can be shipped without special cooling and must then be analysed within 14 days after shipment. Nevertheless it is recommended to perform the assay with freshly collected samples.

PROCEDURE

1-Reagent and 2-Reagent are ready to use.
 For reagent blank 0.9% NaCl is recommended.

REFERENCE VALUES ⁹

serum	mg/l
newborns < 1 month	1.49 – 2.85
children < 5 months	1.01 – 1.92
< 12 months	0.75 – 1.53
12-24 months (F)	0.60 – 1.20
12-24 months (M)	0.77 – 1.85
2-19 years	0.62 – 1.11
adults (F)	0.61 – 1.05
adults (M)	0.71 – 1.21

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY CYSTATIN C CONTROLS (Cat. No 4-460) with each batch of samples.

For the calibration of automatic analysers the CORMAY CYSTATIN C CALIBRATORS (Cat. No 5-185) is recommended. The calibration curve should be prepared every 9 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

These metrological characteristics have been obtained using an automatic analysers ACCENT-300 and Modular P. Results may vary if a different instrument is used.

- Sensitivity:** 0.37 mg/l.
- Linearity:** up to 8.39 mg/l.
- Specificity / Interferences**
 Haemoglobin up to 0.7 g/dl, bilirubin up to 80 mg/dl, triglycerides up to 1400 mg/dl and ascorbate up to 300 mg/l do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean [mg/l]	SD [mg/l]	CV [%]
level 1	1.92	0.021	1.10
Reproducibility (day to day) n = 20	Mean [mg/l]	SD [mg/l]	CV [%]
level 1	0.878	0.024	2.727
level 2	5.245	0.230	4.384

Method comparison

A comparison between cystatin C values determined at ACCENT-300 (y) and at OLYMPUS AU400 (x) using 24 samples gave following results:

$$y = 0.9403x + 0.079 \text{ mg/l}$$

$$R = 0.9886 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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PROGRAM NA ANALIZATOR / APPLICATION / АДАПТАЦИЯ:

Parameters

No.	69	Prim.Wave.	546
Test	CYS C	Sec.Wave.	700
Method	Endpoint	Sample Vol.	5
Direction	Ascend	R1 Vol.	180
Unit	mg/l	R2 Vol.	30
Decimals	1	Line. Limit	20
Incubation	10	Antigen Check	
Reaction	-2 25	Substrat	0
R1 Blank		Mix. R Blank	
Lower	0	Lower	0
Upper	0	Upper	0
Response		Linearity	
Lower	-2.5	Lower	
Upper	2.5	Upper	
Sample	30	Full Name	Cystatin C
Vol.		Print No.	69
Dilution	5		

Calibration

Rule	Spline
K Factor	0
Replicates	1
Interval	0
Sensitivity	0
Correlation	0
Difference	2.5
Blank Response	0 2.5
Coefficient	0
Difference	
Non-linear SD	0

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