ACCENT-300 ASAT

Nr kat. **7-314** (PL)

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Cat. No **7-314** (EN)

INTENDED USE


Diagnostic kit for determination of aspartate aminotransferase activity intended to use in automatic analyser ACCENT-300. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Aspartate aminotransferase (ASAT, AST, GOT) is an enzyme participated in amino acids metabolism. ASAT is found in all tissues but particularly high level of ASAT is observed in heart muscle, skeletal muscle, liver and kidney. This is why elevated ASAT serum level is marker of myocardial infarction and kidney, liver or skeletal muscle injury.

METHOD PRINCIPLE

Optimized, modified method according to International Federation of Clinical Chemistry (IFCC), without pyridoxal phosphate.

The rate of absorbance changing at $\lambda=340$ nm is directly proportional to aspartate aminotransferase activity.

REAGENTS

Package

1-Reagent 4 x 40 ml
 2-Reagent 2 x 20.5 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 12 weeks on board the analyser at 2-10°C.

Concentrations in the test

1-Reagent

Tris buffer, pH 7.7 ≤ 120 mmol/l
 L-aspartate ≤ 360 mmol/l
 MDH ≤ 1.4 U/ml
 LDH ≤ 2.3 U/ml
 stabilizer, preservative

2-Reagent

2-oxoglutarate ≤ 74 mmol/l
 NADH ≤ 1.7 mmol/l
 buffer
 preservatives

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-Reagent meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Warning

H315 Causes skin irritation.
 H319 Causes serious eye irritation.
 P280 Wear protective gloves/protective clothing/eye protection/face protection.
 P302+P352 IF ON SKIN: Wash with plenty of soap and water.
 P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

SPECIMEN

Serum, heparinized or EDTA plasma free from hemolysis. Do not use heparine ammonium salt. Hemolysis should be avoided, since ASAT activity in erythrocytes is 10 times higher than in normal serum. Do not freeze the samples. ASAT activity remains stable in specimen up to 1 day at 15-25°C or up to 4 days at 2-8°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.
 For reagent blank deionized water is recommended.

REFERENCE VALUES ⁶

| serum / plasma | 37°C | |
|----------------|--------------|-------------------------|
| female | up to 31 U/l | up to 0.518 μ kat/l |
| male | up to 37 U/l | up to 0.618 μ kat/l |

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173). For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended. Deionised water should be used as a calibrator 0. The calibration curve should be prepared every 12 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyzer ACCENT-300. Results may vary if a different instrument or a manual procedure is used.

- LOQ (Limit of Quantitation):** 7 U/l (0.12 μ kat/l).
- Linearity:** up to 1000 U/l (16.7 μ kat/l).
- Specificity / Interferences**
 Haemoglobin up to 0.63 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision

| Repeatability (run to run) n = 20 | Mean [U/l] | SD [U/l] | CV [%] |
|--|---------------|-------------|-----------|
| level 1 | 43.8 | 1.5 | 3.5 |
| level 2 | 214 | 2.5 | 1.2 |
| Reproducibility (day to day) n = 80 | Mean [U/l] | SD [U/l] | CV [%] |
| level 1 | 42.9 | 1.4 | 3.4 |
| level 2 | 211 | 4.8 | 2.3 |

Method comparison

A comparison between ASAT values determined at ACCENT-300 (y) and at BECKMAN COULTER AU680 (x) using 120 serum samples gave following results:
 $y = 0,9623 x - 0,3646 \text{ U/l}$;
 $R = 0,994$ (R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

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PROGRAM NA ANALIZATOR / APPLICATION / АДАПТАЦИЯ:

Parameters

| | | | |
|-----------------|-------------|---------------------|------|
| No. | 18 | Prim.Wave. | 340 |
| Test | ASAT II GEN | Sec.Wave. | 450 |
| Method | Kinetic | Sample Vol. | 10 |
| Direction | Descend | R1 Vol. | 200 |
| Unit | U/l | R2 Vol. | 50 |
| Decimals | 1 | Line. Limit | 20 |
| Incubation | 10 | Antigen Check | |
| Reaction | 3 18 | Substrate | 0.55 |
| R1 Blank | | Mix. R Blank | |
| Lower | 0 | Lower | 0 |
| Upper | 0 | Upper | 0 |
| Response | | Linearity | |
| Lower | -2.5 | Lower | 7 |
| Upper | 2.5 | Upper | 1000 |
| Sample Vol. | 45 | Full Name | ASAT |
| Dilution | 5 | Print No. | 18 |

Calibration

| | |
|----------------|--------------------|
| Rule | Multi Point Linear |
| K Factor | 0 |
| Replicates | 2 |
| Interval | 84 |
| Sensitivity | 0 |
| Correlation | 0 |
| Difference | 2.5 |
| Blank Response | 0 2.5 |
| Coefficient | 0 |
| Difference | |
| Non-linear SD | 0 |

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